PRODUCT INFORMATION & MANUAL

Human IL-6 Valukine™ ELISA

Catalog Number: VAL102C

For the quantitative determination of natural and recombinant human Interleukin 6 (IL-6) concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version 202403.1

TABLE OF CONTENTS

I. BACKGROUND.....	2
II. OVERVIEW.....	4
III. ADVANTAGES.....	5
IV. EXPERIMENT.....	8
V. KIT COMPONENTS AND STORAGE.....	9
VI. PREPARATION.....	11
VII. ASSAY PROCEDURE.....	13
VIII. REFERENCES.....	15

I. BACKGROUND

Interleukin 6 (IL-6) is a pleiotropic α -helical 22-28 kDa phosphorylated and variably glycosylated cytokine that plays important roles in the acute phase reaction, inflammation, hematopoiesis, bone metabolism, and cancer progression (1-5). Mature human IL-6 is 183 amino acids (aa) in length and shares 41% aa sequence identity with mouse and rat IL-6 (6). Alternate splicing generates several isoforms with internal deletions, some of which exhibit antagonistic properties (7-10). Cells known to express IL-6 include CD8⁺ T cells, fibroblasts, synoviocytes, adipocytes, osteoblasts, megakaryocytes, endothelial cells (under the influence of endothelins), sympathetic neurons, cerebral cortex neurons, adrenal medulla chromaffin cells, retinal pigment cells, mast cells, keratinocytes, Langerhans cells, fetal and adult astrocytes, neutrophils, monocytes, eosinophils, colonic epithelial cells, B1 B cells, and pancreatic islet beta cells (2, 7, 10-33). IL-6 production is generally correlated with cell activation and is normally kept in control by glucocorticoids, catecholamines, and secondary sex steroids (2). Normal human circulating IL-6 is in the 1 pg/mL range, with slight elevations during the menstrual cycle, modest elevations in certain cancers, and large elevations after surgery (34-38).

IL-6 induces signaling through a cell surface heterodimeric receptor complex composed of a ligand binding subunit (IL-6 R) and a signal transducing subunit (gp130). IL-6 binds to IL-6 R, triggering IL-6 R association with gp130 and gp130 dimerization (39). Gp130 is also a component of the receptors for CLC, CNTF, CT-1, IL-11, IL-27, LIF, and OSM (40). Soluble forms of IL-6 R are generated by both alternative splicing and proteolytic cleavage (3). In a mechanism known as trans-signaling, complexes of soluble IL-6 and IL-6 R elicit responses from gp130-expressing cells that lack cell surface IL-6 R (1, 3). Trans-signaling enables a wider range of cell types to respond to IL-6, as the expression of gp130 is ubiquitous, while that of IL-6 R is predominantly restricted to hepatocytes, monocytes, and resting lymphocytes (1-3). Soluble splice forms of gp130 block trans-signaling from IL-6/IL-6 R but not from other cytokines that use gp130 as a co-receptor (3, 41).

IL-6, along with TNF- α and IL-1, drives the acute inflammatory response, is almost solely responsible for fever and the acute phase response in the liver, and is important in the transition from acute inflammation to either acquired immunity, or chronic inflammatory disease (1-4). It contributes to chronic inflammation in conditions such as obesity, insulin resistance, inflammatory bowel disease, inflammatory arthritis and sepsis when dysregulated, often involving IL-6 trans-signaling (1, 2). It also plays an

important role in the differentiation of naive T cells to Th17 inflammatory cells in the presence of TGF- β . IL-6 modulates bone resorption and is a major effector of inflammatory joint destruction in rheumatoid arthritis through its promotion of Th17 T cell activity (1). It contributes to atherosclerotic plaque development and destabilization (2). However, IL-6 can also have anti-inflammatory effects, such as in skeletal muscle where it is secreted in response to exercise (2). It promotes hematopoiesis by being a growth factor for hematopoietic stem cells, induces B cell maturation to plasma cells and perpetuates multiple myeloma (1, 42). IL-6 also promotes, but probably does not initiate, other types of inflammation-associated carcinogenesis, such as colitis-associated cancer (1).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for human IL-6 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any human IL-6 present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked antibody specific for human IL-6 is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, TMB substrate (Chromogenic agent) is added to the wells and color develops in proportion to the amount of human IL-6 bound in the initial step. The color development is stopped and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- ◆ **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- ◆ This kit is suitable for cell culture supernate, human serum, human EDTA plasma and human heparin plasma.
- ◆ The kit should not be used beyond the expiration date on the kit label.
- ◆ Do not mix or substitute reagents with those from other lots or sources.
- ◆ If samples generate values higher than the highest standard, dilute the samples with Calibrator Diluent (1×) and repeat the assay.
- ◆ Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

Sample	Intra-assay Precision			Inter-assay Precision		
	1	2	3	1	2	3
Mean (pg/mL)	38.8	94.7	205	31.8	97.4	191
Standard Deviation	0.687	2.02	4.69	1.35	3.65	8.25
CV%	1.8	2.1	2.3	4.2	3.8	4.3

B. RECOVERY

The recovery of human IL-6 spiked to levels throughout the range of the assay in various matrices was evaluated.

Sample Type	Average % Recovery	Range (%)
Cell culture supernate (n=4)	102	82 - 122
Human serum (n=4)	96	89 - 103
Human EDTA plasma (n=4)	97	90 - 102
Human heparin plasma (n=4)	96	90 - 99

C. SENSITIVITY

Twenty-four assays were evaluated and the minimum detectable dose (MDD) of human IL-6 ranged from 0.120-0.626 pg/mL. The mean MDD was 0.368 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against a highly purified recombinant human IL-6 produced at R&D Systems.

The NIBSC/WHO Interleukin-6 International Standard 89/548 (IL-6, Human rDNA derived) was evaluated in this kit. The dose response curve of the International Standard (89/548) parallels the Valukine standard curve. To convert sample values obtained with the Valukine Human IL-6 kit to approximate NIBSC 89/548 units, use the equation below.

NIBSC (89/548) approximate value (IU/mL)=0.1006 × Valukine Human IL-6 value (pg/mL)

Additionally, the 2nd International NIBSC/WHO Standard for Interleukin-6 21/308 (IL-6, Human recombinant) was evaluated in this kit. The dose response curve of the International Standard (21/308) parallels the Valukine standard curve. To convert sample values obtained with the Valukine Human IL-6 kit to approximate NIBSC 21/308 units, use the equation below.

NIBSC (21/308) approximate value (IU/mL)=0.0978 × Valukine Human IL-6 value (pg/mL)

Based on data generated from April 2023.

E. LINEARITY

To assess the linearity of the assay, samples were containing or spiked with high concentrations of human IL-6 in various matrices and diluted with Calibrator Diluent (1×) to produce samples with values within the dynamic range of the assay.

Dilution		Cell culture supernate (n=3)	Human serum (n=4)	Human EDTA Plasma (n=4)	Human heparin plasma (n=4)
1:2	Average % of Expected	97	104	103	105
	Range (%)	81-109	94-112	95-111	98-111
1:4	Average % of Expected	100	106	105	106
	Range (%)	88-116	97-115	95-115	98-115
1:8	Average % of Expected	107	108	107	110
	Range (%)	95-121	99-118	99-115	102-120
1:16	Average % of Expected	116	109	108	109
	Range (%)	99-125	103-119	101-116	99-117

F. SAMPLE VALUES

Human serum/plasma- Samples from apparently healthy volunteers were evaluated for the presence of human IL-6 in this assay. No medical histories were available for the donors used in this study.

Sample Type	Mean of Detectable (pg/mL)	%Detectable	Range (pg/mL)
Human serum (n=30)	6.90	40	ND-13.9
Human EDTA plasma (n=30)	6.82	40	ND-13.4
Human heparin plasma (n=30)	6.89	40	ND-12.9

ND=Non-detectable

Cell culture supernates - Human peripheral blood mononuclear cells (1×10^6 cells/mL) were cultured in RPMI supplemented with 10% fetal bovine serum, 50 μ M β -mercaptoethanol, 2 mM L-glutamine, 100 units/mL penicillin, and 100 μ g/mL streptomycin sulfate. PBMCs were left untreated or treated with 10 μ g/mL PHA for 5 days before collecting conditioned media. Aliquots of the cell culture supernates were removed and assayed for levels of human IL-6.

Condition	Day 5 (pg/mL)
Unstimulated	24.0
Stimulated	10,986

G. SPECIFICITY

The following factors were prepared at 3 ng/mL and assayed for cross-reactivity. No significant cross-reactivity or interference was observed.


Recombinant human:	Other recombinants:
gp130	mouse IL-6
IL-6 R α /gp130	rat IL-6
IL-6 sR	porcine IL-6
IL-6 R	

Recombinant cyno macaque IL-6 cross-reacts at approximately 79% and interferes at levels > 50 pg/mL.

IV. EXPERIMENT

EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

(pg/mL)	O.D.	Average	Corrected
0	0.005 0.005	0.005	—
3.13	0.032 0.032	0.032	0.027
6.25	0.063 0.063	0.063	0.058
12.5	0.124 0.125	0.125	0.120
25	0.255 0.263	0.259	0.254
50	0.504 0.504	0.504	0.499
100	0.987 1.010	0.999	0.994
300	2.569 2.575	2.572	2.567

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
Human IL-6 Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with an antibody against human IL-6	1 plate
Human IL-6 Conjugate	Solution of antibody against human IL-6 conjugated to horseradish peroxidase	1 vial
Human IL-6 Standard	Recombinant human IL-6 in a buffered protein base; lyophilized. Refer to the vial label for reconstitution volume.	1 vial
Calibrator Diluent (5×)/RD5P	A 5× concentrated buffered protein base used to dilute standard and samples	1 vial
Wash Buffer Concentrate (25×)	A 25× concentrated solution of buffered surfactant	1 vial
TMB Substrate	TMB ELISA Substrate Solution/ TMB Substrate Solution	2 vials
Stop Solution	2 N sulfuric acid	1 vial
Plate Sealers	Adhesive strip	3 strips

B. STORAGE

Unopened Kit	Store at 2-8°C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Wash Buffer (1×)	May be stored for up to 1 month at 2-8°C.*
	Stop Solution	
	Conjugate	
	TMB Substrate	
	Standard	Aliquot and store for up to 1 month at -20°C in a manual defrost freezer. * Avoid repeated freeze-thaw cycles.
	Calibrator Diluent (5×)/RD5P	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Calibrator Diluent (1×). Prepare fresh for each assay.
Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*	

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- ◆ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ◆ Pipettes and pipette tips.
- ◆ Deionized or distilled water.
- ◆ Squirt bottle, manifold dispenser, or automated microplate washer.
- ◆ 500 mL graduated cylinder.

D. PRECAUTION

- ◆ Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- ◆ The Stop Solution provided with this kit is an acid solution. Wear protective gloves, clothing, eye, and face protection. Wash hands thoroughly after handling.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1 \times).

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at $1000 \times g$. Remove serum and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1 \times).

Plasma - Collect plasma using EDTA or heparin as an anticoagulant. Centrifuge for 15 minutes at $1000 \times g$ within 30 minutes of collection. Assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1 \times).

Note: Citrate plasma has not been validated for this assay.

B. SAMPLE PREPARATION

Human serum samples recommend a 5-fold dilution. A suggested 5-fold dilution is 40 μL of sample + 160 μL of Calibrator Diluent (1 \times). Optimal dilutions should be determined by the end user.

Human plasma samples recommend a 2-fold dilution. A suggested 2-fold dilution is 100 μL of sample + 100 μL of Calibrator Diluent (1 \times). Optimal dilutions should be determined by the end user.

C. REAGENT PREPARATION

Bring all reagents to room temperature before use.

Note: *High concentrations of IL-6 are found in saliva. The use of a face mask and gloves is recommended to protect kit reagents from contamination.*


Wash Buffer (1 \times) - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25 \times) into deionized or distilled water to prepare 500 mL of Wash Buffer (1 \times).

Calibrator Diluent (1 \times) - Use deionized or distilled water to prepare Calibrator Diluent (1 \times).

Human IL-6 Standard - Refer to the vial label for reconstitution volume. This reconstitution produces a stock solution of 3000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*If you have any question, please seek help from our Technical Support.

Pipette 450 μL of Calibrator Diluent ($1\times$) into the 300 pg/mL tube. Pipette 400 μL of Calibrator Diluent ($1\times$) into the 100 pg/mL tube and 300 μL into remaining tubes. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 300 pg/mL standard serves as the high standard. Calibrator Diluent ($1\times$) serves as the zero standard (0 pg/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB Substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB Substrate.

VII. ASSAY PROCEDURE

Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.

Note: *High concentrations of IL-6 are found in saliva. The use of a face mask and gloves is recommended to protect kit reagents from contamination.*

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100 μL of standard and sample per well. Cover with the adhesive strip provided. **Incubate for 2 hours at room temperature.** A plate layout is provided for a record of standards and samples assayed.
4. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Buffer (400 μL) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
5. Add 200 μL of human IL-6 Conjugate to each well. Cover with a new adhesive strip. **Incubate for 2 hours at room temperature.**
6. Repeat the aspiration/wash as in step 4.
7. Add 200 μL of TMB Substrate to each well. **Incubate for 20 minutes at room temperature. Protect from light.**
8. Add 50 μL of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
9. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

10. CALCULATION OF RESULTS

Average the duplicate readings for each standard and sample and subtract the average zero standard optical density (O.D.). Create a standard curve by reducing the data using computer software capable of generating a four-parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis

and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the human IL-6 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Naugler, W.E. and M. Karin (2008) *Trends Mol. Med.* 14:109.
2. Schuett, H. et al. (2009) *Thromb. Haemost.* 102:215.
3. Jones, S.A. (2005) *J. Immunol.* 175:3468.
4. Hodge, D.R. et al. (2005) *Eur. J. Cancer* 41:2502.
5. Rose -John, S. et al. (2006) *J. Leukoc. Biol.* 80:227.
6. Van Snick, J. et al. (1988) *Eur. J. Immunol.* 18:193.
7. Kestler, D.P. et al. (1995) *Blood* 86:4559.
8. Kestler, D.P. et al. (1999) *Am. J. Hematol.* 61:169.
9. Bihl, M.P. et al. (2002) *Am. J. Respir. Cell Mol. Biol.* 27:48.
10. Alberti, L. et al. (2005) *Cancer Res.* 65:2.
11. May, L.T. et al. (1986) *Proc. Natl. Acad. Sci. USA* 83:8957.
12. Sad, S. et al. (1995) *Immunity* 2:271.
13. Cichy, J. et al. (1996) *Biochem. Biophys. Res. Commun.* 227:318.
14. Miyazawa, K. et al. (1998) *Am. J. Pathol.* 152:793.
15. Fried, S.K. et al. (1998) *Endocrinology* 83:847.
16. Ishimi, Y. et al. (1990) *J. Immunol.* 145:3297.
17. Jiang, S. et al. (1994) *Blood* 84:4151.
18. Xin, X. et al. (1995) *Endocrinology* 136:132.
19. Marz, P. et al. (1998) *Proc. Natl. Acad. Sci. USA* 95:3251.
20. Ringheim, G.E. et al. (1995) *J. Neuroimmunol.* 63:113.
21. Gadiant, R.A. et al. (1995) *Neurosci. Lett.* 194:17.
22. Kuppner, M.C. et al. (1995) *Immunology* 84:265.
23. Gagari, E. et al. (1997) *Blood* 89:2654.
24. Cumberbatch, M. et al. (1996) *Immunology* 87:513.
25. Fujisawa, H. et al. (1997) *J. Interferon Cytokine Res.* 17:347.
26. Lee, S.C. et al. (1993) *J. Immunol.* 150:2659.
27. Lafortune, L. et al. (1996) *J. Neuropathol. Exp. Neurol.* 55:515.
28. Ericson, S.G. et al. (1998) *Blood* 91:2099.
29. Melani, C. et al. (1993) *Blood* 81:2744.
30. Lacy, P. et al. (1998) *Blood* 91:2508.
31. Jung, H.C. et al. (1995) *J. Clin. Invest.* 95:55.
32. Spencer, N.F.L. and R.A. Daynes (1997) *Int. Immunol.* 9:745.

33. Campbell, I.L. et al. (1989) *J. Immunol.* 143:1188.
34. D'Auria, L. et al. (1997) *Eur. Cytokine Netw.* 8:383.
35. Yamamura, M. et al. (1998) *Br. J. Haematol.* 100:129.
36. Angstwurm, M.W.A. et al. (1997) *Cytokine* 9:370.
37. Mouawad, R. et al. (1996) *Clin. Cancer Res.* 2:1405.
38. Sakamoto, K. et al. (1994) *Cytokine* 6:181.
39. Murakami, M. et al. (1993) *Science* 260:1808.
40. Muller-Newen, G. (2003) *Sci. STKE* 2003:PE40.
41. Mitsuyama, K. et al. (2006) *Clin. Exp. Immunol.* 143:125.
42. Cerutti, A. et al. (1998) *J. Immunol.* 160:2145.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

1																	
2																	
3																	
4																	
5																	
6																	
7																	
8																	
9																	
10																	
11																	
12																	
	A	B	C	D	E	F	G	H									

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