## **PRODUCT INFORMATION & MANUAL**

**Human HGF Valukine™ ELISA**

**Catalog Number: VAL168**

For the quantitative determination of natural and recombinant  
human HGF concentrations

For research use only.  
Not for diagnostic or therapeutic procedures.

**Bio-Techne China Co. Ltd**

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001  
**info.cn@bio-techne.com**

Please refer to the kit label for expiry date.  
Novus kits are guaranteed for 3 months from date of receipt

Version 202403.2

# TABLE OF CONTENTS

I. BACKGROUND .....	2
II. OVERVIEW .....	4
III. ADVANTAGES .....	5
IV. EXPERIMENT .....	7
V. KIT COMPONENTS AND STORAGE .....	8
VI. PREPARATION .....	11
VII. ASSAY PROCEDURE .....	13
VIII. REFERENCES .....	15

## I. BACKGROUND

Hepatocyte Growth Factor (HGF), also known as scatter factor, hepatopoietin A, and mammary growth factor, is a pleiotropic glycoprotein that regulates the growth and migration of diverse cell types. It is structurally similar to the S1 peptidase Plasminogen. HGF contains an N-terminal PAN/APPLE-like domain, four Kringle domains, and a catalytically inactive serine proteinase-like domain (1, 2). Alternative splicing generates human HGF isoforms that lack the proteinase-like domain and different numbers of the Kringle domains. HGF is secreted as an inactive single chain propeptide that can circulate as a soluble molecule or associate with the extracellular matrix (3, 4). At sites of tissue damage, the propeptide is cleaved after the fourth Kringle domain by serine proteases including HGF Activator and uPA (4-8). The resulting bioactive HGF consists of a disulfidelinked heterodimer of a 60 kDa N-terminal alpha chain and a 30 kDa C-terminal beta chain (4, 5, 9). The serum levels of HGF are elevated in a wide range of pathologies including liver damage (10, 11), acute kidney failure (12), myocardial infarction (13), type 1 diabetes (14), obesity (15), and cancer (16-23), as well as in the synovial fluid of rheumatoid arthritis patients (24). Human HGF shares 91-94% amino acid sequence identity with bovine, canine, feline, mouse, and rat HGF. HGF demonstrates marked species cross-reactivity (25).

HGF exerts its biological activity through the widely expressed receptor tyrosine kinase, HGF R/c-MET (26, 27). This receptor undergoes N-linked glycosylation followed by proteolytic cleavage into 50 kDa N-terminal alpha and 145 kDa C-terminal beta chains (28). The strictly extracellular alpha chain remains disulfide-linked to the beta chain which contains the remaining extracellular, transmembrane, and cytoplasmic domains (26, 27). HGF also binds heparan sulfate proteoglycans, and these interactions enhance the ability of HGF to bind and activate HGF R (29, 30). In the absence of ligand, HGF R forms noncovalent complexes with a variety of membrane proteins including CD44v6, CD151, EGF R, Fas, Integrin  $\alpha 6/\beta 4$ , Plexins B1, B2, B3, and MSP R/Ron (31-38). Ligation of one complex component can trigger activation of the other, followed by cooperative signaling effects (31-38). Formation of some of these heteromeric complexes is a requirement for epithelial cell morphogenesis and tumor cell invasion (34-36). Overexpression and the production of alternate forms of HGF R are implicated in the development of many human cancers (39).

HGF is expressed by fibroblasts, adipocytes, smooth muscle cells, and endothelial cells (1). Expression of HGF R, on the other hand, is found mainly on epithelial cells, suggesting that HGF acts in a paracrine fashion to mediate interactions between stromal and epithelial cells (40). HGF induces the proliferation and migration of epithelial cells as well as multiple other cell types including hepatocytes, chondrocytes, keratinocytes, melanocytes, and endothelial cells (1). It is mitogenic toward most tumor cells but can conversely inhibit their proliferation in some cases (39, 41, 42). During organogenesis, tissue repair, and angiogenesis, HGF promotes epithelial/endothelial morphogenesis by inducing cell scattering and branching tubulogenesis (1, 25, 43, 44). The ability of HGF to regulate angiogenesis and the motility of epithelial cells underlie its importance in the development of solid tumors (39). In addition to its morphogenetic effects, HGF induces a range of responses in diverse tissues (1). It supports the survival, proliferation, and insulin productivity of pancreatic islet cells (45). It functions as a neurotrophic factor during development and in the response to injury (46, 47). It also suppresses inflammation by inducing dendritic cell tolerization, Treg induction, and Th2 bias while inhibiting T cell activation, IL-17 expression, and inflammatory cell infiltration (48-50).

## II. OVERVIEW

### A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for human HGF has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any human HGF present is bound by the immobilized antibody. After washing away any unbound substances, a biotinylated detection antibody specific for human HGF is pipetted into the wells. After washing away any unbound substances, Streptavidin-HRP is pipetted into the wells. Following a wash to remove any unbound reagent, TMB substrate (Chromogenic agent) is added to the wells and color develops in proportion to the amount of human HGF bound in the initial step. The color development is stopped, and the intensity of the color is measured.

### B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- This kit is suitable for cell culture supernate and human serum.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Reagent Diluent (1×) or Calibrator Diluent-S (1×) and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

### III. ADVANTAGES

#### A. PRECISION

**Intra-assay Precision** (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

**Inter-assay Precision** (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

	Intra-assay Precision			Inter-assay Precision		
Sample	1	2	3	1	2	3
Mean (pg/mL)	240.8	981.8	3810.9	250.5	1012.2	3998.7
Standard Deviation	11.1	28.1	119.4	17.3	51.7	321.8
CV%	4.6	2.9	3.1	6.9	5.1	8.0

#### B. RECOVERY

The recovery of human HGF spiked to different levels throughout the range of the assay in cell culture media was evaluated. The recovery ranged from 95.1 to 123.1% with an average of 112.7%.

The recovery of human HGF spiked to different levels throughout the range of the assay in human serum was evaluated. The recovery ranged from 114.0 to 125.9% with an average of 121.0%.

#### C. SENSITIVITY

The minimum detectable dose (MDD) of human HGF is typically less than 9.18 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

#### D. CALIBRATION

This immunoassay is calibrated against a highly purified Sf 21-expressed recombinant human HGF produced at R&D Systems.

## E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with high concentrations of human HGF and diluted with Reagent Diluent 1× (*for cell culture supernate samples*) or Calibrator Diluent-S 1× (*for serum samples*) to produce samples with values within the dynamic range of the assay.

<b>Dilution</b>	<b>Average % of Expected</b>	<b>Range (%)</b>
1:2	90.2	86.2-96.7
1:4	89.6	82.8-101.9
1:8	91.0	80.3-108.1
1:16	89.6	80.2-103.4

## F. SAMPLE VALUES

**Human serum** - Eight human serum samples were evaluated for the presence of human HGF in this assay. All samples measured ranged from 1308.2 to 1846.6 pg/mL with an average of 1632.0 pg/mL.

## G. SPECIFICITY

The following factors prepared at 50 ng/mL were assayed and exhibited no cross-reactivity but do interfere in this assay.


<b>Recombinant human</b>	<b>Recombinant Mouse</b>
HGF R	HGF R
MSP	

## IV. EXPERIMENT

### EXAMPLE STANDARD


The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

CELL CULTURE SUPERNATE ASSAY

pg/mL	OD	Average	Corrected
0	0.030 0.034	0.032	-
125	0.100 0.100	0.100	0.068
250	0.176 0.176	0.176	0.144
500	0.321 0.327	0.324	0.292
1000	0.617 0.621	0.619	0.587
2000	1.161 1.190	1.176	1.144
4000	2.016 2.041	2.029	1.997
8000	2.735 2.739	2.737	2.705

SERUM ASSAY

pg/mL	OD	Average	Corrected
0	0.031 0.032	0.032	-
125	0.087 0.089	0.088	0.057
250	0.142 0.147	0.145	0.113
500	0.259 0.264	0.262	0.230
1000	0.480 0.485	0.483	0.451
2000	0.885 0.900	0.893	0.861
4000	1.488 1.591	1.540	1.508
8000	2.436 2.494	2.465	2.434

## V. KIT COMPONENTS AND STORAGE

### A. MATERIALS PROVIDED

Parts	Description	Size
Human HGF Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with an antibody against human HGF.	1 plate
Human HGF Standard	Recombinant human HGF in a buffered protein base, lyophilized. Refer to the vial label for reconstitution volume.	2 vials
Human HGF Detection Antibody	Biotinylated human HGF antibody, lyophilized. Refer to the vial label for reconstitution volume.	1 vial
Streptavidin-HRP A (200×)	200× Streptavidin conjugated to horseradish peroxidase.	1 vial
Reagent Diluent Concentrate (10×)	A 10× concentrated buffered base used to dilute standard, samples, Detection Antibody and HRP.	1 vial
Calibrator Diluent-S Concentrate (2×)	Concentrated buffer diluent used to dilute standard and samples.	1 vial
Detection Antibody Diluent Concentrate (4×)	A 4× concentrated buffered protein base used to dilute Detection Antibody.	1 vial
Wash Buffer Concentrate (25×)	A 25× concentrated solution of buffered surfactant with preservatives.	1 vial
TMB Substrate	TMB ELISA Substrate Solution/TMB Substrate Solution.	1 vial
Stop Solution	2 N sulfuric acid.	1 vial
Plate Sealers	Adhesive strip.	3 strips

## B. STORAGE

Unopened Kit	Store at 2-8 °C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Streptavidin-HRP A	May be stored for up to 1 month at 2-8 °C.*
	Wash Buffer (1×)	
	TMB Substrate	
	Stop Solution	
	Standard	Prepare fresh for each assay. Standards may be stored for up to 1 month at -20 °C.*
	Detection Antibody	Aliquot and store for up to 1 month at -20 °C in a manual defrost freezer. *
	Reagent Diluent Concentrate (10×)	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Reagent Diluent (1×). Prepare fresh for each assay.
	Calibrator Diluent -S Concentrate (2×)	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Calibrator Diluent-S (1×). Prepare fresh for each assay.
	Detection Antibody Diluent Concentrate (4×)	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Detection Antibody Diluent (1×). Prepare fresh for each assay.
Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8 °C.*	

\* Provided this is within the expiration date of the kit.

### **C. OTHER SUPPLIES REQUIRED**

- ◆ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ◆ Pipettes and pipette tips.
- ◆ Deionized or distilled water.
- ◆ Squirt bottle, manifold dispenser, or automated microplate washer.
- ◆ 500 mL graduated cylinder.
- ◆ Horizontal orbital microplate shaker capable of maintaining a speed of 500±50 rpm.

### **D. PRECAUTION**

- ◆ Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- ◆ The Stop Solution provided with this kit is an acid solution. Wear protective gloves, clothing, eye, and face protection. Wash hands thoroughly after handling.
- ◆ Human HGF is detectable in saliva. Take precautionary measures to prevent contamination of kit reagents while running this assay.

## VI. PREPARATION

### A. SAMPLE COLLECTION AND STORAGE

**Cell Culture Supernates** - Remove particulates by centrifugation and assay immediately or aliquot and store samples at  $\leq -20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles. Samples may require dilution with Reagent Diluent (1 $\times$ ).

**Serum** - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at  $\leq -20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent-S (1 $\times$ ).

### B. SAMPLE PREPARATION

Human serum samples recommend a 2-fold dilution. A suggested 2-fold dilution is 100  $\mu\text{L}$  of sample + 100  $\mu\text{L}$  of **Calibrator Diluent-S (1 $\times$ )**. Optimal dilutions should be determined by the end user.

### C. REAGENT PREPARATION

**Bring all reagents to room temperature before use.**

*Note: human HGF are found in saliva. It is recommended that a face mask and gloves be used to protect kit reagents from contamination.*

**Wash Buffer (1 $\times$ )**- If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25 $\times$ ) into deionized or distilled water to prepare 500 mL of Wash Buffer (1 $\times$ ).

**Calibrator Diluent-S (1 $\times$ )** - Use deionized or distilled water to prepare Calibrator Diluent-S (1 $\times$ ).

**Reagent Diluent (1 $\times$ )** - Use deionized or distilled water to prepare Reagent Diluent (1 $\times$ ).

**Detection Antibody Diluent (1 $\times$ )** - Use deionized or distilled water to prepare Detection Antibody Diluent (1 $\times$ ).

**Detection Antibody (1 $\times$ )**- **Centrifuge briefly before opening. Reconstitution Volume refer to vial label to prepare Detection Antibody (100 $\times$ )**. Allow the Detection Antibody to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions. Aliquot and store if needed. Dilute to Detection Antibody (1 $\times$ ) with Detection Antibody Diluent (1 $\times$ ). Prepare at least 15 minutes prior to use.


**Streptavidin-HRP A (1 $\times$ )** - **Centrifuge briefly before opening.** Dilute to the working

concentration specified on the vial label using Reagent Diluent (1×).

**Human HGF Standard- Centrifuge briefly before opening. Refer to the vial label for the reconstitution volume\***. Reconstitute the human HGF Standard with **Reagent Diluent (1×)** (for cell culture supernate samples) or **Calibrator Diluent-S (1×)** (for serum samples). This reconstitution produces a stock solution of 16000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

\*If you have any question, please seek help from our Technical Support.

Pipette 300 µL of **Reagent Diluent (1×)** or **Calibrator Diluent-S (1×)** into each tubes. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 8000 pg/mL standard serves as the high standard. The **Reagent Diluent (1×)** or **Calibrator Diluent-S (1×)** serves as the zero standard (0 pg/mL).

#### D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB Substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB Substrate.

## VII. ASSAY PROCEDURE

**Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.**

*Note: human HGF are found in saliva. It is recommended that a face mask and gloves be used to protect kit reagents from contamination.*

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100  $\mu$ L of standard and prepared sample per well. Cover with the adhesive strip provided. **Incubate for 2 hours at room temperature on a horizontal orbital microplate shaker set at  $500 \pm 50$  rpm.** A plate layout is provided for a record of standards and samples assayed. (Samples may require dilution. See Sample Preparation section.)
4. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Buffer (400  $\mu$ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
5. Add 100  $\mu$ L of the Detection Antibody (1 $\times$ ) diluted in Detection Antibody Diluent (1 $\times$ ) to each well. **Cover with a new adhesive strip and incubate for 2 hours at room temperature on a horizontal orbital microplate shaker set at  $500 \pm 50$  rpm.**
6. Repeat the aspiration/wash as in step 4.
7. Add 100  $\mu$ L of the working dilution of Streptavidin-HRP A to each well. **Cover the plate and incubate for 20 minutes at room temperature on a horizontal orbital microplate shaker set at  $500 \pm 50$  rpm. Protect from light.**
8. Repeat the aspiration/wash as in step 4.
9. Add 100  $\mu$ L of TMB Substrate to each well. **Incubate for 10 minutes at room temperature on a horizontal orbital microplate shaker set at  $500 \pm 50$  rpm. Protect from light.**

10. Add 50  $\mu\text{L}$  of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
11. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

## 12. **CALCULATION OF RESULTS**

Average the duplicate readings for each standard, control, and sample and subtract the average zero standard optical density (O.D.). Create a standard curve by reducing the data using computer software capable of generating a four-parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the human HGF concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

## VIII. REFERENCES

1. Nakamura, T. and S. Mizuno (2010) Proc. Jpn. Acad., Ser. B. 86:588.
2. Grzelakowska-Sztabert, B. and M. Dudkowska (2011) Growth Factors 4:105.
3. Sasaki, M. et al. (1994) Biochem. Biophys. Res. Commun. 199:772.
4. Naldini, L. et al. (1992) EMBO J. 11:4825.
5. Naka, D. et al. (1992) J. Biol. Chem. 267:20114.
6. Miyazawa, K. et al. (1994) J. Biol. Chem. 269:8966.
7. Miyazawa, K. et al. (1996) J. Biol. Chem. 271:3615.
8. Sisson, T.H. et al. (2009) Blood 114:5052.
9. Nakamura, T. et al. (1989) Nature 342:440.
10. Shiota, G. et al. (1995) Hepatology 21:106.
11. Lindroos, P.M. et al. (1991) Hepatology 13:743.
12. Taman, M. et al. (1997) Clin. Nephrol. 48:241.
13. Seko, Y. et al. (2004) Clin. Sci. (London) 106:439.
14. Kulseng, B. et al. (1998) Acta Diabetol. 35:77.
15. Bell, L.N. et al. (2006) Am. J. Physiol. Endocrinol. Metab. 291:E843.
16. Nakamura, S. et al. (1994) Br. J. Haematol. 87:640.
17. Aune, G. et al. (2011) Gynecol. Oncol. 121:402.
18. Skoldenberg, E.G. et al. (2009) Anticancer Res. 29:3311.
19. Borgiel-Marek, H. et al. (2008) Endokrynol. Pol. 59:467.
20. Osada, S. et al. (2008) Hepatogastroenterology 55:544.
21. Astl, J. et al. (2004) Neuro. Endocrinol. Lett. 25:356.
22. Vessely, D. et al. (2004) Physiol. Res. 53:83.
23. Woodbury, R.L. et al. (2002) J. Proteome Res. 1:233.
24. Feuerherm, A.J. et al. (2001) Scand. J. Rheumatol. 30:229.
25. Grant, D.S. et al. (1993) Proc. Natl. Acad. Sci. USA 90:1937.
26. Bottaro, D.P. et al. (1991) Science 251:802.
27. Giordano, S. et al. (1989) Nature 339:155.

28. Giordano, S. et al. (1989) *Oncogene* 4:1383.
29. Sakata, H. et al. (1997) *J. Biol. Chem.* 272:9457.
30. Zioncheck, T.F. et al. (1995) *J. Biol. Chem.* 270:16871.
31. Klosek, S.K. et al. (2005) *Biochem. Biophys. Res. Commun.* 336:408.
32. Jo, M. et al. (2000) *J. Biol. Chem.* 275:8806.
33. Wang, X. et al. (2002) *Mol. Cell* 9:411.
34. Orian-Rousseau, V. et al. (2002) *Genes Dev.* 16:3074.
35. Trusolino, L. et al. (2001) *Cell* 107:643.
36. Giordano, S. et al. (2002) *Nat. Cell Biol.* 4:720.
37. Conrotto, P. et al. (2004) *Oncogene* 23:5131.
38. Follenzi, A. et al. (2000) *Oncogene* 19:3041.
39. Lesko, E. and M. Majka (2008) *Front. Biosci.* 13:1271.
40. Sonnenberg, E. et al. (1993) *J. Cell Biol.* 123:223.
41. Tajima, H. et al. (1991) *FEBS Lett.* 291:229.
42. Rygaard, K. et al. (1993) *Int. J. Oncol.* 2:991.
43. Montesano, R. et al. (1991) *Cell* 67:901.
44. Maeshima, A. et al. (2000) *Kidney Int.* 58:1511.
45. Dai, C. et al. (2003) *J. Biol. Chem.* 278:27080.
46. Honda, S. et al. (1995) *Brain Res. Mol. Brain Res.* 32:197.
47. Ebens, A. et al. (1996) *Neuron* 17:1157.
48. Benkhoucha, M. et al. (2010) *Proc. Natl. Acad. Sci.* 107:6424.
49. Rutella, S. et al. (2006) *Blood* 108:218.
50. Ito, W. et al. (2005) *Am. J. Respir. Cell. Mol. Biol.* 32:268.

# PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

<b>1</b>								
<b>2</b>								
<b>3</b>								
<b>4</b>								
<b>5</b>								
<b>6</b>								
<b>7</b>								
<b>8</b>								
<b>9</b>								
<b>10</b>								
<b>11</b>								
<b>12</b>								
	<b>A</b>	<b>B</b>	<b>C</b>	<b>D</b>	<b>E</b>	<b>F</b>	<b>G</b>	<b>H</b>

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