## **PRODUCT INFORMATION & MANUAL**

**Human IL-17 Valukine™ ELISA**

**VAL129**

For the quantitative determination of natural and recombinant human  
IL-17 concentrations

For research use only.  
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd  
P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001  
info.cn@bio-techne.com

Please refer to the kit label for expiry date.  
Novus kits are guaranteed for 3 months from date of receipt

Version 201910.1

## TABLE OF CONTENTS

I. BACKGROUND.....	2
II. OVERVIEW.....	3
III. ADVANTAGES.....	4
IV. EXPERIMENT.....	6
V. KIT COMPONENTS AND STORAGE.....	7
VI. PREPARATION.....	10
VII. ASSAY PROCEDURE.....	12
VIII. REFERENCES.....	14

## I. BACKGROUND

Human Interleukin 17 (IL-17), also known as IL-17A and CTLA-8, is a 15-20 kDa, variably glycosylated polypeptide that belongs to the IL-17 family of cytokines (1-5). Human IL-17/17A is synthesized as a 155 amino acid (aa) precursor that contains a 23 aa signal sequence and a 133 aa mature region that possesses a cysteine-knot fold (4-6). In both human and mouse, there is one conserved N-linked glycosylation site that likely contributes 5 kDa to its native molecular weight. IL-17A forms both a 32-38 kDa disulfide-linked homodimer, and a 40-45 kDa covalent heterodimer with IL-17F (7-9). Most secreted IL-17A is in the form of the IL-17A: F heterodimer, however, the IL-17A: A homodimer is the most bioactive of the two forms (8). Mature human IL-17A is 61%, 74%, and 99% aa identical to mouse, porcine, and chimpanzee IL-17A, respectively (10-12). Mammalian cells known to produce IL-17 are the CD4<sup>+</sup> Th17 T cells, Paneth cells, GR1<sup>+</sup>CD11b<sup>+</sup> myeloid suppressor cells, CD27- $\gamma\delta$  T cells, CD1<sup>+</sup>NK1.1<sup>+</sup>iNKT cells, and CD3<sup>-</sup>CD4<sup>+</sup> LTI-like cells (9, 13-17).

A high affinity receptor for human IL-17 has been reported, and appears to be a heteromultimer of IL-17RA and IL-17RC, likely in a 2:1 ratio (1, 18). IL-17RA is a 130 kDa, type I transmembrane glycoprotein that bears no resemblance to members of the cytokine, TNF or immunoglobulin receptor superfamily (2, 10, 15). IL-17RC is also a type I transmembrane protein, approximately 90-95 kDa in size, that shares less than 30% aa identity with IL-17RA (19, 20). Both receptors are needed for IL-17A and IL-17A: F activity. The two receptors appear to form a functional association following ligand binding to IL-17RA (1, 21, 22).

IL-17 is best known for its participation in the recruitment and survival of neutrophils (14, 15, 23, 24, 25). Its induction was initially described to be the result of antigen stimulation of dendritic cells, resulting in IL-23 secretion. In a T cell receptor-independent event, IL-23 induces T cell production of IL-17 (14). Once secreted, IL-17 in the bone marrow would seem to induce stromal/fibroblast expression of both G-CSF and stem cell factor (membrane form), an effect that increases polymorphonuclear neutrophils (PMN) differentiation and production. IL-17 may complement this by directly blocking neutrophil apoptosis, promoting greater circulating PMN numbers (23). In the tissues, IL-17 would also seem to promote neutrophil extravasation, principally through its effects on macrophages and endothelial cells (EC). On macrophages, IL-17 induces TNF- $\alpha$ , IL-1 $\beta$  and IL-6 production (26). TNF- $\alpha$  and IL-1 $\beta$  then act on local ECs to induce G-CSF secretion, an effect that is potentiated by IL-17 (27). IL-17 further contributes to PMN influx by inducing EC CXC chemokine release and NO production, which may increase vascular permeability (14, 28). IL-17 effects are not limited to inflammation. In synovial joints, IL-17 upregulates RANKL expression on osteoblasts. This provides a stimulus for osteoclast formation and subsequent bone resorption (24). In conjunction with IL-4 and CD40L, IL-17A also promotes the generation of IgE secreting cells (29). And in white fat, IL-17A inhibits adipocyte differentiation from preadipocytes, and impairs glucose uptake by mature adipocytes (30).

## II. OVERVIEW

### A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for human IL-17 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any IL-17 present is bound by the immobilized antibody. After washing away any unbound substances, a biotin-linked detect antibody specific for human IL-17 are pipetted into the wells. After washing away any unbound substances, streptavidin-HRP is added. Following a wash to remove any unbound antibody-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of IL-17 bound in the initial step. The color development is stopped, and the intensity of the color is measured.

### B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- This kit is suitable for cell culture supernate and serum.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Calibrator Diluent and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

### III. ADVANTAGES

#### A. PRECISION

**Intra-assay Precision** (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

**Inter-assay Precision** (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

CELL CULTURE SUPERNATE / SERUM

Sample	Intra-assay Precision			Inter-assay Precision		
	1	2	3	1	2	3
Mean (pg/mL)	24.0	102.9	378.4	25.4	103.7	383.1
Standard Deviation	1.1	3.0	6.6	1.1	3.3	20.2
CV%	4.5	2.9	1.7	4.3	3.2	5.3

#### B. RECOVERY

The recovery of human IL-17 spiked to three levels throughout the range of the assay in cell culture media was evaluated. The recovery ranged from 93.4-100.4% with an average of 98.0%.

The recovery of human IL-17 spiked to three levels throughout the range of the assay in human serum was evaluated. The recovery ranged from 85.4-91.5% with an average of 89.5%.

#### C. SENSITIVITY

The minimum detectable dose (MDD) of IL-17 is typically less than 2.93 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

#### D. CALIBRATION

This immunoassay is calibrated against a highly purified *E.coli*-expressed recombinant human IL-17 produced at R&D Systems®.

## E. LINEARITY

To assess the linearity of the assay, different samples containing or spiked with high concentrations of IL-17 and diluted with Calibrator Diluent 1× to produce samples with values within the dynamic range of the assay.

<b>Dilution</b>	<b>Average % of Expected</b>	<b>Range (%)</b>
1:2	106.5	103.6-109.1
1:4	111.6	105.8-117.5
1:8	108.8	107.6-110.7
1:16	91.9	83.9-107.8

## F. SAMPLE VALUES

Five serum samples were evaluated for the presence of IL-17 in this assay. All samples measured less than the lowest human IL-17 standard, 15.6 pg/mL.

## G. SPECIFICITY

This assay recognizes both natural and recombinant human IL-17. The following factors prepared at 50 ng/mL were assayed and exhibited no cross-reactivity or interference.

### **Recombinant Human**

IFN- $\gamma$   
IL-10  
IL-12  
IL-16  
IL-17B  
IL-17C  
IL-17D  
IL-17F

### **Recombinant Mouse**

IL-17

Recombinant human IL-17E does not cross-react in this assay but does interfere at concentrations greater than 6.25 ng/mL.


A sample containing 500 pg/mL of recombinant human IL-17A/F heterodimer reads as 108 pg/mL (21.7% cross-reactivity).

## IV. EXPERIMENT

### EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

CELL CULTURE SUPERNATE / SERUM

## V. KIT COMPONENTS AND STORAGE

### A. MATERIALS PROVIDED

Store the unopened kit at 2-8°C. Do not use past kit expiration date.

<b>Parts</b>	<b>Description</b>	<b>Size</b>
Human IL-17 Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with a mouse antibody against human IL-17.	1 plate
Human IL-17 Standard	Recombinant human IL-17 in a buffered protein base; lyophilized. Refer to the vial label for reconstitution volume.	2 vials
Human IL-17 Detection Antibody	Biotinylated IL-17 polyclonal antibody, lyophilized. Refer to the vial label for reconstitution volume.	1 vial
Streptavidin-HRP B (40×)	40× Streptavidin conjugated to horseradish peroxidase.	1 vial
Reagent Diluent (10×)	A 10× concentrated buffered protein base used to dilute Detection Antibody and HRP.	1 vial
Calibrator Diluent (2×)	A 2× concentrated buffered protein base used to dilute standard and samples.	1 vial
Wash Buffer Concentrate (25×)	A 25× concentrated solution of buffered surfactant with preservatives.	1 vial
Color Reagent A	Stabilized hydrogen peroxide.	1 vial
Color Reagent B	Stabilized chromogen (tetramethylbenzidine).	1 vial
Stop Solution	2 N sulfuric acid.	1 vial
Plate Sealers	Adhesive strips.	3 strips

## B. STORAGE

<b>Unopened Kit</b>	Store at 2-8°C. Do not use past kit expiration date.	
<b>Opened/ Reconstituted Reagents</b>	Streptavidin-HRP B	May be stored for up to 1 month at 2-8°C.*
	Diluted Wash Solution	
	Unmixed Color Reagent A	
	Unmixed Color Reagent B	
	Stop Solution	
	Standard	Prepare fresh for each assay.
	Detection Antibody	Aliquot and store for up to 1 month at -20°C in a manual defrost freezer. *
	Reagent Diluent (10×)	May be stored for up to 1 month at 2-8°C.* Use and discard diluted Reagent Diluent (1×). Prepare fresh for each assay.
Calibrator Diluent (2×)	May be stored for up to 1 month at 2-8°C.* Use and discard diluted Reagent Diluent (1×). Prepare fresh for each assay.	
Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*	

\* Provided this is within the expiration date of the kit.

## C. OTHER SUPPLIES REQUIRED

- Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- Pipettes and pipette tips.
- Deionized or distilled water.
- Squirt bottle, manifold dispenser, or automated microplate washer.
- Test tubes for dilution of standards.
- 100mL and 500 mL graduated cylinder.

#### **D. PRECAUTION**

- Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- Color Reagent B may cause skin, eye, and respiratory irritation. Avoid breathing fumes.
- The Stop Solution provided with this kit is an acid solution. Wear protective gloves, clothing, eye, and face protection. Wash hands thoroughly after handling.

## VI. PREPARATION

### A. SAMPLE COLLECTION AND STORAGE

**Cell Culture Supernates** - Remove particulates by centrifugation and assay immediately or aliquot and store samples at  $\leq -20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles.

**Serum** - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at  $\leq -20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles.

### B. SAMPLE PREPARATION

Cell culture supernate samples require a 2-fold dilution prior to the assay. A suggested 2-fold dilution is 100  $\mu\text{L}$  of sample + 100  $\mu\text{L}$  of Calibrator Diluent (1 $\times$ ).

Serum samples require a 4-fold dilution prior to the assay. A suggested 4-fold dilution is 50  $\mu\text{L}$  of sample + 150  $\mu\text{L}$  of Calibrator Diluent (1 $\times$ ).

### C. REAGENT PREPARATION

**Note:** Bring all reagents to room temperature before use.

**Wash Buffer** - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25 $\times$ ) into deionized or distilled water to prepare 500 mL of Wash Buffer.

**Reagent Diluent (1 $\times$ )** - Add 3 mL of Reagent Diluent (10 $\times$ ) into 27 mL of deionized or distilled water to prepare 30 mL of Reagent Diluent (1 $\times$ ).

**Detection Antibody- Reconstitution Volume refer to vial label with Reagent Diluent (1 $\times$ ).** Aliquot and store if needed. Dilute stock solution in Reagent Diluent (1 $\times$ ) to the working concentration of 20 ng/mL. Prepare at least 15 minutes prior to use.

**Streptavidin-HRP B (1 $\times$ )** - Dilute to the working concentration specified on the vial label using Reagent Diluent (1 $\times$ ).


**Substrate Solution** - Color Reagent A and B should be mixed together in equal volumes within 15 minutes of use. Protect from light. 100  $\mu\text{L}$  of the resultant mixture is required per well.

**IL-17 Standard – Refer to the vial label for the reconstitution volume\* using Calibrator Diluent 1 $\times$ .** This reconstitution produces a stock solution of 2000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

\*if you have any question, please seek help from our Technical Support.

**Pipette 300  $\mu\text{L}$  of Calibrator Diluent 1 $\times$  into each tube.** Use the stock solution to

produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 1000 pg/mL standard serves as the high standard. The Calibrator Diluent 1× serves as the zero standard (0 pg/mL).

#### D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Substrate Solution should remain colorless until added to the plate. Keep Substrate Solution protected from light. Substrate Solution should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the Substrate Solution. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the Substrate Solution.

## VII. ASSAY PROCEDURE

**Note:** Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100  $\mu$ L of Standard, or prepared sample per well. Cover with the adhesive strip provided. Incubate for 2 hours at room temperature. A plate layout is provided for a record of standards and samples assayed. (Samples may require dilution. See Sample Preparation section.)
4. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Solution (400  $\mu$ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
5. Add 100  $\mu$ L of the Detection Antibody diluted in Reagent Diluent, to each well. Cover with a new adhesive strip and incubate 2 hours at room temperature.
6. Repeat the aspiration/wash as in step 4.
7. Add 100  $\mu$ L of the working dilution of Streptavidin-HRP B to each well. Cover the plate and incubate for 20 minutes at room temperature. Avoid placing the plate in direct light.
8. Repeat the aspiration/wash as in step 4.
9. Add 100  $\mu$ L of Substrate Solution to each well. Incubate for 20 minutes at room temperature. Avoid placing the plate in direct light.
10. Add 50  $\mu$ L of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
11. Determine the optical density of each well immediately, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.
12. **CALCULATION OF RESULTS:** Average the duplicate readings for each standard, control, and sample and subtract the average zero standard optical density (O.D.). Create a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a

standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the IL-17 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data. If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

## VIII. REFERENCES

1. Gaffen, S. (2009) *Nat. Rev. Immunol.* 9:556.
2. Iwakura, Y. et al. (2008) *Immunol. Rev.* 226:57.
3. Zhang, X. et al. (2011) *Protein Cell* 2:26.
4. Yao, Z. et al. (1995) *J. Immunol.* 155:5483.
5. Fossiez, F. et al. (1996) *J. Exp. Med.* 183:2593.
6. Hymowitz, S.G. et al. (2001) *EMBO J.* 20:5332.
7. Chang, S. and C. Dong (2007) *Cell Res.* 17:435.
8. Wright, J.F. et al. (2007) *J. Biol. Chem.* 282:13447.
9. Liang, S.C. et al. (2007) *J. Immunol.* 179:7791.
10. Yao, Z. et al. (2005) *Immunity* 3:811.
11. Katoh, S. et al. (2004) *J. Interferon Cytokine Res.* 24:553.
12. GenBank Accession #:XP\_527408.
13. Romagnani, S. et al. (2009) *Mol. Immunol.* 47:3.
14. Kolls, J.K. and A. Linden (2004) *Immunity* 21:467.
15. Witowski, J. et al. (2004) *Cell. Mol. Life Sci.* 61:567.
16. Cua, D.J. and C.M. Tato (2010) *Nat. Rev. Immunol.* 10:479..
17. Shin, H.C. et al. (1998) *Cytokine* 10:841.
18. Yao, Z. et al. (1997) *Cytokine* 9:794.
19. Haudenschild, D. et al. (2002) *J. Biol. Chem.* 277:4309.
20. Toy, D. et al. (2006) *J. Immunol.* 177:36.
21. Hu, Y. et al. (2010) *J. Immunol.* 184:4307.
22. Ely, L.K. et al. (2009) *Nat. Immunol.* 10:1245.
23. Schwarzenberger, P. et al. (2000) *J. Immunol.* 164:4783.
24. Yu, J.J. and S.L. Gaffen (2008) *Front. Biosci.* 13:170.
25. Khader, S.A. and R. Gopal (2010) *Virulence* 1:423.
26. Jovanovic, D.V. et al. (1998) *J. Immunol.* 160:3513.
27. Numasaki, M. et al. (2004) *Immunol. Lett.* 95:97.
28. Miljkovic, D. et al. (2003) *Cell. Mol. Life Sci.* 60:518.
29. Milovanovic, M. et al. (2010) *J. Invest. Dermatol.* 130:2621.
30. Zuniga, L.A. et al. (2010) *J. Immunol.* 185:6947.

**PLATE LAYOUT**

Use this plate layout to record standards and samples assayed.

12									
11									
10									
9									
8									
7									
6									
5									
4									
3									
2									
1									
	<b>A</b>	<b>B</b>	<b>C</b>	<b>D</b>	<b>E</b>	<b>F</b>	<b>G</b>	<b>H</b>	

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