## **PRODUCT INFORMATION & MANUAL**

**Mouse IL-10 Valukine™ ELISA**

**VAL605**

For the quantitative determination of natural and recombinant  
mouse IL-10 concentrations

For research use only.  
Not for diagnostic or therapeutic procedures.

**Bio-Techne China Co. Ltd**

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001  
**info.cn@bio-techne.com**

Please refer to the kit label for expiry date.  
Novus kits are guaranteed for 3 months from date of receipt

Version201912.2

## TABLE OF CONTENTS

I. BACKGROUND.....	2
II. OVERVIEW.....	3
III. ADVANTAGES.....	4
IV. EXPERIMENT.....	7
V. KIT COMPONENTS AND STORAGE.....	8
VI. PREPARATION.....	10
VII. ASSAY PROCEDURE.....	12
VIII. REFERENCES.....	14

## I. BACKGROUND

Interleukin-10 (IL-10), also known as cytokine synthesis inhibitory factor (CSIF), is the charter member of the IL-10  $\alpha$ -helical cytokine family that also includes IL-19, IL-20, IL-22, IL-24, and IL-26/AK155 (1-3). IL-10 is secreted by many activated hematopoietic cell types as well as hepatic stellate cells, keratinocytes, and placental cytotrophoblasts. Whereas human IL-10 is active on mouse cells, mouse IL-10 does not act on human cells (4, 5). Mature mouse IL-10 shares 85% amino acid sequence identity with rat IL-10 and 70-77% with bovine, canine, equine, feline, human, ovine, and porcine IL-10. It contains two intrachain disulfide bridges and is expressed as a 36 kDa noncovalently-associated homodimer (4, 6, 7).

IL-10 mediates its biological activities through a heteromeric receptor complex composed of the type II cytokine receptor subunits IL-10 R $\alpha$  and IL-10 R $\beta$ . IL-10 R $\alpha$  is a 110 kDa transmembrane glycoprotein that is expressed on lymphocytes, NK cells, macrophages, monocytes, astrocytes, intestinal epithelial cells, cytotrophoblasts, and activated hepatic stellate cells (8-13), while the 75 kDa transmembrane IL-10 R $\beta$  is widely expressed (14, 15). The IL-10 dimer binds to two IL-10 R $\alpha$  chains, triggering recruitment of two IL-10 R $\beta$  chains (14, 15). IL-10 R $\beta$  does not bind IL-10 directly but is required for signal transduction. IL-10R $\beta$  also associate with IL-20 R $\alpha$ , IL-22 R  $\alpha$ 1, or IL-28 R $\alpha$  to form the receptor complexes for IL-22, IL-26, IL-28, and IL-29 (16-18).

The involvement of IL-10 in immunoregulation includes both suppressive and stimulatory effects. It functions as an anti-inflammatory cytokine by inhibiting the expansion and activation of Th1 cells and Th17 cells (19-21) and by promoting the development of M2 macrophages (21). Its expression by immunosuppressive regulatory T cells (Treg) and regulatory B cells is important for Treg proliferation (19). Within a tumor microenvironment, however, IL-10 inhibits the expansion of Treg as well as myeloid-derived suppressor cells (22, 23). IL-10 induces the intratumoral accumulation and activation of CD8<sup>+</sup> T cells (24, 25). IL-10 exerts protective effects including limiting tissue damage in arthritic inflammation (19) and promoting muscle regeneration after injury (21), but it also contributes to the persistence of viral infections (26). The levels of IL-10 are elevated in Sjogren's syndrome (saliva), primary CNS lymphoma (30-32).

## II. OVERVIEW

### A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. A monoclonal antibody specific for mouse IL-10 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any mouse IL-10 present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked polyclonal antibody specific for mouse IL-10 is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of mouse IL-10 bound in the initial step. The color development is stopped and the intensity of the color is measured.

### B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- This kit is suitable for cell culture supernate and mouse serum.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Diluent 1× and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

### III. ADVANTAGES

#### A. PRECISION

**Intra-assay Precision** (Precision within an assay)

Two samples were tested twenty times on one plate to assess intra-assay precision.

**Inter-assay Precision** (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

CELL CULTURE SUPERNATE / MOUSE SERUM

Sample	Intra-assay Precision		Inter-assay Precision		
	1	2	1	2	3
n	20	20	20	20	20
Mean (pg/mL)	76.3	123	815	71.4	115
Standard Deviation	2.56	2.85	17.0	5.39	7.64
CV%	3.4	2.3	2.1	7.5	6.6

#### B. RECOVERY

The recovery of mouse IL-10 spiked to different levels throughout the range of the assay in cell culture media was evaluated. The recovery ranged from 95-103% with an average of 99%.

The recovery of mouse IL-10 spiked to different levels throughout the range of the assay in mouse serum was evaluated. The recovery ranged from 75.8-84.0% with an average of 80.3%.

#### C. SENSITIVITY

The minimum detectable dose (MDD) of mouse IL-10 is typically less than 1.97 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

#### D. CALIBRATION

This immunoassay is calibrated against highly purified *E. coli*-expressed recombinant mouse IL-10 produced at R&D Systems.

## E. LINEARITY

To assess the linearity of the assay, samples were spiked with high concentrations of mouse IL-10 and Calibrator Diluent (*for cell culture supernate samples*) or Calibrator Diluent-S (*for mouse serum samples*) to produce samples with values within the dynamic range of the assay.

Dilution	Average % of Expected	Range (%)
1:2	93	90-94
1:4	90	85-94
1:8	89	84-96
1:16	87	84-92

## F. SAMPLE VALUES

**Cell Culture Supernates** - Mouse spleens were removed, rinsed in 1× PBS, and kept on ice in PBS. Organs were then homogenized using a tissue homogenizer. Cells were seeded into media containing RPMI1640 supplemented with 10% fetal bovine serum, 2 mM L-glutamine, 100 U/mL penicillin and 100 ug/mL streptomycin sulfate. Cells were cultured unstimulated or stimulated with 1 pg/mL LPS for 3 days. Aliquots of the cell culture supernates were removed and assayed for levels of natural mouse IL-10.

	Unstimulated (pg/mL)	Stimulated (pg/mL)
Spleen	Non-detectable	77.1

Spleens from individual mice were removed and rinsed in PBS and kept on ice. The tissue was homogenized using a tissue homogenizer and seeded into media containing RPMI1640 supplemented with 10 % fetal bovine serum, 50 µM β-mercaptoethnal, 2 mM L-glutamine, 100 U/mL penicillin and 100 µg/mL streptomycin sulfate. Cells were stimulated with 100 ng/mL recombinant mouse IFN-γ and 1 µg/mL LPS for 4 days. An aliquot of the cell culture supernate was removed, assayed for detectable levels of natural mouse IL-10, and measured 515 pg/mL.

**Serum** - Four mouse serum samples were evaluated for the presence of IL-10 in this assay. All samples measured below the lowest standard, 15.6 pg/mL.

## G. SPECIFICITY

This assay recognizes both natural and recombinant mouse IL-10. The following factors were prepared at 50 ng/mL and assayed for cross-reactivity. Preparations of the following factors at 50 ng/mL in a mid-range recombinant mouse IL-10 control were assayed for interference. No significant cross-reactivity or interference was observed.


<b>Recombinant mouse</b>		<b>Other Recombinants</b>
IFN- $\alpha/\beta$ R2	IL-19	canine IL-10
IFN- $\gamma$ R1	IL-20	equine IL-10
IFN- $\gamma$ R2	IL-22	feline IL-10
IL-10 R	IL-22 R	human IL-10
IL-10 R $\alpha$	IL-24	porcine IL-10
IL-10 R $\beta$		viral IL-10

## IV. EXPERIMENT

### A. EXAMPLE STANDARD


The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

CELL CULTURE SUPERNATE ASSAY

pg/ml	OD	Average	Corrected
0	0.050 0.049	0.050	-
15.6	0.091 0.083	0.087	0.037
31.3	0.134 0.127	0.131	0.081
62.5	0.219 0.209	0.214	0.164
125	0.386 0.386	0.386	0.336
250	0.715 0.697	0.706	0.656
500	1.318 1.312	1.315	1.265
1000	2.266 2.255	2.261	2.211

SERUM ASSAY

pg/ml	OD	Average	Corrected
0	0.102 0.100	0.101	-
15.6	0.144 0.139	0.142	0.041
31.3	0.178 0.171	0.175	0.074
62.5	0.271 0.282	0.277	0.176
125	0.448 0.436	0.442	0.341
250	0.797 0.778	0.788	0.687
500	1.379 1.329	1.354	1.253
1000	2.285 2.120	2.203	2.102

## V. KIT COMPONENTS AND STORAGE

### A. MATERIALS PROVIDED

Parts	Description	Size
Mouse IL-10 Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with monoclonal antibody against mouse IL-10	1 plate
Mouse IL-10 Conjugate	solution of a polyclonal antibody against mouse IL-10 conjugated to horseradish peroxidase	1 vial
Mouse IL-10 Standard	recombinant mouse IL-10 in a buffered protein base; lyophilized	2 vials
Calibrator Diluent (5×)	a 5× concentrated buffered protein base used to dilute standard and cell culture supernate samples	1 vial
Calibrator Diluent-S (2×)	a 2× concentrated buffered protein base used to dilute standard and mouse serum samples	1 vial
Wash Buffer Concentrate (25×)	a 25× concentrated solution of buffered surfactant	1 vial
Color Reagent A	stabilized hydrogen peroxide	1 vial
Color Reagent B	stabilized chromogen (tetramethylbenzidine)	1 vial
Stop Solution	diluted hydrochloric acid solution	1 vial
Plate Covers	adhesive strip	3 strips

## B. STORAGE

Unopened Kit	Store at 2-8°C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Diluted Wash Solution	May be stored for up to 1 month at 2-8°C.*
	Stop Solution	
	Calibrator Dilution (1×)	
	Calibrator Diluent-S (1×)	
	Conjugate	
	Unmixed color reagent A	
	Unmixed Color reagent B	
	Standard	Use a new standard and control for each assay. Discard after use.
Microplate Wells	Return unused wells to the oil ouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8 °C	

\* Provided this is within the expiration date of the kit.

## C. OTHER SUPPLIES REQUIRED

- ◆ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ◆ Pipettes and pipette tips.
- ◆ Deionized or distilled water.
- ◆ Squirt bottle, manifold dispenser, or automated microplate washer.
- ◆ 500 mL graduated cylinder.
- ◆ Horizontal orbital microplate shaker (0.12" orbit) capable of maintaining a speed of 500(±) 50 rpm.
- ◆ Test tubes for dilution of standards and samples.

## D. PRECAUTION

The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material.

## VI. PREPARATION

### A. SAMPLE COLLECTION AND STORAGE

**Cell Culture Supernates** - Remove particulates by centrifugation and assay immediately or aliquot and store samples at  $\leq -20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles.

**Serum** - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at  $\leq -20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles.

### B. SAMPLE PREPARATION

Mouse serum samples require a 2-fold dilution. A suggested 2-fold dilution is 100  $\mu\text{L}$  of mouse serum sample + 100  $\mu\text{L}$  of Calibrator Diluent-S (1 $\times$ ).

### C. REAGENT PREPARATION

**Note:** Bring all reagents to room temperature before use.

**Wash Solution** - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25 $\times$ ) into deionized or distilled water to prepare 500 mL of Wash Buffer.

**Substrate Solution** - Substrates A and B should be mixed together in equal volumes within 15 minutes of use. Protect from light. 100  $\mu\text{L}$  of the resultant mixture is required per well.

**Calibrator Diluent (1 $\times$ )** - Add 20 mL of Calibrator Diluent 5 $\times$  into 80 mL of deionized or distilled water to prepare 100 mL of Calibrator Diluent (1 $\times$ ).


**Calibrator Diluent-S (1 $\times$ )** - Use deionized or distilled water to prepare Calibrator Diluent-S (1 $\times$ ).

**Mouse IL-10 Standard** - Refer to the vial label for reconstitution volume\* using **Calibrator Diluent (1 $\times$ ) (for cell culture supernate samples)** or **Calibrator Diluent-S (1 $\times$ ) (for mouse serum samples)**. This reconstitution produces a stock solution of 10,000 pg/mL. Allow the standard to sit for a minimum of 5 minutes with gentle agitation prior to making dilutions.

\*if you have any question, please seek help from our Technical Support.

**Pipette 900  $\mu\text{L}$  of Calibrator Diluent (1 $\times$ ) (for cell culture supernate samples) or Calibrator Diluent-S (1 $\times$ ) (for mouse serum samples) into the 1000 pg/mL tube.** Pipette 200  $\mu\text{L}$  of Calibrator Diluent (1 $\times$ ) (for cell culture supernate samples) or Calibrator Diluent-S (1 $\times$ ) (for serum samples) into the remaining tubes. Use the stock

solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The standard 1000 pg/mL serves as the high standard. The Calibrator Diluent (1×) (*for cell culture supernate samples*) or Calibrator Diluent-S (1×) (*for serum samples*) serves as the zero standard (0 pg/mL).

#### D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Substrate Solution should remain colorless until added to the plate. Keep Substrate Solution protected from light. Substrate Solution should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the Substrate Solution. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the Substrate Solution.

## VII. ASSAY PROCEDURE

**Note:** Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. A. CELL CULTURE SUPERNATE SAMPLES  
*Add 50  $\mu$ L of Calibrator Diluent (1 $\times$ ) to each well. And then add 50  $\mu$ L of Standard, cell culture supernate sample, or control using Calibrator Diluent (1 $\times$ ) per well.*  
B. MOUSE SERUM SAMPLES  
*Add 100  $\mu$ L of Standard, mouse serum sample, or control using Calibrator Diluent-S (1 $\times$ ) per well.*
4. Cover with the adhesive strip provided. Incubate for 2 hours at room temperature on a horizontal orbital microplate shaker (0.12" orbit) set at 500 $\pm$ 50 rpm. A plate layout is provided for a record of standards and samples assayed.
5. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Solution (400  $\mu$ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
6. Add 100  $\mu$ L of mouse IL-10 conjugate to each well. Cover with a new adhesive strip. Incubate for 2 hours at room temperature on a horizontal orbital microplate shaker (0.12" orbit) set at 500 $\pm$ 50 rpm.
7. Repeat the aspiration/wash as in step 5.
8. Add 100  $\mu$ L of Substrate Solution to each well. Incubate for 30 minutes at room temperature. **Protect from light.**
9. Add 100  $\mu$ L of Stop Solution to each well. The color in the wells should change from blue to yellow. If the color in the wells is green or if the color change does not appear uniform, gently tap the plate to ensure thorough mixing.
10. Determine the optical density of each well within 30 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

11. **CALCULATION OF RESULTS.** Average the duplicate readings for each standard, control, and sample and subtract the average zero standard optical density.

Create a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the IL-10 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

## VIII. REFERENCES

1. Ouyang, W. et al. (2011) *Annu. Rev. Immunol.* 29:71.
2. Sabat, R. et al. (2010) *Cytokine Growth Factor Rev.* 21:331.
3. Saraiva, M. and A. O'Garra (2010) *Nat. Rev. Immunol.* 10:170.
4. Vieira, P. et al. (1991) *Proc. Natl. Acad. Sci. USA* 88:1172.
5. Hsu, D.H. et al. (1990) *Science* 250:830
6. Windsor, W.T. et al. (1993) *Biochemistry* 32:8807
7. Syto, R. et al. (1998) *Biochemistry* 37:16943.
8. Mathurin, P. et al. (2002) *Am. J. Physiol. Gastrointest. Liver Physiol.* 282:G981
9. Szonyi, B.J. et al. (1999) *Mol. Hum. Reprod.* 5:1059.
10. Liu, Y. et al. (1994) *J. Immunol.* 152:1821.
11. Carson, W.E. et al. (1995) *Blood* 85:3577.
12. Denning, T.L. et al. (2000) *Int. Immunol.* 12:133.
13. Ledebroer, A. et al. (2002) *Eur. J. Neurosci.* 16:1175.
14. Kotenko, S.V. et al. (1997) *EMBO* 16:5894.
15. Lutfalla, G. et al. (1993) *Genomics* 16:366.
16. Kotenko, S.V. et al. (2000) *J. Biol. Chem.* 276:2725.
17. Hor, S. et al. (2004) *J. Biol. Chem.* 279:33343.
18. Sheppard, P. et al. (2003) *Nat. Immunol.* 4:63.
19. Carter, N.A. et al. (2011) *J. Immunol.* 186:5569
20. Gu, Y. et al. (2008) *Eur. J. Immunol.* 38:1807.
21. Deng, B. et al. (2012) *J. Immunol.* 189:3669
22. Tanikawa, T. et al. (2012) *Cancer Res.* 72:420.
23. Groux, H. et al. (1997) *Nature* 389:737.
24. Mumm, J.B. et al. (2011) *Cancer Cell* 20:781.
25. Emmerich, J. et al. (2012) *Cancer Res.* 72:3570.
26. Wilson, E.B. and D.G. Brooks (2011) *Curr. Top. Microbiol. Immunol.* 350:39.
27. Bertorello, R. et al. (2004) *Clin. Exp. Med.* 4:148
28. Sasayama, T. et al. (2012) *Neuro. Oncol.* 14:368.
29. Mustea, A. et al. (2006) *Anticancer Res.* 26:1715.
30. George, J. et al. (2012) *Atherosclerosis* 222:519.
31. Borekci, B. et al. (2007) *Am. J. Reprod. Immunol.* 58:56.
32. Camejo, MI (2003) *Arch. Androl.* 49:111.

## PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

<b>1</b>								
<b>2</b>								
<b>3</b>								
<b>4</b>								
<b>5</b>								
<b>6</b>								
<b>7</b>								
<b>8</b>								
<b>9</b>								
<b>10</b>								
<b>11</b>								
<b>12</b>								
	<b>A</b>	<b>B</b>	<b>C</b>	<b>D</b>	<b>E</b>	<b>F</b>	<b>G</b>	<b>H</b>

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