PRODUCT INFORMATION & MANUAL

Mouse VEGF Valukine™ ELISA

VAL608

For the quantitative determination of natural and recombinant
mouse VEGF concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version201912.2

TABLE OF CONTENTS

I. BACKGROUND.....	2
II. OVERVIEW.....	4
III. ADVANTAGES.....	5
IV. EXPERIMENT.....	8
V. KIT COMPONENTS AND STORAGE.....	9
VI. PREPARATION.....	11
VII. ASSAY PROCEDURE.....	13
VIII. RERERENCES.....	15

I. BACKGROUND

Vascular endothelial growth factor (VEGF or VEGF-A), also known as vascular permeability factor (VPF), is a potent mediator of both angiogenesis and vasculogenesis in the fetus and adult (1-3). It is a member of the PDGF family that is characterized by the presence of eight conserved cysteine residues in a cystine knot structure and the formation of antiparallel disulfide-linked dimers (4). Alternately spliced isoforms of 120, 164 and 188 amino acids (aa) have been found in mouse, while 121, 145, 165, 183, 189, and 206 aa isoforms have been identified in human (2,4). In humans, VEGF165 appears to be the most abundant and potent isoform, followed by VEGF121 and VEGF189 (3, 4). The same pattern may exist in mouse. Isoforms other than VEGF120 and 121 contain basic heparin-binding regions and are not freely diffusible (4). Mouse VEGF164 shares 97% aa sequence identity with corresponding regions of rat VEGF. It also shares 89% aa identity with human and porcine VEGF, 88% with bovine VEGF, and 90% with feline, equine and canine VEGF. VEGF is expressed in multiple cells and tissues including skeletal and cardiac muscle (5,6), hepatocytes (7), osteoblasts (8), neutrophils (9), macrophages (10), keratinocytes (11), brown adipose tissue (12), CD34⁺ stem cells (13), endothelial cells (14), fibroblasts, and vascular smooth muscle cells (15). VEGF expression is induced by hypoxia and cytokines such as IL-1, IL-6, IL-8, oncostatin M and TNF- α (3, 4, 9, 16). The isoforms are differentially expressed during development and in the adult (3).

VEGF dimers bind to two related receptor tyrosine kinases, VEGF R1 (also called Flt-1) and VEGF R2 (Flk-1/KDR) and induce their homodimerization and autophosphorylation (3, 4, 7, 17, 18). These receptors have seven extracellular immunoglobulin-like domains and an intracellular split tyrosine kinase domain. They are expressed on vascular endothelial cells and a range of non-endothelial cells. Although VEGF affinity is highest for binding to VEGF R1, VEGF R2 appears to be the primary mediator of VEGF angiogenic activity (3, 4). VEGF165 also binds the semaphorin receptor, neuropilin-1, which promotes complex formation with VEGF R2 (19).

VEGF is best known for its role in vasculogenesis. During embryogenesis, VEGF regulates the proliferation, migration, and survival of endothelial cells (3, 4), thus

regulating blood vessel density and size but playing no role in determining vascular patterns. VEGF promotes bone formation through osteoblast and chondroblast recruitment and is also a monocyte chemoattractant (20-22). In postnatal life, VEGF maintains endothelial cell integrity and is a potent mitogen for micro- and macro-vascular endothelial cells. In adults, VEGF functions mainly in wound healing and the female reproductive cycle (3). In diseased tissues, VEGF promotes vascular permeability. It is thus thought to contribute to tumor metastasis by promoting both extravasation and tumor angiogenesis (23, 24). Various strategies have been employed therapeutically to antagonize VEGF-mediated tumor angiogenesis (25). Circulating VEGF levels correlate with disease activity in autoimmune diseases such as rheumatoid arthritis, multiple sclerosis and systemic lupus erythematosus (26).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. Monoclonal antibody specific for mouse VEGF has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any mouse VEGF present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked polyclonal antibody specific for mouse VEGF is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of mouse VEGF bound in the initial step. The color development is stopped and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- This kit is suitable for cell culture supernate and mouse serum.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Diluent 1× and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Two samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

CELL CULTURE SUPERNATE / MOUSE SERUM

	Intra-assay Precision		Inter-assay Precision		
Sample	1	2	1	2	3
n	20	20	20	20	20
Mean (pg/mL)	52.1	198	19.5	50	252
Standard Deviation	3.8	15	3.0	5.3	16
CV%	7.3	7.5	15.4	10.6	6.3

B. RECOVERY

The recovery of mouse VEGF spiked to different levels throughout the range of the assay in cell culture media was evaluated. The recovery ranged from 71-92% with an average of 80%.

The recovery of mouse VEGF spiked to four levels throughout the range of the assay in mouse serum was evaluated. The recovery ranged from 77.5-88.7% with an average of 83.7%.

C. SENSITIVITY

The minimum detectable dose (MDD) of mouse VEGF is typically less than 4.71 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against highly purified SF21-expressed recombinant mouse VEGF produced at R&D Systems.

E. LINEARITY

To assess the linearity of the assay, samples were spiked with high concentrations of mouse VEGF and Calibrator Diluent (*for cell culture supernate samples*) or Calibrator Diluent-S (*for mouse serum samples*) to produce samples with values within the dynamic range of the assay.

Dilution	Average % of Expected	Range (%)
1:2	104	96 - 111
1:4	108	99 - 117
1:8	108	100 -118
1:16	102	83 -114

F. SAMPLE VALUES

Cell Culture Supernates - LLC1 cells (Lewis lung carcinoma) were seeded at 5×10^5 cells/mL and cultured for 4 days in 25 mL of DMEM supplemented with 10% fetal bovine serum. The cell culture supernate was assayed for mouse VEGF and measured 907 pg/mL.

C2C12 cells (mouse muscle myoblast) were seeded at 2.75×10^6 cells/mL and cultured for 3 days in 100 ml DMEM supplemented with 10% fetal bovine serum and 2mM L-glutamine. The cell culture supernate was assayed for mouse VEGF and measured 885 pg/mL.

WEHI3 cells (Mouse myeloid cells that produced mll-3) were seeded at 2.75×10^6 cells/mL and cultured for 3 days in 100 mL DMEM supplemented with 10% fetal bovine serum and 2mM L-glutamine. The cell culture supernate was assayed for mouse VEGF and measured 889 pg/mL.

Serum - Four mouse serum samples were evaluated for the presence of VEGF in this assay. All samples measured below the lowest standard, 7.8 pg/mL.

G. SPECIFICITY

This assay recognizes both the 164 and 120 amino acid residue forms of mouse VEGF. The factors listed below were prepared at 50 ng/ml in Diluent 1× and assayed for cross-reactivity. Preparations of the following factors prepared at 50 ng/ml in a mid-range mouse VEGF standard were assayed for interference. No significant cross-reactivity or interference was observed.

Recombinant mouse	
VEGF115	VEGF R2
VEGF167	VEGF R3
VEGF186	


At concentrations greater than 78 pg/mL, rmVEGF R1 interferes in this assay.

IV. EXPERIMENT

A. EXAMPLE STANDARD


The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

CELL CULTURE SUPERNATE ASSAY

pg/ml	OD	Average	Corrected
0	0.061	0.066	-
7.8	0.099	0.101	0.035
15.6	0.136	0.138	0.072
31.3	0.211	0.212	0.147
62.5	0.342	0.350	0.285
125	0.600	0.612	0.547
250	1.141	1.149	1.083
500	1.969	2.036	1.971

SERUM ASSAY

pg/ml	OD	Average	Corrected
0	0.120	0.121	-
7.8	0.177	0.183	0.062
15.6	0.238	0.250	0.130
31.3	0.357	0.372	0.251
62.5	0.571	0.600	0.480
125	0.965	0.980	0.860
250	1.610	1.614	1.493
500	2.390	2.396	2.275

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
Mouse VEGF Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with monoclonal antibody against mouse VEGF	1 plate
Mouse VEGF Conjugate	Solution of a polyclonal antibody against mouse VEGF conjugated to horseradish peroxidase with preservatives	1 vial
Mouse VEGF Standard	recombinant mouse VEGF in a buffered protein base with preservatives; lyophilized.	1 vial
Calibrator Diluent (5×)	a 5× concentrated buffered protein base used to dilute standard and cell culture supernate samples	1 vial
Calibrator Diluent-S (2×)	a 2× concentrated buffered protein base used to dilute standard and mouse serum samples	1 vial
Wash Buffer Concentrate (25×)	a 25× concentrated solution of buffered surfactant with preservatives	1 vial
Color Reagent A	stabilized hydrogen peroxide	1 vial
Color Reagent B	stabilized chromogen (tetramethylbenzidine)	1 vial
Stop Solution	diluted hydrochloric acid solution	1 vial
Plate Covers	adhesive strip	3 strips

B. STORAGE

Unopened Kit	Store at 2-8°C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Diluted Wash Solution	May be stored for up to 1 month at 2-8°C.*
	Stop Solution	
	Calibrator Dilution (1×)	
	Calibrator Diluent-S (1×)	
	Conjugate	
	Unmixed Substrate A	
	Unmixed Substrate B	
	Standard	Aliquot and store for up to 1 month at -20°C in a manual defrost freezer.* Avoid repeated freeze-thaw cycles.
Microplate Wells	Return unused wells to the pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*	

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- Pipettes and pipette tips.
- Deionized or distilled water.
- Squirt bottle, manifold dispenser, or automated microplate washer.
- 500 mL graduated cylinder.

D. PRECAUTION

The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles.

B. SAMPLE PREPARATION

Mouse serum samples require a 5-fold dilution. A suggested 5-fold dilution is 40 μ L of mouse serum sample + 160 μ L of Calibrator Diluent-S (1×).

C. REAGENT PREPARATION

Note: Bring all reagents to room temperature before use.

Wash Solution - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25×) into deionized or distilled water to prepare 500 mL of Wash Buffer.

Substrate Solution - Substrates A and B should be mixed together in equal volumes within 15 minutes of use. Protect from light. 100 μ L of the resultant mixture is required per well.

Calibrator Diluent (1×) - Add 20 mL of Calibrator Diluent 5× into 80 mL of deionized or distilled water to prepare 100 mL of Calibrator Diluent (1×).


Calibrator Diluent-S (1×) - Use deionized or distilled water to prepare Calibrator Diluent-S (1×).

Mouse VEGF Standard - Refer to the vial label for reconstitution volume* using Calibrator Diluent (1×) (for cell culture supernate samples) or Calibrator Diluent-S (1×) (for mouse serum samples). This reconstitution produces a stock solution of 500 pg/mL. Allow the standard to sit for a minimum of 5 minutes with gentle

agitation prior to making dilutions.

*if you have any question, please seek help from our Technical Support.

Pipette 200 μ L of Calibrator Diluent (1 \times) (for cell culture supernate samples) or Calibrator Diluent-S (1 \times) (for mouse serum samples) into each tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The undiluted standard 500 pg/mL serves as the high standard. The Calibrator Diluent (1 \times) (for cell culture supernate samples) or Calibrator Diluent-S (1 \times) (for serum samples) serves as the zero standard (0 pg/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Substrate Solution should remain colorless until added to the plate. Keep Substrate Solution protected from light. Substrate Solution should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the Substrate Solution. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the Substrate Solution.

VII. ASSAY PROCEDURE

Note: *Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.*

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 50 μ L of Calibrator Diluent (1 \times) (*for cell culture supernate samples and its standard*) or Calibrator Diluent-S (1 \times) (*for mouse serum samples and its standard*) to each well.
4. Add 50 μ L of Standard, sample, or control per well. Cover with the adhesive strip provided. Incubate for 2 hours at room temperature. A plate layout is provided for a record of standards and samples assayed.
5. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Solution (400 μ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
6. Add 100 μ L of mouse VEGF conjugate to each well. Cover with a new adhesive strip. Incubate for 2 hours at room temperature.
7. Repeat the aspiration/wash as in step 5.
8. Add 100 μ L of Substrate Solution to each well. Incubate for 30 minutes at room temperature. **Protect from light.**
9. Add 100 μ L of Stop Solution to each well. The color in the wells should change from blue to yellow. If the color in the wells is green or if the color change does not appear uniform, gently tap the plate to ensure thorough mixing.
10. Determine the optical density of each well within 30 minutes, using a

microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

11. **CALCULATION OF RESULTS.** Average the duplicate readings for each standard, control, and sample and subtract the average zero standard optical density.

Create a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the VEGF concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. RERERENCES

1. Leung, D.W. et al. (1989) *Science* 246:1306.
2. Shima, D.T. et al. (1996) *J. Biol. Chem.* 271:3877.
3. Byrne, A.M. et al. (2005) *J. Cell. Mol. Med.* 9:777.
4. Robinson, C.J. and Stringer, S.E. (2001) *J. Cell. Sci.* 114:853.
5. Richardson, R.S. et al. (1999) *Am. J. Physiol.* 277:H2247.
6. Sugishita, Y. et al. (2000) *Biochem. Biophys. Res. Commun.* 268:657.
7. Yamane, A. et al. (1994) *Oncogene* 9:2683.
8. Goad, D.L. et al. (1996) *Endocrinology* 137:2262.
9. Gaudry, M. et al. (1997) *Blood* 90:4153.
10. Mdaren, J. et al. (1996) *J. Clin. Invest.* 98:482.
11. Diaz, B.V. et al. (2000) *J. Biol. Chem.* 275:642.
12. Asano, A. et al. (1997) *Biochem. J.* 328:179.
13. Bautz, F. et al. (2000) *Exp. Hematol.* 28:700.
14. Namiki, A. et al. (1995) *J. Biol. Chem.* 270:31189.
15. Nauck, M. et al. (1997) *Am. J. Respir. Cell. Mol. Biol.* 16:398.
16. Angelo, L.S. and R. Kurzrock (2007) *Clin. Cancer Res.* 13:2825.
17. Neufeld, G. et al. (1999) *FASEB. J.* 13:9.
18. Kowalewski, M.P. et al. (2005) Accession #ABB82619.
19. Pan, Q. et al. (2007) *J. Biol. Chem.* 282:24049.
20. Weis, S.M. and D.A. Cheresh (2005) *Nature* 437:497.
21. Breier, G. (2000) *Semin. Thromb. Hemost.* 26:553.
22. Barleon, B. et al. (1996) *Blood* 87:3336.
23. Weis, S.M. and D.A. Cheresh (2005) *Nature* 437:497.
24. Thurston, G. (2002) *J. Anal.* 200:575.
25. Grothey, A. and E. Galanis (2009) *Nat. Rev. Clin. Oncol.* 6:507.
26. Carvalho, J.F. et al. (2007) *J. Clin. Immunol.* 27:246.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								
	A	B	C	D	E	F	G	H

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