


Luminex[®] Multiplex Immunoassays for Classification of Th1, Th2, and Th17 CD4⁺ T Cell Differentiation

Developing CD4⁺ T cells in the thymus differentiate into various distinct subsets. Each subset produces a characteristic combination of cytokines to activate other immune cells, control ongoing immune responses, or manage immune memory.

Bio-Techne has developed R&D Systems[®] Luminex[®] Immunoassays to specifically quantify biomarkers of Th1/Th2 and Th9/Th17/Th22 differentiated CD4⁺ T cells. These two Luminex kits allow for the characterization of T cell lineages present in CD4⁺ T cells which have been induced towards specific T cell subsets.



METHODS

Cell Culture

Using the MagCelect[™] Human CD4⁺T cell Isolation Kit (R&D Systems; Catalog # [MAGH102](#)), CD4⁺T cells were isolated from Human peripheral blood mononuclear cells cultured in RPMI and supplemented with 10% fetal bovine serum, 2 mM L-glutamine, 100 U/mL penicillin, and 100 µg/mL streptomycin sulfate. These cells were either untreated or treated with Th1, Th2 or Th17-inducing reagents.

Th1 Cells: Flasks were coated with 2 µg/mL goat anti-mouse IgG ([AF007](#)) overnight at room temperature. Wash 1X with PBS, and then coat with 1 µg/mL Ms x hCD3ε ([MAB100](#)) 2-3 hours at room temperature. Wash 1X with PBS prior to adding CD4⁺T cells seeded at 0.5 – 1 x 10⁶/mL. To activate the cells the following polarizing cytokines and antibodies were then added to the flask and incubated for 5 days:

20 ng/ml rhIL-2 ([202-IL](#))
40 ng/ml rhIL-12 ([219-IL](#))

5 µg/ml Ms x hIL-4 ([MAB304](#))
5 µg/ml Ms x hCD28 ([MAB342](#))

Further incubate stimulated cells for 24 hours with PMA (10 ng/mL; Tocris; Catalog # [1201/1](#)) and Ionomycin calcium salt (500 ng/mL; Tocris; Catalog # [1704/1](#))

Th17 cells: Flasks were coated with 2 µg/mL goat anti-mouse IgG ([AF007](#)) overnight at room temperature. Wash 1X with PBS, and then coat with 1 µg/mL Ms x hCD3ε ([MAB100](#)) 2-3 hours at room temperature. Wash 1X with PBS prior to adding CD4⁺T cells seeded at 0.5 – 1 x 10⁶/mL. To activate the CD4⁺ T cells the following cytokines and antibodies were then added to the flask and incubated for 5 days:

20 ng/mL rhIL-2 (202-IL)	40 ng/mL rhIL-6 (206-IL)
20 ng/mL rhIL-23 (1290-IL)	5 µg/mL Ms x hCD28 (MAB342)
10 ng/mL rhIL-1β (201-LB)	

Further incubate stimulated cells for 24 hours with PMA (10 ng/mL; Tocris; Catalog # [1201/1](#)) and Ionomycin calcium salt (500 ng/mL; Tocris; Catalog # [1704/1](#))

Th2 cells: Add CD4⁺T cells seeded at 0.5 – 1 x 10⁶/mL to an un-coated flask. To activate the CD4⁺ T cells the following were then added to the flask and incubated for 5 days:

20 ng/ml rhIL-2 (202-IL)	5 ug/ml PHA
--	-------------

Further incubate stimulated cells for 24 hours with PMA (10 ng/mL; Tocris; Catalog # [1201/1](#)) and Ionomycin calcium salt (500 ng/mL; Tocris; Catalog # [1704/1](#))

Luminex Fixed Panel Analysis

Cell culture supernates were analyzed with either the Luminex Th1/Th2 Fixed Panel (Catalog # [LKTM008](#)) or the Luminex Th9/Th17/Th22 Fixed Panel (Catalog # [LKTM009](#)). Samples were assayed according to the procedures outlined in the product inserts. Both kits are available from R&D Systems.

RESULTS

Cells polarized towards a Th1 phenotype secreted increased levels of IFN-α and TNF-β, as expected for Th1 differentiation (Figure **1B**). A hallmark of Th2 differentiation is increased levels of IL-4, IL-5, and IL-13, which we observed (Figure **1A**). Stimulating cells with IL-2, IL-23, IL-1β and IL-6 drove activated CD4⁺ T cell differentiation to a Th17 phenotype, as observed by the secretion of IL-17A (Figure **1D**).

CONCLUSION

The cytokine profiles measured with our Luminex Fixed Panels for Th1/Th2 and Th9/Th17/Th22 were consistent with expected results for each respective T cell phenotype. R&D Systems fixed panels provide a convenient, efficient, and highly-validated immunoassay option to generate biomarker/cytokine profiles in response to intracellular and extracellular pathogenic bacteria, viruses, and fungi. These fixed panels are provided as an off-the-shelf product, making it easier and faster than ever to get the Luminex kits you need.




Figure 1 - Th1, Th2 and Th17 Cytokine Differentiation Profiles

A & B - Cytokine profile for the TH1/Th2 fixed panel analytes.

C & D - Cytokine profile for the Th9/Th17/Th22 fixed panel analytes.

+ = values above the limits of the standard curve

• = stimulating cytokine

aCD3/28 = cells from an uncoated flask incubated in RPMI and 10% FBS for 5 days

NT = Non-treated cells from an antibody coated flask with only Ms x hCD28 added for 5 days

Note: The concentrations represented in this figure have not been corrected for the initial 1:2 sample dilution

REFERENCES

- Suntharalingam, G. *et al.* (2006) *N. Engl. J. Med.* **335**:1018.
- Maude, S.L. *et al.* (2014) *Cancer J.* **20**:119.
- Römer, P.S. *et al.* (2011) *Blood* **118**:6772.
- Findlay, L. *et al.* (2010) *J. Immunol. Methods* **352**:1.
- Stebbing, R. *et al.* (2013) *J. Immunotoxicol.* **10**:75.
- Attarwala, H. (2010). *J. Young Pharm.* **2**:3
- Hünig, T. *et al.* (2011) *Adv Immunol.* **95**:111.

R&D SYSTEMS

NOVUS BIOLOGICALS

TOCRIS

protein simple

ACD™

biotechne®

Global [bio-techne.com](#) [info@bio-techne.com](#) TEL +1 612 379 2956 North America TEL 800 343 7475
Europe | Middle East | Africa TEL +44 (0)1235 529449 China [info.cn@bio-techne.com](#) TEL +86 (21) 52380373

For research use or manufacturing purposes only. Trademarks and registered trademarks are the property of their respective owners.