# Red Systems a biotechne® brand

## Stem Cell-Derived Hepatocyte-Like Cells: An Alternative for High Throughput Toxicity Screening

J. Aho, D. Galitz, J. Sabat, and F. Rinaldi | R&D Systems, Inc., 614 McKinley Place NE, Minneapolis, MN, 55413

#### Abstract

Secondary liver toxicity is a leading cause of commercial pharmaceutical drug recalls and compound failures during drug development. Primary hepatocytes, the most common *in vitro* model for drug-induced liver toxicity testing, are not ideal for high throughput screening due to their limited availability, difficulty to culture, and functional instability in vitro. Hepatocytelike cells, derived from human pluripotent stem cells, are emerging as a stable and renewable model for drug-induced liver toxicity testing. In this study, we hypothesize that hepatocyte-like cells derived from induced pluripotent stem cells using the StemXVivo<sup>®</sup> Hepatocyte Differentiation Kit can be used for high throughput toxicity screens. We have previously demonstrated that the StemXVivo<sup>®</sup> Hepatocyte Differentiation Kit efficiently and consistently generates functional hepatocyte-like cells. To provide further functional validity of these cells we first show that Cyp3A4, an essential hepatic enzyme often analyzed in drug toxicity studies, can be induced in kit-derived hepatocytes. We then used high content imaging to quantitate cell viability against a panel of known hepatotoxic molecules. Our data indicate that the StemXVivo<sup>®</sup> Hepatocyte Differentiation Kit provides a reliable and renewable source of hepatocyte-like cells that can be utilized for high throughput toxicology and drug discovery.

Pluripotent Stem Cell-derived Hepatocyte-like Cells Display Characteristic Hepatocyte Morphology and Cell-specific Protein Expression





### **StemXVivo® Hepatocyte Differentiation Kit**



**Stem Cell-derived Hepatocyte-like Cells have Characteristic Hepatocyte Morphology.** BG01V human embryonic stem cells and iPSK3 induced pluripotent stem cells were differentiated into hepatocyte-like cells using the differentiation kit. Brightfield images of BG01V-derived and iPSK3-derived hepatocytes show the characteristic cobblestone morphology of hepatocytes grown *in vitro*.



Hepatocyte-specific Protein Expression in Pluripotent Stem Cell-derived Hepatocyte-like Cells. iBJ6 human induced pluripotent stem cells were differentiated into hepatocyte-like cells using the differentiation kit. Cells at day 0, 14, and 20 of differentiation were collected, lysed, and analyzed by Simple Western<sup>™</sup> for hepatocyte-specific protein expression. Expression of the stem cell marker, Oct-3/4 (Catalog # AF1759), was observed in undifferentiated iBJ6 cells, but not in cells differentiating into hepatocyte-like cells. Hepatocyte-like cells showed increasing expression of hepatocyte-specific proteins, alpha-Fetoprotein (AFP; Catalog # MAB1368), Serum Albumin (Catalog # AF3329),Cytokeratin 18 (Catalog # MAB7619), and Serpin A1(Catalog # AF1268). HSP60 (Catalog # MAB1800) was used as a housekeeping protein.

#### Pluripotent Stem Cell-derived Hepatocyte-like Cells as an in vitro Liver Toxicity Model



Day 0	Day 5	Day 9	Day 13	Day 20
Starting Pluripotent Stem Cells	Endoderm	Hepatic Endoderm	Hepatoblasts	Hepatocyte-like cells

#### Functional Cyp3A4 Activity in Pluripotent Stem Cell-derived Hepatocytes



Hepatotoxicity to Small Molecules Assessed Using High Content Imaging. iPSK3 (A,B) and iBJ6 (C,D) human pluripotent stem cell-derived hepatocyte-like cells were used to screen small molecules for hepatotoxicity using a cell viability (A,C) and a lipid accumulation (B,D) assay. Hepatocyte-like cells were differentiated in 6-well plates according to the protocol and replated onto Matrigel-coated 96-well plates at 60,000 cells per well for use in high throughput screening experiments. Hepatocyte-like cells were either left untreated (DMSO-alone) or incubated with incremental doses of the known hepatotoxic compounds: Carboplatin (Catalog # 2626), Doxorubicin (Catalog # 2252), Crizotinib (Catalog # 5119), Tacrine (Catalog # 0965), Sunitinib (Catalog # 3768), Staurosporine (Catalog # 1285), and Entacapone (Catalog # 4720). Following 48 hours in culture with compounds, cells were incubated with (A,C) Calcein AM (Catalog # 5119), a fluorescent live-cell marker, or (B,D) LipidTox<sup>\*\*</sup> (ThermoFisher) a fluorescent lipid stain. Cell staining was quantified using the Operetta CLS<sup>\*\*</sup> High Content Analysis System. Statistical significance was determined using an unpaired T-test:  $* = p \le 0.01$ ,  $** = p \le 0.05$ .

## Conclusions

- The StemXVivo® Hepatocyte Differentiation Kit efficiently differentiates human pluripotent stem cells into hepatocyte-like cells.
- Stem cell-derived hepatocyte-like cells demonstrate liver functions, including p450 activity which is important for drug metabolism.
- Stem cell-derived hepatocyte-like cells provide an *in vitro* model to assess hepatotoxicity.
  - Hepatocyte-like cells maintain functional when replated into 96-well plates.
  - 96-well plate format enables high throughput screening of hepatotoxic drugs.
  - Known hepatotoxic compounds are dose-dependently toxic to stem cell-derived hepatocyte-like cells.

• Fluorescent assays are available for high throughput experiments.

For research use or manufacturing purposes only. PS\_Hepatocyte\_ISSCR\_6207

Induction and Activity of Cyp3A4 Expression in Hepatocytes Derived from Human Embryonic Stem Cells. Hepatocyte-like cells differentiated from BG01V human embryonic stem cells using the differentiation kit were treated for 48 hours with 50  $\mu$ M Dexamethasone (Catalog # 1126) and/or 20  $\mu$ M Ketoconazole (Catalog # 1103) to induce and inhibit Cyp3A4 expression, respectively. Dexamethasone treatment resulted in an increase in Cyp3A4 activity compared to untreated DMSO alone controls. Ketoconazole treatment inhibited Cyp3A4 activity in untreated and Dexamethasone-treated samples. Cyp3A4 activity was measured using the P450-Glo CYP3A4 cell-based assay (Promega).

