PRODUCT INFORMATION & MANUAL

Human IL-8/CXCL8 Valukine™ ELISA

VAL103

For the quantitative determination of natural and recombinant
human IL-8/CXCL8 concentrations

For research use only.

Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001

info.cn@bio-techne.com

Novus kits are guaranteed for 3 months from date of receipt

Version 201910.3

TABLE OF CONTENTS

I. BACKGROUND.....	2
II. OVERVIEW.....	4
III. ADVANTAGES.....	5
IV. EXPERIMENT.....	7
V. KIT COMPONENTS AND STORAGE.....	8
VI. PREPARATION.....	10
VII. ASSAY PROCEDURE.....	12
VIII. REFERENCES.....	13

I. BACKGROUND

Interleukin-8 (IL-8; also known as GCP-1, NAP-1, and CXCL8) is a heparin-binding, 8-9 kDa member of the alpha, or CXC family of chemokines. Currently, there are 15 human CXC family members that generally range from 8-12 kDa in size. The majority are found on human chromosome 4, all contain a typical three β -sheet/one α -helix structure, and most show an N-terminal GluLeuArg (ELR) tri-peptide motif (1-3). Human IL-8 is synthesized as a 99 amino acid (aa) precursor that contains a 20 aa signal sequence plus a 79 aa mature region (4-6). It circulates as a monomer, homodimer, and heterodimer with CXCL4/PF4. The monomer is considered the most bioactive, while the heterodimer is proposed to potentiate PF4 activity (7-10). Mature human IL-8 shares 65% and 70% aa sequence identity with porcine and canine IL-8, respectively (11, 12). There is no IL-8 gene counterpart in rodents. Multiple isoforms of IL-8 exist that are generated through both alternative splicing and differential proteolytic processing. Alternative splicing occurs at the C-terminus where there is an

11 aa substitution for amino acids 92-99 (13). Proteolytic processing is likely a cell-specific event and results in IL-8 N-terminal truncation. For example, fibroblasts and endothelial cells cleave amino acids 21 and 22, generating IL-8 (1-7), while monocytes and lymphocytes cleave amino acids 21-25, generating IL-8 (6-7). These short forms generally show increased bioactivity, particularly towards the CXCR1 IL-8 receptor (6, 14). Approximately 15% of IL-8 also undergoes citrullination on Arg27 of the precursor. This creates an IL-8 with an increased half-life and ability to induce leukocytosis (15, 16). A wide variety of cells secrete IL-8, and these include monocytes and neutrophils (17), fibroblasts and keratinocytes (18), mast cells (19), visceral smooth muscle cells (20), dendritic cells (21), type II Greater alveolar cells (22), and endothelial cells (23).

There are two G-protein-coupled receptors for IL-8, termed CXCR1/IL-8RA and CXCR2/IL-8RB, that share 77% aa sequence identity (24). CXCR1 is 45-50 kDa in size and used almost exclusively by IL-8. CXCR2 is 35-40 kDa in size and used by almost all CXC chemokines (25, 26). Both CXCR1 and CXCR2 constitutively homodimerize, and this appears to be their functional configuration. When expressed on the same cell, heterodimerization will also occur, but this configuration undergoes disassembly following IL-8 binding (27). CXCR2 responds to low concentrations of IL-8 and is principally associated with chemotaxis and MMP-9 release. CXCR1, by contrast, responds to high concentrations of IL-8 and is associated with respiratory burst and phospholipase D2 activation (26). Thus, and in the case of neutrophils, CXCR2 is suggested to initiate PMN migration to an inflammatory site followed by CXCR1 mediated priming for anti-microbial activity.

IL-8 is perhaps best known for its proinflammatory effects on immune cells. In essence, IL-8 is secreted by multiple cell types exposed to inflammatory stimuli. For monocytes/macrophages, microbial exposure results in IL-8 release. This is followed by CXCR2-mediated chemotaxis of PMNs to sites of antigenic challenge, and their

subsequent activation/priming for anti-microbial activity. This activity is complemented by IL-8's ability to potentiate M-CSF action in the bone marrow that results in the release and maturation of granulocytes (14). IL-8 is also reported to have an angiogenic effect on tumor-associated endothelial cells. Here, tumor-derived IL-8 can act in a paracrine manner to activate CXCR1 and R2 on endothelial cells. Both receptors have been associated with PI3-K/Akt and RasGTP signaling, resulting in cell survival and proliferation. In addition, IL-8 positively regulates VEGF R2 and EGF R, two TKRs that mediate cell growth and migration (28).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. A monoclonal antibody specific for IL-8 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any IL-8 present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked polyclonal antibody specific for IL-8 is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of IL-8 bound in the initial step. The color development is stopped and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- ◆ **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- ◆ This kit is suitable for cell culture supernate and serum.
- ◆ The kit should not be used beyond the expiration date on the kit label.
- ◆ Do not mix or substitute reagents with those from other lots or sources.
- ◆ If samples generate values higher than the highest standard, dilute the samples with Diluent and repeat the assay.
- ◆ Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

Sample	Intra-assay Precision			Inter-assay Precision		
	1	2	3	1	2	3
Mean (pg/mL)	166	414	708	164	410	751
Standard Deviation	6.3	18.3	33.2	9.4	23.4	33.4
CV%	3.8	4.4	4.7	5.7	5.7	4.4

B. RECOVERY

The recovery of human IL-8 spiked to different levels throughout the range of the assay in cell culture media was evaluated. The recovery ranged from 96-118% with an average of 108%.

The recovery of human IL-8 spiked to different levels throughout the range of the assay in human serum was evaluated. The recovery ranged from 78.7-87.7% with an average of 82.6%.

C. SENSITIVITY

The minimum detectable dose (MDD) of IL-8 is typically less than 7.8 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against highly purified *E.coli*-expressed recombinant human IL-8 produced at R&D Systems®.

E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with high concentrations of IL-8 and diluted with Diluent 1× to produce samples with values within the dynamic range of the assay.

Dilution	Average % of Expected	Range (%)
1:2	114	105 – 125
1:4	114	105 – 122
1:8	112	104 – 124
1:16	108	102 – 107

F. SAMPLE VALUES

Cell Culture Supernates -Human peripheral blood mononuclear cells (1×10^6 cells/mL) were cultured in RPMI supplemented with 10% fetal calf serum, 50 μ M -mercaptoethanol, 2 mM L-glutamine, 100 U/mL penicillin, and 100 μ g/mL streptomycin sulfate and stimulated for 3 days with 10 μ g/mL PHA. An aliquot of the cell culture supernate was removed, assayed for levels of natural IL-8, and measured 403,000 pg/mL.

Serum - Four serum samples were evaluated for the presence of IL-8 in this assay. All samples measured ranged from 92.8 to 100.6 pg/mL with an average of 96.9 pg/mL.

G. SPECIFICITY

This assay recognizes both natural and recombinant human IL-8. The following factors were prepared at 50 ng/mL and assayed for cross-reactivity. Preparations of the following factors at 50 ng/mL in a mid-range rhIL-8 control were assayed for interference. No significant cross-reactivity or interference was observed.

Recombinant human

GRO α
 GRO β
 GRO γ
 I-309
 IP-10
 MCP-1
 MCP-2
 MCP-3
 MIP-1 α
 MIP-1 β


Recombinant mouse

KC
 MIP-1 α
 MIP-1 β

IV. EXPERIMENT

EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

pg/mL	OD	Average	Corrected
0	0.007		
	0.007	0.007	—
	0.114		
15.6	0.11	0.112	0.105
	0.212		
31.2	0.211	0.212	0.205
	0.358		
62.5	0.395	0.392	0.385
	0.743		
125	0.755	0.749	0.742
	1.305		
250	1.274	1.29	1.283
	2.076		
500	2.104	2.091	2.084
	3.065		
1000	2.994	3.03	3.023

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
IL-8 Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with a mouse monoclonal antibody against human IL-8	1 plate
IL-8 Conjugate	solution of polyclonal antibody against IL-8 conjugated to horseradish peroxidase	1 vial
IL-8 Standard	recombinant human IL-8 in a buffered protein base; lyophilized	1 vial
Calibrator Diluent(5×)	A 5× concentrated buffered protein base	1 vial
Wash Buffer Concentrate(25×)	A 25× concentrated solution of buffered surfactant	1 vial
Color Reagent A	Stabilized hydrogen peroxide	1 vial
Color Reagent B	Stabilized chromogen (tetramethylbenzidine)	1 vial
Stop Solution	2 N sulfuric acid	1 vial
Plate Sealers	Adhesive strip	3 strips

B. STORAGE

Unopened Kit	Store at 2-8°C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Diluted Wash Solution	May be stored for up to 1 month at 2-8°C.*
	Stop Solution	
	Diluent 1×	
	Conjugate	
	Unmixed Color Reagent A	
	Unmixed Color Reagent B	
	Standard	Aliquot and store for up to 1 month at -20°C in a manual defrost freezer.* Avoid repeated freeze-thaw cycles.
Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*	

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- ◆ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ◆ Pipettes and pipette tips.
- ◆ Deionized or distilled water.
- ◆ Squirt bottle, manifold dispenser, or automated microplate washer.
- ◆ 500 mL graduated cylinder.

D. PRECAUTION

The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples may require dilution with Diluent 1 \times .

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles.

B. SAMPLE PREPARATION

Serum samples require a 2-fold dilution. A suggested 2-fold dilution is 100 μL of sample + 100 μL of Diluent (1 \times).

C. REAGENT PREPARATION

Note: Bring all reagents to room temperature before use.

Wash Solution - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Solution Concentrate into deionized or distilled water to prepare 500 mL of Wash Buffer.


Substrate Solution - Color Reagent A and B should be mixed together in equal volumes within 15 minutes of use. Protect from light. 200 μL of the resultant mixture is required per well.

Diluent 1 \times - Add 20 mL of Calibrator Diluent 5 \times into 80 mL of deionized or distilled water to prepare 100 mL of Diluent 1 \times .

IL-8 Standard - Refer to the vial label for the reconstitution volume*. This reconstitution produces a stock solution of 2000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*if you have any question, please seek help from our Technical Support.

Pipette 500 μL of Diluent 1 \times into each tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 1000 pg/mL standard serves as the high standard. The Diluent 1 \times serves as the zero standard (0 pg/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Substrate Solution should remain colorless until added to the plate. Keep Substrate Solution protected from light. Substrate Solution should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the Substrate Solution. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the Substrate Solution.

VII. ASSAY PROCEDURE

Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.

***Note:** High concentrations of IL-8 are found in saliva. It is recommended that a face mask and gloves be used to protect kit reagents from contamination.*

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100 μL of Standard, sample, or control per well. Cover with the adhesive strip provided. Incubate for 2 hours at room temperature. A plate layout is provided for a record of standards and samples assayed.
4. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Solution (400 μL) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
5. Add 100 μL of IL-8 Conjugate to each well. Cover with a new adhesive strip. Incubate for 2 hours at room temperature.
6. Repeat the aspiration/wash as in step 4.
7. Add 200 μL of Substrate Solution to each well. Incubate for 30 minutes at room temperature. **Protect from light.**
8. Add 50 μL of Stop Solution to each well. The color in the wells should change from blue to yellow. If the color in the wells is green or if the color change does not appear uniform, gently tap the plate to ensure thorough mixing.
9. Determine the optical density of each well within 30 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.
10. **CALCULATION OF RESULTS :** Average the duplicate readings for each standard, control, and sample and subtract the average zero standard optical density. Create a standard curve by reducing the data using computer software capable of generating a four-parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the IL-8 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data. If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Rosenkilde, M.M. and T.W. Schwartz (2004) *APMIS* 112:481.
2. Colobran, R. et al. (2007) *Clin. Exp. Immunol.* 148:208.
3. Lazennec, G. and A. Richmond (2010) *Trends Mol. Med.* 16:134.
4. Schmid, J. and C. Weissmann (1987) *J. Immunol.* 139:250.
5. Matsushima, K. et al. (1988) *J. Exp. Med.* 167:1883.
6. Van den Steen, P.E. et al. (2000) *Blood* 96:2673.
7. Rajarathnam, K. et al. (2006) *Biochemistry* 45:7882.
8. Burrows, S.D. et al. (1994) *Biochemistry* 33:12741.
9. Nesmelova, I.V. et al. (2005) *J. Biol. Chem.* 280:4948.
10. Nesmelova, I.V. et al. (2008) *J. Biol. Chem.* 283:24155.
11. Goodman, K.B. et al. (1992) *Biochemistry* 31:10483.
12. Ishikawa, J. et al. (1993) *Gene* 131:305.
13. Jang, J. and B. Kim (1998) GenBank Accession #:AAK00048.
14. Mortier, A. et al. (2008) *Pharmacol. Ther.* 120:197.
15. Loos, T. et al. (2009) *Haematologica* 94:1346.
16. Proost, P. et al. (2008) *J. Exp. Med.* 205:2085.
17. Smedman, C. et al. (2009) *J. Immunol. Methods* 346:1.
18. Lim, C.P. et al. (2009) *J. Invest. Dermatol.* 129:851.
19. Moller, A. et al. (1993) *J. Immunol.* 151:3261.
20. Vanaudenaerde, B.M. et al. (2003) *J. Heart Lung Transplant* 22:1280.
21. Sandor, N. et al. (2009) *Mol. Immunol.* 47:438.
22. Brasier, A.R. et al. (1998) *J. Biol. Chem.* 273:3551.
23. Watchorn, T.M. et al. (2002) *Am. J. Physiol. Endocrinol. Metab.* 282:E763.
24. Lee, J. et al. (1992) *J. Biol. Chem.* 267:16283.
25. Nasser, M.W. et al. (2009) *J. Immunol.* 183:3425.
26. Stillie, R. et al. (2009) *J. Leukoc. Biol.* 86:529.
27. Monoz, L.M. et al. (2009) *J. Immunol.* 183:7337.
28. Waugh, D.J.J. and C. Wilson (2008) *Clin. Cancer Res.* 14:6735.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

1													
2													
3													
4													
5													
6													
7													
8													
9													
10													
11													
12													
	A	B	C	D	E	F	G	H					

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