

DESCRIPTION

Formula	$C_{17}H_{24}N_6O_{17}P_2$
Molecular Weight	646.35 Da
Formulation	1 mM provided in 20 mM Tris, pH 8.0
Storage & Stability	Store the unopened product at ≤ -20 °C. Good for 12 months from date of receipt.

APPLICATIONS

- For *in vitro* enzymatic incorporation of azido-sugars into specific, targeted glycans.
- Detecting the presence or absence of terminal GlcNAc residues (including O-GlcNAc).
- Detecting T and Tn antigens.
- Detecting the presence of high mannose glycan.

KEY FEATURES AND BENEFITS

- Can be introduced to proteins and lipids via various GlcNAc transferases.
- Can be conjugated to desired reporter molecules via click chemistry.
- Can be detected via Western Blot, ELISA, and flow cytometry, depending on the type of reporter molecule.
- Contains the smallest possible orthogonal functional group.
- Has minimal side effects on target molecules.
- User-friendly.

For Details: [Wu et al., \(2015\) Carbohydrate Res. 412:1-6.](#)

RELATED REAGENTS

Click Chemistry

- [Biotinylated Alkyne](#)
- [CMP-Azido-Sialic Acid](#)
- [UDP-Azido-GalNAc](#)
- [GDP-Azido-Fucose](#)

Enzymes & Detection Reagents

- GlcNAc specific transferases such as, [Recombinant Human MGAT1](#) , [B3GNT6](#) , [GCNT1](#) and [OGT](#)
- Hexaminodases such as, [Recombinant Human HEXA](#) and [HEXB](#) , [Recombinant B. thetaiotaomicron OGA](#)
- [Streptavidin-HRP](#)

SAMPLE PROTOCOL FOR GLcNAc DETECTION ON CORE-1 O-GLYCAN

Protocols are guidelines. Parameters need to be optimized by end users.

OTHER MATERIALS REQUIRED

- Assay Buffer: 25 mM Tris, 150 mM NaCl, 10 mM MnCl₂, pH 7.5
- Protein Sample
- [Recombinant Human GCNT1 \(R&D Systems[®], Catalog # 7248-GT\)](#)
- [Biotinylated Alkyne \(R&D Systems, Catalog # ES100\)](#)
- CuCl₂, 1 mM in deionized water
- Ascorbic Acid, 20 mM in deionized water
- SDS-PAGE and Western Blot reagents or equivalent
- TBST buffer: 25 mM Tris, 137 mM NaCl, 0.1% Tween[®]-20, pH 7.5
- [Streptavidin-HRP \(R&D Systems, Catalog # DY998\)](#)

ASSAY PROCEDURE

1. Prepare a reaction mixture by combining 5.0 µg of Protein Sample, 1 µg of Recombinant Human GCNT1 in the presence of 1 nmol of UDP-Azido-GlcNAc in the Assay Buffer with the final volume of 25 µL.
2. Prepare negative controls according to step 1 but omit Protein Sample or Recombinant Human GCNT1.
3. Incubate all the reactions and the controls at 37 °C for one hour.
4. To each of the samples, add 5 µL of 1 mM CuCl₂, 5.0 µL of 20 mM Ascorbic Acid, and 5 µL of 1 mM Biotinylated Alkyne. Mix with gentle tapping.
5. Incubate all samples at room temperature for 1 hour.
6. Separate the reactions and controls by SDS-PAGE.
7. Blot the gel to a nitrocellulose membrane.
8. Block the blot with 10% fat-free milk for 5 minutes.
9. Thoroughly wash the membrane with TBST buffer by changing buffer three times for a total of 45 minutes.
10. Incubate the blot with 25 ng/mL Streptavidin-HRP in 30 mL TBST buffer for 30 minutes.
11. Thoroughly wash the membrane with TBST buffer by changing buffer three times for a total of 45 minutes.
12. Detect with commercial ECL (Enhanced Chemiluminescence) reagents.

FINAL ASSAY CONDITIONS PER REACTION

- UDP-Azido-GlcNAc: 1 nmol
- Recombinant Human GCNT1: 1 µg
- Protein Sample: 5 µg
- Reaction volume: 25 µL

CLICK CHEMISTRY REACTION CONDITIONS PER REACTION

- CuCl₂: 5 nmol
- Ascorbic Acid: 100 nmol
- Biotinylated Alkyne: 5 nmol
- Reaction volume: 40 µL