

DESCRIPTION

Formula	C ₅₃ H ₆₇ N ₁₁ O ₂₅ P ₂ S ₃
Molecular Weight	1416.3 Da
Formulation	Lyophilized with Tris, at pH 8.0
Storage & Stability	Store the unopened product at ≤ -20 °C. Good for 12 months from date of receipt.

APPLICATIONS

- Fluorescent labeling with Cy5 of free glycans as well as glycoproteins and glycolipids.
- Fluorescent detection of specific glycan epitopes on glycoproteins as well as cell surface.
- Quantitation of the fucosylation level of specific glycans.
- Together with CMP-Cy3-Sialic Acid, allows for dual labeling and detection of sialoglycans.

KEY FEATURES & BENEFITS

- Excitation at 649 nm and emission at 671 nm, exhibits red fluorescent light under microscope.
- The fluorescent dye is conjugated to the C6 position of the fucose.
- Can be directly introduced into glycoproteins and glycolipids via various fucosyltransferases.
- Can be introduced to live cells for glycan imaging.
- Has minimum side-effects on target molecules.
- Very convenient and user-friendly.

For Details:

[Wu, ZL. *et al.*, \(2020\) Glycobiology 30:970.](#)

RELATED REAGENTS

Click Chemistry

- [CMP-Cy5-Sialic Acid \(ES302\)](#)
- [GDP-Azido-Fucose \(ES101\)](#)
- [CMP-Azido-Sialic Acid \(ES102\)](#)
- [UDP-Azido-GalNAc \(ES103\)](#)
- [UDP-Azido-GlcNAc \(ES104\)](#)
- [CMP-Biotin-Sialic Acid \(ES201\)](#)

Enzymes and Detection Reagents

- [Various fucosyltransferases](#)
- [Various neuraminidase](#)
- [Various fucosidase](#)

SAMPLE PROTOCOL

For direct fluorescent glycan labeling with GDP-Cy5-Fucose. Protocols are guidelines. Parameters need to be optimized by end users.

OTHER MATERIALS REQUIRED

- Assay Buffer: 25 mM Tris, 10 mM MnCl_2 , pH 7.5
- Sample protein
- Recombinant fucosyltransferases such as rhFUT9 ([R&D Systems®, Catalog # 9347-GT](#)) or rhFUT8 ([R&D Systems®, Catalog # 5768-GT](#))
- Protein sample loading dye
- SDS-PAGE and Western Blot reagents or equivalent
- Fluorescent Imager in a far-red fluorescent channel

FINAL ASSAY CONDITIONS PER REACTION

- Sample protein: 0.1 to 5 μg
- GDP-Cy5-Fucose: 0.2 nmol
- Fucosyltransferase: 0.5 μg

ASSAY PROCEDURE

1. Prepare a reaction mixture by combining 0.1 to 5 μg of a sample protein, 0.2 nmol GDP-Cy5-Fucose, 0.5 μg of a fucosyltransferases such as FUT9 or FUT8, add Assay Buffer to the final volume to 30 μL .
2. Prepare a negative control by repeating above but omitting the fucosyltransferases.
3. Incubate all the reactions and controls at 37 °C for 60 minutes.
4. Stop the reactions and controls by adding appropriate volume of protein sample loading dye to each reaction.
5. Separate the reactions and controls by SDS-PAGE.
6. Image the gel with a fluorescent imager in a far-red fluorescent channel.
7. Image the gel with trichloroethanol (TCE) imaging (if TCE is incorporated into the gel) or any other regular protein gel imaging method such as Coomassie® blue staining or silver staining.

Notes:

- FUT9 does not recognize sialylated lactosamine. To increase FUT9 labeling on samples that are highly sialylated, small amount of Recombinant C. perfringens Neuraminidase Protein, CF ([R&D Systems®, Catalog # 5080-NM](#)) may be added into the labeling reaction mixture to increase labeling.
- Recombinant Human Tissue alpha-L-Fucosidase/FUCA1, CF ([R&D Systems®, Catalog #7039-GH](#)) may be used to remove existing fucose residue on samples prior to the reaction to increase labeling. Since the fucosidase is active around pH 4.5 and requires Mg^{2+} , buffer should be conditioned to that of fucosyltransferases after the treatment of fucosidase.