

## PRODUCT DESCRIPTION

This kit contains four conjugated antibodies and four corresponding isotype controls that can be used for single-step staining of human mesenchymal stem/stromal cells (MSCs) (1-4).

## MATERIALS PROVIDED & STORAGE

Store the unopened kit at 2-8 °C **in the dark**. Refer to the kit label for date of expiration.

PART	PART #	DESCRIPTION
Positive Markers	965665	250 µL of CD105/Endoglin-PerCP Mouse IgG <sub>1</sub> ; Clone 166707
	965664	250 µL of CD146/MCAM-CFS Mouse IgG <sub>1</sub> ; Clone 128018
	965663	250 µL of CD90/Thy1-APC Mouse IgG <sub>2A</sub> ; Clone Thy-1A1
Negative Marker	965662	250 µL of CD45-PE Mouse IgG <sub>1</sub> ; Clone 2D1
Isotype Controls	965666	250 µL of Mouse IgG <sub>1</sub> -PE Isotype Control
	965667	250 µL of Mouse IgG <sub>2A</sub> -APC Isotype Control
	965668	250 µL of Mouse IgG <sub>1</sub> -CFS Isotype Control
	965669	250 µL of Mouse IgG <sub>1</sub> -PerCP Isotype Control
Staining Buffer	895027	100 mL of 1X Staining Buffer

## INTENDED USE

This product is designed for the flow cytometric analysis of human MSCs using four fluorochrome-conjugated antibodies.

## PRECAUTION


The Staining Buffer contains 0.09% sodium azide. Sodium azide may react with lead and copper plumbing to form explosive metallic azides. Flush with large volumes of water during disposal.

## SURFACE STAINING PROTOCOL

1. Cell samples should be washed with 2 mL of Staining Buffer, spinning the tube at 300 x g for 5 minutes.
2. Washed cells should be counted and then Fc receptor blocking reagents may be added if desired. If using excess pre-immune IgG to block Fc receptor, use 1 µg of IgG per 1 x 10<sup>5</sup> cells to be stained. The excess IgG does not need to be washed from the cells following the incubation period and can be carried into the staining reaction.
3. Transfer a small volume of cells (1 x 10<sup>6</sup> cells in 100 µL is recommended) into a 5 mL Flow Cytometry tube.
4. Add 10 µL of each antibody, or add 10 µL of each corresponding isotype control antibody to the cells.
5. Incubate the mixture for 30-45 minutes at room temperature **in the dark**.
6. Following the incubation, remove any excess antibody by washing the cells with 2 mL of Staining Buffer. The final cell pellet is resuspended in 200-400 µL of Staining Buffer for flow cytometric analysis.

**Note:** Using multiple fluorochromes requires proper flow cytometric compensation to remove the spillover fluorescence from a particular probe to a certain channel (5).

## DATA EXAMPLES



### Figure 1: Verification of Human Mesenchymal Stem/Stromal Cell Identity by Analysis of MSC Marker Expression.

Undifferentiated (**A**) and osteocyte-differentiated (**B**) human MSCs were stained with the indicated antibodies (filled histograms) or the corresponding isotype control (open histograms), as described in the procedure. Multi-color flow cytometry simultaneously detected cells that were positive for MSC markers (CD105/Endoglin, CD146/MCAM, and CD90/Thy1). Following 21 days of osteocyte-differentiation, the cells showed a characteristic reduction in CD105/Endoglin, CD146/MCAM, and CD90/Thy1 staining. CD45 is a negative control for both cell types.

## REFERENCES

1. Chamberlain, G. *et al.* (2007) *Stem Cells* **25**:2739.
2. Abdallah, B.M. and M. Kassem (2008) *Gene Therapy* **15**:109.
3. Spitkovsky, D. and J. Hescheler (2008) *Min. Inv. Ther.* **17**:79.
4. Delorme, B. *et al.* (2008) *Blood* **111**:2631.
5. Bagwell, B. and E.G. Adams (1993) *Ann. N.Y. Acad. Sci.* **677**:167.