



Mouse Th17 Cell Multi-Color Flow Cytometry Kit

Catalog Number: FMC008

Size: 25 Tests

Product Description

This kit contains four conjugated antibodies (and corresponding isotype controls) that can be used for single-step staining of mouse Th17 cells (1 - 10):

- CCR6-Alexa Fluor[®] 488 (Clone 140706; rat IgG_{2A})
- IL-22-PE (Clone 140301; rat IgG_{2A})
- IL-17-PerCP (goat IgG)
- CD4-APC (Clone GK1.5; rat IgG_{2B})

The kit also contains Fixation/Permeabilization Buffer (30 mL), which contains 1% formaldehyde, saponin, and < 0.05% sodium azide as well as Permeabilization/Wash Buffer (60 mL), which contains saponin and 0.05% sodium azide.

Intended Use

This product is designed for the flow cytometric analysis of mouse Th17 cells using four fluorochrome-conjugated antibodies.

Storage

Store at 2 - 8° C in the dark. Use within 6 months of receipt.

Precautions

- Formaldehyde is a suspected carcinogen. Avoid contact with skin, eyes, and mucous membranes, and avoid inhaling fumes. In case of contact, wash immediately with water and seek medical advice.
- Sodium azide may react with lead and copper plumbing to form explosive metallic azides. Flush with large volumes of water during disposal.

Intracellular Staining Protocol with Simultaneous Fixation/Permeabilization

1. Harvest cells of interest and wash twice in PBS or Hanks' Balanced Salt Solution (HBSS).
2. Approximately 5×10^5 washed cells should be resuspended in 0.5 mL of Fixation/Permeabilization Buffer and incubated at 2 - 8° C for 30 minutes. Cells should be vortexed intermittently in order to maintain a single cell suspension.
3. Centrifuge the cells and resuspend the pellet in 100 - 200 μ L of the Permeabilization/Wash Buffer.
4. Add 10 μ L of each antibody or each corresponding isotype control antibody to the cells.
5. Incubate the mixture for 30 - 45 minutes at 2 - 8° C in the dark.
6. Following the incubation, remove any excess antibody by washing the cells in 2 mL of Permeabilization/Wash Buffer. Resuspend the final cell pellet in 200 - 400 μ L of PBS for flow cytometric analysis.

Notes: *Because saponin-mediated cell permeabilization is a reversible process, it is important to keep the cells in the presence of saponin during intracellular staining. Using multiple fluorochromes requires proper flow cytometric compensation to remove the spillover fluorescence from a particular probe to a certain channel (11).*

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

**R&D Systems, Inc.
1-800-343-7475**

Typical Data




Figure 1: Mouse splenocytes were unstimulated (open histograms) or Th17-stimulated for 5 days (filled histograms), using 5 μ g/mL of anti-CD3 (R&D Systems, Catalog # MAB484), 1 μ g/mL of anti-CD28 (R&D Systems, Catalog # MAB4831), 10 ng/mL of recombinant mouse IL-2 (Cys160Ser) (R&D Systems, Catalog # 1150-ML), 5 ng/mL of human TGF- β (R&D Systems, Catalog # 100-B), 10 ng/mL of recombinant mouse IL-23 (R&D Systems, Catalog # 1887-ML), and 20 ng/mL of recombinant mouse IL-6 (R&D Systems, Catalog # 406-ML). Cells were harvested and stained with the antibodies in this kit or isotype controls following the procedure outlined in this insert and gated on CD4⁺ cells. Histograms show the detected levels of the indicated antibodies in both cell samples.

References

1. Annunziato, F. *et al.* (2009) *Nat. Rev. Rheumatol.* **5**:325.
2. Zhang, Z. *et al.* (2009) *J. Clin. Invest.* **119**:1899.
3. Ghilardi, N. and W. Ouyang (2007) *Seminars in Immunol.* **19**:383.
4. Korn, T. *et al.* (2009) *Annu. Rev. Immunol.* **27**:485.
5. Sundrud, M.S. *et al.* (2009) *Science* **324**:1334.
6. Weaver, C.T. *et al.* (2007) *Ann. Rev. Immunol.* **25**:821.
7. Awasthi, A. and V.K. Kuchroo (2009) *Nature Immunol.* **10**:568.
8. Iwakura, Y. *et al.* (2008) *Immunol. Rev.* **226**:57.
9. Hirota, K. *et al.* (2007) *J. Exp. Med.* **204**:2803.
10. Yamazaki, T. *et al.* (2008) *J. Immunol.* **181**:8391.
11. Bagwell, B. and E.G. Adams (1993) *Ann. N.Y. Acad. Sci.* **677**:167.

Trademarks and registered trademarks are property of their respective owners.

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose, Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300, Fax: (541) 335-0354.