

MATERIAL DATA SHEET

Di-Ub/Ub2 Non-Hydrolyzable Chains (K63-linked), Agarose
Cat. # UCN-302

Linkage specific, non-hydrolyzable di-ubiquitin is resistant to the activity of enzymes (DUB's) that cleave the isopeptide linkage between adjacent ubiquitin molecules. It can be used to investigate binding interactions between di-ubiquitin and proteins that contain elements such as ubiquitin-associated domains (UBAs) or ubiquitin-interacting motifs (UIMs). This product may also be useful in exploring the role of unanchored poly-ubiquitin chains in some signaling pathways.

Product Information

Quantity:	250 µl resin supplied as a 50% slurry in a solution of 20% ethanol.
Stock:	Non-hydrolyzable Ub2 (K63) chains coupled to beads at 1 mg/ml (58 nmol/ml)

Use and Storage

Use:	Prepare resin by rinsing with 10 volumes of water to remove ethanol storage buffer. Equilibrate resin by washing with 10 volumes of desired start buffer. Binding and elution of material is dependent on individual experimental conditions and requirements.
Storage:	Polyubiquitin-agarose can be re-used if properly maintained. After use, clean resin with a wash cycle of 5 volumes 100 mM HEPES pH 8.0, 500 mM NaCl followed by 5 volumes 100 mM NaOAc pH 4.5, 500 mM NaCl. Repeat twice, then rinse resin with a low salt buffer. Store resin at 4°C in neutral aqueous buffer containing 1 mM NaN ₃ or 20 % ethanol as a preservative. DO NOT FREEZE.

Literature

References:	Zeng W., <i>et al.</i> (2010) <i>Cell</i> 141 : 315-330 Komander D. (2009) <i>Biochem Soc Trans</i> 37 : 937-953 Kulathu Y., <i>et al.</i> (2009) <i>Nat. Struct Mol Biol.</i> 16 : 1328-1330
--------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

For Laboratory Research Use Only, Not For Use in Human

Rev: 3/5/2019

840 Memorial Drive, Cambridge, MA 02139 Phone: 617-576-2210 FAX: 617-492-3565
www.bostonbiochem.com

The contents of this datasheet (unless otherwise noted) are Copyright © 2008 Boston Biochem, Inc. All rights reserved. Duplication in whole or in part is strictly prohibited without the express written consent of Boston Biochem, Inc. "Boston Biochem" is a Trademark of Boston Biochem, Inc., 840 Memorial Drive, Cambridge, MA 02139.