PRODUCT INFORMATION & MANUAL

Human IL-6 Valukine™ ELISA

VAL102

For the quantitative determination of natural and recombinant
human Interleukin 6 (IL-6) concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version202006.4

TABLE OF CONTENTS

| | |
|------------------------------------|----|
| I. BACKGROUND..... | 2 |
| II. OVERVIEW..... | 4 |
| III. ADVANTAGES..... | 5 |
| IV. EXPERIMENT..... | 8 |
| V. KIT COMPONENTS AND STORAGE..... | 9 |
| VI. PREPARATION..... | 12 |
| VII. ASSAY PROCEDURE..... | 14 |
| VIII. REFERENCES..... | 15 |

I. BACKGROUND

Interleukin 6 (IL-6) is a pleiotropic α -helical 22-28 kDa phosphorylated and variably glycosylated cytokine that plays important roles in the acute phase reaction, inflammation, hematopoiesis, bone metabolism, and cancer progression (1-5). Mature human IL-6 is 183 amino acids (aa) in length and shares 41% aa sequence identity with mouse and rat IL-6 (6). Alternate splicing generates several isoforms with internal deletions, some of which exhibit antagonistic properties (7-10). Cells known to express IL-6 include CD8⁺ T cells, fibroblasts, synoviocytes, adipocytes, osteoblasts, megakaryocytes, endothelial cells (under the influence of endothelins), sympathetic neurons, cerebral cortex neurons, adrenal medulla chromaffin cells, retinal pigment cells, mast cells, keratinocytes, Langerhans cells, fetal and adult astrocytes, neutrophils, monocytes, eosinophils, colonic epithelial cells, B1 B cells, and pancreatic islet beta cells (2, 7, 10-33). IL-6 production is generally correlated with cell activation and is normally kept in control by glucocorticoids, catecholamines, and secondary sex steroids (2). Normal human circulating IL-6 is in the 1 pg/mL range, with slight elevations during the menstrual cycle, modest elevations in certain cancers, and large elevations after surgery (34-38).

IL-6 induces signaling through a cell surface heterodimeric receptor complex composed of a ligand binding subunit (IL-6 R) and a signal transducing subunit (gp130). IL-6 binds to IL-6 R, triggering IL-6 R association with gp130 and gp130 dimerization (39). Gp130 is also a component of the receptors for CLC, CNTF, CT-1, IL-11, IL-27, LIF, and OSM (40). Soluble forms of IL-6 R are generated by both alternative splicing and proteolytic cleavage (3). In a mechanism known as trans-signaling, complexes of soluble IL-6 and IL-6 R elicit responses from gp130-expressing cells that lack cell surface IL-6 R (1, 3). Trans-signaling enables a wider range of cell types to respond to IL-6, as the expression of gp130 is ubiquitous, while that of IL-6 R is predominantly restricted to hepatocytes, monocytes, and resting lymphocytes (1-3). Soluble splice forms of gp130 block trans-signaling from IL-6/IL-6 R but not from other cytokines that use gp130 as a co-receptor (3, 41).

IL-6, along with TNF- α and IL-1, drives the acute inflammatory response, is almost solely responsible for fever and the acute phase response in the liver, and is important in the transition from acute inflammation to either acquired immunity, or chronic inflammatory disease (1-4). It contributes to chronic inflammation in conditions such as obesity, insulin resistance, inflammatory bowel disease, inflammatory arthritis and sepsis when dysregulated, often involving IL-6 trans-signaling (1, 2). It also plays an important role in the differentiation of naive T cells to Th17 inflammatory cells in the presence of TGF- β . IL-6 modulates bone resorption and is a major effector of inflammatory joint destruction in rheumatoid arthritis through its promotion of Th17 T cell activity (1). It contributes to atherosclerotic plaque development and destabilization (2). However, IL-6 can also have anti-inflammatory effects, such as in skeletal muscle

where it is secreted in response to exercise (2). It promotes hematopoiesis by being a growth factor for hematopoietic stem cells, induces B cell maturation to plasma cells and perpetuates multiple myeloma (1, 42). IL-6 also promotes, but probably does not initiate, other types of inflammation-associated carcinogenesis, such as colitis-associated cancer (1).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. A monoclonal antibody specific for IL-6 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any IL-6 present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked polyclonal antibody specific for IL-6 is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of IL-6 bound in the initial step. The color development is stopped and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- ◆ **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- ◆ This kit is suitable for cell culture supernate, serum and plasma.
- ◆ The kit should not be used beyond the expiration date on the kit label.
- ◆ Do not mix or substitute reagents with those from other lots or sources.
- ◆ If samples generate values higher than the highest standard, dilute the samples with Diluent and repeat the assay.
- ◆ Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

| | Intra-assay Precision | | | Inter-assay Precision | | |
|--------------------|-----------------------|------|------|-----------------------|------|------|
| Sample | 1 | 2 | 3 | 1 | 2 | 3 |
| Mean (pg/mL) | 20.6 | 77.8 | 175 | 23.9 | 83.3 | 177 |
| Standard Deviation | 1.20 | 3.51 | 7.33 | 4.20 | 13.3 | 25.5 |
| CV% | 5.8 | 4.5 | 4.2 | 17.6 | 15.9 | 14.4 |

B. RECOVERY

The recovery of human IL-6 spiked to levels throughout the range of the assay in various matrices was evaluated.

| Sample Type | Average % Recovery | Range |
|--------------------------|--------------------|-----------|
| Cell culture media (n=4) | 95 | 81 - 104% |
| Serum (n=3) | 93 | 80 - 99% |
| Plasma (n=4) | 96 | 81 - 109% |

C. SENSITIVITY

The minimum detectable dose (MDD) of IL-6 is typically less than 1.56 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against highly purified *E. coli*-expressed recombinant human IL-6 produced at R&D Systems. The NIBSC/WHO 1st International Standard for IL-6 (89/548), which was intended as a potency standard, was evaluated in this kit. The NIBSC/WHO standard is a CHO cell-derived recombinant human IL-6.

The dose response curve of the International Standard (89/548) parallels the Valukine standard curve. To convert sample values obtained with the Valukine Human IL-6 kit to approximate NIBSC 89/548 units, use the equation below.

NIBSC (89/548) approximate value (IU/mL) = 0.109 × Valukine Human IL-6 value (pg/mL)

E. LINEARITY

To assess the linearity of the assay, samples were spiked with high concentrations of human IL-6 in various matrices and diluted with Diluent 1× to produce samples with values within the dynamic range of the assay.

| Dilution | | Cell culture media (n=4) | Serum (n=3) | Plasma (n=4) |
|----------|-----------------------|--------------------------|-------------|--------------|
| 1:2 | Average % of Expected | 112 | 102 | 99 |
| | Range (%) | 105 - 117 | 101 - 102 | 89 - 106 |
| 1:4 | Average % of Expected | 111 | 106 | 101 |
| | Range (%) | 104 - 119 | 102 - 111 | 95 - 107 |
| 1:8 | Average % of Expected | 97 | 108 | 99 |
| | Range (%) | 91 - 105 | 103 - 116 | 93 - 104 |
| 1:16 | Average % of Expected | 89 | 109 | 98 |
| | Range (%) | 81 - 98 | 102 - 117 | 90 - 107 |

F. SAMPLE VALUES

Cell Culture Supernates - Human peripheral blood mononuclear cells (1×10^6 cells/mL) were cultured in RPMI supplemented with 10% fetal calf serum, 50 μ M β -mercaptoethanol, 2 mM L-glutamine, 100 U/mL penicillin, and 100 μ g/mL streptomycin sulfate and stimulated for 3 days with 10 μ g/mL PHA. An aliquot of the cell culture supernate was removed, assayed for levels of natural IL-6, and measured 6640 pg/mL.

Serum - Three human serum samples were evaluated for the presence of human IL-6 in this assay. All samples measured ranged from 20.5 to 62.5 pg/mL with an average of 48.0 pg/mL.

Plasma - Four human plasma samples were evaluated for the presence of human IL-6 in this assay. All samples measured ranged from 73.5 to 105 pg/mL with an average of 88.6 pg/mL.

G. SPECIFICITY


This assay recognizes both natural and recombinant human IL-6. The following factors were prepared at 50 ng/mL and assayed for cross-reactivity. Preparations of the following factors at 50 ng/mL in a mid-range rhIL-6 control were assayed for interference. No significant cross-reactivity or interference was observed.

| Recombinant human | Recombinant mouse |
|-------------------|-------------------|
| sgp130 | IL-6 |
| IL-6 sR | |
| IL-6 sR/sgp130 | |

IV. EXPERIMENT

EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

| pg/mL | OD | Average | Corrected |
|-------|----------------|---------|-----------|
| 0 | 0.024 0.022 | 0.023 | - |
| 3.12 | 0.067 0.065 | 0.066 | 0.043 |
| 6.25 | 0.108 0.105 | 0.107 | 0.084 |
| 12.5 | 0.190 0.180 | 0.185 | 0.162 |
| 25 | 0.355 0.343 | 0.349 | 0.326 |
| 50 | 0.714 0.682 | 0.698 | 0.675 |
| 100 | 1.110 1.141 | 1.126 | 1.103 |
| 300 | 2.629 2.808 | 2.719 | 2.696 |

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

| Parts | Description | Size |
|-------------------------------|--|----------|
| Human IL-6 Microplate | 96 well polystyrene microplate (12 strips of 8 wells) coated with a mouse monoclonal antibody against human IL-6 | 1 plate |
| Human IL-6 Conjugate | Solution of polyclonal antibody against human IL-6 conjugated to horseradish peroxidase | 1 vial |
| Human IL-6 Standard | recombinant human IL-6 in a buffered protein base; lyophilized | 1 vial |
| Calibrator Diluent (5×) | a 5× concentrated buffered protein base | 1 vial |
| Wash Buffer Concentrate (25×) | a 25× concentrated solution of buffered surfactant | 1 vial |
| Color Reagent A | stabilized hydrogen peroxide | 1 vial |
| Color Reagent B | stabilized chromogen (tetramethylbenzidine) | 1 vial |
| Stop Solution | 2 N sulfuric acid | 1 vial |
| Plate Sealers | adhesive strip | 3 strips |

B. STORAGE

| | | |
|--------------------------------------|---|---|
| Unopened Kit | Store at 2-8°C. Do not use past kit expiration date. | |
| Opened/ Reconstituted Reagents | Diluted Wash Buffer | May be stored for up to 1 month at 2-8°C.* |
| | Stop Solution | |
| | Diluent 1× | |
| | Conjugate | |
| | Unmixed Color Reagent A | |
| | Unmixed Color Reagent B | |
| | Standard | Aliquot and store for up to 1 month at -20°C in a manual defrost freezer. * Avoid repeated freeze-thaw cycles. |
| Microplate Wells | Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.* | |

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- ◆ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ◆ Pipettes and pipette tips.
- ◆ Deionized or distilled water.
- ◆ Squirt bottle, manifold dispenser, or automated microplate washer.
- ◆ 500 mL graduated cylinder.

D. PRECAUTION

The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent 1 \times .

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles.

Plasma - Collect plasma using EDTA, heparin, or citrate as an anticoagulant. Centrifuge for 15 minutes at 1000 x g within 30 minutes of collection. Assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles.

B. SAMPLE PREPARATION

Serum samples require a 5-fold dilution. A suggested 5-fold dilution is 40 μL of sample + 160 μL of Diluent (1 \times).

Plasma samples require a 2-fold dilution. A suggested 2-fold dilution is 100 μL of sample + 100 μL of Diluent (1 \times).

C. REAGENT PREPARATION

Note: Bring all reagents to room temperature before use.

Wash Buffer - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25 \times) into deionized or distilled water to prepare 500 mL of Wash Buffer.


Substrate Solution - Color Reagent A and B should be mixed together in equal volumes within 15 minutes of use. Protect from light. 200 μL of the resultant mixture is required per well.

Diluent 1 \times - Add 20 mL of Calibrator Diluent Concentrate 5 \times into 80 mL of deionized or distilled water to prepare 100 mL of Diluent 1 \times .

IL-6 Standard - Refer to the vial label for reconstitution volume*. This reconstitution produces a stock solution of 300 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*if you have any question, please seek help from our Technical Support.

Pipette 667 μL of Diluent 1 \times into the 100 pg/mL tube. Pipette 500 μL of Diluent 1 \times into each remaining tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The undiluted standard serves as the high standard (300 pg/mL). The Diluent 1 \times serves as the zero standard (0 pg/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Substrate Solution should remain colorless until added to the plate. Keep Substrate Solution protected from light. Substrate Solution should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the Substrate Solution. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the Substrate Solution.

VII. ASSAY PROCEDURE

Note: Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100 μ L of Standard, sample, or control per well. Cover with the adhesive strip provided. Incubate for 2 hours at room temperature. A plate layout is provided for a record of standards and samples assayed.
4. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Buffer (400 μ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
5. Add 200 μ L of human IL-6 Conjugate to each well. Cover with a new adhesive strip. Incubate for 2 hours at room temperature.
6. Repeat the aspiration/wash as in step 4.
7. Add 200 μ L of Substrate Solution to each well. Incubate for 20 minutes at room temperature. **Protect from light.**
8. Add 50 μ L of Stop Solution to each well. The color in the wells should change from blue to yellow. If the color in the wells is green or if the color change does not appear uniform, gently tap the plate to ensure thorough mixing.
9. Determine the optical density of each well within 30 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.
10. **CALCULATION OF RESULTS:** Average the duplicate readings for each standard, control, and sample and subtract the average zero standard optical density. Create a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the IL-6 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data. If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Naugler, W.E. and M. Karin (2008) *Trends Mol. Med.* 14:109.
2. Schuett, H. et al. (2009) *Thromb. Haemost.* 102:215.
3. Jones, S.A. (2005) *J. Immunol.* 175:3468.
4. Hodge, D.R. et al. (2005) *Eur. J. Cancer* 41:2502.
5. Rose -John, S. et al. (2006) *J. Leukoc. Biol.* 80:227.
6. Van Snick, J. et al. (1988) *Eur. J. Immunol.* 18:193.
7. Kestler, D.P. et al. (1995) *Blood* 86:4559.
8. Kestler, D.P. et al. (1999) *Am. J. Hematol.* 61:169.
9. Bihl, M.P. et al. (2002) *Am. J. Respir. Cell Mol. Biol.* 27:48.
10. Alberti, L. et al. (2005) *Cancer Res.* 65:2.
11. May, L.T. et al. (1986) *Proc. Natl. Acad. Sci. USA* 83:8957.
12. Sad, S. et al. (1995) *Immunity* 2:271.
13. Cichy, J. et al. (1996) *Biochem. Biophys. Res. Commun.* 227:318.
14. Miyazawa, K. et al. (1998) *Am. J. Pathol.* 152:793.
15. Fried, S.K. et al. (1998) *Endocrinology* 83:847.
16. Ishimi, Y. et al. (1990) *J. Immunol.* 145:3297.
17. Jiang, S. et al. (1994) *Blood* 84:4151.
18. Xin, X. et al. (1995) *Endocrinology* 136:132.
19. Marz, P. et al. (1998) *Proc. Natl. Acad. Sci. USA* 95:3251.
20. Ringheim, G.E. et al. (1995) *J. Neuroimmunol.* 63:113.
21. Gadiant, R.A. et al. (1995) *Neurosci. Lett.* 194:17.
22. Kuppner, M.C. et al. (1995) *Immunology* 84:265.
23. Gagari, E. et al. (1997) *Blood* 89:2654.
24. Cumberbatch, M. et al. (1996) *Immunology* 87:513.
25. Fujisawa, H. et al. (1997) *J. Interferon Cytokine Res.* 17:347.
26. Lee, S.C. et al. (1993) *J. Immunol.* 150:2659.
27. Lafortune, L. et al. (1996) *J. Neuropathol. Exp. Neurol.* 55:515.
28. Ericson, S.G. et al. (1998) *Blood* 91:2099.
29. Melani, C. et al. (1993) *Blood* 81:2744.
30. Lacy, P. et al. (1998) *Blood* 91:2508.
31. Jung, H.C. et al. (1995) *J. Clin. Invest.* 95:55.
32. Spencer, N.F.L. and R.A. Daynes (1997) *Int. Immunol.* 9:745.

33. Campbell, I.L. et al. (1989) *J. Immunol.* 143:1188.
34. D'Auria, L. et al. (1997) *Eur. Cytokine Netw.* 8:383.
35. Yamamura, M. et al. (1998) *Br. J. Haematol.* 100:129.
36. Angstwurm, M.W.A. et al. (1997) *Cytokine* 9:370.
37. Mouawad, R. et al. (1996) *Clin. Cancer Res.* 2:1405.
38. Sakamoto, K. et al. (1994) *Cytokine* 6:181.
39. Murakami, M. et al. (1993) *Science* 260:1808.
40. Muller-Newen, G. (2003) *Sci. STKE* 2003:PE40.
41. Mitsuyama, K. et al. (2006) *Clin. Exp. Immunol.* 143:125.
42. Cerutti, A. et al. (1998) *J. Immunol.* 160:2145.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

| | | | | | | | | | | | | | | | | | |
|-----------|----------|----------|----------|----------|----------|----------|----------|----------|--|--|--|--|--|--|--|--|--|
| 1 | | | | | | | | | | | | | | | | | |
| 2 | | | | | | | | | | | | | | | | | |
| 3 | | | | | | | | | | | | | | | | | |
| 4 | | | | | | | | | | | | | | | | | |
| 5 | | | | | | | | | | | | | | | | | |
| 6 | | | | | | | | | | | | | | | | | |
| 7 | | | | | | | | | | | | | | | | | |
| 8 | | | | | | | | | | | | | | | | | |
| 9 | | | | | | | | | | | | | | | | | |
| 10 | | | | | | | | | | | | | | | | | |
| 11 | | | | | | | | | | | | | | | | | |
| 12 | | | | | | | | | | | | | | | | | |
| | A | B | C | D | E | F | G | H | | | | | | | | | |

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