PRODUCT INFORMATION & MANUAL

Human IFN- γ Valukine™ ELISA

VAL104

For the quantitative determination of natural and
recombinant human IFN- γ concentrations

For research use only.

Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001

info.cn@bio-techne.com

Novus kits are guaranteed for 3 months from date of receipt

Version 201911.3

TABLE OF CONTENTS

I. BACKGROUND.....	2
II. OVERVIEW.....	3
III. ADVANTAGES.....	4
IV. EXPERIMENT.....	6
V. KIT COMPONENTS AND STORAGE.....	7
VI. PREPARATION.....	9
VII. ASSAY PROCEDURE.....	11
VIII. REFERENCES.....	12

I. BACKGROUND

Interferon-gamma (IFN- γ , also known as type II interferon) is an important immunoregulatory cytokine that was originally identified through its anti-viral activity (1, 2). It plays key roles in host defense by exerting anti-viral, anti-proliferative, and immunoregulatory activities (3, 4). On many cell types, IFN- γ induces the production of cytokines and upregulates the expression of various membrane proteins including class I and II MHC antigens, Fc receptors, leukocyte adhesion molecules, and B7 family antigens. IFN- γ is a potent activator of macrophage effector functions. It directs the synthesis, class switching, and secretion of immunoglobulins by B cells. IFN- γ also influences T-helper cell phenotype development by inhibiting Th2 differentiation and stimulating Th1 development (3, 4). IFN- γ plays a central role in the progression of inflammatory diseases such as autoimmunity and atherosclerosis (5, 6).

Biologically active IFN- γ consists of a noncovalently linked homodimer of 20-25 kDa variably glycosylated subunits (7). Mature human IFN- γ shares 90% amino acid (aa) sequence identity with rhesus macaque IFN- γ , 59%-64% with bovine, canine, equine, feline and porcine IFN- γ , and 37%-43% aa identity with cotton rat, mouse, and rat IFN- γ . IFN- γ dimers bind to transmembrane IFN- γ RI (alpha subunits) which then interact with transmembrane IFN- γ RII (beta subunits) to form a functional receptor complex of two and two subunits (8, 9). Inclusion of IFN- γ RII in the receptor complex increases the ligand binding affinity as well as the efficiency of signal transduction (9, 10). Whereas the β -chain is expressed constitutively on many cell types, the cellular regulation of the β -chain correlates with an IFN- γ responsive state and is tightly regulated (8).

IFN- γ is produced by a number of cell types under inflammatory conditions, including dendritic epidermal/ $\gamma\delta$ T cells (11), keratinocytes (12), peripheral blood $\gamma\delta$ T cells (13), mast cells (14), neurons (15), CD8⁺ T cells (16), macrophages (17), B cells (18), neutrophils (19), NK cells (20), CD4⁺ T cells (21), and testicular spermatids (22).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. A polyclonal antibody specific for IFN- γ has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any IFN- γ present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked polyclonal antibody specific for IFN- γ is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of IFN- γ bound in the initial step. The color development is stopped and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- ◆ **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- ◆ This kit is suitable for cell culture supernate and serum.
- ◆ The kit should not be used beyond the expiration date on the kit label.
- ◆ Do not mix or substitute reagents with those from other lots or sources.
- ◆ If samples generate values higher than the highest standard, dilute the samples with Diluent and repeat the assay.
- ◆ Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

	Intra-assay Precision			Inter-assay Precision		
Sample	1	2	3	1	2	3
Mean (pg/mL)	47.8	148	321	61.5	164	342
Standard Deviation	3.67	7.57	14.2	5.55	10.4	28.4
CV%	7.7	5.1	4.4	9.0	6.4	8.3

B. RECOVERY

The recovery of human IFN- γ spiked to different levels throughout the range of the assay in cell culture media was evaluated. The recovery ranged from 90-116% with an average of 103%.

The recovery of human IFN- γ spiked to different levels throughout the range of the assay in human serum was evaluated. The recovery ranged from 88.8-100.5% with an average of 96.3%.

C. SENSITIVITY

The minimum detectable dose (MDD) of IFN- γ is typically less than 7.8 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against highly purified *E. coli*-expressed recombinant human IFN- γ produced at R&D Systems. The NIBSC 1st British Standard for human leukocyte IFN- γ (82/587), which was intended as a potency standard, was evaluated in this kit.

The dose response curve of the NIBSC Standard (82/587) parallels the Valukine standard curve. To convert sample values obtained with the Valukine Human IFN- γ kit to approximate NIBSC 82/587 units, use the equation below.

NIBSC (82/587) approximate value (IU/mL) = 0.010 × Valukine Human IFN- γ value (pg/mL)

E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with high concentrations of IFN- γ and diluted with Diluent 1 \times to produce samples with values within the dynamic range of the assay.

Dilution	Average % of Expected	Range (%)
1:2	107	96 – 109
1:4	111	107 – 115
1:8	111	104 – 121
1:16	106	91 – 114

F. SAMPLE VALUES

Cell Culture Supernates -Human peripheral blood mononuclear cells (1×10^6 cells/mL) were cultured in RPMI supplemented with 10% fetal calf serum, 50 μ M β -mercaptoethanol, 2 mM L-glutamine, 100 U/mL penicillin, and 100 μ g/mL streptomycin sulfate. Cells were cultured unstimulated or stimulated with 10 μ g/mL PHA. Aliquots of the cell culture supernates was removed on days 1 and 5 and assayed for levels of natural IFN- γ .

Condition	Day 1 (pg/mL)	Day 5 (pg/mL)
Unstimulated	453	2651
Stimulated	8360	7851

Serum - Four serum samples were evaluated for the presence of IFN- γ in this assay. All samples measured below the lowest standard, 15.6pg/mL.

G. SPECIFICITY

This assay recognizes both natural and recombinant human IFN- γ . The following factors were prepared at 50 ng/mL and assayed for cross-reactivity. Preparations of the following factors at 50 ng/mL in a mid-range rhIFN- γ control were assayed for interference. No significant cross-reactivity or interference was observed.


Recombinant human	Recombinant mouse	Recombinant rat
IL-1 α	IFN- γ	IFN- γ
IL-1 β	IL-1 α	
IL-1 α	IL-1 β	
IL-6	IL-6	
SCF		

Recombinant rhesus macaque IFN- γ was found to cross-react approximately 20%.

IV. EXPERIMENT

EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

pg/mL	OD	Average	Corrected
0	0.085 0.087 0.161	0.086	—
15.6	0.170 0.228	0.166	0.080
31.3	0.232 0.343	0.230	0.144
62.5	0.375 0.574	0.359	0.273
125	0.624 1.013	0.599	0.513
250	1.069 1.762	1.041	0.955
500	1.914 2.969	1.838	1.752
1000	3.000	2.985	2.899

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
IFN- γ Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with a mouse polyclonal antibody against human IFN- γ	1 plate
IFN- γ Conjugate	a solution of polyclonal antibody against IFN- γ conjugated to horseradish peroxidase	1 vial
IFN- γ Standard	recombinant human IFN- γ in a buffered protein base; lyophilized	1 vial
Calibrator Diluent(5 \times)	a 5 \times concentrated buffered protein base	1 vial
Wash Buffer Concentrate(25 \times)	a 25 \times concentrated solution of buffered surfactant	1 vial
Color Reagent A	stabilized hydrogen peroxide	1 vial
Color Reagent B	stabilized chromogen (tetramethylbenzidine)	1 vial
Stop Solution	2 N sulfuric acid	1 vial
Plate Covers	adhesive strip	3 strips

B. STORAGE

Unopened Kit	Store at 2-8°C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Diluted Wash Solution	May be stored for up to 1 month at 2-8°C.*
	Stop Solution	
	Diluent 1×	
	Conjugate	
	Unmixed Color Reagent A	
	Unmixed Color Reagent B	
	Standard	Aliquot and store for up to 1 month at <-20°C. in a manual defrost freezer. * Avoid repeated freeze-thaw cycles.
Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*	

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- ◆ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ◆ Pipettes and pipette tips.
- ◆ Deionized or distilled water.
- ◆ Squirt bottle, manifold dispenser, or automated microplate washer.
- ◆ 500 mL graduated cylinder.

D. PRECAUTION

The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples may require dilution with Diluent 1 \times .

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles.

B. SAMPLE PREPARATION

Serum samples require a 2-fold dilution. A suggested 2-fold dilution is 100 μL of sample + 100 μL of Diluent (1 \times).

C. REAGENT PREPARATION

Note: Bring all reagents to room temperature before use.

Wash Solution - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25 \times) into deionized or distilled water to prepare 500 mL of Wash Buffer.


Substrate Solution - Color Reagent A and B should be mixed together in equal volumes within 15 minutes of use. Protect from light. 200 μL of the resultant mixture is required per well.

Diluent 1 \times - Add 20 mL of Calibrator Diluent 5 \times into 80 mL of deionized or distilled water to prepare 100 mL of Diluent 1 \times .

IFN- γ Standard - Refer to the vial label for reconstitution volume for reconstitution*. This reconstitution produces a stock solution of 1000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*if you have any question, please seek help from our Technical Support.

Use polypropylene tubes. Pipette 500 μL of Diluent 1 \times into each tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 1000 pg/mL standard serves as the high standard. The Diluent 1 \times serves as the zero standard (0pg/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Substrate Solution should remain colorless until added to the plate. Keep Substrate Solution protected from light. Substrate Solution should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the Substrate Solution. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the Substrate Solution.

VII. ASSAY PROCEDURE

Note: Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100 μL of Standard, sample, or control per well. Cover with the adhesive strip provided. Incubate for 2 hours at room temperature. A plate layout is provided for a record of standards and samples assayed.
4. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Solution (400 μL) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
5. Add 200 μL of IFN- γ Conjugate to each well. Cover with a new adhesive strip. Incubate for 2 hours at room temperature.
6. Repeat the aspiration/wash as in step 4.
7. Add 200 μL of Substrate Solution to each well. Incubate for 30 minutes at room temperature. **Protect from light.**
8. Add 50 μL of Stop Solution to each well. The color in the wells should change from blue to yellow. If the color in the wells is green or if the color change does not appear uniform, gently tap the plate to ensure thorough mixing.
9. Determine the optical density of each well within 30 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.
10. **CALCULATION OF RESULTS :** Average the duplicate readings for each standard, control, and sample and subtract the average zero standard optical density. Create a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the IFN- γ concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data. If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Billiau, A. and P. Matthys (2009) Cytokine Growth Factor Rev. 20:97.
2. Wheelock, E.F. (1965) Science 146:310.
3. Schoenborn, J.R & C.B. Wilson (2007) Adv. Immunol. 96:41.
4. Pestka, S. et al. (2004) Immunol. Rev. 205:8.
5. Kelchtermans, H. et al. (2008) Trends Immunol. 29:479.
6. McLaren, J.E. and D.P. Ramji (20 09) Cytokine Growth Factor Rev. 20:125.
7. Gray, P.W. and D.V. Goeddel (1982) Nature 298:859.
8. Bach, E.E. et al. (1997) Annu. Rev. Immunol. 15:563.
9. Marsters, S.A. et al. (1995) Proc. Natl. Acad. Sci. 92:5401.
10. Krause, C.D. et al. (2000) J. Biol. Chem. 275:22995.
11. Sugaya, M. et al. (1999) J. Invest. Dermatol. 113:350.
12. Howie, S.E.M. et al. (1996) J. Invest. Dermatol. 106:1218.
13. Battistini, L. et al. (1997) J. Immunol. 159:3723.
14. Gupta, A.A. et al. (1996) J. Immunol. 157:2123.
15. Neumann, H. et al. (1997) J. Exp. Med. 186:2023.
16. Hoiden, I. & G. Moller (1996) Scand. J. Immunol. 44:501.
17. Puddu, P. et al. (1997) J. Immunol. 159:3490.
18. Yoshimoto, T. et al. (1997) Proc. Natl. Acad. Sci. USA 94:3948.
19. Yeaman, G.R. et al. (1998) J. Immunol. 160:5145.
20. Asea, A. et al. (1996) Clin. Exp. Immunol. 105:376.
21. Briscoe, D.M. et al. (1997) J. Immunol. 159:3247.
22. Dejuco, N. et al. (1995) Endocrinology 136:4925.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

1																	
2																	
3																	
4																	
5																	
6																	
7																	
8																	
9																	
10																	
11																	
12																	
	A	B	C	D	E	F	G	H									

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