PRODUCT INFORMATION & MANUAL

Human VEGF Valukine™ ELISA

VAL106

For the quantitative determination of natural and recombinant
human VEGF concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version202006.3

TABLE OF CONTENTS

I. BACKGROUND.....	2
II. OVERVIEW.....	4
III. ADVANTAGES.....	5
IV. EXPERIMENT.....	8
V. KIT COMPONENTS AND STORAGE.....	9
VI. PREPARATION.....	11
VII. ASSAY PROCEDURE.....	13
VIII. REFERENCES.....	15

I. BACKGROUND

Vascular endothelial growth factor (VEGF or VEGF-A), also known as vascular permeability factor (VPF), is a potent mediator of both angiogenesis and vasculogenesis in the fetus and adult (1-3). It is a member of the PDGF family that is characterized by the presence of eight conserved cysteine residues in a cystine knot structure and the formation of antiparallel disulfide-linked dimers (4). Humans express alternately spliced isoforms of 121, 145, 165, 183, 189, and 206 amino acids (aa) in length (4). VEGF165 appears to be the most abundant and potent isoform, followed by VEGF121 and VEGF189 (3, 4). Isoforms other than VEGF121 contain basic heparin-binding regions and are not freely diffusible (4). Human VEGF165 shares 88% aa sequence identity with corresponding regions of mouse and rat VEGF. VEGF is expressed in multiple cells and tissues including skeletal and cardiac muscle (5, 6), hepatocytes (7), osteoblasts (8), neutrophils (9), macrophages (10), keratinocytes (11), brown adipose tissue (12), CD34⁺ stem cells (13), endothelial cells (14), fibroblasts, and vascular smooth muscle cells (15). VEGF expression is induced by hypoxia and cytokines such as IL-1, IL-6, IL-8, oncostatin M and TNF- α (3, 4, 9, 16). VEGF isoforms are differentially expressed during development and in the adult (3).

VEGF dimers bind to two related receptor tyrosine kinases, VEGF R1 (also called Flt-1) and VEGF R2 (Flk-1/KDR), and induce their homodimerization and autophosphorylation (3, 4, 7, 17, 18). These receptors have seven extracellular immunoglobulin-like domains and an intracellular split tyrosine kinase domain. They are expressed on vascular endothelial cells and a range of non-endothelial cells. Although VEGF affinity is highest for binding to VEGF R1, VEGF R2 appears to be the primary mediator of VEGF angiogenic activity (3, 4). VEGF165 also binds the semaphorin receptor, neuropilin-1, which promotes complex formation with VEGF R2 (19).

VEGF is best known for its role in vasculogenesis. During embryogenesis, VEGF regulates the proliferation, migration, and survival of endothelial cells (3, 4), thus regulating blood vessel density and size, but playing no role in determining vascular patterns. VEGF promotes bone formation through osteoblast and chondroblast recruitment and is also a monocyte chemoattractant (20-22). After birth, VEGF maintains endothelial cell integrity and is a potent mitogen for micro- and macro-vascular endothelial cells. In adults, VEGF functions mainly in wound healing

and the female reproductive cycle (3). In diseased tissues, VEGF promotes vascular permeability. It is thus thought to contribute to tumor metastasis by promoting both extravasation and tumor angiogenesis (23, 24). Various strategies have been employed therapeutically to antagonize VEGF-mediated tumor angiogenesis (25). Circulating VEGF levels correlate with disease activity in autoimmune diseases such as rheumatoid arthritis, multiple sclerosis and systemic lupus erythematosus (26).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. A monoclonal antibody specific for VEGF has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any VEGF present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked polyclonal antibody specific for VEGF is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of VEGF bound in the initial step. The color development is stopped and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- ◆ **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- ◆ This kit is suitable for cell culture supernate, serum, EDTA Plasma and citrate plasma.
- ◆ The kit should not be used beyond the expiration date on the kit label.
- ◆ Do not mix or substitute reagents with those from other lots or sources.
- ◆ If samples generate values higher than the highest standard, dilute the samples with Diluent and repeat the assay.
- ◆ Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

Sample	Intra-assay Precision			Inter-assay Precision		
	1	2	3	1	2	3
Mean (pg/mL)	128	401	810	129	388	749
Standard Deviation	8.24	26.0	40.2	12.0	31.1	70.2
CV%	6.4	6.5	5.0	9.3	8.0	9.4

B. RECOVERY

The recovery of human VEGF spiked to levels throughout the range of the assay in various matrices was evaluated.

Sample Type	Average % Recovery	Range
Cell culture media (n=4)	121	113 - 135%
Serum (n=3)	124	116 - 127%
Plasma (n=6)	100	78 - 119%

C. SENSITIVITY

The minimum detectable dose (MDD) of human VEGF is typically less than 7.8 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against highly purified *Sf-21*-expressed recombinant human VEGF produced at R&D Systems.

E. LINEARITY

To assess the linearity of the assay, samples were spiked with high concentrations of human VEGF in various matrices and diluted with Diluent 1× to produce samples with values within the dynamic range of the assay.

Dilution		Cell culture media (n=4)	Serum (n=3)	Plasma (n=6)
1:2	Average % of Expected	100	88	102
	Range (%)	96 - 104	87 - 89	87 - 110
1:4	Average % of Expected	101	91	108
	Range (%)	95 - 108	90 - 91	81 - 124
1:8	Average % of Expected	102	89	117
	Range (%)	94 - 111	88 - 91	106 - 125
1:16	Average % of Expected	101	80	102
	Range (%)	88 - 110	74 - 86	88 - 119

F. SAMPLE VALUES

Cell Culture Supernates - Human peripheral blood mononuclear cells (1×10^6 cells/mL) were cultured in RPMI supplemented with 5% fetal calf serum, 50 μ M β -mercaptoethanol, 2 mM L-glutamine, 100 U/mL penicillin, and 100 μ g/mL streptomycin sulfate. Cells were cultured unstimulated or stimulated 10 μ g/mL PHA. Aliquots of the cell culture supernates was removed and assayed for levels of natural VEGF and measured 332 pg/ml and 1440 pg/ml respectively.

Serum - Three serum samples were evaluated for the presence of human VEGF in this assay. All samples measured ranged from 217.5 to 375.5 pg/mL with an average of 282.4 pg/mL.

Plasma - Six human plasma samples were evaluated for the presence of human VEGF in this assay. All samples measured below the lowest standard, 15.6 pg/mL.

G. SPECIFICITY

This assay recognizes both natural and recombinant human VEGF. The following factors were prepared at 50 ng/mL and assayed for cross-reactivity. Preparations of the following factors at 50 ng/mL in a mid-range rhVEGF control were assayed for interference. No significant cross-reactivity or interference was observed.

Recombinant human				Recombinant mouse	
ANG	sgp130	IL-8	PDGF-AA	IL-1 α	IL-13
AR	GRO α	IL-10	PDGF-AB	IL-1 β	MIP-1 α
CNTF	GRO β	IL-11	PDGF-BB	IL-3	MIP-1 β
β -ECGF EGF	GRO γ	IL-12	P/GF	IL-4	P/GF-2
Epo	HB-EGF	IL-13	TGF- β 1	IL-5	TNF- α
FGF acidic	HGF	LAP (TGF- β 1)	TNF- α	IL-6	VEGF120
FGF basic	IL-1 α	M-CSF	TNF R	IL-9	VEGF164
FGF-4	IL-1 β	MCP-1	sTNF RII	IL-10	VEGF-R3
FGF-5	IL-2	MIP-1 α	VEGF165/P/GF		
FGF-6	IL-3	MIP-1 β	VEGF-B167		
FGF-7	IL-4	β -NGF	VEGF-C	Recombinant rat	
G-CSF	IL-5	OSM	VEGF-D	VEGF	
GM-CSF	IL-6	PD-ECGF	VEGF-R3		


VEGF-related factors showing cross-reactivity or interference:

- ◆ rhVEGF R1/Fc Chimera interfered at levels >500 pg/mL
- ◆ rhVEGF R2/Fc Chimera interfered at levels >2000 pg/mL
- ◆ rmVEGF R1 interfered at levels >500 pg/mL
- ◆ rmVEGF R2 interfered at levels >4000 pg/mL
- ◆ Canine VEGF showed 67% cross-reactivity

IV. EXPERIMENT

EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

pg/mL	OD	Average	Corrected
0	0.013 0.016	0.015	-
15.63	0.036 0.035	0.036	0.021
31.25	0.059 0.062	0.061	0.046
62.5	0.131 0.115	0.123	0.109
125	0.283 0.222	0.253	0.238
250	0.471 0.531	0.501	0.487
500	0.995 1.024	1.010	0.995
1000	2.041 2.006	2.024	2.009

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
Human VEGF Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with a mouse monoclonal antibody against human VEGF	1 plate
Human VEGF Conjugate	solution of polyclonal antibody against human VEGF conjugated to horseradish peroxidase	1 vial
Human VEGF Standard	recombinant human VEGF in a buffered protein base; lyophilized	2 vials
Calibrator Diluent (5×)	a 5× concentrated buffered protein base	1 vial
Wash Buffer Concentrate (25×)	a 25× concentrated solution of buffered surfactant	1 vial
Color Reagent A	stabilized hydrogen peroxide	1 vial
Color Reagent B	stabilized chromogen (tetramethylbenzidine)	1 vial
Stop Solution	2 N sulfuric acid	1 vial
Plate Covers	adhesive strip	3 strips

B. STORAGE

Unopened Kit	Store at 2-8°C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Diluted Wash Solution	May be stored for up to 1 month at 2-8°C.*
	Stop Solution	
	Diluent 1×	
	Conjugate	
	Unmixed Color Reagent A	
	Unmixed Color Reagent B	
Reagents	Standard	Discard the VEGF stock solution and dilutions after 4 hours. Use a fresh standard for each assay.
	Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- ◆ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ◆ Pipettes and pipette tips.
- ◆ Deionized or distilled water.
- ◆ Squirt bottle, manifold dispenser, or automated microplate washer.
- ◆ 500 mL graduated cylinder.

D. PRECAUTION

- ◆ The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material.
- ◆ VEGF is detectable in saliva. Take the necessary precautions to prevent contamination of the kit reagents.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Diluent 1 \times .

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles.

Plasma - Collect plasma using EDTA, or citrate as an anticoagulant. Centrifuge for 15 minutes at 1000 x g within 30 minutes of collection. Assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles.

B. SAMPLE PREPARATION

Serum samples require a 5-fold dilution. A suggested 5-fold dilution is 40 μ L of sample + 160 μ L of Diluent (1 \times).

Plasma samples require a 2-fold dilution. A suggested 2-fold dilution is 100 μ L of sample + 100 μ L of Diluent (1 \times).

C. REAGENT PREPARATION

Note: Bring all reagents to room temperature before use.

Wash Solution - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25 \times) into deionized or distilled water to prepare 500 mL of Wash Buffer.

Substrate Solution - Color Reagent A and B should be mixed together in equal volumes within 15 minutes of use. Protect from light. 200 μ L of the resultant mixture is required per well.


Diluent 1 \times - Add 20 mL of Calibrator Diluent 5 \times into 80 mL of deionized or distilled water to prepare 100 mL of Diluent 1 \times .

VEGF Standard - Refer to the vial label for the reconstitution volume*. This reconstitution produces a stock solution of 2000 pg/mL. Allow the standard to sit for a

minimum of 15 minutes with gentle agitation prior to making dilutions.

*if you have any question, please seek help from our Technical Support.

Pipette 500 μL of Diluent 1 \times into each tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 1000 pg/mL standard serves as the high standard. Diluent 1 \times serves as the zero standard (0 pg/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Substrate Solution should remain colorless until added to the plate. Keep Substrate Solution protected from light. Substrate Solution should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the Substrate Solution. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the Substrate Solution.

VII. ASSAY PROCEDURE

Note: Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100 μ L of Standard, sample, or control per well. Cover with the adhesive strip provided. Incubate for 2 hours at room temperature. A plate layout is provided for a record of standards and samples assayed.
4. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Solution (400 μ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
5. Add 200 μ L of human VEGF Conjugate to each well. Cover with a new adhesive strip. Incubate for 2 hours at room temperature.
6. Repeat the aspiration/wash as in step 4.
7. Add 200 μ L of Substrate Solution to each well. Incubate for 30 minutes at room temperature. **Protect from light.**
8. Add 50 μ L of Stop Solution to each well. The color in the wells should change from blue to yellow. If the color in the wells is green or if the color change does not appear uniform, gently tap the plate to ensure thorough mixing.
9. Determine the optical density of each well within 30 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.
10. **CALCULATION OF RESULTS:** Average the duplicate readings for each standard, control, and sample and subtract the average zero standard optical density. Create

a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the VEGF concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data. If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Leung, D.W. et al. (1989) *Science* 246:1306.
2. Keck, P.J. et al. (1989) *Science* 246:1309.
3. Byrne, A.M. et al. (2005) *J. Cell. Mol. Med.* 9:777.
4. Robinson, C.J. and Stringer, S.E. (2001) *J. Cell. Sci.* 114:853.
5. Richardson, R.S. et al. (1999) *Am. J. Physiol.* 277:H2247.
6. Sugishita, Y. et al. (2000) *Biochem. Biophys. Res. Commun.* 268:657.
7. Yamane, A. et al. (1994) *Oncogene* 9:2683.
8. Goad, D.L. et al. (1996) *Endocrinology*137:2262.
9. Gaudry, M. et al. (1997) *Blood* 90:4153.
10. McLaren, J. et al. (1996) *J. Clin. Invest.* 98:482.
11. Diaz, B.V. et al. (2000) *J. Biol. Chem.* 275:642.
12. Asano, A. et al. (1997) *Biochem. J.* 328:179.
13. Bautz, F. et al. (2000) *Exp. Hematol.* 28:700.
14. Namiki, A. et al. (1995) *J. Biol. Chem.* 270:31189.
15. Nauck, M. et al. (1997) *Am. J. Respir. Cell. Mol. Biol.* 16:398.
16. Angelo, L.S. and R. Kurzrock (2007) *Clin. Cancer Res.* 13:2825.
17. Neufeld, G. et al. (1999) *FASEB. J.* 13:9.
18. Kowalewski, M.P. et al. (2005) Accession #ABB82619.
19. Pan, Q. et al. (2007) *J. Biol. Chem.* 282:24049.
20. Weis, S.M. and D.A. Cheresh (2005) *Nature* 437:497.
21. Breier, G. (2000) *Semin. Thromb. Hemost.* 26:553.
22. Barleon, B. et al. (1996) *Blood* 87:3336.
23. Weis, S.M. and D.A. Cheresh (2005) *Nature* 437:497.
24. Thurston, G. (2002) *J. Anat.* 200:575.
25. Grothey, A. and E. Galanis (2009) *Nat. Rev. Clin. Oncol.* 6:507.
26. Carvalho, J.F. et al. (2007) *J. Clin. Immunol.* 27:246.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								
	A	B	C	D	E	F	G	H

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