## **PRODUCT INFORMATION & MANUAL**

**Human MMP-9 Valukine™ ELISA**

**VAL113**

For the quantitative determination of natural and recombinant  
human MMP-9 concentrations

For research use only.  
Not for diagnostic or therapeutic procedures.

**Bio-Techne China Co. Ltd**

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001  
**info.cn@bio-techne.com**

Please refer to the kit label for expiry date.  
Novus kits are guaranteed for 3 months from date of receipt

Version202209.3

# TABLE OF CONTENTS

- I. BACKGROUND ..... 2
- II. OVERVIEW ..... 4
- III. ADVANTAGES ..... 5
- IV. EXPERIMENT ..... 7
- V. KIT COMPONENTS AND STORAGE ..... 8
- VI. PREPARATION ..... 10
- VII. ASSAY PROCEDURE ..... 12
- VIII. REFERENCES ..... 13

## I. BACKGROUND

Matrix metalloproteinases (MMPs), also called matrixins, constitute a family of zinc and calcium dependent endopeptidases that function in the breakdown of the extracellular matrix (ECM) and in the processing of a variety of molecules in different subcellular environments. They play an important role in many normal physiological processes such as embryonic development, morphogenesis, reproduction, and tissue remodeling (1, 2). They also participate in inflammatory and autoimmune disorders such as arthritis, cancer, and cardiovascular disease (3-5). While the amounts of newly synthesized MMPs are regulated mainly at the levels of transcription, the proteolytic activities of existing MMPs are controlled through both the activation of proenzymes or zymogens and the inhibition of active enzymes by endogenous inhibitors,  $\alpha$ 2-Macroglobulin, and tissue inhibitors of metalloproteinases (TIMPs) (6).

MMP-9 (also referred to as gelatinase B, 92 kDa type IV collagenase, 92 kDa gelatinase, and type V collagenase) is secreted as a glycosylated proenzyme (6-8). Activation of the proenzyme involves proteolytic removal of the N-terminal pro region, resulting in the 82 kDa active enzyme (9, 10). Active human MMP-9 shares 72% and 74% amino acid sequence identity with mouse and rat MMP-9, respectively. In addition to the zinc-binding site, the catalytic domain also contains three contiguous fibronectin type II homology units responsible for binding gelatin (11). A proline-rich hinge region links the catalytic domain to the C-terminal hemopexin-like domain. In vitro treatment of the proenzyme with 4-aminophenylmercuric acetate (APMA) produces not only the active enzyme but also a C-terminal truncated form with activity comparable to that of the active form (12). MMP-9 degrades components of the ECM with high specific activity for denatured collagens (gelatin). It can cleave native collagens of type III, IV, V, and XI, as well as Elastin, Nidogen-1, and Vitronectin (2, 3). MMP-9 can also cleave a variety of chemokines and growth factors (e.g. IL-1 $\beta$ , CXCL8/IL-8, CXCL7, CXCL4, CXCL1, Latent TGF- $\beta$ , membrane bound TNF- $\alpha$ , VEGF, and FGF basic), Amyloid  $\beta$  peptide, Substance P, and Myelin Basic Protein (3, 13-15). This action can increase or decrease the biological activity of soluble factors and can also liberate them from association with the ECM (16, 17). MMP-9 can also trigger signaling through various transmembrane proteins or inhibit signaling by inducing their shedding from the cell surface (e.g. CD44, E-Cadherin, Integrins, ICAM-1, and

IL-2 R $\alpha$ ) (3, 18-20).

MMP-9 is produced by a variety of normal and transformed cells including neutrophils, monocytes, macrophages, astrocytes, fibroblasts, osteoclasts, chondrocytes, keratinocytes, endothelial and epithelial cells. It exerts physiological and pathological angiogenic and remodeling effects on the vasculature (21-25). Activated neutrophils release proMMP-9 which is free of TIMP-1, allowing the liberation of pro-angiogenic FGF-2 from the ECM (17). MMP-9 in complex with TIMP-1 does not induce FGF-2 release (17). Neutrophil-derived MMP-9 exacerbates the inflammatory response, in part by generating collagen-derived peptides that induce the release of additional neutrophil MMP-9 (26). MMP-9 also plays a role in bone formation and remodeling (1, 21, 27), methamphetamine induced behavioral sensitization and reward (28), the regulation of neuronal synapse remodeling (29), trophoblast invasion during implantation (30), and the inactivation of Serpin  $\alpha$ 1-Proteinase Inhibitor (31). The shedding of adhesion proteins by MMP-9 has a direct effect on tumor cell invasiveness (18-20). Circulating levels of MMP-9 are increased in many inflammatory disorders including intraluminal thrombus formation (32), atherosclerosis (33), Crohn's disease (34), hepatitis C virus infection (35), colorectal cancer (36), and Duchenne muscular dystrophy (37). The ratio of MMP-9 to TIMP-1 is also increased in multiple sclerosis serum (38) and cystic fibrosis sputum (39), but it is decreased in the serum during cytomegalovirus infection (40). Levels of free MMP-9 and complexes of MMP-9 with Lipocalin-2/NGAL are elevated in the urine of ovarian cancer and uterine tract infection patients, respectively (41, 42).

## II. OVERVIEW

### A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. A monoclonal antibody specific for MMP-9 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any MMP-9 present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked polyclonal antibody specific for MMP-9 is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of MMP-9 bound in the initial step. The color development is stopped and the intensity of the color is measured.

### B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- This kit is suitable for cell culture supernate and serum.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Calibrator Diluent and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

### III. ADVANTAGES

#### A. PRECISION

##### **Intra-assay Precision** (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

##### **Inter-assay Precision** (Precision between assays)

Three samples were tested in forty separate assays to assess inter-assay precision.

Sample	Intra-assay Precision			Inter-assay Precision		
	1	2	3	1	2	3
n	20	20	20	40	40	40
Mean (ng/mL)	0.833	2.04	11.0	0.972	2.35	12.2
Standard Deviation	0.017	0.039	0.316	0.077	0.184	0.845
CV%	2.0	1.9	2.9	7.9	7.8	6.9

#### B. RECOVERY

The recovery of human MMP-9 spiked to different levels throughout the range of the assay in cell culture media was evaluated. The recovery ranged from 85-104% with an average of 97%.

#### C. SENSITIVITY

The minimum detectable dose (MDD) of MMP-9 is typically less than 0.156 ng/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

#### D. CALIBRATION

This immunoassay is calibrated against highly purified *CHO cell*-expressed recombinant human MMP-9 produced at R&D Systems®.

#### E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with high concentrations of MMP-9 and diluted with Calibrator Diluent 1x to produce samples with values within the dynamic range of the assay.

Dilution	Average % of Expected	Range (%)
1:2	104	101 – 107
1:4	104	101 – 106
1:8	102	97 – 110
1:16	102	96 – 110

## F. SAMPLE VALUES

**Cell Culture Supernates** - Human peripheral blood mononuclear cells ( $5 \times 10^6$  cells/mL) were cultured in RPMI supplemented with 5% fetal calf serum, 50  $\mu$ M  $\beta$ -mercaptoethanol, 2 mM L-glutamine, 100 U/mL penicillin, and 100  $\mu$ g/mL streptomycin sulfate. Cells were cultured unstimulated or stimulated with 10  $\mu$ g/mL PHA. Aliquots of the cell culture supernate was removed on day 1 and 5 and assayed for human MMP-9.

Condition	Day 1 (ng/mL)	Day 3 (ng/mL)
Unstimulated	132	33.8
Stimulated	522	210

**Serum** - Four serum samples were evaluated for the presence of MMP-9 in this assay. All samples measured ranged from 557.9 to 610.8 ng/mL with an average of 581.6 ng/mL.

## G. SPECIFICITY

This assay recognizes both natural and recombinant 92 kDa Pro-MMP-9 and 82 kDa active MMP-9. It does not measure the 65 kDa form. The following factors were prepared at 100 or 200 ng/mL in Diluent 1x and assayed for cross-reactivity. Preparations of the following factors at 200 ng/mL in a mid-range rhMMP-9 control were assayed for interference. No significant cross-reactivity or interference was observed.

### Recombinant Human

ADAM8	MMP-8
ADAM10	MMP-10
ADAM15	MMP-12
Lipocalin-2/NGAL	MMP-13
MMP-1	TACE (ADAM17)
MMP-2	TIMP-2
MMP-3	TIMP-3
MMP-7	TIMP-4

No cross-reactivity was observed with recombinant human TIMP-1, but interference was observed at concentrations  $\geq 6.25$  ng/mL.

### Recombinant mouse


ADAM9	MMP-2
ADAM10	MMP-3
ADAM15	MMP-9
Lipocalin-2/NGAL	TIMP-1

Interference was also observed with recombinant rat TIMP-1 at concentrations  $\geq 200$  ng/mL.

## IV. EXPERIMENT

### EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

(ng/mL)	O.D.	Average	Corrected
0	0.040 0.041	0.040	—
0.313	0.111 0.110	0.110	0.070
0.625	0.172 0.170	0.172	0.132
1.25	0.292 0.288	0.290	0.250
2.5	0.532 0.508	0.520	0.480
5	0.910 0.899	0.904	0.864
10	1.612 1.668	1.640	1.600
20	2.607 2.638	2.622	2.582

## V. KIT COMPONENTS AND STORAGE

### A. MATERIALS PROVIDED

Parts	Description	Size
MMP-9 Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with a mouse monoclonal antibody against human MMP-9	1 plate
MMP-9 Conjugate	monoclonal antibody against MMP-9 conjugated to horseradish peroxidase with preservatives	1 vial
Assay Diluent	a buffered protein base with preservatives	2 vials
MMP-9 Standard	recombinant human MMP-9 in a buffered protein base with preservatives; lyophilized	1 vial
Calibrator Diluent (5×)	a solution of 5× buffered protein base with preservatives	2 vials
Wash Buffer Concentrate (25×)	a 25× concentrated solution of buffered surfactant with preservatives	1 vial
TMB Substrate	TMB ELISA Substrate Solution	2 vials
Stop Solution	2 N sulfuric acid	1 vial
Plate Covers	adhesive strip	3 strips

### B. STORAGE

<b>Unopened Kit</b>	Store at 2-8°C. Do not use past kit expiration date.	
<b>Opened/ Reconstituted Reagents</b>	Diluted Wash Solution	May be stored for up to 1 month at 2-8°C.*
	Stop Solution	
	Assay Diluent	
	Calibrator Diluent 1×	
	Conjugate	
	TMB Substrate	Use a new standard for each assay. Discard after use.
	Standard	
Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*	

\* Provided this is within the expiration date of the kit.

### C. OTHER SUPPLIES REQUIRED

- Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- Pipettes and pipette tips.
- Deionized or distilled water.
- Squirt bottle, manifold dispenser, or automated microplate washer.
- **Horizontal orbital microplate shaker (0.12" orbit) capable of maintaining a speed of 500 ± 50 rpm.**
- 500 mL graduated cylinder.
- Polypropylene tubes for dilution of standards and samples.

### D. PRECAUTION

- MMP-9 is detectable in saliva. Use a face mask and gloves to protect kit reagents from contamination.
- Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- Wear protective gloves, clothing, eye, and face protection. Wash hands thoroughly after handling.

## VI. PREPARATION

### A. SAMPLE COLLECTION AND STORAGE

**Cell Culture Supernates** - Remove particulates by centrifugation and assay immediately or aliquot and store samples at  $\leq -20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles.

**Serum** - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at  $\leq -20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles.

### B. SAMPLE PREPARATION

Cell culture supernate samples require at least a 100-fold dilution prior to the assay. A suggested 100-fold dilution is 10  $\mu\text{L}$  of Activated sample + 990  $\mu\text{L}$  of Calibrator Diluent (1 $\times$ ).

Serum samples require a 200-fold dilution. For example, add 10  $\mu\text{L}$  of serum into a tube with 190  $\mu\text{L}$  Calibrator Diluent (1 $\times$ ) to prepare a 20-fold diluted sample. Mix through and then pipette 20  $\mu\text{L}$  of prepared 20-fold diluted sample into a tube with 180  $\mu\text{L}$  Calibrator Diluent (1 $\times$ ) to prepare a final 200-fold diluted sample.

### C. REAGENT PREPARATION

**Note:** *Bring all reagents to room temperature before use.*


**Wash Solution** - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25 $\times$ ) into deionized or distilled water to prepare 500 mL of Wash Buffer.

**Calibrator Diluent 1 $\times$**  - Add 20 mL of Calibrator Diluent 5 $\times$  into 80 mL of deionized or distilled water to prepare 100 mL of Calibrator Diluent 1 $\times$ .

**MMP-9 Standard-** Refer to the vial label for reconstitution volume\*. This reconstitution produces a stock solution of 20 ng/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

\*if you have any question, please seek help from our Technical Support.

**Use polypropylene tubes.** Pipette 500  $\mu\text{L}$  of the appropriate Calibrator Diluent 1 $\times$  into each tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 20 ng/mL standard serves as the high standard. The appropriate Calibrator Diluent 1 $\times$  serves as the zero standard (0 pg/mL).

#### D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep Substrate Solution protected from light. Substrate Solution should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the Substrate Solution. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the Substrate Solution.

## VII. ASSAY PROCEDURE

**Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.**

*Note: High concentrations of MMP-9 are found in saliva. It is recommended that a face mask and gloves be used to protect kit reagents from contamination.*

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plat frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100  $\mu$ L of Assay Diluent to each well.
4. Add 100  $\mu$ L of Standard, sample, or control per well. Cover with the adhesive strip provided. Incubate for 2 hours at room temperature on a horizontal orbital microplate shaker (0.12" orbit) set at  $500 \pm 50$  rpm. A plate layout is provided for a record of standards and samples assayed.
5. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Solution (400  $\mu$ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
6. Add 200  $\mu$ L of MMP-9 Conjugate to each well. Cover with a new adhesive strip. Incubate for 1 hours at room temperature.
7. Repeat the aspiration/wash as in step 5.
8. Add 200  $\mu$ L of TMB Substrate to each well. Incubate for 30 minutes at room temperature. **Protect from light.**
9. Add 50  $\mu$ L of Stop Solution to each well. The color in the wells should change from blue to yellow. If the color in the wells is green or if the color change does not appear uniform, gently tap the plate to ensure thorough mixing.
10. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.
11. **CALCULATION OF RESULTS:** Average the duplicate readings for each standard, control, and sample and subtract the average zero standard optical density. Create a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the MMP-9 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data. If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

## VIII. REFERENCES

1. Page-McCaw, A. et al. (2007) *Nat. Rev. Mol. Cell. Biol.* 9:221.
2. Rodriguez, D. et al. (2010) *Biochim. Biophys. Acta* 1803:39.
3. Bauvois, B. (2012) *Biochim. Biophys. Acta* 1825:29.
4. Heissig, B. et al. (2010) *Histol. Histopathol.* 25:765.
5. Ram, M. et al. (2006) *J. Clin. Immunol.* 26:299.
6. Nagase, H. et al. (2006) *Cardiovasc. Res.* 69:562.
7. Okada, A. et al. (1995) *Gene* 164:317.
8. Xia, Y. et al. (1996) *FEBS Lett.* 382:285.
9. Van Wart, H.E. and H. Birkedal-Hansen (1990) *Proc. Natl. Acad. Sci. USA* 87:5578.
10. Okada, Y. et al. (1992) *J. Biol. Chem.* 267:21712.
11. Collier, I.E. et al. (1992) *J. Biol. Chem.* 267:6776.
12. O'Connell, J.P. et al. (1994) *J. Biol. Chem.* 269:14967.
13. Backstrom, J.R. et al. (1996) *J. Neurosci.* 16:7910.
14. Backstrom, J.R. and Z.A. Tokes (1995) *J. Neurochem.* 64:1312.
15. Shiryayev, S.A. et al. (2009) *PLoS ONE* 4:e4952.
16. Van den Steen, P.E. et al. (2000) *Blood* 96:2673.
17. Ardi, V.C. et al. (2009) *J. Biol. Chem.* 284:25854.
18. Chetty, C. et al. (2012) *Cell Signal.* 24:549.
19. Zuo, J.H. et al. (2011) *J. Cell. Biochem.* 112:2508.
20. Pal-Ghosh, S. et al. (2011) *J. Cell Sci.* 124:2666.
21. Vu, T.H. et al. (1998) *Cell* 93:411.
22. Pyo, R. et al. (2000) *J. Clin. Invest.* 105:1641.
23. Bergers, G.R. et al. (2000) *Nat. Cell Biol.* 2:737.
24. Svedin, P. et al. (2007) *J. Neurosci.* 27:1511.
25. Johnson, J.L. et al. (2011) *Arterioscler. Throm. Vasc. Biol.* 31:e35.
26. Xu, X. et al. (2011) *PLoS ONE* 6:e15781.

27. Engsig, M.T. et al. (2000) *J. Cell Biol.* 151:879.
28. Mizoguchi, H. et al. (2005) *J. Neurochem.* 100:1579.
29. Michaluk, P. et al. (2011) *J. Cell Sci.* 124:3369.
30. Bischof, P. et al. (2002) *J. Reprod. Immunol.* 55:3.
31. Liu, Z. et al. (2000) *Cell* 102:647.
32. Wiernicki, I. et al. (2011) *Dis. Markers* 31:67.
33. Tanindi, A. et al. (2011) *Open Cardiovasc. Med. J.* 5:110.
34. Kofla-Dlubacz, A. et al. (2011) *Dig. Dis. Sci.* Oct. 14 Epub.
35. Helaly, G.F. (2011) *Br. J. Biomed. Sci.* 68:38.
36. Dragutinovic, V.V. et al. (2011) *Mol. Cell. Biochem.* 355:173.
37. Nadarajah, V.D. et al. (2011) *Neuromuscul. Disord.* 21:569.
38. Avolio, C. et al. (2005) *Mult. Scler.* 11:441.
39. Jackson, P.L. et al. (2010) *Mol. Med.* 16:159.
40. Straat, K. et al. (2009) *J. Virol.* 83:830.
41. Coticchia, C.M. et al. (2011) *Gynecol. Oncol.* 123:295.
42. Hatipoglu, S. et al. (2011) *Pediatr. Nephrol.* 26:1263.

## PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

<b>1</b>									
<b>2</b>									
<b>3</b>									
<b>4</b>									
<b>5</b>									
<b>6</b>									
<b>7</b>									
<b>8</b>									
<b>9</b>									
<b>10</b>									
<b>11</b>									
<b>12</b>									
	<b>A</b>	<b>B</b>	<b>C</b>	<b>D</b>	<b>E</b>	<b>F</b>	<b>G</b>	<b>H</b>	

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