PRODUCT INFORMATION & MANUAL

Human GM-CSF Valukine™ ELISA

VAL124

For the quantitative determination of natural and recombinant human
GM-CSF concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd
P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version 201907.01

TABLE OF CONTENTS

I. BACKGROUND.....	2
II. OVERVIEW.....	3
III. ADVANTAGES.....	4
IV. EXPERIMENT.....	6
V. KIT COMPONENTS AND STORAGE.....	7
VI. PREPARATION.....	10
VII. ASSAY PROCEDURE.....	12
VIII. REFERENCES.....	14

I. BACKGROUND

Granulocyte-Macrophage Colony Stimulating Factor (GM-CSF), also known as CSF-2, is a pleiotrophic 30 kDa member of the Common beta Chain (β c) cytokine family that also includes IL-3 and IL-5. GM-CSF is secreted by a wide variety of activated immune, mesenchymal, and epithelial cell types and circulates as a variably glycosylated monomer (1, 2). It is upregulated in multiple cell types during inflammation including encephalitogenic T cells (3-5), allergen exposed lung endothelial cells (6), and IgE activated mast cells (7). Mature human GM-CSF shares 54% and 63% amino acid sequence identity with mouse and rat GM-CSF, respectively (8).

Number of functions have been attributed to GM-CSF. It induces monocyte, neutrophil, and eosinophil production from CD34+ stem cell precursors (9, 10). It can act togetherwith IL-4 or Flt-3 Ligand to induce the development and maintenance of myeloid and dermal dendritic cells (10-14). It also acts as a neutrophil and dendritic cell chemoattractant (6, 15). GM-CSF promotes Th1 and Th17 cell mediated autoimmune inflammation as well as the inflammatory activation of dendritic cells, microglia, alveolar macrophages, and eosinophils (3-5, 16-20). In addition, it cooperates with G-CSF in promoting tumor cell proliferation and invasion (21).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for human GM-CSF has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any GM-CSF present is bound by the immobilized antibody. After washing away any unbound substances, a biotin-linked detect antibody specific for human GM-CSF are pipetted into the wells. After washing away any unbound substances, streptavidin-HRP is added. Following a wash to remove any unbound antibody-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of GM-CSF bound in the initial step. The color development is stopped, and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- This kit is suitable for cell culture supernate and serum.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Reagent Diluent or Calibrator Diluent and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

CELL CULTURE SUPERNATE / SERUM

Sample	Intra-assay Precision			Inter-assay Precision		
	1	2	3	1	2	3
Mean (pg/mL)	31.3	123.2	531.9	31.6	123.9	505.6
Standard Deviation	1.1	3.4	17.5	1.2	2.9	34.9
CV%	3.5	2.8	3.3	3.9	2.4	6.9

B. RECOVERY

The recovery of human GM-CSF spiked to three levels throughout the range of the assay in cell culture media was evaluated. The recovery ranged from 84.3-96.2% with an average of 89.7%.

The recovery of human GM-CSF spiked to three levels throughout the range of the assay in serum samples was evaluated. The recovery ranged from 111.5-118.9% with an average of 116.2%.

C. SENSITIVITY

The minimum detectable dose (MDD) of GM-CSF is typically less than 3.8 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against a highly purified *E. coli*-expressed recombinant human GM-CSF produced at R&D Systems®.

E. LINEARITY

To assess the linearity of the assay, samples spiked with high concentrations of GM-CSF and diluted with Reagent Diluent 1× (*for cell culture supernate samples*) or Calibrator Diluent 1× (*for serum samples*) to produce samples with values within the dynamic range of the assay.

Dilution	Average % of Expected	Range (%)
1:2	105.2%	94.1-113.6%
1:4	106.8%	90.0-119.6%
1:8	101.9%	86.2-116.6%
1:16	107.8%	87.0-119.6%

F. SAMPLE VALUES

Five Serum samples were evaluated for the presence of GM-CSF in this assay. All samples measured less than the lowest human GM-CSF standard, 15.6pg/mL.

G. SPECIFICITY

This assay recognizes both natural and recombinant human GM-CSF. The following factors prepared at 50 ng/mL were assayed and exhibited no cross-reactivity or interference.

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
IV. EXPERIMENT

EXAMPLE STANDARD


The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

CELL CULTURE SUPERNATE ASSAY

pg/mL	OD	Average	Corrected
0	0.014	0.014	-
15.6	0.114	0.115	0.101
31.3	0.217	0.215	0.202
62.5	0.408	0.409	0.396
125	0.750	0.758	0.745
250	1.299	1.294	1.281
500	1.980	1.952	1.939
1000	2.484	2.455	2.442

SERUM ASSAY

pg/mL	OD	Average	Corrected
0	0.020	0.021	-
15.6	0.086	0.086	0.066
31.3	0.140	0.144	0.123
62.5	0.275	0.277	0.257
125	0.516	0.520	0.500
250	0.927	0.958	0.937
500	1.539	1.570	1.549
1000	2.241	2.252	2.232

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Store the unopened kit at 2-8°C. Do not use past kit expiration date.

Parts	Description	Size
Human GM-CSF Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with a mouse antibody against human GM-CSF.	1 plate
Human GM-CSF Standard	Recombinant human GM-CSF in a buffered protein base; lyophilized. Refer to the vial label for reconstitution volume.	2 vials
Human GM-CSF Detection Antibody	Biotinylated GM-CSF polyclonal antibody, lyophilized. Refer to the vial label for reconstitution volume.	1 vial
Streptavidin-HRP B (40×)	40× Streptavidin conjugated to horseradish peroxidase.	1 vial
Reagent Diluent (10×)	A 10× concentrated buffered protein base used to dilute Detection Antibody, HRP, standard and cell culture supernate samples.	1 vial
Calibrator Diluent (1×)	A buffered protein base used to dilute standard and serum samples.	2 vial
Wash Buffer Concentrate (25×)	A 25× concentrated solution of buffered surfactant with preservatives.	1 vial
Color Reagent A	Stabilized hydrogen peroxide.	1 vial
Color Reagent B	Stabilized chromogen (tetramethylbenzidine).	1 vial
Stop Solution	2 N sulfuric acid.	1 vial
Plate Sealers	Adhesive strips.	3 strips

B. STORAGE

Unopened Kit	Store at 2-8 °C. Do not use past kit expiration date.
Opened/ Reconstituted Reagents	Streptavidin-HRP B
	Diluted Wash Solution
	Unmixed Color Reagent A
	Unmixed Color Reagent B
	Stop Solution
	Calibrator Diluent (1×)
	Standard
	Detection Antibody
	Reagent Diluent (10×)
Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8 °C.*

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- ◆ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ◆ Pipettes and pipette tips.
- ◆ Deionized or distilled water.
- ◆ Squirt bottle, manifold dispenser, or automated microplate washer.
- ◆ Test tubes for dilution of standards.
- ◆ 100mL and 500 mL graduated cylinder.

D. PRECAUTION

- ◆ Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- ◆ Color Reagent B may cause skin, eye, and respiratory irritation. Avoid breathing fumes.
- ◆ The Stop Solution provided with this kit is an acid solution.
- ◆ Wear protective gloves, clothing, eye, and face protection. Wash hands thoroughly after handling.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles.

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles.

B. SAMPLE PREPARATION

Cell culture supernate samples require a 4-fold dilution prior to the assay. A suggested 4-fold dilution is 50 μL of sample + 150 μL of **Reagent Diluent (1 \times)**.

Serum and plasma samples require a 2-fold dilution prior to the assay. A suggested 2-fold dilution is 100 μL of sample + 100 μL of **Calibrator Diluent (1 \times)**.

C. REAGENT PREPARATION

Note: Bring all reagents to room temperature before use.

Wash Buffer - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25 \times) into deionized or distilled water to prepare 500 mL of Wash Buffer.

Reagent Diluent (1 \times) - Add 7 mL of Reagent Diluent (10 \times) into 63 mL of deionized or distilled water to prepare 70 mL of Reagent Diluent (1 \times).

Detection Antibody- Reconstitution Volume refer to vial label with Reagent Diluent (1 \times). Aliquot and store if needed. Dilute stock solution in Reagent Diluent (1 \times) to the working concentration of 10 ng/mL. Prepare at least 15 minutes prior to use.

Streptavidin-HRP B (1 \times) - Dilute to the working concentration specified on the vial label using Reagent Diluent (1 \times).


Substrate Solution - Color Reagent A and B should be mixed together in equal volumes within 15 minutes of use. Protect from light. 100 μL of the resultant mixture is required per well.

GM-CSF Standard – Refer to the vial label for the reconstitution volume* using Reagent Diluent 1 \times (for cell culture supernate samples) or Calibrator Diluent 1 \times (for

serum samples). This reconstitution produces a stock solution of 2000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*if you have any question, please seek help from our Technical Support.

Pipette 300 µL of Reagent Diluent 1× (for cell culture supernate samples) or Calibrator Diluent 1× (for serum samples) into each tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 1000 pg/mL standard serves as the high standard. The **Reagent Diluent 1× (for cell culture supernate samples) or Calibrator Diluent 1× (for serum samples)** serves as the zero standard (0 pg/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Substrate Solution should remain colorless until added to the plate. Keep Substrate Solution protected from light. Substrate Solution should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the Substrate Solution. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the Substrate Solution.

VII. ASSAY PROCEDURE

Note: Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100 µL of Standard, or prepared sample per well. Cover with the adhesive strip provided. Incubate for 2 hours at room temperature. A plate layout is provided for a record of standards and samples assayed. (Samples may require dilution. See Sample Preparation section.)
4. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Solution (400 µL) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
5. Add 100 µL of the Detection Antibody diluted in Reagent Diluent, to each well. Cover with a new adhesive strip and incubate 2 hours at room temperature.
6. Repeat the aspiration/wash as in step 4.
7. Add 100 µL of the working dilution of Streptavidin-HRP B to each well. Cover the plate and incubate for 20 minutes at room temperature. Avoid placing the plate in direct light.
8. Repeat the aspiration/wash as in step 4.
9. Add 100 µL of Substrate Solution to each well. Incubate for 30 minutes at room temperature. Avoid placing the plate in direct light.
10. Add 50 µL of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
11. Determine the optical density of each well immediately, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

12. CALCULATION OF RESULTS

Average the duplicate readings for each standard, control, and sample and subtract the average zero standard optical density (O.D.).

Create a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a

standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the GM-CSF concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.


If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Shiomi, A. and T. Usui (2015) *Mediators Inflamm.* 2015:568543.
2. Broughton, S.E. et al. (2012) *Immunol. Rev.* 250:277.
3. Ponomarev, E.D. et al. (2007) *J. Immunol.* 178:39.
4. El-Behi, M. et al. (2011) *Nat. Immunol.* 12:568.
5. Codarri, L. et al. (2011) *Nat. Immunol.* 12:560.
6. Willart, M.A. et al. (2012) *J. Exp. Med.* 209:1505.
7. Wodnar-Filipowicz, A. et al. (1989) *Nature* 339:150.
8. Wong, G.G. et al. (1985) *Science* 228:810.
9. Wognum, A.W. (1994) *Blood* 84:764.
10. Inaba, K. et al. (1992) *J. Exp. Med.* 176:1693.
11. Herbst, B. et al. (1998) *Br. J. Haematol.* 101:231.
12. Kingston, D. et al. (2009) *Blood* 114:835.
13. Gilliet, M. et al. (2002) *J. Exp. Med.* 195:953.
14. Xu, Y. et al. (2007) *J. Immunol.* 179:7577.
15. Gomez-Cambronero, J. et al. (2003) *J. Immunol.* 171:6846.
16. Sonderegger, I. et al. (2008) *J. Exp. Med.* 205:2281.
17. Berclaz, P.Y. et al. (2002) *Blood* 100:4193.
18. Parajuli, B. et al. (2012) *J. Neuroinflammation* 9:268.
19. Sorgi, C.A. et al. (2012) *PLoS One* 7:e40523.
20. Horie, S. and H. Kita (1994) *J. Immunol.* 512:5457.
21. Obermueller, E. et al. (2004) *Cancer Res.* 64:7801.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

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