## **PRODUCT INFORMATION & MANUAL**

**Hyaluronan Valukine™ ELISA**

**Catalog Number: VAL152**

For the quantitative determination of natural and recombinant  
Hyaluronan concentrations

For research use only.  
Not for diagnostic or therapeutic procedures.

**Bio-Techne China Co. Ltd**

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001  
**info.cn@bio-techne.com**

Please refer to the kit label for expiry date.  
Novus kits are guaranteed for 3 months from date of receipt

Version 202308.2

# TABLE OF CONTENTS

|                                     |    |
|-------------------------------------|----|
| I. BACKGROUND .....                 | 2  |
| II. OVERVIEW .....                  | 4  |
| III. ADVANTAGES .....               | 5  |
| IV. EXPERIMENT .....                | 7  |
| V. KIT COMPONENTS AND STORAGE ..... | 8  |
| VI. PREPARATION .....               | 11 |
| VII. ASSAY PROCEDURE .....          | 13 |
| VIII. REFERENCES .....              | 15 |

## I. BACKGROUND

Hyaluronan (HA), also called hyaluronic acid or sodium hyaluronate, is a naturally occurring, highly hydrated linear polymer composed of the repeating disaccharide D-glucuronic acid-beta-1, 3-N-acetylglucosamine-beta-1,4 (1-6). It is a glycosaminoglycan (GAG) that is ubiquitously present in the extracellular matrix of all vertebrates. Unlike other GAGs such as chondroitin sulfate, dermatan sulfate, and heparan sulfate, HA lacks sulfate groups and is almost never covalently associated with other proteins (3, 6, 7). Mammalian HA is synthesized by one of three distinct multipass transmembrane HA synthases, HAS1, 2, and 3 (1-4, 7, 8). During HA synthesis, the growing polymer is translocated across the membrane to the extracellular surface as UDP-sugar substrates are added (7). The three HA synthases differ in their expression patterns and rates of synthesis, and may produce HA polymers with different chain lengths (8). In general, HA synthases produce high molecular weight HA (HMW, >500 kDa), while degradation is the main source of low molecular weight HA fragments (LMW, <500 kDa). HA turns over rapidly in mammals via the concerted actions of a family of hyaluronidases, and is cleared by the liver and kidneys (1, 3, 5, 7, 9).

HA is important in the tissues for hydration of the extracellular matrix, cell adhesion and migration (1, 6, 7). It interacts with a large number of hyaluronan-binding proteins, or hyaladherins. In the extracellular matrix, HA binds Link proteins such as aggrecan and HAPLN1 in cartilage, versican in smooth muscle, and brevican in the brain, and neurocan in nervous tissue (1, 6, 7, 10). In the circulation, it binds plasma binding proteins such as ITIH1/2, TSG-6, and HABP1/C1qBP (1, 3, 5, 7, 10). HA binds many cell surface receptors, such as the widely expressed CD44, RHAMM on vascular endothelia and other cells, LYVE-1 on lymphatic endothelia, and the endocytic receptor Stabilin-2 in endothelial sinuses such as in lymph nodes and liver (5-7, 10). Myeloid cell TLR2 and TLR4 also bind LMW HA fragments (11-13).

Cell signaling pathways activated through HA receptors are in part dependent on the size of HA (4, 6, 8, 10-13). While HMW HA is considered anti-angiogenic, anti-inflammatory, and immunosuppressive, LMW HA can be highly angiogenic and pro-inflammatory, and promotes dendritic cell maturation (1, 4-6, 12, 14). HA is mainly produced by fibroblasts, and contributes to fibroblast production of inflammatory

mediators and migration to areas of fibrosis (1, 6). Its production by aortic smooth muscle and fibroblasts is increased in the aging arterial wall (6, 15). HA promotes mesodermal and cardiac differentiation from embryonic stem cells, and hematopoietic differentiation from hematopoietic stem cells in vitro (1). Plasma HMW HA mainly comes from the lymphatic circulation (7). It is elevated in angiopathy, liver fibrosis, autoimmune diseases such as advanced scleroderma, rheumatoid arthritis and lupus nephritis, end-stage renal failure, idiopathic pulmonary arterial hypertension, adult respiratory distress syndrome, and in diabetics with poor glucose control (1, 15, 16). HA is used therapeutically as a drug carrier, and injected into diseased joints to increase hydration and suppress inflammation (5, 17, 18). It protects tissues by scavenging reactive oxygen species; however, reactive oxygen species can also cause HA fragmentation, contributing to inflammation (1). HMW HA binding to CD44 promotes activity of regulatory T cells, which control inflammation (14, 19). HA promotes epithelial to mesenchymal transition, both during development and during oncogenic transformation (1, 7, 20). Cancer cells frequently produce hyaluronidases which degrade HA, promoting angiogenesis and inflammation (4, 20).

## II. OVERVIEW

### A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. Recombinant human Aggrecan has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any Hyaluronan present is bound by the immobilized recombinant human Aggrecan. After washing away any unbound substances, a biotinylated recombinant human Aggrecan specific for Hyaluronan is pipetted into the wells. After washing away any unbound substances, Streptavidin-HRP is pipetted into the wells. Following a wash to remove any unbound reagent, TMB substrate solution (Chromogenic agent) is added to the wells and color develops in proportion to the amount of Hyaluronan bound in the initial step. The color development is stopped, and the intensity of the color is measured.

### B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- This kit is suitable for cell culture supernate and serum.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Reagent Diluent (1×) and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

### III. ADVANTAGES

#### A. PRECISION

##### **Intra-assay Precision** (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

##### **Inter-assay Precision** (Precision between assays)

Three samples were tested in twenty-eight separate assays to assess inter-assay precision.

| Sample             | Intra-assay Precision |     |      | Inter-assay Precision |     |      |
|--------------------|-----------------------|-----|------|-----------------------|-----|------|
|                    | 1                     | 2   | 3    | 1                     | 2   | 3    |
| Mean (ng/mL)       | 1.2                   | 4.3 | 17.9 | 1.2                   | 4.5 | 18.4 |
| Standard Deviation | 0.1                   | 0.3 | 1.1  | 0.1                   | 0.4 | 1.5  |
| CV%                | 7.3                   | 6.1 | 6.2  | 7.2                   | 8.0 | 8.1  |

#### B. RECOVERY

The recovery of Hyaluronan spiked to different levels throughout the range of the assay in cell culture media was evaluated. The recovery ranged from 84-106.9% with an average of 95.0%.

The recovery of Hyaluronan spiked to different levels throughout the range of the assay in human serum was evaluated. The recovery ranged from 88.1-118.8% with an average of 102.7%.

#### C. SENSITIVITY

The minimum detectable dose (MDD) of Hyaluronan is typically less than 0.081 ng/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

## D. CALIBRATION

This immunoassay is calibrated against a highly purified form of Hyaluronan (Catalog #GLR004) from R&D Systems®.

## E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with high concentrations of Hyaluronan and diluted with Reagent Diluent (1×) to produce samples with values within the dynamic range of the assay.

| Dilution | Average % of Expected | Range (%)   |
|----------|-----------------------|-------------|
| 1:2      | 100.1                 | 90.3-110.3  |
| 1:4      | 105.2                 | 86.1-114.9  |
| 1:8      | 107.2                 | 106.2-120.9 |
| 1:16     | 100.8                 | 87.9-115.1  |

## F. SAMPLE VALUES

**Cell Culture Supernates** - TF-1 human blood leukemia cells ( $1 \times 10^5$  cells/mL) and SNU-1 gastric cancer cells ( $2 \times 10^5$  cells/mL) were cultured in RPMI 1640 supplemented with 10% fetal bovine serum, 2 mM L-glutamine, 100 U/mL penicillin, and 100 µg/mL streptomycin sulfate for 3 days. MCF-7 human breast cancer cells ( $4 \times 10^5$  cells/mL) were cultured in DMEM supplemented with 10% fetal bovine serum, 2 mM L-glutamine, 100 U/mL penicillin, and 100 µg/mL streptomycin sulfate for 3 days. An aliquot of the cell culture supernate was removed, assayed for human Hyaluronan, and measured ranged from 35.8 to 44.4 ng/mL with an average of 38.8 ng/mL.

**Serum** - Six human serum samples were evaluated for the presence of Hyaluronan in this assay. All samples measured ranged from 44 to 126 ng/mL with an average of 77 ng/mL.

## G. SPECIFICITY


Billirubin, Hemoglobin, and Heparin were assayed at 10, 50, and 1000 µg/mL, respectively, and exhibited no cross-reactivity or interference in this assay.

The low molecular weight (15-40 kDa), medium molecular weight (75-350 kDa), and high molecular weight (>950 kDa) forms of Hyaluronan are all detected in this assay.

## IV. EXPERIMENT

### EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

| ng/mL | OD    | Average | Corrected |
|-------|-------|---------|-----------|
| 0     | 0.052 | 0.052   | -         |
|       | 0.052 |         |           |
| 0.625 | 0.136 | 0.141   | 0.089     |
|       | 0.145 |         |           |
| 1.25  | 0.235 | 0.236   | 0.184     |
|       | 0.236 |         |           |
| 2.5   | 0.411 | 0.419   | 0.367     |
|       | 0.426 |         |           |
| 5     | 0.651 | 0.692   | 0.640     |
|       | 0.733 |         |           |
| 10    | 1.111 | 1.197   | 1.145     |
|       | 1.283 |         |           |
| 20    | 1.665 | 1.685   | 1.633     |
|       | 1.705 |         |           |
| 40    | 2.163 | 2.166   | 2.114     |
|       | 2.168 |         |           |

## V. KIT COMPONENTS AND STORAGE

### A. MATERIALS PROVIDED

| Parts                         | Description   | Size     |
|-------------------------------|---|----------|
| Hyaluronan Microplate         | 96 well polystyrene microplate (12 strips of 8 wells) coated with recombinant human Aggrecan.           | 1 plate  |
| Hyaluronan Standard           | Hyaluronan in a buffered protein base; lyophilized. Refer to the vial label for reconstitution volume.  | 2 vials  |
| Hyaluronan Detection Antibody | Biotinylated recombinant human Aggrecan lyophilized. Refer to the vial label for reconstitution volume. | 1 vial   |
| PBS (10×)                     | A 10× concentrated PBS  | 1 vial   |
| Streptavidin-HRP B (40×)      | 40× Streptavidin conjugated to horseradish peroxidase.  | 1 vial   |
| Reagent Diluent (5×)          | A 5× concentrated buffered protein base used to dilute Standard, Samples, Detection Antibody and HRP.   | 1 vial   |
| Wash Buffer Concentrate (25×) | A 25× concentrated solution of buffered surfactant with preservatives.                                  | 1 vial   |
| TMB Substrate                 | TMB ELISA Substrate Solution  | 1 vial   |
| Stop Solution                 | 2 N sulfuric acid.  | 1 vial   |
| Plate Sealers                 | Adhesive strip.   | 3 strips |

## B. STORAGE

|                                      |  |   |
|--------------------------------------|--|---|
| Unopened Kit                         | Store at 2-8 °C. Do not use past kit expiration date.  |   |
| Opened/<br>Reconstituted<br>Reagents | Streptavidin-HRP<br>B  | May be stored for up to 1 month at 2-8 °C.*   |
|                                      | Diluted Wash<br>Solution   |   |
|                                      | TMB Substrate  |   |
|                                      | PBS (10×)  |   |
|                                      | Stop Solution  |   |
|                                      | Standard   | Prepare fresh for each assay. Standards may be stored for up to 1 month at -20°C*.  |
|                                      | Detection Antibody   | Aliquot and store for up to 1 month at -20 °C in a manual defrost freezer.*   |
|                                      | Reagent Diluent<br>(5×)  | May be stored for up to 1 month at 2-8 °C.*<br>Use and discard diluted Reagent Diluent (1×).<br>Prepare fresh for each assay. |
| Microplate Wells                     | Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8 °C.* |   |

\* Provided this is within the expiration date of the kit.

### **C. OTHER SUPPLIES REQUIRED**

- ◆ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ◆ Pipettes and pipette tips.
- ◆ Deionized or distilled water.
- ◆ Squirt bottle, manifold dispenser, or automated microplate washer.
- ◆ 500 mL graduated cylinder.
- ◆ An incubator which can provide stable incubation conditions up to 25°C ±0.5°C.

### **D. PRECAUTION**

- ◆ Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- ◆ The Stop Solution provided with this kit is an acid solution. Wear protective gloves, clothing, eye, and face protection. Wash hands thoroughly after handling.

## VI. PREPARATION

### A. SAMPLE COLLECTION AND STORAGE

**Cell Culture Supernates** - Remove particulates by centrifugation and assay immediately or aliquot and store samples at  $\leq -20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles. Samples may require dilution with Reagent Diluent (1 $\times$ ).

**Serum** - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000  $\times$  g. Remove serum and assay immediately or aliquot and store samples at  $\leq -20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles. Samples may require dilution with Reagent Diluent (1 $\times$ ).

### B. SAMPLE PREPARATION

Cell Culture Supernates samples recommend a 10-fold dilution. A suggested 10-fold dilution is 20  $\mu\text{L}$  of sample + 180  $\mu\text{L}$  of Reagent Diluent (1 $\times$ ). Optimal dilutions should be determined by the end user.

Serum samples recommend a 20-fold dilution. A suggested 20-fold dilution is 10  $\mu\text{L}$  of sample + 190  $\mu\text{L}$  of Reagent Diluent (1 $\times$ ). Optimal dilutions should be determined by the end user.

### C. REAGENT PREPARATION

**Note:** *Bring all reagents to room temperature before use.*

**Wash Buffer** - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25 $\times$ ) into deionized or distilled water to prepare 500 mL of Wash Buffer.

**PBS (1 $\times$ )**- Use deionized or distilled water to prepare PBS (1 $\times$ ).

**Reagent Diluent (1 $\times$ )** - Dilute Reagent Diluent (5 $\times$ ) with PBS (1 $\times$ ) to prepare Reagent Diluent (1 $\times$ ).

**Detection Antibody (1 $\times$ )** - **Centrifuge briefly before opening. Reconstitution volume refer to vial label to prepare Detection Antibody (100 $\times$ ).** Allow the Detection Antibody to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions. Aliquot and store if needed. Dilute to Detection Antibody (1 $\times$ ) with Reagent Diluent (1 $\times$ ). Prepare at least 15 minutes prior to use.


**Streptavidin-HRP B (1 $\times$ )** - **Centrifuge briefly before opening.** Dilute to the working concentration specified on the vial label using Reagent Diluent (1 $\times$ ).

**Hyaluronan Standard** - **Centrifuge briefly before opening. Refer to the vial label**

**for the reconstitution volume\***. This reconstitution produces a stock solution of 80 ng/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

\*If you have any question, please seek help from our Technical Support.

**Pipette 300  $\mu$ L of the appropriate Reagent Diluent (1 $\times$ ) into each tube.** Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 40 ng/mL standard serves as the high standard. The Reagent Diluent (1 $\times$ ) serves as the zero standard (0 ng/mL).

#### D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep Substrate Solution protected from light. Substrate Solution should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB substrate.

## VII. ASSAY PROCEDURE

**Note:** *Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.*

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100  $\mu$ L of standard, or prepared sample per well. Cover with the adhesive strip provided. **Incubate for 2 hours in an incubator at 25°C.** A plate layout is provided for a record of standards and samples assayed. (Samples may require dilution. See Sample Preparation section.)
4. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Solution (400  $\mu$ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
5. Add 100  $\mu$ L of the Detection Antibody (1 $\times$ ) in Reagent Diluent (1 $\times$ ), to each well. Cover with a new adhesive strip and **incubate 2 hours in an incubator at 25°C.**
6. Repeat the aspiration/wash as in step 4.
7. Add 100  $\mu$ L of the working dilution of Streptavidin-HRP B to each well. Cover the plate and **incubate for 20 minutes in an incubator at 25°C. Avoid placing the plate in direct light.**
8. Repeat the aspiration/wash as in step 4.
9. Add 100  $\mu$ L of TMB Substrate to each well. **Incubate for 20 minutes in an incubator at 25°C. Avoid placing the plate in direct light.**
10. Add 50  $\mu$ L of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
11. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical

imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

## 12. **CALCULATION OF RESULTS**

Average the duplicate readings for each standard, control, and sample and subtract the average zero standard optical density (O.D.). Create a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the Hyaluronan concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

## VIII. REFERENCES

1. Noble, P.W. et al. (2011) *Physiol. Rev.* 91:221.
2. Weigel, P.H. et al. (1997) *J. Biol. Chem.* 272:13997.
3. Bost, F. et al. (1998) *Eur. J. Biochem.* 252:339.
4. Stern, R. (2004) *Eur. J. Cell Biol.* 83:317.
5. Gaffney, J. et al. (2010) *Mol. BioSyst.* 6:437.
6. Karangelis, D.E. et al. (2010) *Curr. Med. Chem.* 17:4018.
7. Jackson, D.G. (2009) *Immunol. Rev.* 230:216.
8. Itano, N. and K. Kimata (2002) *IUBMB Life* 54:195.
9. Csoka, A.B. et al. (2001) *Matrix Biol.* 20:499.
10. Day, A.J. (2001)  
<http://www.glycoforum.gr.jp/science/hyaluronan/HA16/HA16J.html>.
11. Gariboldi, S. et al. (2008) *J. Immunol.* 181:2103.
12. Scheibner, K.A. et al. (2006) *J. Immunol.* 177:1272.
13. Powell, J.D. and M.R. Horton (2005) *Immunol. Res.* 32:207.
14. Bollyky, P.L. et al. (2009) *J. Leukoc. Biol.* 86:567.
15. Aytekin, M. et al. (2008) *Am. J. Physiol. Lung Cell Mol. Physiol.* 295:L789.
16. Mine, S. et al. (2006) *Endocrine J.* 53:761.
17. Ferguson, E.L. et al. (2001) *Int. J. Pharm.* 420:84.
18. Masuko, K. et al. (2009) *Int. J. Gen. Med.* 2:77.
19. Bollyky, P.L. et al. (2007) *J. Immunol.* 179:744.
20. Sironen, R.K. et al. (2011) *Exp. Cell Res.* 317:3

## PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

|           |          |          |          |          |          |          |          |          |
|-----------|----------|----------|----------|----------|----------|----------|----------|----------|
| <b>1</b>  |          |          |          |          |          |          |          |          |
| <b>2</b>  |          |          |          |          |          |          |          |          |
| <b>3</b>  |          |          |          |          |          |          |          |          |
| <b>4</b>  |          |          |          |          |          |          |          |          |
| <b>5</b>  |          |          |          |          |          |          |          |          |
| <b>6</b>  |          |          |          |          |          |          |          |          |
| <b>7</b>  |          |          |          |          |          |          |          |          |
| <b>8</b>  |          |          |          |          |          |          |          |          |
| <b>9</b>  |          |          |          |          |          |          |          |          |
| <b>10</b> |          |          |          |          |          |          |          |          |
| <b>11</b> |          |          |          |          |          |          |          |          |
| <b>12</b> |          |          |          |          |          |          |          |          |
|           | <b>A</b> | <b>B</b> | <b>C</b> | <b>D</b> | <b>E</b> | <b>F</b> | <b>G</b> | <b>H</b> |

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