PRODUCT INFORMATION & MANUAL

Human IFN- α 2 Valukine™ ELISA

VAL169

For the quantitative determination of natural and recombinant
human IFN- α 2 concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version202305.1

TABLE OF CONTENTS

I. BACKGROUND	2
II. OVERVIEW	3
III. ADVANTAGES	4
IV. EXPERIMENT	7
V. KIT COMPONENTS AND STORAGE	8
VI. PREPARATION	11
VII. ASSAY PROCEDURE	13
VIII. REFERENCES	15

I. BACKGROUND

There are 3 major classes of interferons (IFNs): Type I, Type II and Type III. Interferon alpha (IFN- α), along with IFN- β , IFN- δ , IFN- ϵ , IFN- κ , IFN- ω and IFN- τ are all Type I IFNs (1). The sole type II IFN is IFN- γ . Type III IFNs include IFN- λ 1, IFN- λ 2, IFN- λ 3 and IFN- λ 4 (2). As a part of the innate immune response, Type I IFNs are rapidly induced in response to viral nucleic acids such as double stranded DNA or RNA (dsDNA, dsRNA) and single stranded RNA (ssRNA), viral glycoproteins, microbial cytosine-phosphate-guanosine (CpG) DNA, DNA damage, and chromosomal instability (3,4).

IFN- α subtypes are well described (5-9). There are 16 human IFN- α subtypes with 80% amino acid identity (10). The number of IFN- α subtypes varies by species with 6 equine subtypes, 17 porcine subtypes, 14 bovine subtypes, and 9 canine subtypes known currently. Human IFN- α subtypes include: IFN- α 1a, IFN- α 1b, IFN- α 2a, IFN- α 2b, IFN- α 2c, IFN- α 4a, IFN- α 4b, IFN- α 5, IFN- α 6, IFN- α 7, IFN- α 8, IFN- α 10, IFN- α 14, IFN- α 16, IFN- α 17, and IFN- α 21. Although there is one known heterodimeric IFN- α receptor (IFN- α R, described below), each IFN- α subtype has been correlated with differing biological activities (7). Variability in biological responses can be attributed to differences in binding affinity and duration, receptor density, feedback responses and intracellular characteristics (11). IFN- α responses have been described as robust, especially in the context of viral infection responsiveness by all cells or tunable in a cell type specific manner.

IFN- α signaling is well characterized (2, 5,11,12,13). IFN- α is a ligand for IFN α R, which includes two subunits IFN α R1 and IFN α R2. IFN- α ligand binding to the ubiquitously expressed IFN α R1 triggers a conformational change which allows for the heterodimerization of IFN α R1 and IFN α R2 (10). Heterodimerization results in the cross phosphorylation of the Janus-activated Kinase 1 (JAK1) on IFN α R2 and tyrosine kinase 2 (TYK2) on IFN α R1 respectively, as well as the intracellular domain of IFN α R1 and IFN α R2 (9). The transcription factors Signal transducer and activator of transcription (STAT) 1 and 2 are subsequently recruited to IFN α R via their Src homology 2 (SH2) domain and phosphorylated. In the canonical IFN- α signaling pathway, phosphoSTAT1/STAT2 heterodimers associate with Interferon Regulatory Factor 9 (IRF9) to form Interferon Stimulated Gene Factor 3 (ISGF3), which translocates into the nucleus to bind to gamma-activated sequences (GASs) or interferon-stimulated response elements (ISREs). This stimulates the transcription of interferon-stimulated genes (ISGs). STAT 2 homodimers and monomers can also associate with IRF-9 to form ISGF3-like complex, which also binds to ISRE to stimulate the transcription of ISGs (9).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for human IFN- α 2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any human IFN- α 2 present is bound by the immobilized antibody. A biotinylated detection antibody specific for human IFN- α 2 is pipetted into the wells. After washing away any unbound substances, Streptavidin-HRP is pipetted into the wells. Following a wash to remove any unbound reagent, TMB substrate solution (Chromogenic agent) is added to the wells and color develops in proportion to the amount of human IFN- α 2 bound in the initial step. The color development is stopped, and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- This kit is suitable for cell culture supernates and human serum.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Calibrator Diluent (1 \times) or Calibrator Diluent-S (1 \times) and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

	Intra-assay Precision			Inter-assay Precision		
Sample	1	2	3	1	2	3
Mean (pg/mL)	5.9	22.8	92.1	6.0	24.1	96.7
Standard Deviation	0.2	0.9	3.6	0.2	2.4	7.6
CV%	3.9	3.8	3.9	4.0	9.9	7.8

B. RECOVERY

The recovery of human IFN- α 2 spiked to different levels throughout the range of the assay in cell culture media was evaluated. The recovery ranged from 84.7 to 119.3% with an average of 101.2%.

The recovery of human IFN- α 2 spiked to different levels throughout the range of the assay in human serum was evaluated. The recovery ranged from 84.4 to 122.5% with an average of 101.6%.

C. SENSITIVITY

The minimum detectable dose (MDD) of human IFN- α 2 is typically less than 0.18 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against a highly purified HEK293 expressed recombinant human IFN- α 2 produced at R&D Systems.

E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with high concentrations of human IFN- α 2 and diluted with Calibrator Diluent (1 \times) (*for cell culture supernate samples*) or Calibrator Diluent-S (1 \times) (*for serum samples*) to produce samples with values within the dynamic range of the assay.

Dilution	Average % of Expected	Range (%)
1:2	100.3	92.7-112.1
1:4	102.3	88.0-105.8
1:8	102.9	89.0-111.6
1:16	102.1	87.7-110.4

F. SAMPLE VALUES

Serum - Six human serum samples were evaluated for the presence of human IFN- α 2 in this assay. All samples measured less than the lowest standard, 3.1 pg/mL.

G. SPECIFICITY

The following factors prepared at 50 ng/mL were assayed and exhibited no cross-reactivity or interference.

Recombinant human:		Recombinant mouse:
IFN- α 1A	IFN- α 21	IFN- α 2
IFN- α 1B	IFN- α / β R1	
IFN- α 4B	IFN- α / β R2/Fc Chimera	
IFN- α 5	IFN- β	


Recombinant Human Factor	Concentration Tested (pg/mL)	Observed Value (pg/mL)	% Cross-reactivity
IFN- α 2A	156.3	78	50.2
IFN- α 4A	312.5	44	13.9
IFN- α 6	78.1	32	40.5
IFN- α 7	2500	26	1.0
IFN- α 8	312.5	30	9.6
IFN- α 10	312.5	52	16.8
IFN- α 14	1250	36	2.9
IFN- α 16	312.5	21	6.8
IFN- α 17	312.5	34	10.7

IV. EXPERIMENT

EXAMPLE STANDARD


The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

CELL CULTURE SUPERNATE ASSAY

pg/mL	OD	Average	Corrected
0	0.020 0.022	0.021	-
3.1	0.079 0.081	0.080	0.059
6.3	0.142 0.144	0.143	0.122
12.5	0.269 0.273	0.271	0.250
25	0.531 0.533	0.532	0.511
50	1.049 1.054	1.052	1.031
100	1.942 1.947	1.945	1.924
200	2.764 2.777	2.771	2.750

SERUM ASSAY

pg/mL	OD	Average	Corrected
0	0.019 0.020	0.020	-
3.1	0.059 0.061	0.060	0.041
6.3	0.102 0.102	0.102	0.083
12.5	0.185 0.196	0.191	0.171
25	0.359 0.363	0.361	0.342
50	0.697 0.718	0.708	0.688
100	1.341 1.367	1.354	1.335
200	2.332 2.349	2.341	2.321

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
Human IFN- α 2 Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with a mouse antibody against human IFN- α 2.	1 plate
Human IFN- α 2 Standard	Recombinant human IFN- α 2 in a buffered protein base; lyophilized. Refer to the vial label for reconstitution volume.	2 vials
Human IFN- α 2 Detection Antibody	Biotinylated human IFN- α 2 antibody, lyophilized. Refer to the vial label for reconstitution volume.	1 vial
Calibrator Diluent Concentrate (4 \times)	A 4 \times concentrated buffered diluent used to dilute standard and cell culture supernate samples.	1 vial
Calibrator Diluent-S (1 \times)	A buffered diluent used to dilute standard and serum samples.	2 vials
Detection Antibody Diluent Concentrate (10 \times)	A 10 \times concentrated buffered diluent used to dilute Detection Antibody.	1 vial
Streptavidin-HRP A (200 \times)	200 \times Streptavidin conjugated to horseradish peroxidase.	1 vial
Reagent Diluent Concentrate (10 \times)	A 10 \times concentrated buffered protein base used to dilute HRP.	1 vial
Wash Buffer Concentrate (25 \times)	A 25 \times concentrated solution of buffered surfactant with preservatives.	1 vial
TMB Substrate	TMB ELISA Substrate Solution.	1 vial
Stop Solution	2 N sulfuric acid.	1 vial
Plate Sealers	Adhesive strip.	3 strips

B. STORAGE

Unopened Kit	Store at 2-8 °C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Streptavidin-HRP A	May be stored for up to 1 month at 2-8 °C.*
	Diluted Wash Solution	
	TMB Substrate	
	Stop Solution	
	Standard	Prepare fresh for each assay. Standards may be stored for up to 1 month at -20 °C.*
	Detection Antibody	Aliquot and store for up to 1 month at -20 °C in a manual defrost freezer. *
	Reagent Diluent Concentrate (10×)	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Reagent Diluent (1×). Prepare fresh for each assay.
	Calibrator Diluent Concentrate (4×)	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Calibrator Diluent (1×). Prepare fresh for each assay.
	Calibrator Diluent-S (1×)	May be stored for up to 1 month at 2-8 °C.*
	Detection Antibody Diluent Concentrate (10×)	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Detection Antibody Diluent (1×). Prepare fresh for each assay.
Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8 °C.*	

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- ◆ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ◆ Pipettes and pipette tips.
- ◆ Deionized or distilled water.
- ◆ Squirt bottle, manifold dispenser, or automated microplate washer.
- ◆ 500 mL graduated cylinder.
- ◆ Horizontal orbital microplate shaker capable of maintaining a speed of 500±50 rpm.

D. PRECAUTION

- ◆ Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- ◆ The Stop Solution provided with this kit is an acid solution. Wear protective gloves, clothing, eye, and face protection. Wash hands thoroughly after handling.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples require dilution with Calibrator Diluent (1 \times).

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples require dilution with Calibrator Diluent-S (1 \times).

B. SAMPLE PREPARATION

Serum samples recommend a 4-fold dilution. A suggested 4-fold dilution is 50 μL of sample + 150 μL of Calibrator Diluent-S (1 \times). Optimal dilutions should be determined by the end user.

C. REAGENT PREPARATION

Note: Bring all reagents to room temperature before use.

Wash Buffer (1 \times) - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25 \times) into deionized or distilled water to prepare 500 mL of Wash Buffer (1 \times).

Calibrator Diluent (1 \times) - Use deionized or distilled water to prepare Calibrator Diluent (1 \times).

Detection Antibody Diluent (1 \times) - Use deionized or distilled water to prepare Detection Antibody Diluent (1 \times).

Reagent Diluent (1 \times) - Use deionized or distilled water to prepare Reagent Diluent (1 \times).

Detection Antibody - **Centrifuge briefly before opening. Reconstitution volume refer to vial label to prepare Detection Antibody (100 \times).** Allow the detection antibody to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions. Aliquot and store if needed. Dilute to Detection Antibody (1 \times) with Detection Antibody Diluent (1 \times). Prepare at least 15 minutes prior to use.


Streptavidin-HRP A (1 \times) - **Centrifuge briefly before opening.** Dilute to the working concentration specified on the vial label using Reagent Diluent (1 \times).

Human IFN- α 2 Standard - **Centrifuge briefly before opening. Refer to the vial label**

for the reconstitution volume*. This reconstitution produces a stock solution of 400 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*If you have any question, please seek help from our Technical Support.

Pipette 300 μ L of the appropriate Calibrator Diluent (1 \times) (for cell culture supernate samples) or Calibrator Diluent-S (1 \times) (for serum samples) into each tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 200 pg/mL standard serves as the high standard. The Calibrator Diluent (1 \times) (for cell culture supernate samples) or Calibrator Diluent-S (1 \times) (for serum samples) serves as the zero standard (0 pg/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB Substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB Substrate.

VII. ASSAY PROCEDURE

Note: Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100 μ L of standard, or prepared sample per well. And add 50 μ L of the Detection Antibody (1 \times) diluted in Detection Antibody Diluent, to each well. Cover with the adhesive strip provided. **Incubate for 3 hours at room temperature on a horizontal orbital microplate shaker set at 500 ± 50 rpm.** A plate layout is provided for a record of standards and samples assayed. (Samples may require dilution. See Sample Preparation section.)
4. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Solution (400 μ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
5. Add 100 μ L of the working dilution of Streptavidin-HRP A to each well. Cover the plate and **incubate for 30 minutes at room temperature on a horizontal orbital microplate shaker set at 500 ± 50 rpm. Protect from light.**
6. Repeat the aspiration/wash as in step 4.
7. Add 100 μ L of TMB Substrate to each well. **Incubate for 30 minutes at room temperature on a horizontal orbital microplate shaker set at 500 ± 50 rpm. Protect from light.**
8. Add 50 μ L of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
9. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in

the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

10. **CALCULATION OF RESULTS**

Average the duplicate readings for each standard, control, and sample and subtract the average zero standard optical density (O.D.). Create a standard curve by reducing the data using computer software capable of generating a four-parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the human IFN- α 2 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Murira, A. & A. Lamarre (2016) *Front. Immunol.* 7:1
2. Lazear, H. M. et al. (2019) *Immunity.* 50:909
3. Tomasello, E. et al. (2014) *Front. Immunol.* 5:1
4. Fuertes, M.B. et al. (2013) *Trends Immunol.* 34:67
5. Ivashkiv, L. B. et al. (2014) *Nat. Rev. Immunol.* 14:36
6. Li, S. et al. (2018) *Cell. Phys. And Biochem.* 51:2377
7. Gibbert, K. et al. (2010) *Br. J of Pharm.* 168:1048
8. Moll, H. P. et al. (2011) *Cytokine.* 53:52
9. Provance, O. K. & J. Lewis-Wambi (2019) *Breast Cancer Res.* 21:59
10. Paul, F. et al. (2015) *Gene.* 567:132
11. Schreiber, G. (2017) *JBC.* 292:7285
12. Heise, R. et al. (2015) *PLOS One.* 11:1
13. Blaszczyk, K. et al. (2016) *Cytokine & Growth Factors Rev.* 29:71

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

1											
2											
3											
4											
5											
6											
7											
8											
9											
10											
11											
12											
	A	B	C	D	E	F	G	H			

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