PRODUCT INFORMATION & MANUAL

Human Fibronectin Valukine™ ELISA

VAL179

For the quantitative determination of natural and recombinant human Fibronectin concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version202211.1

TABLE OF CONTENTS

I. BACKGROUND	2
II. OVERVIEW	3
III. ADVANTAGES	4
IV. EXPERIMENT	6
V. KIT COMPONENTS AND STORAGE	7
VI. PREPARATION	10
VII. ASSAY PROCEDURE	12
VIII. REFERENCES	14

I. BACKGROUND

Fibronectin (FN) is a large modular glycoprotein that is found as a polymeric fibrillar network in the extracellular matrix (ECM) and as soluble disulfide-linked dimeric protomers in plasma and other body fluids (1-4). Fibronectin is a ligand for many molecules including fibrin, heparin, chondroitin sulfate, collagen/gelatin, various integrins and growth factors, myocilin, and Apolipoprotein A (4-10). It is involved in multiple cellular processes such as cell adhesion/ migration, blood clotting, morphogenesis, tissue repair, and cell signaling. Fibronectin functions are mediated by the insoluble polymeric fibrils. Conversion of soluble Fibronectin to Fibronectin fibrils in the ECM is initiated by binding to cell surface integrins, resulting in exposure of cryptic epitopes necessary for polymerization (4). Mature human Fibronectin shares 88% aa sequence identity with mouse and rat Fibronectin.

Fibronectin is made up of three types of homologous structural motifs termed FN type I, type II, and type III repeats (11-13). Alternative splicing generates multiple isoforms of Fibronectin which may have insertions of extra type III domains (EDA and EDB) or alteration of the type III connecting segment (IIICS) (3, 13). Differential splicing within the IIICS domain determines the presence of CS1 and CS2 sequences and the sensitivity to proteases (14, 15). The IIICS domain contains two sites (CS1 and CS2) that interact with Integrin $\alpha4\beta1$ (14, 16, 17). The CS1 sequence is not accessible in full length Fibronectin but is exposed by protease digestion, thereby enabling cell adhesion via Integrin $\alpha4\beta1$ (17). A fragment from the first type III repeat, known as Anastellin, binds to Fibronectin and induces the formation of high molecular weight disulfide-linked superfibronectin (sFN) (18-20). sFN exhibits increased protease sensitivity and cellular adhesiveness compared to Fibronectin (18, 20). Anastellin acts in vivo to inhibit angiogenesis, tumor growth, and metastasis (19).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for human Fibronectin has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any human Fibronectin present is bound by the immobilized antibody. After washing away any unbound substances, a biotinylated detection antibody specific for human Fibronectin is pipetted into the wells. After washing away any unbound substances, Streptavidin-HRP is pipetted into the wells. Following a wash to remove any unbound reagent, TMB substrate is added to the wells and color develops in proportion to the amount of human Fibronectin bound in the initial step. The color development is stopped, and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- This kit is suitable for cell culture supernate and human serum.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Calibrator Diluent (1×) and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

Sample	Intra-assay Precision			Inter-assay Precision		
	1	2	3	1	2	3
Mean (ng/mL)	1.6	5.8	24.1	1.5	5.8	24.1
Standard Deviation	0.1	0.2	1.5	0.1	0.2	1.4
CV%	3.3	3.4	6.2	4.7	3.4	5.7

B. RECOVERY

The recovery of human Fibronectin spiked to different levels throughout the range of the assay in cell culture media was evaluated. The recovery ranged from 82.8-121.1% with an average of 101.8%.

C. SENSITIVITY

The minimum detectable dose (MDD) of human Fibronectin is typically less than 0.011 ng/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against human Fibronectin purified from human plasma at R&D Systems.

E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with high concentrations of human Fibronectin and diluted with Calibrator Diluent (1×) to produce samples with values within the dynamic range of the assay.

Dilution	Average % of Expected	Range (%)
1:2	99.9	91.6-109.2
1:4	103.7	85.9-119.7
1:8	104.4	85.9-123.9
1:16	106.3	87.7-121.8

F. SAMPLE VALUES

Cell Culture Supernates - MCF-7 cells were seeded in 100 mL of DMEM/F12 supplemented with 10% fetal bovine serum, 100 U/mL penicillin, and 100 µg/mL streptomycin sulfate for 2 days. The cell culture supernate was assayed for human Fibronectin and measured 18.5 ng/mL.

TF-1 cells were seeded in 100 mL of 1640 supplemented with 10% fetal bovine serum, 100 U/mL penicillin, and 100 µg/mL streptomycin sulfate, 20ng/mL GM-CSF for 2 days. The cell culture supernate was assayed for human Fibronectin and measured 4 ng/mL.

Serum - Six human serum samples were evaluated for the presence of human Fibronectin in this assay. All samples measured ranged from 14.6 to 117.0 µg/mL with an average of 54.9 µg/mL.

G. SPECIFICITY

The following factors prepared at 2000 ng/mL were assayed and exhibited no cross-reactivity or interference.

Recombinant human:	
COMP	Integrin αVβ1
Fibulin 1C	Integrin αVβ3
Galectin-3BP	Integrin αVβ6
Integrin α4β7	Integrin αVβ8
Integrin α5β1	Tenascin R


A sample containing 250 ng/mL of recombinant human Fibronectin (aa 631-705) reads as 4 ng/mL (1.6% cross-reactivity).

A sample containing 1000 ng/mL of bovine Fibronectin reads as 4 ng/mL (0.4% cross-reactivity). Bovine Fibronectin may be present in Bovine Sera.

IV. EXPERIMENT

EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

ng/mL	OD	Average	Corrected
0	0.029 0.029	0.029	-
0.8	0.126 0.129	0.128	0.099
1.6	0.227 0.232	0.230	0.201
3.1	0.411 0.414	0.413	0.384
6.3	0.733 0.736	0.735	0.706
12.5	1.205 1.227	1.216	1.187
25	1.856 1.872	1.864	1.835
50	2.445 2.458	2.452	2.423

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
Human Fibronectin Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with a mouse antibody against human Fibronectin.	1 plate
Human Fibronectin Standard	Human Fibronectin in a buffered protein base; lyophilized. Refer to the vial label for reconstitution volume.	2 vials
Human Fibronectin Detection Antibody	Biotinylated human Fibronectin antibody, lyophilized. Refer to the vial label for reconstitution volume.	1 vial
Calibrator Diluent Concentrate (2×)	Concentrated buffered diluent used to dilute standard and samples.	2 vials
Streptavidin-HRP A (200×)	200× Streptavidin conjugated to horseradish peroxidase.	1 vial
Reagent Diluent Concentrate (10×)	A 10× concentrated buffered protein base used to dilute HRP.	1 vial
Wash Buffer Concentrate (25×)	A 25× concentrated solution of buffered surfactant with preservatives.	1 vial
TMB Substrate	TMB ELISA Substrate Solution.	1 vial
Stop Solution	2 N sulfuric acid.	1 vial
Plate Sealers	Adhesive strip.	3 strips

B. STORAGE

Unopened Kit	Store at 2-8 °C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Streptavidin-HRP A	May be stored for up to 1 month at 2-8 °C.*
	Diluted Wash Solution	
	TMB Substrate	
	Stop Solution	
	Standard	Prepare fresh for each assay. Standards may be stored for up 1 month at -20°C.*
	Detection Antibody	Aliquot and store for up to 1 month at -20 °C in a manual defrost freezer. *
	Reagent Diluent Concentrate (10×)	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Reagent Diluent (1×). Prepare fresh for each assay.
	Calibrator Diluent Concentrate (2×)	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Calibrator Diluent (1×). Prepare fresh for each assay.
Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8 °C.*	

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- ◆ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ◆ Pipettes and pipette tips.
- ◆ Deionized or distilled water.
- ◆ Squirt bottle, manifold dispenser, or automated microplate washer.
- ◆ 500 mL graduated cylinder.

D. PRECAUTION

- ◆ Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- ◆ The Stop Solution provided with this kit is an acid solution. Wear protective gloves, clothing, eye, and face protection. Wash hands thoroughly after handling.
- ◆ Human Fibronectin is detectable in saliva. Take precautionary measures to prevent contamination of kit reagents while running this assay.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1 \times).

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1 \times).

B. SAMPLE PREPARATION

Cell culture supernate samples require a 2-fold dilution prior to the assay. A suggested 2-fold dilution is 100 μL of sample + 100 μL of Calibrator Diluent (1 \times).

Human serum samples require a 2000-fold dilution. For example, add 10 μL of serum into a tube with 90 μL Calibrator Diluent (1 \times) to prepare a 10-fold diluted sample. Mix through and then pipette 10 μL of prepared 10-fold diluted sample into a tube with 90 μL Calibrator Diluent (1 \times) to prepare a 100-fold diluted sample. Mix through and then pipette 10 μL of prepared 100-fold diluted sample into a tube with 190 μL Calibrator Diluent (1 \times) to prepare a 2000-fold diluted sample.

C. REAGENT PREPARATION

Bring all reagents to room temperature before use.

Note: Human Fibronectin are found in saliva. It is recommended that a face mask and gloves be used to protect kit reagents from contamination.

Wash Buffer (1 \times) - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25 \times) into deionized or distilled water to prepare 500 mL of Wash Buffer (1 \times).

Calibrator Diluent (1 \times) - Use deionized or distilled water to prepare Calibrator Diluent (1 \times).

Reagent Diluent (1 \times) - Use deionized or distilled water to prepare Reagent Diluent (1 \times).

Detection Antibody - **Centrifuge briefly before opening. Reconstitution Volume refer to vial label with Reagent Diluent (1 \times).** Aliquot and store if needed. Dilute stock solution in Reagent Diluent (1 \times) to the working concentration of 12.5 ng/mL. Prepare at least 15 minutes prior to use.


Streptavidin-HRP A (1 \times) - **Centrifuge briefly before opening.** Dilute to the working

concentration specified on the vial label using Reagent Diluent (1×).

Human Fibronectin Standard - Centrifuge briefly before opening. Refer to the vial label for the reconstitution volume*. This reconstitution produces a stock solution of 100 ng/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*If you have any question, please seek help from our Technical Support.

Pipette 300 µL of the appropriate Calibrator Diluent (1×) into each tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 50 ng/mL standard serves as the high standard. The Calibrator Diluent (1×) serves as the zero standard (0 pg/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB substrate should remain colorless until added to the plate. Keep TMB substrate protected from light. TMB substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB substrate.

VII. ASSAY PROCEDURE

Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.

Note: Human Fibronectin are found in saliva. It is recommended that a face mask and gloves be used to protect kit reagents from contamination.

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100 μ L of standard, or prepared sample per well. Cover with the adhesive strip provided. Incubate for 2 hours at room temperature. A plate layout is provided for a record of standards and samples assayed. (Samples may require dilution. See Sample Preparation section.)
4. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Solution (400 μ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
5. Add 100 μ L of the Detection Antibody diluted in Reagent Diluent, to each well. Cover with a new adhesive strip and incubate for 2 hours at room temperature.
6. Repeat the aspiration/wash as in step 4.
7. Add 100 μ L of the working dilution of Streptavidin-HRP A to each well. Cover the plate and incubate for 20 minutes at room temperature. **Protect from light.**
8. Repeat the aspiration/wash as in step 4.
9. Add 100 μ L of TMB Substrate to each well. Incubate for 20 minutes at room temperature. **Protect from light.**
10. Add 50 μ L of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
11. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm.

If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

12. **CALCULATION OF RESULTS**

Average the duplicate readings for each standard, control, and sample and subtract the average zero standard optical density (O.D.). Create a standard curve by reducing the data using computer software capable of generating a four-parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the human Fibronectin concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Potts, J.R. and I.D. Campbell (1996) *Matrix Biol.* 15:313.
2. Singh, P. et al. (2010) *Annu. Rev. Cell Dev. Biol.* 26:397.
3. White, E.S. and A.F. Muro (2011) *IUBMB Life* 63:538.
4. Mao, Y. and J.E. Schwarzbauer (2005) *Matrix Biol.* 24:389.
5. Peters, D.M. et al. (2005) *Exp. Cell Res.* 303:218.
6. Edelstein, C. et al. (2005) *J. Lipid Res.* 46:2673.
7. Pytela, R. et al. (1985) *Cell* 40:191.
8. Wayner, E.A. et al. (1989) *J. Cell Biol.* 109:1321.
9. Altroff, H. et al. (2004) *J. Biol. Chem.* 279:55995.
10. Zhu, J. and R.A.F. Clark (2014) *J. Invest. Dermatol.* 134:895.
11. Bernard, M.P. et al. (1985) *Biochemistry* 24:2698.
12. Kornblihtt, A.R. et al. (1983) *Proc. Natl. Acad. Sci. USA* 80:3218.
13. Kornblihtt, A.R. et al. (1985) *EMBO J.* 4:1755.
14. Mould, A.P. et al. (1991) *J. Biol. Chem.* 266:3579.
15. Abe, Y. et al. (2005) *Biochem. Biophys. Res. Commun.* 338:1640.
16. Komoriya, A. et al. (1991) *J. Biol. Chem.* 266:15075.
17. Ugarova, T.P. et al. (1996) *Biochemistry* 35:10913.
18. Morla, A. et al. (1994) *Nature* 367:193.
19. Yi, M. and E. Rouslahti (2001) *Proc. Natl. Acad. Sci. USA* 98:620.
20. Ohashi, T. and H.P. Erickson (2005) *J. Biol. Chem.* 280:39143.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								
	A	B	C	D	E	F	G	H

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