PRODUCT INFORMATION & MANUAL

Human Angiopoietin-2 Valukine™ ELISA

Catalog Number: VAL204

For the quantitative determination of natural and recombinant human
Angiopoietin-2 concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version 202411.1

TABLE OF CONTENTS

I.	BACKGROUND	2
II.	OVERVIEW	3
III.	ADVANTAGES	4
IV.	EXPERIMENT	7
V.	KIT COMPONENTS AND STORAGE	8
VI.	PREPARATION	11
VII.	ASSAY PROCEDURE	13
VIII.	REFERENCES	15

I. BACKGROUND

Angiopoietin-2 (Ang2) is an approximately 70 kDa secreted glycoprotein that plays a complex role in angiogenesis, inflammation, and vascular development. Angiopoietin-2 contains a coiled-coil domain that mediates multimerization and a C-terminal Fibrinogen-like domain that mediates receptor binding (1, 2). It forms disulfide-linked dimers, trimers, tetramers, and pentamers (2-4). Mature human Angiopoietin-2 shares 86% amino acid (aa) sequence identity with mouse and rat Angiopoietin-2. Angiopoietin-2 is expressed in vascular endothelial (EC) and smooth muscle cells in highly angiogenic tissues (e.g. placenta, ovaries, uterus, and tumor vasculature) (5-9), lung epithelial cells, differentiating myotubes, and neural progenitor cells (10-12). It is upregulated by cellular stress (7, 10, 13, 14) and circulates at elevated levels in sepsis and acute lung injury (10, 15).

Both Angiopoietin-2 and the related Angiopoietin-1 (Ang-1) are ligands for the receptor tyrosine kinase Tie-2 (5). Angiopoietin-2 acts as a Tie-2 agonist, although it can also partially antagonize the activation of Tie-2 by Angiopoietin-1 (5, 7, 15-17). Angiopoietin-2 promotes EC survival, proliferation, and migration, and also promotes sprouting angiogenesis (7, 17, 18). It induces the loss of pericytes from vessel walls and an increase of vascular permeability (14, 15, 19-21). Angiopoietin-2 is required for the development of lymphatic vessels as well as the postnatal remodeling of both lymphatic and vascular vessels (6, 19, 22). In cancer, Angiopoietin-2 production is induced by tumor-derived VEGF and augments VEGF-induced angiogenesis and tumor growth (8, 19, 23). In addition, Angiopoietin-2 modulates cell adhesion and promotes tumor cell migration through interactions with Integrins containing $\alpha 3$, $\alpha 5$, αV , $\beta 1$, or $\beta 3$ chains (24-26). It promotes leukocyte adhesion to the vascular endothelium and extravasation to inflammatory sites (10, 13, 19, 27). Aside from the circulatory system, Angiopoietin-2 promotes the differentiation of myoblasts, regulatory T cells, and neurons (11, 12, 28).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for human Angiopoietin-2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any human Angiopoietin-2 present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked antibody specific for human Angiopoietin-2 is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, TMB substrate (Chromogenic agent) is added to the wells and color develops in proportion to the amount of human Angiopoietin-2 bound in the initial step. The color development is stopped, and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- ♦ **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- ♦ This kit is suitable for cell culture supernates, human serum, human plasma and human saliva.
- ♦ The kit should not be used beyond the expiration date on the kit label.
- ♦ Do not mix or substitute reagents with those from other lots or sources.
- ♦ If samples generate values higher than the highest standard, dilute the samples with Calibrator Diluent RD5-5 and repeat the assay.
- ♦ Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in forty separate assays to assess inter-assay precision.

	Intra-assay Precision			Inter-assay Precision		
Sample	1	2	3	1	2	3
Mean (pg/mL)	237	703	1301	276	805	1494
Standard Deviation	16.4	45.6	54.4	28.8	73.0	111
CV%	6.9	6.5	4.2	10.4	9.1	7.4

B. RECOVERY

The recovery of human Angiopoietin-2 spiked to different levels throughout the range of the assay in various matrices was evaluated.

Sample Type	Average % Recovery	Range (%)
Cell culture media (n=4)	103	98-113
Human serum* (n=4)	100	90-107
Human EDTA plasma* (n=4)	95	85-104
Human heparin plasma* (n=4)	94	87-104

*Samples were diluted prior to assay as direct in the Sample Preparation section.

C. SENSITIVITY

Thirty-eight assays were evaluated and the minimum detectable dose (MDD) of human Angiopoietin-2 ranged from 1.20-21.3 pg/mL. The mean MDD was 8.29 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against a highly purified NS0-expressed recombinant human Angiopoietin-2 produced at R&D Systems.

E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with high concentrations of human Angiopoietin-2 and diluted with Calibrator Diluent RD5-5 to produce samples with values within the dynamic range of the assay.

Dilution		Cell culture media (n=4)	Human serum* (n=4)	Human EDTA plasma* (n=4)	Human heparin plasma* (n=4)	Human saliva (n=4)
1:2	Average % of Expected	100	95	97	95	105
	Range (%)	98-102	87-103	93-103	89-101	101-110
1:4	Average % of Expected	97	95	95	94	107
	Range (%)	90-103	86-112	93-98	88-101	100-112
1:8	Average % of Expected	96	95	100	91	107
	Range (%)	87-105	85-114	89-113	87-100	104-112

*Samples were diluted prior to assay as directed in the Sample Preparation section.

F. SAMPLE VALUES

Human Serum/Plasma - Samples from apparently healthy volunteers were evaluated for the presence of human Angiopoietin-2 in this assay. No medical histories were available for the donors used in this study.

Sample Type	Mean (pg/mL)	Range (pg/mL)	Standard Deviation (pg/mL)
Human serum (n=60)	2494	1065-8907	1341
Human EDTA plasma (n=35)	1964	1071-4389	808
Human heparin plasma (n=35)	2049	1009-4973	913

Human saliva - Nine samples were evaluated for the presence of human Angiopoietin-2 in this assay and ranged from 247-822 pg/mL.

Cell Culture Supernates - Two cultures of HUVEC human umbilical vein endothelial cells were grown to confluency in EGM supplemented with 2% fetal bovine serum, 2 mM L-glutamine, 100 U/mL penicillin, 100 µg/mL streptomycin sulfate and bovine brain extract. Aliquots of the cell culture supernates were removed and assayed for levels of human Angiopoietin-2.

Stimulant	Values (pg/mL)
HUVEC Culture 1	2326
HUVEC Culture 2	7611

Note: Significant levels of Angiopoietin-2 are found in fetal bovine, bovine, porcine, equine, and rabbit serum. The background level of Angiopoietin-2 in control medium should be determined and subtracted from samples of conditioned medium.

G. SPECIFICITY

This assay recognizes natural and recombinant human Angiopoietin-2.


The factors listed below were prepared at 50 ng/mL in Calibrator Diluent RD5-5 and assayed for cross-reactivity. Preparations of the following factors at 50 ng/mL in a mid-range recombinant human Angiopoietin-2 standard were also assayed for interference. No significant cross-reactivity or interference was observed.

Recombinant human:	
ANG	VEGF ₁₂₁
Ang-1	VEGF ₁₆₅
Tie-2	VEGF/PIGF

IV. EXPERIMENT

EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
Human Angiopoietin-2 Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with an antibody against human Angiopoietin-2	1 plate
Human Angiopoietin-2 Conjugate	Solution of antibody against human Angiopoietin-2 conjugated to horseradish peroxidase	1 vial
Human Angiopoietin-2 Standard	Recombinant human Angiopoietin-2 in a buffered protein base; lyophilized. Refer to the vial label for reconstitution volume	2 vials
Assay Diluent RD1-76	A buffered protein base	1 vial
Calibrator Diluent RD5-5	A buffered protein base used to dilute standard and samples	1 vial
Wash Buffer Concentrate (25×)	A 25× concentrated solution of buffered surfactant	1 vial
TMB Substrate	TMB ELISA Substrate Solution/TMB Substrate Solution	2 vials
Stop Solution	2 N sulfuric acid	1 vial
Plate Sealers	Adhesive strip	3 strips

B. STORAGE

Unopened Kit	Store at 2-8°C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Wash Buffer (1×)	May be stored for up to 1 month at 2-8°C.*
	Stop Solution	
	Conjugate	
	Assay Diluent RD1-76	
	Calibrator Diluent RD5-5	
	TMB Substrate	
	Standard	Use a fresh standard for each assay. Discard after use.
	Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- ♦ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ♦ Pipettes and pipette tips.
- ♦ Deionized or distilled water.
- ♦ Squir bottle, manifold dispenser, or automated microplate washer.
- ♦ 500 mL graduated cylinder.
- ♦ Horizontal orbital microplate shaker (0.12" orbit) capable of maintaining a speed of 500 ± 50 rpm
- ♦ Collection device for saliva samples that has no protein binding or filtering capabilities such as a Salivette® or equivalent
- ♦ Polypropylene test tubes for dilution of standards and samples.

D. PRECAUTION

- ♦ Angiopoietin-2 is detectable in saliva. Take precautionary measures to prevent contamination of kit reagents while running this assay.
- ♦ Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- ♦ The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

The sample collection and storage conditions listed below are intended as general guidelines. Sample stability has not been evaluated.

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent RD5-5.

Caution: *Significant levels of Angiopoietin-2 are found in fetal bovine, bovine, porcine, equine, and rabbit sera. The background level of Angiopoietin-2 in control medium should be determined and subtracted from samples of conditioned medium.*

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at $1000 \times g$. Remove serum and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent RD5-5.

Plasma - Collect plasma using EDTA or heparin as an anticoagulant. Centrifuge for 15 minutes at $1000 \times g$ within 30 minutes of collection. Assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent RD5-5.

Note: *Citrate plasma has not been validated for use in this assay.*

Saliva - Collect saliva using a collection device such as a Salivette® or equivalent. Assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent RD5-5.

Note: *Saliva collector must not have any protein binding or filtering capabilities.*

B. SAMPLE PREPARATION

Use polypropylene tubes.

Human serum and plasma samples recommend a 5-fold dilution. A suggested 5-fold dilution is 50 μ L of sample + 200 μ L of Calibrator Diluent RD5-5. Optimal dilutions should be determined by the end user.

C. REAGENT PREPARATION

Bring all reagents to room temperature before use.

Note: *High concentrations of Angiopoietin-2 are found in saliva. We recommend using a face mask and gloves to protect kit reagents from contamination.*


Wash Buffer (1 $\times$) - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25 $\times$) into deionized or distilled water to prepare 500 mL

of Wash Buffer (1×).

Human Angiopoietin-2 Standard - Refer to the vial label for the reconstitution volume* Reconstitute the Human Angiopoietin-2 Standard with deionized or distilled water. This reconstitution produces a stock solution of 30000 pg/mL. Mix the standard to ensure complete reconstitution and allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*If you have any question, please seek help from our Technical Support.

Use polypropylene tubes. Pipette 450 µL of the Calibrator Diluent RD5-5 into the 3000 pg/mL tube. Pipette 200 µL into the remaining tubes. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 3000 pg/mL standard serves as the high standard. The appropriate Calibrator Diluent RD5-5 serves as the zero standard (0 pg/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB Substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB Substrate.

VII. ASSAY PROCEDURE

Bring all other reagents and samples to room temperature before use. It is recommended that all standards and samples be assayed in duplicate.

NOTE: *Angiopoietin-2 is detectable in saliva. Take precautionary measures to prevent contamination of kit reagents while running this assay.*

1. Prepare all reagents, working standards and samples as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100 μ L of Assay Diluent RD1-76 to each well.
4. Add 50 μ L of standard and prepared sample per well. Cover with the adhesive strip provided. **Incubate for 2 hours at room temperature on a horizontal orbital microplate shaker (0.12" orbit) set at 500 rpm $\pm$ 50 rpm.** A plate layout is provided for a record of standards and samples assayed.
5. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Buffer (400 μ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
6. Add 200 μ L of Human Angiopoietin-2 Conjugate to each well. Cover with a new adhesive strip. **Incubate for 2 hours at room temperature on a horizontal orbital microplate shaker (0.12" orbit) set at 500 rpm $\pm$ 50 rpm.**
7. Repeat the aspiration/wash as in step 5.
8. Add 200 μ L of TMB Substrate to each well. **Incubate for 30 minutes at room temperature on the benchtop. Protect from light.**
9. Add 50 μ L of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
10. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

11. CALCULATION OF RESULTS

Average the duplicate readings for each standard and sample and subtract the average zero standard optical density (O.D.). Create a standard curve by reducing the data using computer software capable of generating a four-parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the human Angiopoietin-2 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Augustin, H.G. et al. (2009) *Nat. Rev. Mol. Cell Biol.* 10:165.
2. Procopio, W.N. et al. (1999) *J. Biol. Chem.* 274:30196.
3. Kim, K-T. et al. (2005) *J. Biol. Chem.* 280:20126.
4. Xu, Y. and Q. Yu (2001) *J. Biol. Chem.* 276:34990.
5. Maisonpierre, P.C. et al. (1997) *Science* 27:55.
6. Gale, N.W. et al. (2002) *Dev. Cell* 3:411.
7. Daly, C. et al. (2006) *Proc. Natl. Acad. Sci. USA* 103:15491.
8. Zhang, L. et al. (2003) *Cancer Res.* 63:3403.
9. Holash, J. et al. (1999) *Science* 284:1994.
10. Bhandari, V. et al. (2006) *Nat. Med.* 12:1286.
11. Mofarrahi, M. and S.N.A. Hussain (2011) *PLoS ONE* 7:e22882.
12. Liu, X.S. et al. (2009) *J. Biol. Chem.* 284:22680.
13. Tressel, S.L. et al. (2008) *Arterioscler. Thromb. Vasc. Biol.* 28:1989.
14. Hammes, H-P. et al. (2004) *Diabetes* 53:1104.
15. Parikh, S.M. et al. (2006) *PLoS Med.* 3:e46.
16. Teichert-Kuliszewska, K. et al. (2001) *Cardiovasc. Res.* 49:659.
17. Yuan, H.T. et al. (2009) *Mol. Cell. Biol.* 29:2011.
18. Lobov, I.B. et al. (2002) *Proc. Natl. Acad. Sci. USA* 99:11205.
19. Fagiani, E. et al. (2011) *Cancer Res.* 71:5717.
20. Cao, Y. et al. (2007) *Cancer Res.* 67:3835.
21. Nasarre, P. et al. (2009) *Cancer Res.* 69:1324.
22. Dellinger, M. et al. (2008) *Dev. Biol.* 319:309.
23. Yoshiji, H. et al. (2005) *Gut* 54:1768.
24. Carlson, T.R. et al. (2001) *J. Biol. Chem.* 276:26516.
25. Hu, B. et al. (2006) *Cancer Res.* 66:775.
26. Imanishi, Y. et al. (2007) *Cancer Res.* 67:4254.
27. Fiedler, U. et al. (2006) *Nat. Med.* 12:235.
28. Coffelt, S.B. et al. (2011) *J. Immunol.* 186:4183.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

12								
11								
10								
9								
8								
7								
6								
5								
4								
3								
2								
1								
	A	B	C	D	E	F	G	H

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