PRODUCT INFORMATION & MANUAL

Human TIMP-1 Valukine™ ELISA

Catalog Number: VAL213

For the quantitative determination of natural and recombinant human
Tissue Inhibitor of Metalloproteinases 1 (TIMP-1) concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version 202410.1

TABLE OF CONTENTS

I.	BACKGROUND	2
II.	OVERVIEW	3
III.	ADVANTAGES	4
IV.	EXPERIMENT	8
V.	KIT COMPONENTS AND STORAGE	9
VI.	PREPARATION	12
VII.	ASSAY PROCEDURE	14
VIII.	REFERENCES	16

I. BACKGROUND

Matrix Metalloproteinases (MMPs) are zinc-dependent endopeptidases that catalyze degradation of extracellular matrix proteins, thereby controlling such processes as development, tissue remodeling, wound healing and tumor metastasis (1-3). The activity of MMPs is controlled by regulation of expression and secretion, by proteolytic activation of pro-enzymes and by the Tissue Inhibitors of Metalloproteinases (TIMPs) (4, 5). TIMPs form 1:1, non-covalent complexes with MMPs, blocking access of substrates to the MMP catalytic site. TIMPs are highly specific for MMPs in general but not for any particular MMP. Functional specificity is conferred by other characteristics. TIMP-1 is an inducible protein and TIMP-2 is a constitutive protein and both are soluble and widely distributed. TIMP-3 is restricted to the extracellular matrix and TIMP-4 is largely restricted to cardiac tissue. For reviews on MMPs and TIMPs, see references 1-5.

TIMP-1 is a 184 amino acid residue glycosylated protein, though glycosylation is not necessary for activity (6). It has 12 cysteines (conserved among all TIMPs) that form disulfide bonds in a pattern that gives distinct N- and C-terminal domains (7). The N-terminal domain contains sites that bind to the MMP substrate-binding site (8). Binding of TIMP-1 does not leave a peptide bond in position for proteolysis and is not cleaved (5). The TIMP/MMP complex can dissociate to yield enzyme and active TIMP-1 (9). The C-terminal domain binds to an external site on MMPs, increasing overall affinity (5). TIMP-1 binds with high affinity to the inactive pro-MMP-9, forming a complex in which TIMP-1 retains its ability to inhibit the activity of another active MMP via its N-terminal domain (10).

TIMP-1 is widely synthesized by many cells and tissues (4). Transcription of the TIMP-1 gene is induced by pro-inflammatory cytokines (IL-1, IL-6, OSM, LIF and TNF- α), TGF- β 1 and phorbol esters (4, 11). Many physiological functions of TIMP-1 are closely tied to the functions of MMPs, and an improper balance of MMP and TIMP production correlates with pathological conditions such as arthritis, tumor growth and metastasis (4). On the other hand, TIMP-1 was independently discovered as an erythroid potentiating activity (12, 13), an activity that appears to be functionally distinct from MMP inhibitory activity (14). TIMP-1 binds to certain cell lines and is translocated to the nucleus (15). It inhibits apoptosis in B cells (16), further suggesting that it independently functions in multiple ways to support survival and growth of cells in contrast to its function of inhibition of MMPs.

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for human TIMP-1 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any human TIMP-1 present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked antibody specific for human TIMP-1 is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, TMB substrate (Chromogenic agent) is added to the wells and color develops in proportion to the amount of human TIMP-1 bound in the initial step. The color development is stopped, and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- ◆ **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- ◆ This kit is suitable for cell culture supernates, human serum, human plasma and human saliva.
- ◆ The kit should not be used beyond the expiration date on the kit label.
- ◆ Do not mix or substitute reagents with those from other lots or sources.
- ◆ If samples generate values higher than the highest standard, dilute the samples with Calibrator Diluent (1×) and repeat the assay.
- ◆ Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in forty separate assays to assess inter-assay precision.

Sample	Intra-assay Precision			Inter-assay Precision		
	1	2	3	1	2	3
Mean (ng/mL)	0.48	1.27	6.95	0.51	1.28	6.90
Standard Deviation	0.02	0.05	0.35	0.02	0.05	0.34
CV%	4.2	3.9	5.0	3.9	3.9	4.9

B. RECOVERY

The recovery of human TIMP-1 spiked to levels throughout the range of the assay in various matrices was evaluated.

Sample Type	Average % Recovery	Range (%)
Cell culture media (n=5)	102	94-113
Human Serum* (n=5)	98	87-102
Human Heparin Plasma* (n=5)	99	89-108
Human EDTA Plasma* (n=5)	99	93-108
Human Saliva* (n = 4)	105	89-121

*Samples were initially diluted prior to assay as directed in the Sample Preparation section.

C. SENSITIVITY

The minimum detectable dose (MDD) of human TIMP-1 is typically less than 0.08 ng/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against a highly purified NS0-expressed recombinant human TIMP-1 produced at R&D Systems.

E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with high concentrations of human TIMP-1 and diluted with Calibrator Diluent (1×) to produce samples with values within the dynamic range of the assay.

Dilution		Cell culture media (n=10)	Human Serum* (n=9)	Human EDTA plasma* (n=9)	Human Heparin plasma* (n=9)	Human Saliva* (n=8)
1:2	Average % of Expected	100	103	102	104	101
	Range (%)	97-105	100-110	98-108	100-108	94-105
1:4	Average % of Expected	100	103	103	103	100
	Range (%)	96-106	98-111	96-110	95-110	96-109
1:8	Average % of Expected	103	101	100	103	100
	Range (%)	89-110	93-110	92-110	94-111	95-113
1:16	Average % of Expected	104	100	100	100	101
	Range (%)	91-112	93-108	93-108	85-111	87-115

*Samples were diluted prior to assay as directed in the Sample Preparation section.

F. SAMPLE VALUES

Human Serum/Plasma/Saliva - Samples from apparently healthy volunteers were evaluated for the presence of human TIMP-1 in this assay. No medical histories were available for the donors used in this study.

Sample Type	Mean (ng/mL)	Range (ng/mL)	Standard Deviation (ng/mL)
Human Serum (n=60)	190	87-524	72.1
Human Heparin plasma (n=60)	84	39-279	33.1
Human EDTA plasma (n=60)	98	44-304	35.4
Human Saliva (n=4)	121	46-208	82.3

Cell Culture Supernates - Human peripheral blood mononuclear cells (5×10^6 cells/mL) were cultured in RPMI supplemented with 5% fetal bovine serum, 50 μ M β -mercaptoethanol, 2 mM L-glutamine, 100 U/mL penicillin, and 100 μ g/mL streptomycin sulfate. The cells were cultured unstimulated or stimulated with 10 μ g/mL PHA. Aliquots of the culture supernate were removed on days 1 and 5 and assayed for levels of natural human TIMP-1.

Condition	Day 1 (ng/mL)	Day 5 (ng/mL)
Unstimulated	76	163
Stimulated	201	274

G. SPECIFICITY

This assay recognizes natural and recombinant human (rh) TIMP-1.

The factors listed below were prepared at 200 ng/mL in Calibrator Diluent (1 \times) and assayed for cross-reactivity. Preparations of the following factors at 200 ng/mL in a mid-range recombinant human TIMP-1 control were assayed for interference. No significant cross-reactivity or interference was observed.

Recombinant human:
MMP-1
MMP-2
MMP-3
TIMP-2
TIMP-3
TIMP-4


No cross-reactivity was observed with rhMMP-9, but interference was observed at concentrations ≥ 100 ng/mL.

rhMMP-9 Concentration (ng/mL)	Observed TIMP-1 Value (ng/mL)
200	227.6
100	249.8
50	263.4
25	269.8
12.5	287.0
0	276.9

IV. EXPERIMENT

EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

(ng/mL)	O.D.	Average	Corrected
0	0.033 0.034	0.034	—
0.156	0.086 0.087	0.086	0.052
0.313	0.133 0.135	0.134	0.100
0.625	0.235 0.237	0.236	0.202
1.25	0.419 0.423	0.421	0.387
2.5	0.740 0.774	0.757	0.723
5	1.324 1.347	1.336	1.302
10	2.142 2.190	2.166	2.132

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
Human TIMP-1 Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with an antibody against human TIMP-1	1 plate
Human TIMP-1 Conjugate	Solution of antibody against human TIMP-1 conjugated to horseradish peroxidase	1 vial
Human TIMP-1 Standard	Recombinant human TIMP-1 in a buffered protein base; lyophilized. Refer to the vial label for reconstitution volume	1 vial
Assay Diluent RD1X	A buffered protein base	1 vial
Calibrator Diluent Concentrate (5×)/ RD5P	A 5× concentrated buffered protein base used to dilute standard and samples	1 vial
Wash Buffer Concentrate (25×)	A 25× concentrated solution of buffered surfactant	1 vial
TMB Substrate	TMB ELISA Substrate Solution/TMB Substrate Solution	2 vials
Stop Solution	2 N sulfuric acid	1 vial
Plate Sealers	Adhesive strip	3 strips

B. STORAGE

Unopened Kit	Store at 2-8°C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Wash Buffer (1×)	May be stored for up to 1 month at 2-8°C.*
	Stop Solution	
	Conjugate	
	TMB Substrate	
	Standard	May be stored for up to 1 month at 2-8 °C.*
	Assay Diluent RD1X	May be stored for up to 1 month at 2-8 °C.*
	Calibrator Diluent Concentrate (5×)/ RD5P	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Calibrator Diluent (1×). Prepare fresh for each assay.
Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*	

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- ◆ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ◆ Pipettes and pipette tips.
- ◆ Deionized or distilled water.
- ◆ Squir bottle, manifold dispenser, or automated microplate washer.
- ◆ 100 mL and 500 mL graduated cylinder.
- ◆ Horizontal orbital shaker (0.12" orbit) capable of maintaining a speed of 500±50 rpm

D. PRECAUTION

- ◆ Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- ◆ The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

The sample collection and storage conditions listed below are intended as general guidelines. Sample stability has not been evaluated.

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 × g. Remove serum and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

Plasma - Collect plasma using heparin or EDTA as an anticoagulant. Centrifuge for 15 minutes at 1000 × g within 30 minutes of collection. Assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

Note: *Citrate plasma has not been validated for use in this assay.*

Grossly lipemic samples should not be used in this assay.

Saliva - Collect saliva using a collection device such as a Salivette or equivalent. Assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

Note: *Saliva collector cannot have any enzyme binding or filtering capabilities.*

B. SAMPLE PREPARATION

All samples recommend a 100-fold dilution. A suggested 100-fold dilution is 10 µL of sample + 990 µL of Calibrator Diluent (1×). Optimal dilutions should be determined by the end user.

C. REAGENT PREPARATION

Bring all reagents to room temperature before use.

Note: *High levels of TIMP-1 are found in saliva. Take necessary precautions (e.g. mask and gloves) to protect kit reagents.*


Wash Buffer (1×) - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25×) into deionized or distilled water to prepare 500 mL of Wash Buffer (1×).

Calibrator Diluent (1×) - Use deionized or distilled water to prepare Calibrator Diluent (1×).

Human TIMP-1 Standard - Refer to the vial label for the reconstitution volume*. Reconstitute the Human TIMP-1 Standard with deionized or distilled water. This reconstitution produces a stock solution of 10 ng/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*If you have any question, please seek help from our Technical Support.

Pipette 200 μ L of Calibrator Diluent (1 \times) into each tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The undiluted standard stock serves as the high standard (10 ng/mL). Calibrator Diluent (1 \times) serves as the zero standard (0 ng/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB Substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB Substrate.

VII. ASSAY PROCEDURE

Bring all other reagents and samples to room temperature before use. It is recommended that all standards and samples be assayed in duplicate.

Note: *High levels of TIMP-1 are found in saliva. Take necessary precautions (e.g. mask and gloves) to protect kit reagents.*

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100 μ L of Assay Diluent RD1X to each well. *Assay Diluent RD1X may contain a crystalline precipitate. Mix well before and during use.*
4. Add 50 μ L of standard and prepared sample per well. Cover with the adhesive strip provided. **Incubate for 2 hours at room temperature on a horizontal orbital shaker (0.12" orbit) set at 500 \pm 50 rpm.** A plate layout is provided for a record of standards and samples assayed.
5. Aspirate each well and wash, repeating the process two times for a total of three washes. Wash by filling each well with Wash Buffer (400 μ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
6. Add 200 μ L of Human TIMP-1 Conjugate to each well. Cover with a new adhesive strip. **Incubate for 1 hour at room temperature on a horizontal orbital shaker (0.12" orbit) set at 500 \pm 50 rpm.**
7. Repeat the aspiration/wash as in step 5.
8. Add 200 μ L of TMB Substrate to each well. **Incubate for 30 minutes at room temperature. Protect from light.**
9. Add 50 μ L of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
10. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

11. CALCULATION OF RESULTS

Average the duplicate readings for each standard and sample and subtract the average zero standard optical density (O.D.). Create a standard curve by

reducing the data using computer software capable of generating a four-parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the human TIMP-1 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Borkakoti, N. (1998) *Prog. Biophys. Mol. Biol.* 70:73.
2. Massova, I. et al. (1998) *FASEB J.* 12:1075.
3. Bode, W. et al. (1999) *Cell. Mol. Life Sci.* 55:639.
4. Gomez, D.E. et al. (1997) *Eur. J. Cell Biol.* 74:111.
5. Bode, W. et al. (1999) *APMIS* 107:3.
6. Schultz, R.M. et al. (1988) *Cancer Res.* 48:5539.
7. Williamson, R.A. et al. (1990) *Biochem. J.* 268:267.
8. Murphy, G. et al. (1991) *Biochemistry* 30:10362.
9. Murphy, G. et al. (1989) *Biochem. J.* 26:1031.
10. Kolkenbrok, H. et al. (1995) *Biol. Chem. Hoppe-Seyler* 376:495.
11. Borden, P. and R.A. Heller (1997) *Crit. Rev. Eukaryot. Gene Expr.* 7:159.
12. Westbrook, C.A. et al. (1984) *J. Biol. Chem.* 259:9992.
13. Docherty, A.J. et al. (1985) *Nature* 318:66.
14. Chesler, L. et al. (1995) *Blood* 86:4506.
15. Ritter, L.M. et al. (1999) *Biochem. Biophys. Res. Commun.* 257:494.
16. Guedez, L. et al. (1998) *J. Clin. Invest.* 102:2002.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

1									
2									
3									
4									
5									
6									
7									
8									
9									
10									
11									
12									
	A	B	C	D	E	F	G	H	

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