PRODUCT INFORMATION & MANUAL

Human Thrombomodulin/BDCA-3 Valukine™ ELISA

Catalog Number: VAL252

For the quantitative determination of natural and recombinant human
Thrombomodulin concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version 202411.1

TABLE OF CONTENTS

I.	BACKGROUND	2
II.	OVERVIEW	3
III.	ADVANTAGES	4
IV.	EXPERIMENT	8
V.	KIT COMPONENTS AND STORAGE	9
VI.	PREPARATION	12
VII.	ASSAY PROCEDURE	15
VIII.	REFERENCES	17

I. BACKGROUND

Human Thrombomodulin, also known as BDCA-3 and CD141, is a transmembrane protein mainly expressed by vascular endothelial cells. Thrombomodulin is an important component in the anti-coagulation and fibrinolysis system. When the coagulation cascade is activated, prothrombin is converted to thrombin by coagulation factors Va and VIIIa, ultimately leading to fibrin clot formation. Thrombomodulin functions as a cell surface receptor for thrombin. The thrombomodulin-thrombin complex activates protein C to degrade coagulation factors Va and VIIIa, thereby reducing the amount of thrombin generated and inhibiting coagulation. Thrombomodulin-thrombin complex also activates thrombin-activatable fibrinolysis inhibitor (TAFI), creating a carboxypeptidase that inhibits fibrinolysis (1-3).

Thrombomodulin is a type I transmembrane protein, consisting of an N-terminal C-type lectin domain, six contiguous epidermal growth factor (EGF)-like domains, a highly glycosylated serine/threonine-rich domain, a transmembrane segment, and a short C-terminal cytoplasmic tail. The EGF-like domains are critical for thrombin and protein C binding (4). When vascular endothelial cell injury occurs, the extracellular domain of thrombomodulin is released from the cell surface by proteolytic cleavage, giving rise to multiple fragments in the circulation (5, 6). Quantification of soluble thrombomodulin in the circulation therefore provides insight into the extent of vascular endothelial cell damage. It has been documented that patients suffering from cardiovascular diseases, such as acute coronary syndrome, pulmonary thromboembolism, and severe hemorrhage, have increased circulating thrombomodulin concentrations (7-9). Soluble thrombomodulin retains anticoagulant effect and its therapeutic potential has been explored. Preliminary data has shown that it might have value in preventing deep vein thrombosis and in treatment of disseminated intravascular coagulation (10-12).

Thrombomodulin also contributes to other biological processes beyond anticoagulation. It is well established that the protein C system as well as the thrombin-protease activated receptor system are important coagulation factors that play important roles in inflammation. By modulating the activity of these mediators, thrombomodulin exhibits an anti-inflammatory effect (13-14). Thrombomodulin also regulates cell adhesion and proliferation through its lectin-binding domain (15, 16). Additionally, since mouse embryos lacking functional thrombomodulin cannot survive, thrombomodulin might have unique functions during development (17, 18). Furthermore, accumulating experimental evidence has demonstrated that thrombomodulin may be able to suppress the metastatic capacity of tumor cells to confer a protective role in many types of malignancies (19).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for human Thrombomodulin has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any human Thrombomodulin present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked antibody specific for human Thrombomodulin is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, TMB substrate (Chromogenic agent) is added to the wells and color develops in proportion to the amount of human Thrombomodulin bound in the initial step. The color development is stopped, and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- ♦ **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- ♦ This kit is suitable for cell culture supernates, cell lysates, human serum, human plasma and human urine.
- ♦ The kit should not be used beyond the expiration date on the kit label.
- ♦ Do not mix or substitute reagents with those from other lots or sources.
- ♦ If samples generate values higher than the highest standard, dilute the samples with Calibrator Diluent (1×) and repeat the assay.
- ♦ Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in forty separate assays to assess inter-assay precision.

	Intra-assay Precision			Inter-assay Precision		
Sample	1	2	3	1	2	3
Mean (pg/mL)	571	1123	2318	518	1073	2159
Standard Deviation	13.2	31.1	83.9	41.2	76.1	123
CV%	2.3	2.8	3.6	8.0	7.1	5.7

B. RECOVERY

The recovery of human Thrombomodulin spiked to levels throughout the range of the assay in various matrices was evaluated.

Sample Type	Average % Recovery	Range (%)
Cell culture media* (n=4)	96	88-106
Human Serum* (n=4)	89	85-96
Human EDTA plasma* (n=4)	95	86-105
Human Heparin plasma* (n=4)	91	85-96
Human Citrate plasma* (n=4)	94	85-102
Human Urine* (n=4)	99	85-115

*Samples were initially diluted prior to assay as directed in the Sample Preparation section.

C. SENSITIVITY

Forty assays were evaluated and the minimum detectable dose (MDD) of human

Thrombomodulin ranged from 2.91-27.0 pg/mL. The mean MDD was 7.82 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against a highly purified NS0-expressed recombinant human Thrombomodulin produced at R&D Systems.

E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with high concentrations of human Thrombomodulin and diluted with Calibrator Diluent (1×) to produce samples with values within the dynamic range of the assay.

Dilution		Cell culture media* (n=4)	Human Serum* (n=4)	Human EDTA plasma* (n=4)	Human Heparin plasma* (n=4)	Human Citrate plasma* (n=4)	Human Urine* (n=4)
1:2	Average % of Expected	106	98	98	97	99	102
	Range (%)	103-110	93-102	93-101	94-98	98-100	96-112
1:4	Average % of Expected	107	95	100	96	102	107
	Range (%)	102-111	93-98	94-107	93-101	101-102	103-114
1:8	Average % of Expected	111	94	105	97	105	110
	Range (%)	104-115	85-98	93-113	96-99	104-107	103-115
1:16	Average % of Expected	111	88	107	98	106	110
	Range (%)	108-114	85-92	90-113	95-101	102-111	100-119

*Samples were diluted prior to assay as directed in the Sample Preparation section.

F. SAMPLE VALUES

Human Serum/Plasma/Urine - Samples from apparently healthy volunteers were evaluated for the presence of human Thrombomodulin in this assay. No medical histories were available for the donors used in this study.

Sample Type	Mean (pg/mL)	Range (pg/mL)	Standard Deviation (pg/mL)
Human Serum (n=35)	4032	2866-5318	637
Human EDTA plasma (n=35)	4226	2815-5407	671
Human Heparin plasma (n=35)	4092	2805-5479	645
Human Citrate plasma (n=35)	3511	2353-4541	535
Human Urine (n=11)	15,800	3540-29,100	8340

Cell Culture Supernates/Cell Lysates:

HUVEC human umbilical vein endothelial cells were cultured in EGM-2 until confluent.

MG-63 Human osteosarcoma cells were cultured in DMEM supplemented with 10% fetal bovine serum, 2 mM L-glutamine, 100 U/mL penicillin, and 100 µg/mL streptomycin sulfate for 4 days.

U-87 MG human glioblastoma/astrocytoma cells were cultured in MEM supplemented with 10% fetal bovine serum, 2 mM L-glutamine, 100 U/mL penicillin, 1.0 mM sodium pyruvate, and 100 µg/mL streptomycin sulfate until confluent.

Cell Line	Cell culture supernates (pg/mL)	Cell lysates (pg/mg cell lysate)
HUVEC	ND	5022
MG-63	ND	670
U-87 MG	186	817

ND=Non detectable

G. SPECIFICITY

This assay recognizes natural and recombinant human Thrombomodulin. This assay also recognizes Thrombomodulin complexed with Thrombin, Protein C, and Thrombin/Protein C.


The factors listed below were prepared at 50 ng/mL in Calibrator Diluent (1×) and assayed for cross-reactivity. Preparations of the following factors at 50 ng/mL in a mid-range recombinant human Thrombomodulin control were assayed for interference. No significant cross-reactivity or interference was observed.

Recombinant human:	Recombinant mouse:	Natural proteins:
Coagulation Factor VII	Coagulation Factor VII	human LDL
Coagulation Factor Xa	Coagulation Factor XI	human HDL
Coagulation Factor XI	Thrombomodulin	human VLDL
Protein C		
Protein S		
Thrombin		

IV. EXPERIMENT

EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

(pg/mL)	O.D.	Average	Corrected
0	0.030 0.030	0.030	—
62.5	0.075 0.076	0.076	0.046
125	0.125 0.127	0.126	0.096
250	0.222 0.222	0.222	0.192
500	0.423 0.429	0.426	0.396
1000	0.801 0.822	0.812	0.782
2000	1.555 1.561	1.558	1.528
4000	2.733 2.832	2.783	2.753

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
Human Thrombomodulin Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with an antibody against human Thrombomodulin	1 plate
Human Thrombomodulin Conjugate	Solution of antibody against human Thrombomodulin conjugated to horseradish peroxidase	1 vial
Human Thrombomodulin Standard	Recombinant human Thrombomodulin in a buffered protein; lyophilized. Refer to the vial label for reconstitution volume	1 vial
Assay Diluent RD1X	A buffered protein base	1 vial
Calibrator Diluent Concentrate (5×)/ RD5P	A 5× concentrated buffered protein base used to dilute standard and samples	1 vial
Wash Buffer Concentrate (25×)	A 25× concentrated solution of buffered surfactant	1 vial
TMB Substrate	TMB ELISA Substrate Solution/TMB Substrate Solution	2 vials
Stop Solution	2 N sulfuric acid	1 vial
Plate Sealers	Adhesive strip	3 strips

B. STORAGE

Unopened Kit	Store at 2-8°C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Wash Buffer (1×)	May be stored for up to 1 month at 2-8°C.*
	Stop Solution	
	Conjugate	
	TMB Substrate	
	Standard	May be stored for up to 1 month at 2-8 °C.*
	Assay Diluent RD1X	May be stored for up to 1 month at 2-8 °C.*
	Calibrator Diluent Concentrate (5×)/ RD5P	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Calibrator Diluent (1×). Prepare fresh for each assay.
	Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- ♦ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ♦ Pipettes and pipette tips.
- ♦ Deionized or distilled water.
- ♦ Squirt bottle, manifold dispenser, or automated microplate washer.
- ♦ 500 mL graduated cylinder.
- ♦ Test tubes for dilution of standards and samples.
- ♦ If using cell lysate samples, the following is also required:
 - PBS
 - Cell Lysis Buffer 1 (R&D Systems, Catalog # 890713).

D. PRECAUTION

- ♦ Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- ♦ The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material.
- ♦ Human Thrombomodulin is detectable in human saliva. Take precautionary measures to prevent contamination of kit reagents while running this assay.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

The sample collection and storage conditions listed below are intended as general guidelines. Sample stability has not been evaluated.

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

Cell Lysates - Cells must be lysed prior to assay. Refer to the Cell Lysis Procedure. Samples may require dilution with Calibrator Diluent (1×).

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 × g. Remove serum and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

Plasma - Collect plasma using EDTA, heparin, or citrate as an anticoagulant. Centrifuge for 15 minutes at 1000 × g within 30 minutes of collection. Assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

Urine - Aseptically collect the first urine of the day (mid-stream), voided directly into a sterile container. Centrifuge to remove particulate matter, and assay immediately or aliquot and store at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

B. CELL LYSIS PROCEDURE

Use the following procedure for the preparation of cell lysate samples.

1. Gently wash cells with cold PBS. Pour off and discard the PBS.
2. Add sufficient Cell Lysis Buffer 1 to cover the cells (for example, 10 mL if using a T75 flask) and incubate at room temperature for 30 minutes with gentle agitation.
3. Collect the cell lysates and centrifuge at 12,000 rpm for 10 minutes to remove cell debris.
4. Collect the supernate and determine the total protein concentration of the lysate using the Bradford method (20).
5. Aliquot the lysis supernatant and store at ≤ -80 °C until ready for use.

C. SAMPLE PREPARATION

Cell culture supernate and cell lysate samples recommend a 2-fold dilution. A suggested 2-fold dilution is 100 μ L of sample + 100 μ L of Calibrator Diluent (1 $\times$). Optimal dilutions should be determined by the end user.

Human serum, plasma, and urine samples recommend a 10-fold dilution. A suggested 10-fold dilution is 20 μ L of sample + 180 μ L of Calibrator Diluent (1 $\times$). Optimal dilutions should be determined by the end user.

D. REAGENT PREPARATION

Bring all reagents to room temperature before use.


Wash Buffer (1 $\times$) - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25 $\times$) into deionized or distilled water to prepare 500 mL of Wash Buffer (1 $\times$).

Calibrator Diluent (1 $\times$) - Use deionized or distilled water to prepare Calibrator Diluent (1 $\times$).

Human Thrombomodulin Standard - Refer to the vial label for the reconstitution volume*. Reconstitute the Human Thrombomodulin Standard with deionized or distilled water. This reconstitution produces a stock solution of 40000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*If you have any question, please seek help from our Technical Support.

Pipette 900 μ L of Calibrator Diluent (1 $\times$) into the 4000 pg/mL tube. Pipette 500 μ L into the remaining tubes. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 4000 pg/mL standard serves as the high standard. Calibrator Diluent (1 $\times$) serves as the zero standard (0 pg/mL).

E. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB Substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB Substrate.

VII. ASSAY PROCEDURE

Bring all other reagents and samples to room temperature before use. It is recommended that all standards and samples be assayed in duplicate.

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100 μ L of Assay Diluent RD1X to each well. *May contain crystals. Warm to room temperature and mix well to dissolve.*
4. Add 50 μ L of standard and prepared sample per well. Cover with the adhesive strip provided. **Incubate for 2 hours at room temperature.** A plate layout is provided for a record of standards and samples assayed.
5. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Buffer (400 μ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
6. Add 200 μ L of Human Thrombomodulin Conjugate to each well. Cover with a new adhesive strip. **Incubate for 2 hours at room temperature.**
7. Repeat the aspiration/wash as in step 5.
8. Add 200 μ L of TMB Substrate to each well. **Incubate for 30 minutes at room temperature. Protect from light.**
9. Add 50 μ L of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
10. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

11. CALCULATION OF RESULTS

Average the duplicate readings for each standard and sample and subtract the average zero standard optical density (O.D.).

Create a standard curve by reducing the data using computer software capable of generating a four-parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the human Thrombomodulin concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Davie, E.W. et al. (1991) *Biochemistry* 30:10363.
2. Esmon, C.T. (1987) *Science* 235:1348.
3. Binette, T.M. et al. (2007) *Blood* 110:3168.
4. Dittman, W.A. and P.W. Majerus (1990) *Blood* 75:329.
5. Ishii, H. and P.W. Majerus (1985) *J. Clin. Invest.* 76:2178.
6. Ohlin, A.K. et al. (2005) *J. Thromb. Haemost.* 3:976.
7. Takano, S. et al. (1990) *Blood* 76:2024.
8. Salomaa, V. et al. (1999) *Lancet* 353:1729.
9. Jansson, J.H. et al. (1997) *Circulation* 96:2938.
10. Kurosawa, S. et al. (1987) *J. Biol. Chem.* 262:2206.
11. Kearon, C. et al. (2005) *J. Thromb. Haemost.* 3:962.
12. Saito, H. et al. (2007) *J. Thromb. Haemost.* 5:31.
13. Bernard, G.R. et al. (2001) *N. Engl. J. Med.* 344:699.
14. Ossovskaya, V.S. and N.W. Bunnett (2004) *Physiol. Rev.* 84:579.
15. Huang, H.C. et al. (2003) *J. Biol. Chem.* 278:46750.
16. Zhang, Y. et al. (1998) *J. Clin. Invest.* 101:1301.
17. Healy, A.M. et al. (1995) *Proc. Natl. Acad. Sci. USA* 92:850.
18. Conway, E.M. et al. (1999) *Blood* 93:3442.
19. Hanly, A.M. and D.C. Winter (2007) *Semin. Thromb. Hemost.* 33:673.
20. Bradford, M.M. (1976) *Anal. Biochem.* 72:248.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

12								
11								
10								
9								
8								
7								
6								
5								
4								
3								
2								
1								
	A	B	C	D	E	F	G	H

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