## **PRODUCT INFORMATION & MANUAL**

**Mouse IL-6 Valukine™ ELISA**

**VAL604G**

For the quantitative determination of natural and recombinant  
Mouse IL-6 concentrations

For research use only.  
Not for diagnostic or therapeutic procedures.

**Bio-Techne China Co. Ltd**

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001  
**info.cn@bio-techne.com**

Please refer to the kit label for expiry date.  
Novus kits are guaranteed for 3 months from date of receipt

Version202209.1

# TABLE OF CONTENTS

I. BACKGROUND .....	2
II. OVERVIEW .....	4
III. ADVANTAGES .....	5
IV. EXPERIMENT .....	8
V. KIT COMPONENTS AND STORAGE .....	9
VI. PREPARATION .....	11
VII. ASSAY PROCEDURE .....	13
VIII. REFERENCES .....	15

## I. BACKGROUND

Interleukin 6 (IL-6) is a pleiotropic  $\alpha$ -helical 22-28 kDa phosphorylated and variably glycosylated cytokine that plays important roles in acute phase reactions, inflammation, hematopoiesis, bone metabolism, and cancer progression (1-5). Mature Mouse IL-6 is 187 amino acids (aa) in length and shares 42% and 85% as sequence identity with human and rat IL-6, respectively (6). Alternate splicing generates several isoforms with internal deletions some of which exhibit antagonistic properties (7-10). Cells known to express IL-6 include CD8<sup>+</sup> T cells, fibroblasts, synoviocytes, adipocytes, osteoblasts, megakaryocytes, endothelial cells (under the influence of endothelin), sympathetic neurons, cerebral cortex neurons, adrenal medulla chromaffin cells, retinal pigment cells, mast cells, keratinocytes, Langerhans cells, fetal and adult astrocytes, neutrophils, monocytes, eosinophils, colonic epithelial cells, B1 B cells and pancreatic islet beta cells (2,7,10-33). IL-6 production is generally correlated with cell activation and is normally kept in control by glucocorticoids, catecholamines and secondary sex steroids (2). Normal human circulating IL-6 is in the 1 pg/mL range, with slight elevations during the menstrual cycle, modest elevations in certain cancers, and large elevations after surgery (34-38).

IL-6 induces signaling through a cell surface heterodimeric receptor complex composed of a ligand binding subunit (IL-6R) and a signal transducing subunit (gp130). IL-6 binds to IL-6R, triggering IL-6R association with gp130 and gp130 dimerization (39). Gp130 is also a component of the receptors for CLC, CNTF, CT-1, IL-11, IL-27, LIF, and OSM (40). Soluble forms of IL-6R are generated by both alternate splicing and proteolytic cleavage (3). In a mechanism known as trans-signaling, complexes of soluble IL-6 and IL-6R elicit responses from gp130-expressing cells that lack cell surface IL-6R (1, 3). Trans-signaling enables a wider range of cell types to respond to IL-6, as the expression of gp130 is ubiquitous, while that of IL-6R is predominantly restricted to hepatocytes, monocytes, and resting lymphocytes (1-3). Soluble splice forms of gp130 block trans-signaling from IL-6/IL-6R but not from other cytokines that use gp130 as a co-receptor (3, 41).

IL-6, along with TNF- $\alpha$  and IL-1, drives the acute inflammatory response, is almost solely responsible for fever and the acute phase response in the liver, and is important

in the transition from acute inflammation to either acquired immunity or chronic inflammatory disease (1-4). It contributes to chronic inflammation in conditions such as obesity, insulin resistance, inflammatory bowel disease, inflammatory arthritis and sepsis when dysregulated, often involving IL-6 trans-signaling (1, 2). It also plays an important role in the differentiation of naive T cells to Th17 inflammatory cells in the presence of TGF- $\beta$ . IL-6 modulates bone resorption and is a major effector of inflammatory joint destruction in rheumatoid arthritis through its promotion of Th17 T cell activity (1). It contributes to atherosclerotic plaque development and destabilization (2). However, IL-6 can also have anti-inflammatory effects, such as in skeletal muscle where it is secreted in response to exercise (2). It promotes hematopoiesis by being a growth factor for hematopoietic stem cells, induces B cell maturation to plasma cells and perpetuates multiple myeloma (1, 42). IL-6 also promotes, but probably does not initiate, other types of inflammation-associated carcinogenesis, such as colitis-associated cancer (1).

## II. OVERVIEW

### A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. A monoclonal antibody specific for mouse IL-6 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any mouse IL-6 present is bound by the immobilized antibody. After washing away any unbound substances, a biotinylated detection antibody specific for mouse IL-6 is pipetted into the wells. After washing away any unbound substances, streptavidin-HRP is pipetted into the wells. Following a wash to remove any unbound reagent, TMB substrate is added to the wells and color develops in proportion to the amount of mouse IL-6 bound in the initial step. The color development is stopped, and the intensity of the color is measured.

### B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- This kit is suitable for cell culture supernate, mouse serum and plasma.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Calibrator Diluent (1×) and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

### III. ADVANTAGES

#### A. PRECISION

##### **Intra-assay Precision** (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

##### **Inter-assay Precision** (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

	Intra-assay Precision			Inter-assay Precision		
Sample	1	2	3	1	2	3
Mean (pg/mL)	245.8	62.4	15.2	247.7	61.9	15.1
Standard Deviation	10.4	2.8	0.7	9.7	2.6	0.6
CV%	4.2	4.5	4.5	3.9	4.2	4.0

#### B. RECOVERY

The recovery of mouse IL-6 spiked to levels throughout the range of the assay in various matrices was evaluated.

Sample Type	Average % Recovery	Range %
Cell culture media (n=5)	102.7	95.0-109.6
Serum (n=5)	99.8	94.9-106.8
Plasma (n=3)	89.8	82.8-95.6

#### C. SENSITIVITY

The minimum detectable dose (MDD) of mouse IL-6 is typically less than 0.62 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

## D. CALIBRATION

This immunoassay is calibrated against highly purified *E. coli*-expressed recombinant mouse IL-6 produced at R&D Systems. The NIBSC Standard for Murine IL-6 (93/730), which was intended as a potency standard, was evaluated in this kit.

The dose response curve of the NIBSC Standard (93/730) parallels the Valukine standard curve. To convert sample values obtained with the Valukine Mouse IL-6 kit to approximate NIBSC 93/730 units, use the equation below.

NIBSC (93/730) approximate value (IU/mL) = 1.30 × Valukine Mouse IL-6 value (pg/mL)

## E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with high concentrations of mouse IL-6 and diluted with Calibrator Diluent 1× to produce samples with values within the dynamic range of the assay.

Dilution	Average % of Expected	Range (%)
1:2	101.1	92.0-110.2
1:4	101.5	90.9-113.9
1:8	99.4	90.6-113.0
1:16	98.3	88.5-113.5

## F. SAMPLE VALUES

**Cell Culture Supernates** - Two spleen organ tissues from a mouse were homogenized and seeded in 100 mL of RPMI1640 supplemented with 10% fetal bovine serum, 2 mM L-glutamine, 100 U/mL penicillin, 100 µg/mL streptomycin sulfate, and 10 µg/mL ConA for 2 days. The cell culture supernate was assayed for mouse IL-6 and measured 194 pg/mL.

EL-4 cells (mouse thymoma) were seeded at  $2 \times 10^5$  cell/mL and cultured for 4 days in 100 mL of DMEM supplemented with 10% horse serum, 10 µg/mL PHA and 10 ng/mL PMA. The cell culture supernate was assayed for mouse IL-6 and measured 317 pg/mL.

BALB/3T3 cells were seeded in 100 mL of DMEM supplemented with 10% fetal bovine serum for 2 days. The cell culture supernate was assayed for mouse IL-6 and measured 157.1 pg/mL.

**Serum** - Five serum samples were evaluated for the presence of IL-6 in this assay. All samples measured ranged from 45.7 to 94.2 pg/mL with an average of 65.7 pg/mL.

**Plasma** - Three plasma samples were evaluated for the presence of IL-6 in this assay. All samples measured ranged from 35.3 to 109.1 pg/mL with an average of 81.6 pg/mL.

## G. SPECIFICITY

This assay recognizes both natural and recombinant mouse IL-6. The following factors were prepared at 50 ng/mL and assayed for cross-reactivity. Preparations of the following factors at 50 ng/mL in a mid-range recombinant mouse IL-6 control were assayed for interference. No significant cross-reactivity or interference was observed.

<b>Recombinant Mouse</b>	<b>Recombinant Human</b>
CT-1	IL-6
GP130	IL-6R
IL11	
LIF	
OSM	
IL-6R	


Recombinant Rat IL-6 cross-reacts approximately 0.24% in this assay.

## IV. EXPERIMENT

### EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

pg/mL	OD	Average	Corrected
0	0.053	0.053	-
7.8	0.152	0.149	0.097
15.6	0.242	0.245	0.192
31.3	0.435	0.434	0.381
62.5	0.826	0.822	0.769
125	1.522	1.524	1.471
250	2.500	2.478	2.425
500	2.987	2.987	2.934

## V. KIT COMPONENTS AND STORAGE

### A. MATERIALS PROVIDED

Parts	Description	Size
Mouse IL-6 Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with monoclonal antibody against mouse IL-6.	1 plate
Mouse IL-6 Detection Antibody	Biotinylated mouse IL-6 antibody, lyophilized. Refer to the vial label for reconstitution volume.	1 vial
Mouse IL-6 Standard	Recombinant mouse IL-6 in a buffered protein base, lyophilized. Refer to the vial label for reconstitution volume.	2 vials
Calibrator Diluent Concentrate (4×)	A 4× concentrated buffered protein base used to dilute standard, samples and Detection Antibody.	1 vial
Reagent Diluent Concentrate (10×)	A 10× concentrated buffered protein base used to dilute HRP.	1 vial
Streptavidin-HRP A (200×)	200× Streptavidin conjugated to horseradish peroxidase.	1 vial
Wash Buffer Concentrate (25×)	A 25× concentrated solution of buffered surfactant.	1 vial
TMB Substrate	TMB ELISA Substrate Solution.	1 vial
Stop Solution	Diluted hydrochloric acid solution.	1 vial
Plate Covers	Adhesive strip.	3 strips

## B. STORAGE

Unopened Kit	Store at 2-8 °C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Streptavidin-HRP A	May be stored for up to 1 month at 2-8 °C.*
	Diluted Wash Solution	
	TMB Substrate	
	Stop Solution	
	Standard	Prepare fresh for each assay.
	Detection Antibody	Aliquot and store for up to 1 month at -20 °C in a manual defrost freezer. *
	Calibrator Diluent Concentrate (4×)	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Calibrator Diluent (1×). Prepare fresh for each assay.
	Reagent Diluent Concentrate (10×)	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Reagent Diluent (1×). Prepare fresh for each assay.
Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8 °C.*	

\* Provided this is within the expiration date of the kit.

## C. OTHER SUPPLIES REQUIRED

- Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- Pipettes and pipette tips.
- Deionized or distilled water.
- Squirt bottle, manifold dispenser, or automated microplate washer.
- 500 mL graduated cylinder.

## D. PRECAUTION

- The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material.

## VI. PREPARATION

### A. SAMPLE COLLECTION AND STORAGE

**Cell Culture Supernates** - Remove particulates by centrifugation and assay immediately or aliquot and store samples at  $\leq -20$  °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

**Serum** - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at  $\leq -20$  °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

**Plasma** - Collect plasma using EDTA, heparin, or citrate as an anticoagulant. Centrifuge for 15 minutes at 1000 x g within 30 minutes of collection. Assay immediately or aliquot and store samples at  $\leq -20$  °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

### B. SAMPLE PREPARATION

Cell Culture Supernate, serum and plasma samples require a 2-fold dilution. A suggested 2-fold dilution is 100  $\mu$ L of sample + 100  $\mu$ L of Calibrator Diluent (1×).

### C. REAGENT PREPARATION

**Note:** *Bring all reagents to room temperature before use.*

**Wash Solution** - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25×) into deionized or distilled water to prepare 500 mL of Wash Buffer.

**Calibrator Diluent (1×)** - Use deionized or distilled water to prepare Calibrator Diluent (1×).

**Reagent Diluent (1×)** - Use deionized or distilled water to prepare Reagent Diluent (1×).

**Detection Antibody** - **Centrifuge briefly before opening. Reconstitution volume refer to vial label with Calibrator Diluent (1×).** Aliquot and store if needed. Dilute stock solution in Calibrator Diluent (1×) to the working concentration of 0.5  $\mu$ g/mL. Prepare at least 15 minutes prior to use.


**Streptavidin-HRP A (1×)** - **Centrifuge briefly before opening.** Dilute to the working concentration specified on the vial label using Reagent Diluent (1×).

**Mouse IL-6 Standard** - **Centrifuge briefly before opening. Refer to the vial label**

**for reconstitution volume\***. This reconstitution produces a stock solution of 5000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

\*If you have any question, please seek help from our Technical Support.

Pipette 900  $\mu\text{L}$  of Calibrator Diluent ( $1\times$ ) into the 500 pg/mL tube. Pipette 500  $\mu\text{L}$  of the appropriate Calibrator Diluent ( $1\times$ ) into the remaining tubes. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 500 pg/mL standard serves as the high standard. The Calibrator Diluent ( $1\times$ ) serves as the zero standard (0 pg/mL).

#### D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB Substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB Substrate.

## VII. ASSAY PROCEDURE

**Note:** *Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.*

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100  $\mu\text{L}$  of standard, or prepared sample per well. Cover with the adhesive strip provided. Incubate for 2 hours at room temperature. A plate layout is provided for a record of standards and samples assayed. (Samples may require dilution. See Sample Preparation section.)
4. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Solution (400  $\mu\text{L}$ ) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
5. Add 100  $\mu\text{L}$  of the Detection Antibody diluted in Calibrator Diluent (1x), to each well. Cover with a new adhesive strip and incubate 2 hours at room temperature.
6. Repeat the aspiration/wash as in step 4.
7. Add 100  $\mu\text{L}$  of the working dilution of Streptavidin-HRP A to each well. Cover the plate and incubate for 30 minutes at room temperature. **Protect from light.**
8. Repeat the aspiration/wash as in step 4.
9. Add 100  $\mu\text{L}$  of TMB Substrate to each well. Incubate for 20 minutes at room temperature. **Protect from light.**
10. Add 50  $\mu\text{L}$  of Stop Solution to each well. The color in the wells should change from blue to yellow. If the color in the wells is green or if the color change does not appear uniform, gently tap the plate to ensure thorough mixing.
11. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical

imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

12. **CALCULATION OF RESULTS.** Average the duplicate readings for each standard, control, and sample and subtract the average zero standard optical density. Create a standard curve by reducing the data using computer software capable of generating a four-parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the Mouse IL-6 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data. If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

## VIII. REFERENCES

1. Naugler, W.E. and M. Karin (2008) *Trends Mol. Med.* 14:109.
2. Schuett, H. et al. (2009) *Thromb. Haemost.* 102:215.
3. Jones, SA (2005) *J. Immunol.* 175:3468.
4. Hodge, D.R. et al. (2005) *Eur. J. Cancer*41:2502.
5. Rose-John, S. et al. (2006) *J. leukoc. Biol.* 80:227.
6. Van Snick, J. et al. (1988) *Eur. J. Immunol.*18:193.
7. Kestler, D.P. et al. (1995) *Blood* 86:4559.
8. Kestler, D.P. et al. (1999) *Am. J. Hematol.* 61:169.
9. Bihl, M.P. at al. (2002) *Am. J. Respir. Cell Mol. Biol.* 27:48.
10. Alberti, L. et al. (2005) *Cancer Res.* 65:2.
11. May, L.T. et al. (1986) *Proc. Natl. Acad. Sci. USA* 83:8957.
12. Sad, S. et al. (1995) *Immunity* 2:271.
13. Cichy, J. et Bf. (1996) *Biochem. Biophys. Res. Commun.* 227:318
14. Miyazawa, K. et al. (1998) *Am. J. Pathol.* 152:793.
15. Fried, S.K. et al. (1998) *Endocrinology* 83:847.
16. Ishimi, Y. et al. (1990) *J. Immunol.*145:3297.
17. Jiang, S. et al. (1994) *Blood* 84:4151.
18. Xin, X. et al. (1995) *Endocrinology* 136:132.
19. Marz, P. et al. (1998) *Proc. Natl. Acad. Sci. USA*95:3251.
20. Ringheim, G.E. et al. (1995) *J. Neuroimmunol.* 63:113.
21. Gadiant, RA et al. (1995) *Neurosci. Lett.*194:17.
22. Kuppner, M.C. et al. (1995} *Immunology* 84:265.
23. Gagari, E, et al. (1997) *Blood* 89:2654.
24. Cumberbatch, M. et al. (1996) *Immunology* 87:513.
25. Fujisawa, H. et al. (1997) *J. Interferon Cytokine Res.* 17:347.
26. Lee, S.C. et al. (1993) *J. Immunol.* 150:2659.
27. Lafortune, L. et al. (1996) *J. Neuropathol. Exp. Neural.* 55:515.

28. Ericson, S.G. et al. (1998) *Blood* 91:2099.
29. Melani, C. et al. (1993) *Blood* 81:2744.
30. Lacy, P. et al. (1998) *Blood* 91:2508.
31. Jung, H.C. et al. (1995) *J. Clin. Invest.* 95:55.
32. Spencer, N.F.L. and RA Daynes (1997) *Int. Immunol.* 9:745.
33. Campbell, I.L. et al. (1989) *J. Immunol.* 143:1188.
34. D'Auria, L. et al. (1997) *Eur. Cytokine Netw.* 8:383.
35. Yamamura, M. et al. (1998) *Br. J. Haematol.* 100:129.
36. Angstwurm, MW.A. et al. (1997) *Cytokine* 9:370.
37. Mouawad, R. et al. (1996) *Clin. Cancer Res.* 2:1405.
38. Sakamoto, K. et al. (1994) *Cytokine* 6:181.
39. Murakami, M. et al. (1993) *Science* 260:1808.
40. Muller-Newen, G. (2003) *Sci. STKE* 2003:PE40.
41. Mitsuyama, K. et al. (2006) *Clin. Exp. Immunol.* 143:125.
42. Cerutti, A et al. (1998) *J. Immunol.* 160:2145.

# PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

<b>1</b>									
<b>2</b>									
<b>3</b>									
<b>4</b>									
<b>5</b>									
<b>6</b>									
<b>7</b>									
<b>8</b>									
<b>9</b>									
<b>10</b>									
<b>11</b>									
<b>12</b>									
	<b>A</b>	<b>B</b>	<b>C</b>	<b>D</b>	<b>E</b>	<b>F</b>	<b>G</b>	<b>H</b>	

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