PRODUCT INFORMATION & MANUAL

Mouse Total MMP-9 Valukine™ ELISA

VAL621

For the quantitative determination of natural and recombinant Mouse active and Pro-Matrix Metalloproteinase 9 (Total MMP-9) concentration

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version202210.1

TABLE OF CONTENTS

I. BACKGROUND	2
II. OVERVIEW	3
III. ADVANTAGES	4
IV. EXPERIMENT	6
V. KIT COMPONENTS AND STORAGE	7
VI. PREPARATION	9
VII. ASSAY PROCEDURE	11
VIII. REFERENCES	13

I. BACKGROUND

Matrix MetalloProteinases (MMPs), also called matrixins, constitute a family of zinc and calcium dependent endopeptidases that function in the breakdown of extracellular matrix (ECM). They play an important role in many normal physiological processes such as embryonic development, morphogenesis, reproduction, and tissue remodeling (1, 2). They also participate in many pathological processes such as arthritis, cancer, and cardiovascular disease (3). While the amounts of newly synthesized MMPs are regulated mainly at the levels of transcription, the proteolytic activities of existing MMPs are controlled through both the activation of pro-enzymes or zymogens and the inhibition of active enzymes by endogenous inhibitors, α 2-macroglobulin, and tissue inhibitors of metalloproteinases (TIMPs).

MMP-9 (also referred to as gelatinase B, 92 kDa type IV collagenase, 92 kDa gelatinase, and type V collagenase) is secreted as a glycosylated pro-enzyme (4). Activation of the pro-enzyme involves a proteolytic removal of the N-terminal pro region containing the cysteine switch motif conserved in MMPs (5). The resulting 82 kDa active enzyme consists of a catalytic domain with a zinc-binding motif conserved in metzincins (6, 7). The catalytic domain also contains three contiguous fibronectin type II homology units responsible for binding gelatin (8). A proline-rich hinge region links the catalytic domain to the C-terminal hemopexin-like domain. In vitro treatment of the pro-enzyme with 4-aminophenylmercuric acetate (APMA) produces not only the active enzyme but also a C-terminal truncated form with activity comparable to that of the active form (9).

MMP-9 degrades components of the ECM with high specific activity for denatured collagens (gelatin). It can cleave native collagens of type IV, V, and XI, as well as elastin. MMP-9 can also cleave a variety of non-ECM molecules such as IL-1 β , IL-8, connective tissue-activating peptide-III, platelet factor-4, GRO α , substance P, myelin basic protein, and amyloid β peptide. MMP-9 can increase or decrease the biological activity of these molecules, depending upon the site of cleavage (10, 11).

MMP-9 is produced by a variety of normal and transformed cells including neutrophils, monocytes, macrophages, astrocytes, fibroblasts, osteoclasts, chondrocytes, keratinocytes, endothelial and epithelial cells, and is regulated by various agents. Transgenic mouse models report that MMP-9 regulates growth plate angiogenesis and apoptosis of hypertrophic chondrocytes (12), mediates early bone development through osteoclast recruitment (13), suppresses development of experimental abdominal aortic aneurysms (14), influences sensitization and reward (15), allows cytotrophoblast invasion (16), contributes to skin carcinogenesis (17), triggers the angiogenic switch during carcinogenesis (18), inactivates the serpin α 1-proteinase inhibitor (19), and associates with ischemia-induced blood brain barrier permeability (20).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for mouse MMP-9 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any MMP-9 present is bound by the immobilized antibody. After washing away any unbound substances, a biotinylated detection antibody specific for mouse MMP-9 is pipetted into the wells. After washing away any unbound substances, Streptavidin-HRP is pipetted into the wells. Following a wash to remove any unbound reagent, TMB Substrate is added to the wells and color develops in proportion to the amount of mouse MMP-9 bound in the initial step. The color development is stopped, and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- This kit is suitable for cell culture supernate and mouse serum.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Calibrator Diluent (1×) and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

Sample	Intra-assay Precision			Inter-assay Precision		
	1	2	3	1	2	3
Mean (pg/mL)	169.9	640.4	2662.3	169.9	649.0	2679.3
Standard Deviation	6.8	16.5	116.6	6.6	21.2	116.5
CV%	4.0	2.6	4.4	3.9	3.3	4.3

B. RECOVERY

The recovery of mouse MMP-9 spiked to different levels throughout the range of the assay in cell culture media was evaluated. The recovery ranged from 84.4-92.6% with an average of 90.1%.

The recovery of mouse MMP-9 spiked to different levels throughout the range of the assay in mouse serum was evaluated. The recovery ranged from 94.6-101.5% with an average of 98.1%.

C. SENSITIVITY

The minimum detectable dose (MDD) of mouse MMP-9 is typically less than 2.56 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against a highly purified NS0-expressed recombinant mouse MMP-9 produced at R&D Systems.

E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with high concentrations of mouse MMP-9 and diluted with Calibrator Diluent (1×) to produce samples with values within the dynamic range of the assay.

Dilution	Average % of Expected	Range (%)
1:2	103.0	98.7-113.5
1:4	107.3	100.3-120.2
1:8	109.4	97.5-127.9
1:16	110.4	100.2-122.2

F. SAMPLE VALUES

Serum - Three mouse serum samples were evaluated for the presence of mouse MMP-9 in this assay. All samples measured ranged from 52.1 to 69.4 ng/mL with an average of 58.3 ng/mL.

G. SPECIFICITY

The following factors prepared at 50 ng/mL were assayed and exhibited no cross-reactivity or interference.

Recombinant mouse	Other recombinants
MMP-2	human MMP-9
MMP-3	rat MMP-9
MMP-7	
MMP-8	
MMP-12	
TIMP-1	


This kit recognizes the Pro form of recombinant mouse (rm) MMP-9 and the mature form of rmMMP-9 by itself or when complexed to rmTIMP-1 or rmTIMP-2.

Recombinant mouse TIMP-2 does not cross-react in this assay but does interfere at concentrations > 6.25 ng/mL.

IV. EXPERIMENT

EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
Mouse Total MMP-9 Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with a Rat antibody against Mouse MMP-9.	1 plate
Mouse Total MMP-9 Standard	Recombinant Mouse MMP-9 in a buffered base; lyophilized. Refer to the vial label for reconstitution volume.	2 vials
Mouse Total MMP-9 Detection Antibody	Biotinylated Mouse MMP-9 antibody, lyophilized. Refer to the vial label for reconstitution volume.	1 vial
Calibrator Diluent (1×)	Buffered diluent used to dilute standard and samples.	2 vials
Reagent Diluent (1×)	Buffered protein base used to dilute Detection Antibody and Streptavidin-HRP.	2 vials
Streptavidin-HRP A (200×)	200× Streptavidin conjugated to horseradish peroxidase.	1 vial
Wash Buffer Concentrate (25×)	A 25× concentrated solution of buffered surfactant with preservatives.	1 vial
TMB Substrate	TMB ELISA Substrate Solution.	1 vial
Stop Solution	2 N sulfuric acid.	1 vial
Plate Sealers	Adhesive strip.	3 strips

B. STORAGE

Unopened Kit	Store at 2-8 °C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Streptavidin-HRP A	May be stored for up to 1 month at 2-8 °C.*
	Diluted Wash Solution	
	TMB Substrate	
	Stop Solution	
	Standard	Prepare fresh for each assay.
	Detection Antibody	Aliquot and store for up to 1 month at -20 °C in a manual defrost freezer. *
	Reagent Diluent (1×)	May be stored for up to 1 month at 2-8 °C.*
	Calibrator Diluent (1×)	May be stored for up to 1 month at 2-8 °C.*
Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8 °C.*	

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- ◆ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ◆ Pipettes and pipette tips.
- ◆ Deionized or distilled water.
- ◆ Squirt bottle, manifold dispenser, or automated microplate washer.
- ◆ 500 mL graduated cylinder.
- ◆ Horizontal orbital microplate shaker capable of maintaining a speed of 500±50 rpm.

D. PRECAUTION

- ◆ Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- ◆ The Stop Solution provided with this kit is an acid solution. Wear protective gloves, clothing, eye, and face protection. Wash hands thoroughly after handling.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples require dilution with Calibrator Diluent (1 \times).

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples require dilution with Calibrator Diluent (1 \times).

B. SAMPLE PREPARATION

Cell culture supernate samples require a 4-fold dilution. A suggested 4-fold dilution is 50 μL of sample + 150 μL of Calibrator Diluent (1 \times).

Serum samples require a 50-fold dilution. A suggested 50-fold dilution is 10 μL of sample + 490 μL of Calibrator Diluent (1 \times).

C. REAGENT PREPARATION

Note: Bring all reagents to room temperature before use.

Wash Buffer (1 \times)- If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25 \times) into deionized or distilled water to prepare 500 mL of Wash Buffer (1 \times).

Detection Antibody- **Centrifuge briefly before opening. Reconstitution Volume refer to vial label with Reagent Diluent (1 \times)**. Aliquot and store if needed. Allow the detection antibody to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions. Dilute stock solution in Reagent Diluent (1 \times) to the working concentration of 100 ng/mL. Prepare at least 15 minutes prior to use.


Streptavidin-HRP A (1 \times) - **Centrifuge briefly before opening**. Dilute to the working concentration specified on the vial label using Reagent Diluent (1 \times).

Mouse Total MMP-9 Standard - **Centrifuge briefly before opening. Refer to the vial label for the reconstitution volume***. This reconstitution produces a stock

solution of 10,000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*If you have any question, please seek help from our Technical Support.

Pipette 300 μ L of the appropriate Calibrator Diluent (1 \times) into each tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 5000 pg/mL standard serves as the high standard. The Calibrator Diluent (1 \times) serves as the zero standard (0 pg/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB Substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB Substrate.

VII. ASSAY PROCEDURE

Note: Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100 μ L of standard, or prepared sample per well. Cover with the adhesive strip provided. Incubate for 2 hours at room temperature **on a horizontal orbital microplate shaker set at 500 \pm 50rpm**. A plate layout is provided for a record of standards and samples assayed. (Samples may require dilution. See Sample Preparation section.)
4. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Solution (400 μ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
5. Add 100 μ L of the Detection Antibody diluted in Reagent Diluent, to each well. Cover with a new adhesive strip and incubate 2 hours at room temperature **on a horizontal orbital microplate shaker set at 500 \pm 50rpm**.
6. Repeat the aspiration/wash as in step 4.
7. Add 100 μ L of the working dilution of Streptavidin-HRP A to each well. Cover the plate and incubate for 20 minutes at room temperature. **Avoid placing the plate in direct light.**
8. Repeat the aspiration/wash as in step 4.
9. Add 100 μ L of TMB Substrate to each well. Incubate for 20 minutes at room temperature. **Avoid placing the plate in direct light.**
10. Add 50 μ L of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
11. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570

nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

12. CALCULATION OF RESULTS

Average the duplicate readings for each standard, sample and subtract the average zero standard optical density (O.D.). Create a standard curve by reducing the data using computer software capable of generating a four-parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the mouse MMP-9 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Nagase, H. and J.F. Woessner, Jr. (1999) *J. Biol. Chem.* 274:2191.
2. Page-McCaw, A. et al. (2007) *Nat. Rev. Mol. Cell. Biol.* 9:221.
3. Parks, W.C. and R.P. Mecham (1998) in *Matrix Metalloproteinases*, Parks, W.C. and R.P. Mecham, eds., Academic Press, San Diego.
4. Collier, I.E. and G.I. Goldberg (1998) in *Handbook of Proteolytic Enzymes*, A.J. Barrett, et al. eds., Academic Press, San Diego, pp. 1205-1210.
5. Van Wart, H.E. and H. Birkedal-Hansen (1990) *Proc. Natl. Acad. Sci. USA* 87:5578.
6. Jiang, W. and J.S. Bond (1992) *FEBS Lett.* 312:110.
7. Bode, W. et al. (1993) *FEBS Lett.* 331:134.
8. Collier, I.E. et al. (1992) *J. Biol. Chem.* 267:6776.
9. O'Connell, J.P. et al. (1994) *J. Biol. Chem.* 269:14967.
10. Vu, T.H. and Z. Werb (1998) in *Matrix Metalloproteinases*, Parks, W.C. and R.P. Mecham, eds., Academic Press, San Diego, pp. 115-148.
11. Van den Steen, P.E. et al. (2000) *Blood* 96:2673.
12. Vu, T.H. et al. (1998) *Cell* 93:411.
13. Engsig, M.T. et al. (2000) *J. Cell Biol.* 151:879.
14. Pyo, R. et al. (2000) *J. Clin. Invest.* 105:1641.
15. Mizoguchi, H. et al. (2005) *J. Neurochem.* 100:1579.
16. Bischof, P. et al. (2002) *J. Reprod. Immunol.* 55:3.
17. Coussens, L.M. et al. (2000) *Cell* 103:481.
18. Bergers, G.R. et al. (2000) *Nat. Cell Biol.* 2:737.
19. Liu, Z. et al. (2000) *Cell* 102:647.
20. Svedin, P. et al. (2007) *J. Neurosci.* 27:1511.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

1																	
2																	
3																	
4																	
5																	
6																	
7																	
8																	
9																	
10																	
11																	
12																	
	A	B	C	D	E	F	G	H									

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