PRODUCT INFORMATION & MANUAL

Human IL-10 Valukine™ ELISA

Catalog Number: VAL112

For the quantitative determination of natural and recombinant human IL-10 concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd
P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version 202606.6

TABLE OF CONTENTS

I. BACKGROUND 2

II. OVERVIEW 3

III. ADVANTAGES 4

IV. EXPERIMENT 8

V. KIT COMPONENTS AND STORAGE 9

VI. PREPARATION 11

VII. ASSAY PROCEDURE 13

VIII. REFERENCES 15

I. BACKGROUND

Interleukin-10 (IL-10), also known as cytokine synthesis inhibitory factor (CSIF), is the charter member of the IL-10 α-helical cytokine family that also includes IL-19, IL-20, IL-22, IL-24, and IL-26/AK155 (1-3). IL-10 is secreted by many activated hematopoietic cell types as well as hepatic stellate cells, keratinocytes, and placental cytotrophoblasts. Whereas human IL-10 is active on mouse cells, mouse IL-10 does not act on human cells (4, 5). Mature human IL-10 shares 86% amino acid sequence identity with equine IL-10 and 72%-80% with bovine, canine, feline, guinea pig, mouse, ovine, porcine, and rat IL-10. It contains two intrachain disulfide bridges and is expressed as a 36 kDa noncovalently-associated homodimer (4, 6, 7).

IL-10 mediates its biological activities through a heteromeric receptor complex composed of the type II cytokine receptor subunits IL-10 Rα and IL-10 Rβ. IL-10 Rα is a 110 kDa transmembrane glycoprotein that is expressed on lymphocytes, NK cells, macrophages, monocytes, astrocytes, intestinal epithelial cells, cytotrophoblasts, and activated hepatic stellate cells (8-13), while the 75 kDa transmembrane IL-10 Rβ is widely expressed (14, 15). The IL-10 dimer binds to two IL-10 Rα chains, triggering recruitment of two IL-10 Rβ chains (14, 15). IL-10 Rβ does not bind IL-10 directly but is required for signal transduction. IL-10 Rβ also associates with IL-20 Rα, IL-22 Rα1, or IL-28 Rα to form the receptor complexes for IL-22, IL-26, IL-28, and IL-29 (16-18).

The involvement of IL-10 in immunoregulation includes both suppressive and stimulatory effects. It functions as an anti-inflammatory cytokine by inhibiting the expansion and activation of Th1 cells and Th17 cells (19-21) and by promoting the development of M2 macrophages (21). Its expression by immunosuppressive regulatory T cells (Treg) and regulatory B cells is important for Treg proliferation (19). Within a tumor microenvironment, however, IL-10 inhibits the expansion of Treg as well as myeloid-derived suppressor cells (22, 23). IL-10 induces the intratumoral accumulation and activation of $CD8^+$ T cells (24, 25). IL-10 exerts protective effects including limiting tissue damage in arthritic inflammation (19) and promoting muscle regeneration after injury (21), but it also contributes to the persistence of viral infections (26). The levels of IL-10 are elevated in Sjogren's syndrome (saliva), primary CNS lymphoma (cerebrospinal fluid), and ovarian cancer (serum and ascites) (27-29). Its levels are decreased in the serum in patients with recurrent heart attacks or during preeclampsia and also in the seminal fluid of infertile men (30-32).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for IL-10 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any IL-10 present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked antibody specific for IL-10 is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, TMB substrate solution (Chromogenic agent) is added to the wells and color develops in proportion to the amount of IL-10 bound in the initial step. The color development is stopped and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with the appropriate Calibrator Diluent and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in forty separate assays to assess inter-assay precision.

CELL CULTURE SUPERNATE ASSAY

	Intra-assay Precision			**Inter-assay Precision**		
Sample	1	2	3	1	2	3
n	20	20	20	40	40	40
Mean (pg/mL)	24.4	79.6	223	23.6	74.2	227
Standard Deviation	1.6	3.9	5.5	1.8	4.9	12.7
CV%	6.6	4.9	2.5	7.6	6.6	5.6

SERUM/PLASMA ASSAY

	Intra-assay Precision			**Inter-assay Precision**		
Sample	1	2	3	1	2	3
n	20	20	20	40	40	40
Mean (pg/mL)	23.9	76.9	231	23.2	75.1	228
Standard Deviation	1.2	3.3	3.9	1.7	5.6	13.4
CV%	5.0	4.3	1.7	7.3	7.5	5.9

B. RECOVERY

The recovery of human IL-10 spiked to three different levels in samples throughout the range of the assay in various matrices was evaluated.

Sample Type	Average % Recovery	Range (%)
Cell culture media (n=5)	98	88 - 110
Serum (n=5)	100	92 - 110
EDTA plasma (n=5)	92	75 - 105
Heparin plasma (n=5)	94	84 - 105
Citrate plasma (n=5)	96	87 - 108

C. SENSITIVITY

The minimum detectable dose (MDD) of IL-10 is typically less than 3.9 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against highly purified *Sf* 21-expressed recombinant human IL-10 produced at R&D Systems®.

The NIBSC/WHO Reference Reagent recombinant human IL-10 93/722 was evaluated in this kit. The dose response curve of the NIBSC/WHO Reference Reagent 93/722 parallels the Valukine standard curve. To convert sample values obtained with the Valukine Human IL-10 kit to approximate NIBSC 93/722 units, use the equation below.

NIBSC (93/722) approximate value (U/mL) = 0.0103 x Valukine Human IL-10 value (pg/mL)

Note: *Based on data generated in May 2010.*

E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with high concentrations of IL-10 and diluted with the appropriate Calibrator Diluent to produce samples with values within the dynamic range of the assay.

		Cell culture media (n=5)	Serum (n=5)	EDTA plasma (n=5)	Heparin plasma (n=5)	Citrate plasma (n=5)
1:2	Average % of Expected	98	101	100	99	97
	Range (%)	92 - 108	96 - 106	97 - 102	91 - 104	93 - 100
1:4	Average % of Expected	96	102	100	100	100
	Range (%)	86 - 109	97 - 108	91 - 106	92 - 106	93 - 105
1:8	Average % of Expected	94	102	105	100	97
	Range (%)	84 - 104	97 - 109	96 - 112	93 - 107	89 - 104
1:16	Average % of Expected	90	102	108	95	89
	Range (%)	82 - 98	91 - 114	99 - 114	85 - 110	82 - 96

F. SAMPLE VALUES

Serum/Plasma - Forty serum and plasma samples from apparently healthy volunteers were evaluated for the presence of human IL-10 in this assay. All samples measured less than the lowest standard, 7.8 pg/mL. No medical histories were available for the donors used in this study.

Cell Culture Supernates -Human peripheral blood mononuclear cells (1×10^6 cells/mL) were cultured in RPMI supplemented with 10% fetal calf serum, 50 µM β-mercaptoethanol, 2 mM L-glutamine, 100 U/mL penicillin, and 100 µg/mL streptomycin sulfate. Cells were cultured unstimulated or stimulated with 10 µg/mL PHA for 1, 3 and 5 days. Aliquots of the cell culture supernate was removed and assayed for IL-10.

Condition	Day 1 (pg/mL)	Day 3 (pg/mL)	Day 5 (pg/mL)
Unstimulated	98	134	127
Stimulated	2756	3563	2255

G. SPECIFICITY

This assay recognizes both natural and recombinant IL-10. The following factors were prepared at 50 ng/mL in the appropriate Calibrator Diluent and assayed for cross-reactivity. Preparations of the following factors at 50 ng/mL in a mid-range rhIL-10 control were assayed for interference. No significant cross-reactivity or interference was observed.


Recombinant Human	**Recombinant Mouse**
IL-10 (cytomegalovirus)	IL-10
IL-10 Rα	IL-10 Rα
IL-10 Rβ/Fc Chimera	
Other Recombinants	
Canine IL-10	Porcine IL-10
Equine IL-10	Rat IL-10
Feline IL-10	Viral IL-10 (Epstein-Barr)

IV. EXPERIMENT

EXAMPLE STANDARD


The standard curves are provided for demonstration only. A standard curve should be generated for each set of samples assayed.

CELL CULTURE SUPERNATE ASSAY

(pg/mL)	O.D.	Average	Corrected
0	0.028 0.032	0.030	—
7.8	0.071 0.088	0.080	0.050
15.6	0.108 0.110	0.109	0.079
31.3	0.203 0.211	0.207	0.177
62.5	0.418 0.432	0.425	0.395
125	0.718 0.755	0.737	0.707
250	1.440 1.506	1.473	1.443
500	2.729 2.779	2.754	2.724

SERUM/PLASMA ASSAY

(pg/mL)	O.D.	Average	Corrected
0	0.040 0.044	0.042	—
7.8	0.091 0.106	0.098	0.056
15.6	0.135 0.137	0.136	0.094
31.3	0.213 0.216	0.214	0.172
62.5	0.378 0.391	0.384	0.342
125	0.705 0.729	0.717	0.675
250	1.305 1.322	1.313	1.271
500	2.330 2.332	2.331	2.289

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
IL-10 Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with an antibody against IL-10	1 plate
IL-10 Conjugate	A solution of antibody against IL-10 conjugated to horseradish peroxidase with preservatives	1 vial
IL-10 Standard	Recombinant IL-10 in a buffered protein base with preservatives; lyophilized. *Refer to the vial label for reconstitution volume.*	1 vial
Assay Diluent RD1W	A buffered protein base with preservatives. *For serum/plasma samples.*	1 vial
Calibrator Diluent (5×)/RD5C	A 5× concentrated buffered protein base with preservatives. *For cell culture supernate samples.*	1 vial
Calibrator Diluent RD6P	Animal serum with preservatives. *For serum/plasma samples.*	1 vial
Wash Buffer Concentrate (25×)	A 25× concentrated solution of buffered surfactant	1 vial
TMB Substrate	TMB ELISA Substrate Solution/TMB Substrate Solution	2 vials
Stop Solution	2 N sulfuric acid	1 vial
Plate Sealers	Adhesive strip	3 strips

B. STORAGE

<table>
<tr><td>Unopened Kit</td><td colspan="2">Store at 2-8°C. Do not use past kit expiration date.</td></tr>
<tr><td rowspan="7">Opened/ Reconstituted Reagents</td><td>Wash Buffer (1×)</td><td rowspan="4">May be stored for up to 1 month at 2-8°C.*</td></tr>
<tr><td>Stop Solution</td></tr>
<tr><td>Conjugate</td></tr>
<tr><td>TMB Substrate</td></tr>
<tr><td>Standard</td><td>Aliquot and store for up to 1 month at <-20°C in a manual defrost freezer. * Avoid repeated freeze-thaw cycles.</td></tr>
<tr><td>Assay Diluent RD1W</td><td rowspan="2">May be stored for up to 1 month at 2-8 °C.*</td></tr>
<tr><td>Calibrator Diluent RD6P</td></tr>
<tr><td>Calibrator Diluent (5×)/RD5C</td><td>May be stored for up to 1 month at 2-8 °C.*
Use and discard diluted Calibrator Diluent (1×). Prepare fresh for each assay.</td></tr>
<tr><td>Microplate Wells</td><td>Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*</td></tr>
</table>

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- Pipettes and pipette tips.
- Deionized or distilled water.
- Squirt bottle, manifold dispenser, or automated microplate washer.
- 500 mL graduated cylinder.

D. PRECAUTION

The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

The sample collection and storage conditions listed below are intended as general guidelines. Sample stability has not been evaluated.

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at ≤-20°C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent RD5C (1×).

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at ≤ -20°C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent RD6P.

Plasma - Collect plasma using EDTA, heparin, or citrate as an anticoagulant. Centrifuge for 15 minutes at 1000 x g within 30 minutes of collection. Assay immediately or aliquot and store samples at ≤ -20°C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent RD6P.

Note: *Hemolyzed samples are not suitable for use this assay.*

B. SAMPLE PREPARATION

Polypropylene tubes must be used. Do not use glass.

Cell culture supernate samples require at least a 10-fold dilution. A suggested 10-fold dilution is 50 µL of sample + 450 µL of Calibrator Diluent RD5C (1×)*. Optimal dilutions should be determined by the end user.

*See Reagent Preparation section.

C. REAGENT PREPARATION

Note: *Bring all reagents to room temperature before use.*


Wash Buffer (1×) - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25×) into deionized or distilled water to prepare 500 mL of Wash Buffer (1×).

Calibrator Diluent RD5C (1×) - Use deionized or distilled water to prepare Calibrator Diluent RD5C (1×).

IL-10 Standard - Refer to the vial label for reconstitution volume*. This reconstitution produces a stock solution of 5000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*If you have any question, please seek help from our Technical Support.

Use polypropylene tubes. Pipette 900 µL of Calibrator Diluent RD5C (1×) (*for cell culture supernate samples*) or Calibrator Diluent RD6P (*for serum/plasma samples*) into 500 pg/mL tube. Pipette 500 µL of the appropriate Calibrator Diluent into each remaining tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 500 pg/mL standard serves as the high standard. The appropriate Calibrator Diluent serves as the zero standard (0 pg/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB Substrate Solution. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB Substrate.

VII. ASSAY PROCEDURE

Note: *Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.*

1. Prepare all reagents, samples and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. **For Serum/Plasma Samples Only:** Add 50 µL of Assay Diluent RD1W to each well.
4. Add 200 µL of standard or sample per well. Cover with the adhesive strip provided. **Incubate for 2 hours at room temperature.** A plate layout is provided for a record of standards and samples assayed.
5. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Buffer (400 µL) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
6. Add 200 µL of Human IL-10 Conjugate to each well. Cover with a new adhesive strip.

 For Cell Culture Supernate Samples: Incubate for 1 hour at room temperature.

 For Serum/Plasma Samples: Incubate for 2 hours at room temperature.
7. Repeat the aspiration/wash as in step 5.
8. Add 200 µL of TMB Substrate to each well. **Protect from light.**

 For Cell Culture Supernate Samples: Incubate for 20 minutes at room temperature.

 For Serum/Plasma Samples: Incubate for 30 minutes at room temperature.
9. Add 50 µL of Stop Solution to each well. The color in the wells should change from blue to yellow. If the color in the wells is green or the color change does not appear uniform, gently tap the plate to ensure thorough mixing.

10. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

11. **CALCULATION OF RESULTS**

 Average the duplicate readings for each standard and sample and subtract the average zero standard optical density (O.D.). Create a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the IL-10 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.


 If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Ouyang, W. et al. (2011) Annu. Rev. Immunol. 29:71.
2. Sabat, R. et al. (2010) Cytokine Growth Factor Rev. 21:331.
3. Saraiva, M. and A. O'Garra (2010) Nat. Rev. Immunol. 10:170.
4. Vieira, P. et al. (1991) Proc. Natl. Acad. Sci. USA 88:1172.
5. Hsu, D.-H. et al. (1990) Science 250:830.
6. Windsor, W.T. et al. (1993) Biochemistry 32:8807.
7. Syto, R. et al. (1998) Biochemistry 37:16943.
8. Mathurin, P. et al. (2002) Am. J. Physiol. Gastrointest. Liver Physiol. 282:G981.
9. Szony, B.J. et al. (1999) Mol. Hum. Reprod. 5:1059.
10. Liu, Y. et al. (1994) J. Immunol. 152:1821.
11. Carson, W.E. et al. (1995) Blood 85:3577.
12. Denning, T.L. et al. (2000) Int. Immunol. 12:133.
13. Ledeboer, A. et al. (2002) Eur. J. Neurosci. 16:1175.
14. Kotenko, S.V. et al. (1997) EMBO J. 16:5894.
15. Lutfalla, G. et al. (1993) Genomics 16:366.
16. Kotenko, S.V. et al. (2000) J. Biol. Chem. 276:2725.
17. Hor, S. et al. (2004) J. Biol. Chem. 279:33343.
18. Sheppard, P. et al. (2003) Nat. Immunol. 4:63.
19. Carter, N.A. et al. (2011) J. Immunol. 186:5569.
20. Gu, Y. et al. (2008) Eur. J. Immunol. 38:1807.
21. Deng, B. et al. (2012) J. Immunol. 189:3669.
22. Tanikawa, T. et al. (2012) Cancer Res. 72:420.
23. Groux, H. et al. (1997) Nature 389:737.
24. Mumm, J.B. et al. (2011) Cancer Cell 20:781.
25. Emmerich, J. et al. (2012) Cancer Res. 72:3570.
26. Wilson, E.B. and D.G. Brooks (2011) Curr. Top. Microbiol. Immunol. 350:39.
27. Bertorello, R. et al. (2004) Clin. Exp. Med. 4:148.
28. Sasayama, T. et al. (2012) Neuro. Oncol. 14:368.
29. Mustea, A. et al. (2006) Anticancer Res. 26:1715.
30. George, J. et al. (2012) Atherosclerosis 222:519.
31. Borekci, B. et al. (2007) Am. J. Reprod. Immunol. 58:56.
32. Camejo, M.I. (2003) Arch. Androl. 49:111.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

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