PRODUCT INFORMATION & MANUAL

Mouse IFN-γ Valukine™ ELISA

Catalog Number: VAL607

For the quantitative determination of natural and recombinant mouse IFN-γ concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd
P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version 202609.5

TABLE OF CONTENTS

I. BACKGROUND 2

II. OVERVIEW 3

III. ADVANTAGES 4

IV. EXPERIMENT 7

V. KIT COMPONENTS AND STORAGE 8

VI. PREPARATION 11

VII. ASSAY PROCEDURE 13

VIII. REFERENCES 15

I. BACKGROUND

Interferon-gamma (IFN-γ, also known as type II interferon) is an important immunoregulatory cytokine that was originally identified through its anti-viral activity (1, 2). It plays key roles in host defense by exerting antiviral, antiproliferative, and immunoregulatory activities (3, 4). On many cell types, IFN-γ induces the production of cytokines and upregulates the expression of various membrane proteins including class I and II MHC antigens, Fc receptors, leukocyte adhesion molecules, and B7 family antigens. IFN-γ is a potent activator of macrophage effector functions. It directs the synthesis, class switching, and secretion of immunoglobulins by B cells. IFN-γ also influences T-helper cell phenotype development by inhibiting Th2 differentiation and stimulating Th1 development (3, 4). IFN-γ plays a central role in the progression of inflammatory diseases such as autoimmunity and atherosclerosis (5, 6).

Biologically active IFN-γ consists of a noncovalently linked homodimer of 20-25 kDa variably glycosylated subunits (7). Mature mouse IFN-γ shares 86% amino acid (aa) sequence identity with rat IFN-γ and 38-44% aa identity with bovine, canine, cotton rat, equine, feline, human, porcine, and rhesus IFN-γ. IFN-γ dimers bind to transmembrane IFN-γ RI (alpha subunits) which then interact with transmembrane IFN-γ RII (beta subunits) to form the functional receptor complex of two α and two β subunits (8, 9). Inclusion of IFN-γ RII in the receptor complex increases the ligand binding affinity as well as the efficiency of signal transduction (9, 10). Whereas the α-chain is expressed constitutively on many cell types, the cellular regulation of the β-chain correlates with an IFN-γ responsive state and is tightly regulated (8).

IFN-γ is produced by a number of cell types under inflammatory conditions, including dendritic epidermal/γδ T cells (11), keratinocytes (12), peripheral blood γδ T cells (13), mast cells (14), neurons (15), $CD8^+$ T cells (16), macrophages (17), B cells (18), neutrophils (19), NK cells (20), $CD4^+$ T cells (21), and testicular spermatids (22).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for mouse IFN-γ has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any mouse IFN-γ present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked antibody specific for mouse IFN-γ is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, TMB substrate (Chromogenic agent) is added to the wells and color develops in proportion to the amount of mouse IFN-γ bound in the initial step. The color development is stopped and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- This kit is suitable for cell culture supernate and mouse serum.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with appropriate calibrator diluent and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

CELL CULTURE SUPERNATE ASSAY

	Intra-assay Precision			**Inter-assay Precision**		
Sample	1	2	3	1	2	3
Mean (pg/mL)	44.8	156	518	43.3	144	509
Standard Deviation	2.2	5.0	14.2	3.6	14.7	48.6
CV%	4.9	3.2	2.7	8.3	10.2	9.5

SERUM ASSAY

	Intra-assay Precision			**Inter-assay Precision**		
Sample	1	2	3	1	2	3
Mean (pg/mL)	45.1	143	528	42.8	139	518
Standard Deviation	2.1	4.4	11.8	3.6	13.6	46.5
CV%	4.7	3.1	2.2	8.4	9.8	9.0

B. RECOVERY

The recovery of mouse IFN-γ spiked to three levels throughout the range of the assay in various matrices was evaluated.

Sample Type	**Average % Recovery**	**Range%**
Cell culture supernates (n=7)	105	98-115
Serum (n=5)	97	91-106

C. SENSITIVITY

The minimum detectable dose (MDD) of mouse IFN-γ is typically less than 2 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against highly purified *E. coli*-expressed recombinant mouse IFN-γ produced at R&D Systems.

The NIH reference preparation mouse IFN-γ Gg02-901-533 which was intended as a bioassay standard, was evaluated in this kit. Each ampule contains a nominal 1.0 µg of natural mouse IFN-γ and was assigned an arbitrary unitage of 1000 Units/ampule.

NIH Gg02-901-533: 1 Unit of standard = 230 pg of Valukine™ Mouse IFN-γ

E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with high concentrations of mouse IFN-γ and diluted with appropriate calibrator diluent to produce samples with values within the dynamic range of the assay.

		Cell culture supernates (n=4)	**Serum (n=4)**
1:2	Average % of Expected	101	97
	Range (%)	99-104	92-101
1:4	Average % of Expected	101	95
	Range (%)	98-104	92-98
1:8	Average % of Expected	104	94
	Range (%)	100-110	91-96
1:16	Average % of Expected	100	92
	Range (%)	95-105	90-93

F. SAMPLE VALUES

Cell Culture Supernates:

Mouse splenocytes (2 x 10^6 cells/mL) were cultured for 3 days in RPMI plus 10% fetal bovine serum and stimulated with 10 µg/mL PHA. An aliquot of the cell culture

supernate was removed, assayed for mouse IFN-γ , and measured 32 ng/mL.

EL-4 mouse lymphoblast cells (9 x 10^5 cells/mL) were cultured for 2 days in DMEM plus 10% fetal bovine serum and stimulated with 10 µg/mL PHA and 10 ng/mL PMA. An aliquot of the cell culture supernate was removed, assayed for mouse IFN-γ, and measured 4 ng/mL.

Mouse lung conditioned media (1 lung, 1-2 mm pieces in 10 mL of medium) was collected after culturing for 5 days in RPMI plus 10% fetal bovine serum. An aliquot of the cell culture supernate was removed, assayed for mouse IFN-γ, and measured 198 pg/mL.

Mouse serum - Forty mouse serum samples were evaluated for the presence of mouse IFN-γ in this assay. Thirty-eight samples measured less than the lowest Mouse IFN-g Standard, 9.4 pg/mL. Two samples read 24 pg/mL and 60 pg/mL, respectively.

G. SPECIFICITY

This assay recognizes both natural and recombinant mouse IFN-γ. The following factors were prepared at 50 ng/mL in calibrator diluent RD5Y/RD6-12 and assayed for cross-reactivity. Preparations of the following factors at 50 ng/mL in a mid-range recombinant mouse IFN-γ control were assayed for interference. No significant cross-reactivity or interference was observed.

Recombinant mouse			**Recombinant human**
C10	IL-4	KC	IFN-γ
G-CSF	IL-5	LIF	
GM-CSF	IL-6	M-CSF	
IFN-γ RI	IL-7	MIP-1 α	
IFN-γ RII	IL-9	MIP-1 β	
IFN-κ	IL-10	MIP-2	
IL-1 α	IL-10 R	SCF	
IL-1 β	IL-12	TNF-α	
IL-2	IL-13	Tpo	
IL-3	JE	VEGF	


A sample containing 66 ng/mL of recombinant rat IFN-γ measured 27 pg/mL in this assay (0.04% cross-reactivity). Upon dilution, the dose-curve of the recombinant rat IFN-γ was parallel to the mouse IFN-γ standard curve.

IV. EXPERIMENT

EXAMPLE STANDARD


The standard curves are provided for demonstration only. A standard curve should be generated for each set of samples assayed.

CELL CULTURE SUPERNATE ASSAY

(pg/mL)	O.D.	Average	Corrected
0	0.042 0.043	0.042	—
9.4	0.077 0.077	0.077	0.035
18.8	0.113 0.112	0.112	0.070
37.5	0.183 0.186	0.184	0.142
75	0.327 0.329	0.328	0.286
150	0.580 0.594	0.587	0.545
300	1.130 1.124	1.127	1.085
600	2.020 2.042	2.031	1.989

SERUM ASSAY

(pg/mL)	O.D.	Average	Corrected
0	0.035 0.037	0.036	—
9.4	0.069 0.071	0.070	0.034
18.8	0.099 0.107	0.103	0.067
37.5	0.170 0.179	0.174	0.138
75	0.309 0.318	0.314	0.278
150	0.551 0.576	0.564	0.528
300	1.053 1.081	1.067	1.031
600	1.950 2.005	1.978	1.942

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
Mouse IFN-γ Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with an antibody against mouse IFN-γ	1 plate
Mouse IFN-γ Conjugate	Antibody against mouse IFN-γ conjugated to horseradish peroxidase	1 vial
Mouse IFN-γ Standard	Recombinant mouse IFN-γ in a buffered protein base; lyophilized. *Refer to the vial label for reconstitution volume.*	2 vials
Assay Diluent RD1-21	A buffered protein base with preservatives. May contain a precipitate. *Mix well before and during use.*	1 vial
Calibrator Diluent RD5Y	A buffered protein base used to dilute standard and samples. *For cell culture supernate samples*	1 vial
Calibrator Diluent RD6-12	A buffered protein base used to dilute standard and samples. *For serum samples*	1 vial
Wash Buffer Concentrate (25×)	A 25× concentrated solution of buffered surfactant	1 vial
TMB Substrate	TMB ELISA Substrate Solution/TMB Substrate Solution	1 vial
Stop Solution	Diluted hydrochloric acid solution	1 vial
Plate Sealers	Adhesive strip	3 strips

B. STORAGE

<table>
<tr><td>Unopened Kit</td><td colspan="2">Store at 2-8 °C. Do not use past kit expiration date.</td></tr>
<tr><td rowspan="8">Opened/ Reconstituted Reagents</td><td>Wash Buffer (1×)</td><td rowspan="6">May be stored for up to 1 month at 2-8 °C.*</td></tr>
<tr><td>Assay Diluent
RD1-21</td></tr>
<tr><td>Calibrator
Diluent RD5Y</td></tr>
<tr><td>Calibrator Diluent RD6-12</td></tr>
<tr><td>Stop Solution</td></tr>
<tr><td>Conjugate</td></tr>
<tr><td>TMB Substrate</td><td rowspan="2">Prepare fresh for each assay.
Discard after use.
Standards may be stored for up to 1 month at -20°C *.</td></tr>
<tr><td>Standard</td></tr>
<tr><td>Microplate Wells</td><td>Return unused wells to the oil ouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*</td></tr>
</table>

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- Pipettes and pipette tips.
- Deionized or distilled water.
- Squirt bottle, manifold dispenser, or automated microplate washer.
- 500 mL graduated cylinder.

D. PRECAUTION

- The Stop Solution provided with this kit is an acid solution.
- Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- TMB substrate may cause skin, eye, and respiratory irritation. Avoid breathing fumes.
- Wear protective gloves, clothing, eye, and face protection. Wash hands thoroughly after handling

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

The sample collection and storage conditions listed below are intended as general guidelines. Sample stability has not been evaluated.

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles.

Serum - Use a serum separator tube (SST) and allow samples to clot for 2 hours at room temperature before centrifugation for 20 minutes at 2000 × g. Remove serum and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles.

Note: *Grossly hemolyzed or lipemic samples may not be suitable for use in this assay.*

B. SAMPLE PREPARATION

Cell culture supernate samples may require a dilution.

Optimal dilutions should be determined by the end user.

C. REAGENT PREPARATION


Note: *Bring all reagents to room temperature before use.*

Wash Buffer (1×) - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25×) into deionized or distilled water to prepare 500 mL of Wash Buffer (1×).

Mouse IFN-γ Standard - Centrifuge briefly before opening. Refer to the vial label for reconstitution volume*. Reconstitute the Mouse IFN-γ Standard with Calibrator Diluent RD5Y (for cell culture supernate samples) or Calibrator Diluent RD6-12 (for serum samples). Do not substitute other diluents. This reconstitution produces a stock solution of 3000 pg/mL. Allow the standard to sit for a minimum of 5 minutes with gentle mixing prior to making dilutions.

*If you have any question, please seek help from our Technical Support.

Use polypropylene tubes. Pipette 400 µL of Calibrator Diluent RD5Y (for cell culture supernate samples) or Calibrator Diluent RD6-12 (for serum samples) into the 600 pg/mL tube. Pipette 200 µL of appropriate calibrator diluent into each remaining tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 600 pg/mL standard serves as the high standard. The appropriate calibrator diluent serves as the zero standard (0 pg/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB Substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB Substrate.

VII.ASSAY PROCEDURE

Note: *Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.*

1. Prepare all reagents, samples and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 50 µL of Assay Diluent RD1-21 to each well.
4. Add 50 µL of standard and sample per well. Mix by gently tapping the plate frame for 1 minute. Cover with the adhesive strip provided. **Incubate for 2 hours at room temperature.** A plate layout is provided for a record of standards and samples assayed.
5. Aspirate each well and wash, repeating the process four times for a total of five washes. Wash by filling each well with Wash Buffer (400 µL) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
6. Add 100 µL of mouse IFN-γ conjugate to each well. Cover with a new adhesive strip. **Incubate for 2 hours at room temperature.**
7. Repeat the aspiration/wash as in step 5.
8. Add 100 µL of TMB Substrate to each well. **Incubate for 30 minutes at room temperature. Protect from light**.
9. Add 100 µL of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
10. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction

may be higher and less accurate.

11. **CALCULATION OF RESULTS.** Average the duplicate readings for each standard and sample and subtract the average zero standard optical density. Create a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the mouse IFN-γ concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.


 If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Billiau, A. and P. Matthys (2009) Cytokine Growth Factor Rev. 20:97.
2. Wheelock, E.F. (1965) Science 146:310.
3. Schoenborn, J.R & C.B. Wilson (2007) Adv. Immunol. 96:41.
4. Pestka, S. et al. (2004) Immunol. Rev. 202:8.
5. Kelchtermans, H. et al. (2008) Trends Immunol. 29:479.
6. McLaren, J.E. and D.P. Ramji (2009) Cytokine Growth Factor Rev. 20:125.
7. Gray, P.W. and D.V. Goeddel (1983) Proc. Natl. Acad. Sci. 80:5842.
8. Bach, E.A. et al. (1997) Annu. Rev. Immunol. 15:563.
9. Marsters, S.A. et al. (1995) Proc. Natl. Acad. Sci. 92:5401.
10. Krause, C.D. et al. (2000) J. Biol. Chem. 275:22995.
11. Sugaya, M. et al. (1999) J. Invest. Dermatol. 113:350.
12. Howie, S.E.M. et al. (1996) J. Invest. Dermatol. 106:1218.
13. Battistini, L. et al. (1997) J. Immunol. 159:3723.
14. Gup ta, A.A. et al. (1996) J. Immunol. 157:2123.
15. Neumann, H. et al. (1997) J. Exp. Med. 186:2023.
16. Hoiden, I. & G. Moller (1996) Scand. J. Immunol. 44:501.
17. Puddu, P. et al. (1997) J. Immunol. 159:3490.
18. Yoshimoto, T. et al. (1997) Proc. Natl. Acad. Sci. 94:3948.
19. Yeaman, G.R. et al. (1998) J. Immunol. 160:5145.
20. Asea, A. et al. (1996) Clin. Exp. Immunol. 105:376.
21. Briscoe, D.M. et al. (1997) J. Immunol. 159:3247.
22. Dejuco, N. et al. (1995) Endocrinology 136:4925.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

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