PRODUCT INFORMATION & MANUAL

Mouse TNF-α Valukine™ ELISA

Catalog Number: VAL609

For the quantitative determination of natural and recombinant mouse TNF-α concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd
P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version 202606.6

TABLE OF CONTENTS

I. BACKGROUND 2

II. OVERVIEW 3

III. ADVANTAGES 4

IV. EXPERIMENT 9

V. KIT COMPONENTS AND STORAGE 10

VI. PREPARATION 13

VII. ASSAY PROCEDURE 15

VIII. REFERENCES 17

I. BACKGROUND

Tumor necrosis factor alpha (TNF-α), also known as cachectin and TNFSF1A, is the prototypic ligand of the TNF superfamily (1). It is a pleiotropic molecule that plays a central role in inflammation, immune system development, apoptosis, and lipid metabolism (2-5). TNF-α is also involved in a number of pathological conditions including asthma, Crohn's disease, rheumatoid arthritis, neuropathic pain, obesity, type 2 diabetes, septic shock, autoimmunity, and cancer (5-11).

Mouse TNF-α is synthesized as a 26 kDa type II transmembrane protein that consists of a 35 amino acid (aa) cytoplasmic domain, a 21 aa transmembrane segment, and a 179 aa extracellular domain (ECD) (12). Within the ECD mouse TNF-α shares 95% aa sequence identity with rat, and 80% aa identity with canine, equine, feline, human, rabbit, and porcine TNF-α. It is produced by a wide variety of immune, epithelial, endothelial, and tumor cells. TNF-α is assembled intracellularly to form a nonconvalently linked homotrimer which is expressed on the cell surface (13). Cell surface TNF-α can both induce the lysis of tumor cells and virus infected cells, and generate its own downstream cell signaling following ligation by soluble TNF RI (14, 15). Shedding of membrane bound TNF-α by TACE/ADAM17 releases the bioactive cytokine, a 55 kDa soluble trimer containing the TNF-α extracellular domain (16-18).

TNF-α binds the ubiquitous 55-60kDa TNF RI (19, 20) and the hematopoietic cell restricted 78-80 kDa TNF RII (21, 22), both of which are also expressed as homotrimers (1, 23). Both type I and type II receptors bind TNF-α with comparable affinity and can promote NFκB activation (24-27). Only TNF RI, however, contains a cytoplasmic death domain which triggers the activation of apoptosis (3, 28). Soluble forms of both types of receptors are released into human serum and urine and can neutralize the biological activity of TNF (29-31).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for mouse TNF-α has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any mouse TNF-α present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked antibody specific for mouse TNF-α is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, TMB substrate (Chromogenic agent) is added to the wells and color develops in proportion to the amount of mouse TNF-α bound in the initial step. The color development is stopped and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- This kit is suitable for cell culture supernate and mouse serum, mouse plasma.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Calibrator Diluent RD5K/RD6-12 and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

CELL CULTURE SUPERNATE ASSAY

	Intra-assay Precision			**Inter-assay Precision**		
Sample	1	2	3	1	2	3
Mean (pg/mL)	54.6	116	405	52.6	138	373
Standard Deviation	1.50	4.51	12.5	4.61	8.49	29.7
CV%	2.7	3.9	3.1	8.8	6.2	8.0

SERUM/PLASMA ASSAY

	Intra-assay Precision			**Inter-assay Precision**		
Sample	1	2	3	1	2	3
Mean (pg/mL)	50.5	134	375	53.0	146	388
Standard Deviation	1.91	4.95	11.5	3.69	8.42	29.7
CV%	3.8	3.7	3.1	7.0	5.8	7.7

B. RECOVERY

The recovery of mouse TNF-α spiked to three levels throughout the range of the assay in various matrices was evaluated.

Sample Type	Average % Recovery	Range (%)
Cell culture supernates (n=6)	104	94 - 111
Mouse serum (n=4)	94	88 - 106
EDTA plasma (n=4)	95	85 - 103
Heparin plasma (n=4)	92	86 - 98

C. SENSITIVITY

One hundred seventeen assays were evaluated and the minimum detectable dose (MDD) of mouse TNF-α ranged from 0.36-7.21 pg/mL. The mean MDD was 1.88 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against a recombinant mouse TNF-α produced at R&D Systems®.

The Non WHO Reference Material 88/532 was evaluated in this kit. Each ampule contains a nominal 1 µg of recombinant mouse TNF-α, and was assigned an arbitrary unitage of 200,000 U/ampule. To convert sample values obtained with the Valukine™ Mouse TNF-α kit to approximate NIBSC units, use the equation below:

NIBSC/WHO (88/532) approximate value (U/mL) = 1.8553 x Valukine Mouse TNF-α (pg/mL)

Note: *Based on data generated in March 2017.*

E. LINEARITY

To assess the linearity of the assay, samples were spiked with various concentrations of mouse TNF-α in each matrix, diluted with the appropriate Calibrator Diluent RD5K/RD6-12 and then assayed.

		Cell culture supernates (n=4)	Serum (n=4)	EDTA plasma (n=4)	Heparin plasma (n=4)
1:2	Average % of Expected	106	103	104	105
	Range (%)	102 - 110	90 - 111	102 - 105	98 - 108
1:4	Average % of Expected	105	105	105	110
	Range (%)	98 - 112	92 - 110	101 - 109	105 - 115
1:8	Average % of Expected	100	107	107	112
	Range (%)	89 - 113	99 - 111	101 - 113	110 - 115
1:16	Average % of Expected	105	107	109	109
	Range (%)	99 - 110	97 - 112	101 - 115	101 - 117

F. SAMPLE VALUES

Serum/Plasma - Twenty individual mouse serum and plasma samples were evaluated for levels of mouse TNF-α in this assay. All samples measured less than the lowest Mouse TNF-α Standard, 10.9 pg/mL.

Cell Culture Supernates:

EL-4 mouse lymphoblast cells (2 x 10^5 cells/mL) were cultured in DMEM supplemented with 10% fetal bovine serum, 2 mM L-glutamine, 100 U/mL penicillin, and 100 µg/mL streptomycin sulfate. Cells were cultured unstimulated or stimulated with 10 µg/mL lipopolysaccharide (LPS) and 100 ng/mL recombinant mouse (rm) IL-10 for 4 days. Aliquots of the cell culture supernates were removed and assayed for levels of mouse TNF-α.

Condition	(pg/mL)
Unstimulated	ND
Stimulated	1220

ND=Non-detectable

RAW 264.7 mouse monocyte/macrophage cells (3.5 x 10^6 cells/mL) were cultured in DMEM supplemented with 10% fetal bovine serum, 2 mM L-glutamine, 100 U/mL penicillin, and 100 µg/mL streptomycin sulfate. Cells were cultured unstimulated or stimulated with 5.0 µg/mL LPS and 100 ng/mL rmIL-10 for 7 days. Aliquots of the cell culture supernates were removed and assayed for levels of mouse TNF-α.

Condition	(pg/mL)
Unstimulated	151
Stimulated	9605

Organs from individual mice were removed, rinsed in 1X PBS and kept on ice in tubes containing 1X PBS. Organs were cut into 1-2 mm pieces and homogenized using a tissue homogenizer. Cells were cultured in RPMI 1640 supplemented with 10% fetal bovine serum, 2 mM L-glutamine, 100 U/mL penicillin, and 100 µg/mL streptomycin sulfate. Cells were cultured unstimulated or stimulated with 1.0 µg/mL LPS for 1 day. Aliquots of the cell culture supernates were removed and assayed for levels of mouse TNF-α.

Tissue Type	Unstimulated (pg/mL)	Stimulated (pg/mL)
Mouse Liver	ND	25.7
Mouse Lung	55.9	97.3
Mouse spleen	ND	34.6

ND=Non-detectable

G. SPECIFICITY

This assay recognizes both natural and recombinant mouse TNF-α.

The following factors were prepared at 100 ng/mL in each calibrator diluent and assayed for cross-reactivity. Preparations of the following factors prepared at 100 ng/mL in a mid-range mouse TNF-α control were assayed for interference. No significant cross-reactivity or interference was observed.

Recombinant mouse	**Recombinant human**
CD40	TNF-α
CD40 ligand	**Recombinant canine**
Fas	TNF-α
Fas ligand	**Recombinant porcine**
LIF	TNF-α
OPG	
RANK	
RANK-L	
TNF-β	
TRAIL	
TROY	

47% cross-reactivity was observed with recombinant rat TNF-α.


Recombinant mouse TNF RI and TNF RII were found to interfere in this assay at concentrations above 1.25 ng/mL and 12.5 ng/mL, respectively.

IV. EXPERIMENT

EXAMPLE STANDARD


The standard curves are provided for demonstration only. A standard curve should be generated for each set of samples assayed.

CELL CULTURE SUPERNATE ASSAY

(pg/mL)	O.D.	Average	Corrected
0	0.020 0.021	0.021	—
10.9	0.060 0.061	0.061	0.040
21.9	0.099 0.101	0.100	0.079
43.8	0.179 0.181	0.180	0.159
87.5	0.336 0.344	0.340	0.319
175	0.632 0.637	0.635	0.614
350	1.203 1.242	1.223	1.202
700	2.265 2.281	2.273	2.252

SERUM/PLASMA ASSAY

(pg/mL)	O.D.	Average	Corrected
0	0.013 0.016	0.015	—
10.9	0.051 0.051	0.051	0.036
21.9	0.091 0.099	0.095	0.080
43.8	0.167 0.171	0.169	0.154
87.5	0.285 0.302	0.294	0.279
175	0.571 0.574	0.573	0.558
350	1.057 1.090	1.074	1.059
700	2.111 2.139	2.125	2.110

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
Mouse TNF-α Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with antibody against mouse TNF-α	1 plate
Mouse TNF-α Conjugate	Solution of an antibody against mouse TNF-α conjugated to horseradish peroxidase	1 vial
Mouse TNF-α Standard	Recombinant mouse TNF-α in a buffered protein base with preservatives; lyophilized. *Refer to the vial label for the reconstitution volume*	2 vials
Assay Diluent RD1-63	A buffered protein base with preservatives.	1 vial
Calibrator Diluent RD5K	A buffered protein base with preservatives. *For cell culture supernate samples.*	1 vial
Calibrator Diluent RD6-12	A buffered protein base with preservatives. *For serum/plasma samples.*	1 vial
Wash Buffer Concentrate (25×)	A 25× concentrated solution of buffered surfactant	1 vial
TMB Substrate	TMB ELISA Substrate Solution/TMB Substrate Solution	1 vial
Stop Solution	Diluted hydrochloric acid solution	1 vial
Plate Sealers	Adhesive strip	3 strips

B. STORAGE

<table>
<tr><td>Unopened Kit</td><td colspan="2">Store at 2-8 °C. Do not use past kit expiration date.</td></tr>
<tr><td rowspan="8">Opened/ Reconstituted Reagents</td><td>Wash Buffer (1×)</td><td rowspan="4">May be stored for up to 1 month at 2-8 °C.*</td></tr>
<tr><td>Stop Solution</td></tr>
<tr><td>Conjugate</td></tr>
<tr><td>TMB Substrate</td></tr>
<tr><td>Standard</td><td>Use a new standard for each assay. Discard after use.</td></tr>
<tr><td>Assay Diluent RD1-63</td><td rowspan="3">May be stored for up to 1 month at 2-8 °C.*</td></tr>
<tr><td>Calibrator Diluent RD5K</td></tr>
<tr><td>Calibrator Diluent RD6-12</td></tr>
<tr><td>Microplate Wells</td><td>Return unused wells to the foil ouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*</td></tr>
</table>

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- Pipettes and pipette tips.
- Deionized or distilled water.
- Squirt bottle, manifold dispenser, or automated microplate washer.
- 500 mL graduated cylinder.

D. PRECAUTION

- Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

The sample collection and storage conditions listed below are intended as general guidelines. Sample stability has not been evaluated.

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at ≤-20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent RD5K.

Serum - Allow blood samples to clot for 2 hours at room temperature before centrifuging for 20 minutes at 2000 x g. Remove serum and assay immediately or aliquot and store samples at ≤ -20 °C Avoid repeated freeze-thaw cycles.Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent RD6-12.

Plasma - Collect plasma using EDTA or heparin as an anticoagulant. Centrifuge for 20 minutes at 2000 x g within 30 minutes of collection. Assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent RD6-12.

Note: *Grossly hemolyzed samples are not suitable for use in this assay.*

Citrate plasma has not been validated for use in this assay.

B. REAGENT PREPARATION

Note: *Bring all reagents to room temperature before use.*

Wash Buffer (1×) - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25×) into deionized or distilled water to prepare 500 mL of Wash Buffer (1×).


Mouse TNF-α Standard - Refer to the vial label for the reconstitution volume*.

Reconstitute the Mouse TNF-α Standard with deionized or distilled water. This reconstitution produces a stock solution of 7000 pg/mL. Allow the standard to sit for a minimum of 5 minutes with gentle mixing prior to making dilutions.

*If you have any question, please contact our Technical Support.

Pipette 900 µL of Calibrator Diluent RD5K (*for cell culture supernate samples*) or

Calibrator Diluent RD6-12 (*for serum/plasma samples*) into the 700 pg/mL tube. Pipette 200 µL of appropriate Calibrator Diluent into each remaining tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The standard 700 pg/mL serves as the high standard. The appropriate Calibrator Diluent serves as the zero standard (0 pg/mL).

C. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB Substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB Substrate.

VII. ASSAY PROCEDURE

Note: *Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.*

1. Prepare all reagents, samples and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 50 µL of Assay Diluent RD1-63 to each well.
4. Add 50 µL of Standard and sample per well. Mix by gently tapping the plate frame for 1 minute. Cover with the adhesive strip provided. **Incubate for 2 hours at room temperature**. A plate layout is provided for a record of standards and samples assayed.
5. Aspirate each well and wash, repeating the process four times for a total of five washes. Wash by filling each well with Wash Buffer (400 µL) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
6. Add 100 µL of mouse TNF-α conjugate to each well. Cover with a new adhesive strip. **Incubate for 2 hours at room temperature**.
7. Repeat the aspiration/wash as in step 5.
8. Add 100 µL of TMB Substrate to each well. **Incubate for 30 minutes at room temperature**. **Protect from light**.
9. Add 100 µL of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
10. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

11. CALCULATION OF RESULTS.

Average the duplicate readings for each standard and sample and subtract the average zero standard optical density (O.D.).Create a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the mouse TNF-α concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.


If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Idriss, H.T. and J.H. Naismith (2000) Microsc. Res. Tech. 50:184.
2. Salek-Ardakani, S. and M. Croft (2010) J. Interferon Cytokine Res. 30:205.
3. Van Herreweghe, F. et al. (2010) Cell. Mol. Life Sci. 67:1567.
4. Chen, X. et al. (2009) Cell Biochem. Funct. 27:407.
5. Hehlgans, T. and K. Pfeffer (2005) Immunology 115:1.
6. Berry, M. et al. (2007) Curr. Opin. Pharmacol. 7:279.
7. D'Haens, G. (2003) Curr. Pharm. Des. 9:289.
8. Feldmann, M. and R.N. Maini (2001) Annu. Rev. Immunol. 19:163.
9. Leung, L. and C.M. Cahill (2010) J. Neuroinflammation 7:27.
10. Tzanavari, T. et al. (2010) Curr. Dir. Autoimmun. 11:145.
11. Wu, Y. and B.P. Zhou (2010) Br. J. Cancer 102:639.
12. Pennica, D. et al. (1984) Nature 312:724.
13. Tang, P. et al. (1996) Biochemistry 35:8216.
14. Perez, C. et al. (1990) Cell 63:251.
15. Watts, A.D. et al. (1999) EMBO J. 18:2119.
16. Black, R.A. et al. (1997) Nature 385:729.
17. Moss, M.L. et al. (1997) Nature 385:733.
18. Gearing, A.J.H. et al. (1994) Nature 370:555.
19. Schall, T.J. et al. (1990) Cell 61:361.
20. Loetscher, H. et al. (1990) Cell 61:351.
21. Dembic, Z. et al. (1990) Cytokine 2:231.
22. Smith, C.A. et al. (1990) Science 248:1019.
23. Loetscher, H. et al. (1991) J. Biol. Chem. 266:18324.
24. Rothe, M. et al. (1995) Science 269:1424.
25. Ruby, J. et al. (1997) J. Exp. Med. 186:1591.
26. Pinckard, J.K. et al. (1997) J. Biol. Chem. 272:10784.
27. Mukhopadhyay, A. et al. (2001) J. Biol. Chem. 276:31906.
28. Hsu, H. et al. (1995) Cell 81:495.
29. Seckinger, P. et al. (1989) J. Biol. Chem. 264:11966.
30. Olsson, I. et al. (1989) Eur. J. Haematol. 42:270.
31. Engelmann, H. et al. (1990) J. Biol. Chem. 265:1531.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

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