### **R** SYSTEMS a biotechne brand


# Recombinant SARS-CoV-2 B.1.1.7 Spike S1 Subunit His-tag

Catalog Number: 11137-CV

DESCRIPTION	
Source	Human embryonic kidney cell, HEK293-derived sars-cov-2 Spike S1 Subunit protein Val16-Pro681 (His69del, Val70del, Tyr145del, Asn501Tyr, Ala570Asp, Asp614Gly, Pro681His), with a C-terminal 6-His tag Accession # YP_009724390.1
N-terminal Sequence Analysis	Val16
Predicted Molecular Mass	75 kDa

SPECIFICATIONS	
SDS-PAGE	106-123 kDa, under reducing conditions.
Activity	Measured by its binding ability in a functional ELISA with Recombinant Human ACE-2 His-tag (Catalog # 933-ZN).
Endotoxin Level	<0.10 EU per 1 µg of the protein by the LAL method.
Purity	>95%, by SDS-PAGE visualized with Silver Staining and quantitative densitometry by Coomassie® Blue Staining.
Formulation	Lyophilized from a 0.2 μm filtered solution in PBS with Trehalose. See Certificate of Analysis for details.

PREPARATION AND STORAGE	
Reconstitution	Reconstitute at 500 μg/mL in PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.
	<ul> <li>12 months from date of receipt, -20 to -70 °C as supplied.</li> </ul>
	<ul> <li>1 month, 2 to 8 °C under sterile conditions after reconstitution.</li> </ul>
	<ul> <li>3 months, -20 to -70 °C under sterile conditions after reconstitution.</li> </ul>



µg/lane of Recombinant SARS-CoV-2 B.1.1.7 Spike S1 Subunit His-tag Protein (Catalog # 11137-CV) was resolved with SDS-PAGE under reducing (R) and non-reducing (NR) conditions and visualized by Coomassie® Blue staining, showing bands at 106 -

Rev. 5/18/2022 Page 1 of 2



Global bio-techne.com info@bio-techne.com techsupport@bio-techne.com TEL +1 612 379 2956 USA TEL 800 343 7475 Canada TEL 855 668 8722 China TEL +86 (21) 52380373 Europe | Middle East | Africa TEL +44 (0)1235 529449



## Recombinant SARS-CoV-2 B.1.1.7 Spike S1 Subunit His-tag

Catalog Number: 11137-CV

#### BACKGROUND

SARS-CoV-2, which causes the global pandemic coronavirus disease 2019 (Covid-19), belongs to a family of viruses known as coronaviruses that also include MERS-CoV and SARS-CoV-1. Coronaviruses are commonly comprised of four structural proteins: Spike protein (S), Envelope protein (E), Membrane protein (M) and Nucleocapsid protein (N) (1). The SARS-CoV-2 S protein is a glycoprotein that mediates membrane fusion and viral entry. The S protein is homotrimeric, with each ~180-kDa monomer consisting of two subunits, S1 and S2 (2). In SARS-CoV-2, as with most coronaviruses, proteolytic cleavage of the S protein into S1 and S2 subunits is required for activation. The S1 subunit is focused on attachment of the protein to the host receptor while the S2 subunit is involved with cell fusion (35). The S1 subunit of SARS-CoV-2 shares 65% and 22% amino acid sequence identity with that of SARS-CoV-1 and MERS, respectively. The S Protein of the SARS-CoV-2 virus, like the SARS-CoV-1 counterpart, binds a metallopeptidase, Angiotensin-Converting Enzyme 2 (ACE-2), but with much higher affinity and faster binding kinetics through the receptor binding domain (RBD) located in the C-terminal region of S1 subunit (6). It has been demonstrated that the S Protein can invade host cells through the CD147/EMMPRIN receptor and mediate membrane fusion (7, 8). Polyclonal antibodies to the RBD of the SARS-CoV-2 protein have been shown to inhibit interaction with the ACE-2 receptor, confirming RBD as an attractive target for vaccinations or antiviral therapy (9). There is also promising work showing that the RBD may be used to detect presence of neutralizing antibodies present in a patient's bloodstream, consistent with developed immunity after exposure to the SARS-CoV-2 (10). Several emerging SARS-CoV-2 genomes have been identified including the B 1.1.7 (United Kingdom) variant (11). The B 1.1.7 variant appears to more easily transmissible, exhibit increased viral loads and, potentially, be associated with higher mortality rates compared to pre

#### References:

- 1. Wu, F. et al. (2020) Nature 579:265
- 2. Tortorici, M.A. and D. Veesler (2019) Adv. Virus Res. 105:93.
- 3. Bosch, B.J. et al. (2003). J. Virol. 77:8801.
- 4. Belouzard, S. et al. (2009) Proc. Natl. Acad. Sci. 106:5871.
- 5. Millet, J.K. and G.R. Whittaker (2015) Virus Res. 202:120.
- 6. Ortega, J.T. *et al*. (2020) EXCLI J. **19**:410.
- 7. Wang, K. et al. (2020) Sig. Transduct. Target Ther. 5:283.
- 8. Isabel, et al. (2020) Sci. Rep. 10:14031.
- 9. Tai, W. et al. (2020) Cell. Mol. Immunol. 17:613.
- 10. Okba, N.M.A. et al. (2020) Emerg. Infect. Dis. https://doi.org/10.3201/eid2607.200841.
- 11. Kidd, M. et al. (2021) The Journal of Infectious Diseases 223:1666.
- 12. Zahradník, J. *et al.* (2021) Nat. Microbiol. **6**:1188.
- 13. Davies, N.G. (2020) Science 372:eabg3055.

Rev. 5/18/2022 Page 2 of 2



Global bio-techne.com info@bio-techne.com techsupport@bio-techne.com TEL +1 612 379 2956 USA TEL 800 343 7475 Canada TEL 855 668 8722 China TEL +86 (21) 52380373 Europe | Middle East | Africa TEL +44 (0)1235 529449