bio-techne[®]

Recombinant SARS-CoV-2 XBB.1 Spike RBD His-tag

Catalog Number: 11330-CV

RDsystems

DESCRIPTION				
Source	Human embryonic kidney cell, HEK293-derived sars-cov-2 Spike RBD protein			
	SARS-CoV-2 XBB.1 (Arg319-Phe541) (Gly339His, Arg346Thr, Leu368IIe, Ser371Phe, Ser373Pro, Ser375Phe, Thr376Ala, Asp405Asn, Arg408Ser, Lys417Asn, Asn440Lys, Val445Pro, Gly446Ser, Asn460Lys, Ser477Asn, Thr478Lys, Glu484Ala, Phe486Ser, Phe490Ser, Gln498Arg, Asn501Tyr, Tyr505His) Accession # YP_009724390.1	6-His tag		
	N-terminus	C-terminus		
N-terminal Sequence Analysis	Arg319			
Predicted Molecular Mass	26 kDa			
SPECIFICATIONS				

OF LOIT IOATIONS			
SDS-PAGE	30-37 kDa, under reducing conditions.		
Activity	Measured by its binding ability in a functional ELISA with Recombinant Human ACE-2 Fc Chimera (Catalog # 10544-ZN).		
Endotoxin Level	<0.10 EU per 1 µg of the protein by the LAL method.		
Purity	>95%, by SDS-PAGE under reducing conditions and visualized by Colloidal Coomassie® Blue stain at 5 µg per lane.		
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details.		

PREPARATION AND STORAGE		
ion Reconstitute at 500 μg/mL in PBS.		
The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.		
Use a manual defrost freezer and avoid repeated freeze-thaw cycles.		
 12 months from date of receipt, -20 to -70 °C as supplied. 		
 1 month, 2 to 8 °C under sterile conditions after reconstitution. 		

3 months, -20 to -70 °C under sterile conditions after reconstitution.



Rev. 1/23/2023 Page 1 of 2



Global bio-techne.com info@bio-techne.com techsupport@bio-techne.com TEL +1 612 379 2956 USA TEL 800 343 7475 Canada TEL 855 668 8722 China TEL +86 (21) 52380373 Europe | Middle East | Africa TEL +44 (0)1235 529449

biotechne

Recombinant SARS-CoV-2 XBB.1 Spike RBD His-tag

Catalog Number: 11330-CV

RDsystems

BACKGROUND

SARS-CoV-2, which causes the global pandemic coronavirus disease 2019 (Covid-19), belongs to a family of viruses known as coronaviruses that also include MERS-CoV and SARS-CoV-1. Coronaviruses are commonly comprised of four structural proteins: Spike protein (S), Envelope protein (E), Membrane protein (M) and Nucleocapsid protein (N) (1). The SARS-CoV-2 S protein is a glycoprotein that mediates membrane fusion and viral entry. The S protein is homotrimeric, with each ~180-kDa monomer consisting of two subunits, S1 and S2 (2). In SARS-CoV-2, as with most coronaviruses, proteolytic cleavage of the S protein into S1 and S2 subunits is required for activation. The S1 subunit is focused on attachment of the protein to the host receptor while the S2 subunit is involved with cell fusion (3-5). The S Protein of the SARS-CoV-2 virus, like the SARS-CoV-1 counterpart, binds a metallopeptidase, Angiotensin-Converting Enzyme 2 (ACE-2), but with much higher affinity and faster binding kinetics through the receptor binding domain (RBD) located in the C-terminal region of S1 subunit (6). It has been demonstrated that the S Protein can invade host cells through the CD147/EMMPRIN receptor and mediate membrane fusion (7, 8). Polyclonal antibodies to the RBD of the SARS-CoV-2 protein have been shown to inhibit interaction with the ACE-2 receptor, confirming RBD as an attractive target for vaccinations or antiviral therapy (9). There is also promising work showing that the RBD may be used to detect presence of neutralizing antibodies present in a patient's bloodstream, consistent with developed immunity after exposure to the SARS-CoV-2 (10). Based on amino acid (aa) sequence homology, the SARS-CoV-2 RBD has 73% identity with SARS-CoV-1 RBD, but only 22% homology with the MERS RBD. Several emerging SARS-CoV-2 genomes have been identified including the Omicron, or B.1.1.529, variant. First identified in November 2021 in South Africa, the Omicron variant quickly became the predominant SARS-CoV-2 variant and is considered a variant of concern (VOC). The Omicron variant contains 32 mutations in the S protein, 3 to 4 times more than in other SARS-CoV-2 variants, that potentially affect viral fitness and transmissibility (11). Of these mutations, 15 are located in the RBD domain and allow the Omicron variant to bind ACE-2 with greater affinity and, potentially, increased transmissibility (11, 12). Several additional mutations throughout the S protein have been shown or are predicted to enhance spike cleavage and could aid transmission (13-15). The study of the Omicron variant's impact on immune escape and reduced neutralization activity to monoclonal antibodies along with an increased risk of reinfection, even among vaccinated individuals, remains ongoing (16). The XBB.1 variant is a combination of BA.2.10.1 and BA.2.75 sublineage, and is predicted to resist some COVID treatment (17).

References:

- 1. Wu, F. et al. (2020) Nature 579:265.
- 2. Tortorici, M.A. and D. Veesler (2019) Adv. Virus Res. 105:93
- 3. Bosch, B.J. et al. (2003) J. Virol. 77:8801.
- 4. Belouzard, S. et al. (2009) Proc. Natl. Acad. Sci. 106:5871.
- 5. Millet, J.K. and G.R. Whittaker (2015) Virus Res. 202:120.
- 6. Ortega, J.T. et al. (2020) EXCLI J. 19:410.
- 7. Wang, K. et al. (2020) bioRxiv https://www.biorxiv.org/content/10.1101/2020.03.14.988345v1.
- 8. Isabel, S. et al. (2020) Sci Rep. 10, 14031. <u>https://doi.org/10.1038/s41598-020-70827-z</u>.
- 9. Tai, W. et al. (2020) Cell. Mol. Immunol. 17:613.
- 10. Okba, N. M. A. et al. (2020). Emerg. Infect. Dis. https://doi.org/10.3201/eid2607.200841.
- 11. Shah, M. and Woo, H.G. (2021) bioRxiv https://doi.org/10.1101/2021.12.04.471200.
- 12. Lupala, C.S. *et al.* (2021) bioRxiv <u>https://doi.org/10.1101/2021.12.10.472102</u>.
- 13. Zhang, L. et al. (2020) Nat Commun. 11:6013.
- 14. Lasek-Nesselquist, E. et al. (2021) medRxiv https://doi.org/10.1101/2021.03.10.21253285.
- 15. Scheepers, C. et al. (2021) medRxiv https://doi.org/10.1101/2021.08.20.21262342.
- 16. Callaway, E. and Ledford, H. (2021) Nature 600:197.
- 17. Callaway E. et al. (2022) Nature 611:213.

Rev. 1/23/2023 Page 2 of 2



Global bio-techne.com info@bio-techne.com techsupport@bio-techne.com TEL +1 612 379 2956 USA TEL 800 343 7475 Canada TEL 855 668 8722 China TEL +86 (21) 52380373 Europe | Middle East | Africa TEL +44 (0)1235 529449