bio-techne[®] RDSYSTEMS

Catalog Number: 5439-DK/CF

DESCRIPTION	
Source	<i>Spodoptera frugiperda, Sf</i> 21 (baculovirus)-derived human Dkk-1 protein Thr32-His266 Accession # O94907
N-terminal Sequence Analysis	Thr32
Predicted Molecular Mass	25.8 kDa

SPECIFICATIONS	
SDS-PAGE	33-38 kDa, reducing conditions
Activity	Measured by its ability to inhibit Wnt induced TCF reporter activity in HEK293 human embryonic kidney cells. The ED ₅₀ for this effect is 10.0-250 ng/mL.
Endotoxin Level	<1.0 EU per 1 μ g of the protein by the LAL method.
Purity	>95%, by SDS-PAGE with silver staining.
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS. See Certificate of Analysis for details.


PREPARATION AND STORAGE	
Reconstitution	Reconstitute at 100 μg/mL in PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	 Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 3 months, -20 to -70 °C under sterile conditions after reconstitution.

DATA


bio-techne® RDsystems

Recombinant Human Dkk-1

Catalog Number: 5439-DK/CF



Cell Culture



iPSC-derived Cardiomyocytes **Differentiated on Cultrex Stem** Cell Qualified RGF BME Express Stage-specific Markers. JOY6 human iPSCs were cultured in media containing Cultrex[™] Stem Cell Qualified RGF BME (Catalog # 3434-010-02) and N21-MAX Insulin Free Media Supplement (Catalog # AR010) along with Recombinant Human Activin A on day 0, Recombinant Human FGF-basic (Catalog # 4114-TC), Recombinant Human BMP-4, and CHIR99021 (Catalog # 4423) on days 1-5, and Recombinant Human Dkk-1 (Catalog # 5439-DK/CF) on days 5-7 to induce cardiomyocyte differentiation. Cells were then cultured in media supplemented with N21-MAX Insulin Free Media Supplement (Catalog # AR010) on days 7-12 and media supplemented with N21-MAX Media Supplement (Catalog # AR008) from day 12 and beyond. Cells were fixed and stained for stage-specific markers at select time points during the procedure. The pluripotency marker Oct-4 was detected using a Mouse Anti-Human Oct-4A Monoclonal Antibody (Catalog # MAB17591). The mesoderm marker, Snail is expressed intermediately during differentiation (Day 1). It was detected using a Goat Anti-Human Snail Polyclonal Antibody (Catalog #AF3639). The cardiomyocyte markers NKX2.5 and Troponin T are not present in cells during early (Day 0) and intermediate (Day1) differentiation and become more highly expressed during the later stages of differentiation (Day 7, Day 30). NKX2.5 was detected using a Goat Anti-Human NKX2.5 Polyclonal Antibody (Catalog # AF2444) and Troponin T was detected using a Mouse Anti-Human Cardiac Troponin T Monoclonal Antibody (Catalog # MAB1874). Snail and NKX2.5 primary antibodies were visualized with a NorthernLights™ 557conjugated Donkey Anti-Goat IgG Secondary Antibody (Catalog # NL001). Oct-4 and Troponin T were visualized with a NorthernLights 557-conjugated Donkey Anti-Mouse IgG Secondary Antibody (Catalog # NL007).


Rev. 5/14/2024 Page 2 of 3

Global | bio-techne.com info@bio-techne.com techsupport@bio-techne.com TEL: 1.612.379.2956 USA | TEL: 800.343.7475 Canada | TEL: 855.668.8722 Europe | Middle East | Africa TEL: +44.0.1235.529449 China | info.cn@bio-techne.com TEL: 400.821.3475



Catalog Number: 5439-DK/CF

Cell Culture



iPSC-derived Cardiomyocytes Express Atrial and Ventricular Markers. JOY6 human induced pluripotent stem cells (iPSCs) were differentiated into cardiomyocytes in media containing Cultrex™ Stem Cell Qualified Reduced Growth Factor Basement Membrane Extract (Catalog # 3434-010-02) and N21-MAX Insulin Free Media Supplement (Catalog # AR010) along with Recombinant Human Activin A on day 0, Recombinant Human FGF-basic (Catalog # 4114-TC), Recombinant Human BMP-4, and CHIR99021 (Catalog # 4423) on days 1-5, and Recombinant Human Dkk-1 (Catalog # 5439-DK/CF) on days 5-7. Following culture in media supplemented with N21-MAX Insulin Free Media Supplement (Catalog # AR010) on days 7-12 and media supplemented with N21-MAX Media Supplement (Catalog # AR008) from day 12 and beyond, cells were fixed and stained for the atrial-specific marker, MLC2a, and the ventriclespecific marker, MYH7 using a Rabbit Anti-Human MYH7 Monoclonal Antibody (Catalog # MAB90961).

BACKGROUND

Dickkopf related protein 1 (Dkk-1) is the founding member of the Dickkopf family of proteins that includes Dkk-1, -2, -3, -4, and a related protein, Soggy (1, 2). Dkk proteins are secreted proteins that contain two conserved cysteine-rich domains separated by a linker region. Each domain contains ten cysteine residues (1-3). Mature human Dkk-1 is a 40 kDa glycosylated protein that shares 86%, 87%, 90% and 91% aa sequence identity with mouse, rat, rabbit and bovine Dkk-1, respectively. It also shares 42% and 36% aa identity with human Dkk-2 and Dkk-4, respectively. Dkk-1 and Dkk-4 are well documented antagonists of the canonical Wnt signaling pathway (1, 2). This pathway is activated by Wnt engagement of a receptor complex composed of the Frizzled proteins and one of two low-density lipoprotein receptor-related proteins, LRP5 or LRP6 (4). Dkk-1 antagonizes Wnt by forming ternary complexes of LRP5/6 with Kremen1 or Kremen2 (4, 5). Dkk-1/LRP6/Krm2 complex internalization has been shown to down-regulate Wnt signaling (4, 5). Dkk-1 is expressed throughout development and antagonizes Wnt-7a during limb development (6, 7). Other sites of expression include developing neurons, hair follicles and the retina of the eye (8, 9). The balance between Wnt signaling and Dkk-1 includes, Dkk-1 is expressed in osteoblasts and osteocytes, and neurons. Cerebral ischemia induces Dkk-1 expression, which contributes to neuronal cell death (12).

References:

- 1. Krupnik, V.E. et al. (1999) Gene 238:301.
- 2. Niehrs, C. (2006) Oncogene 25:7469.
- 3. Bullock, C.M. et al. (2004) Mol. Pharmacol. 65:582.
- 4. Mao, B. et al. (2001) Nature 411:321.
- 5. Mao, B. et al. (2002) Nature 417:664.
- 6. Kemp, C. et al. (2005) Dev. Dyn. 233:1064.
- 7. Adamska, M. et al. (2004) Dev. Biol. 272:134.
- 8. Li, J. et al. (2006) Bone 36:754.
- 9. Verani, R. et al. (2006) J. Neurochem. 101:242.
- 10. Pinzone, J.J. et al. (2009) Blood 113:517.
- 11. Morvan, F. et al. (2006) J. Bone Miner. Res. 21:934.
- 12. Cappuccio, I. et al. (2005) J. Neurosci. 25:2647.

Rev. 5/14/2024 Page 3 of 3
Bio-Techne®
Global | bio-techne.com info@bio-techne.com techsupport@bio-techne.com TEL: 1.612.379.2956
USA | TEL: 800.343.7475 Canada | TEL: 855.668.8722 Europe | Middle East | Africa TEL: +44.0.1235.529449
China | info.cn@bio-techne.com TEL: 400.821.3475