

## PRODUCT DESCRIPTION

GMP N-2 MAX Media Supplement is manufactured using animal-free raw materials in an animal-free laboratory following cGMP guidelines. This media supplement can be added to base media formulations to support the expansion of neural stem cells (1) under cGMP-grade culture conditions. GMP N-2 MAX Media Supplement has been modified from the Bottenstein's formulation (2) to be compliant in *ex vivo* cell and tissue manufacturing work flows. It is supplied as a 100X concentrate.

## STABILITY & STORAGE

Product is shipped on dry ice. Store **in the dark** at  $\leq -20$  °C in a manual defrost freezer. Do not use beyond the expiration date.

## MEDIA SUPPLEMENT COMPONENTS

COMPONENT	AMOUNT (µg/mL)
Recombinant Human Insulin	2500
Recombinant Human Transferrin	10,000
Putrescine	1611
Selenite	0.52
Progesterone	0.63

## PRODUCT TESTING

TEST	SPECIFICATION
Endotoxin Level	$\leq 2$ EU/mL for a 1X solution
Mycoplasma	Negative when tested in a ribosomal RNA hybridization assay
Sterility	Tested to USP <71> sterility guidelines
Activity	Supports expansion, maintenance, and differentiation of rat cortical stem cells

## SAFETY INFORMATION

- For use in preclinical research or for *ex vivo* cell and tissue manufacturing.
- Not intended for direct administration into humans or animals. Not for parenteral use.
- The safety and efficacy of this product in diagnostic or other clinical uses has not been established.

## PRECAUTION

When handling bio-hazardous materials such as human cells, safe laboratory procedures should be followed and protective clothing should be worn.

## LIMITATIONS

- This reagent should not be used beyond the expiration date indicated on the label.
- Do not use if package is damaged. Use undamaged and sealed bottles only.
- Results may vary due to variations among cells derived from different donors.

## MEDIA PREPARATION

Dilute GMP N-2 MAX Media Supplement 1:100 in basal media before use. The diluted medium may be stored **in the dark** at 2-8 °C for up to 2 weeks.

**Note:** For *ex vivo* cell or tissue manufacturing applications, this supplement should be used in combination with GMP-compliant base media and raw materials.

## MANUFACTURING SPECIFICATIONS

R&D Systems®, a Bio-Techne® Brand's GMP media are produced according to relevant sections of the following documents: USP Chapter 1043, Ancillary Materials for Cell, Gene and Tissue-Engineered Products and Eu. Ph. 5.2.12, Raw Materials of Biological Origin for the Production of Cell-based and Gene Therapy Medicinal Products.

R&D Systems quality focus includes:

- Manufacturing and testing under an ISO 9001:2015 and ISO 13485:2016 certified quality system.
- Documented processes and QA control of documentation and process changes.
- Personnel training programs.
- Raw material testing and vendor qualification/monitoring.
- Fully validated equipment, processes and test methods.
- Equipment calibration schedules using a computerized calibration program.
- Facility maintenance, safety programs and pest control.
- Material review process for variances.
- Monitoring of stability over product shelf life.

Additional testing and documentation requested by the customer can be arranged at an additional cost.

Production records and facilities are available for examination by appropriate personnel on site at R&D System in Minneapolis, Minnesota USA.

R&D System sells GMP products for preclinical or clinical *ex vivo* use. They are not for *in vivo* use. Read the following End of Terms prior to using this product.

## END USER TERMS OF USE OF PRODUCT

The following terms are offered to you upon your acceptance of these End User Terms of Use of Product. By using this product, you indicate your acknowledgment and agreement to these End User Terms of Use of Product. If you do not agree to be bound by and comply with all of the provisions of these End User Terms of Use of Product, you should contact your supplier of the product and make arrangements to return the product.

We suggest you print and retain a copy of these End User Terms of Use of Product for your records.

The End User is aware that R&D Systems, Inc. sells GMP products for preclinical or clinical *ex vivo* use and not for *in vivo* use. The End User further agrees, as a condition of the sale of R&D Systems GMP products that: a) the End User will not use this GMP Product in any procedure wherein the product may be administered to humans, unless the End User has obtained, or prior to their use will have obtained, an Investigational New Drug (IND) exemption from the FDA and will use the product only in accordance with the protocols of such IND and of the Institutional Review Board overseeing the proposed research, or b) the End User will use the products outside of the United States in accordance with the protocols of research approved by the Institutional Review Board or authorized ethics committee and regulatory agencies to which the End User is subject to in their territory.

R&D Systems, Inc. has the right, at its sole discretion, to modify, add or remove any terms or conditions of these End User Terms of Use without notice or liability to you. Any changes to these End User Terms of Use are effective immediately following the printing of such changes on this product insert. The most recent version of these End User Terms of Use of Product may be found at: [RnDSystems.com/Legal](http://RnDSystems.com/Legal).

You agree to review these End User Terms of Use of Product to ensure any subsequent use by you of R&D Systems GMP Products following changes to these End User Terms of Use of Product constitutes your acceptance of all such changes.

## REFERENCES

1. Johe, K.K. *et al.* (1996) *Genes & Development* **10**:3129.
2. Bottenstein, J.E. (1985) *Cell Culture in the Neurosciences*, Plenum Press: New York and London.