

## Human GluR7/GRIK3 Alexa Fluor® 700-conjugated Antibody

Recombinant Monoclonal Rabbit IgG Clone # 2914A Catalog Number: FAB113281N

100 µg

DESCRIPTION		
Species Reactivity	Human	
Specificity	Detects human GluR7/GRIK3 in direct ELISA.	
Source	Recombinant Monoclonal Rabbit IgG Clone # 2914A	
Purification	Protein A or G purified from hybridoma culture supernatant	
Immunogen	Human embryonic kidney cell, HEK293-derived human GluR7/GRIK3 Gly249-Met1 Accession # Q13003	
Conjugate	Alexa Fluor 700 Excitation Wavelength: 675-700 nm Emission Wavelength: 723 nm	
Formulation	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide.  *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet	
	(SDS) for additional information and handling instructions.	

### APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. General Protocols are available in the Technical Information section on our website.

Flow Cytometry	Titration recommended for optimal concentration with starting range of 0.1-1 μg/1 million cells. Sample used for this
	experiment was CHO cells transfected with hGRIK3 vs irrelevant CHO transfectant cells.

# PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.	
Stability & Storage Protect from light. Do not freeze.		
	<ul> <li>12 months from date of receipt, 2 to 8 °C as supplied.</li> </ul>	

### BACKGROUND

GluR7 is a receptor for glutamate. Glutamate receptors are the predominant excitatory neurotransmitter receptors in the mammalian brain. This gene can co-assemble with GRIK4 or GRIK5 to form heteromeric receptors. GluR7 is one of the five subtypes of kainate receptors. Kainate receptors are widely expressed in the CNS and are involved in the regulation of transmitter release, synapse formation, and in the pathophysiology of brain diseases. Genetic variants in the GluR7 gene are associated with schizophrenia, major depression, and bipolar disorder.

- 1. Entrez Gene: GRIK3 glutamate receptor, ionotropic, kainate 3.
- 2. Matute C. "Therapeutic Potential of Kainate Receptors". CNS Neuroscience & Therapeutics. 2011 Dec; 17(6):661.
- 3. Begni S, Popoli M, Moraschi S, Bignotti S, Tura G.B, Gennarelli M., "Association Between the Ionotropic Glutamate Receptor Kainate 3 (GRI3) Ser310Ala Polymorphism and Schizophrenia". Mol Psychiatry. 2002; 7(4):416.

## PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Rev. 8/3/2023 Page 1 of 1