bio-techne® RD SYSTEMS

Catalog Number: IC1195N 100 µg

	DESCRIPTION	
Purification Protein A or G purified from hybridoma culture supernatant Immunogen Rat brain-derived microtubules Conjugate Alexa Fluor 700 Excitation Wavelength: 675-700 nm Emission Wavelength: 723 nm Formulation Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data S	Specificity	Detects mammalian and chicken neuron-specific β-III tubulin but not other β-tubulin isotypes in Western blots.
Immunogen Rat brain-derived microtubules Conjugate Alexa Fluor 700 Excitation Wavelength: 675-700 nm Emission Wavelength: 723 nm Formulation Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data S	Source	Monoclonal Mouse IgG _{2A} Clone # TuJ-1
Conjugate Alexa Fluor 700 Excitation Wavelength: 675-700 nm Emission Wavelength: 723 nm Formulation Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data S	Purification	Protein A or G purified from hybridoma culture supernatant
Excitation Wavelength: 675-700 nm Emission Wavelength: 723 nm Formulation Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data S	Immunogen	Rat brain-derived microtubules
*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data S	Conjugate	Excitation Wavelength: 675-700 nm
	Formulation	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide.
(SDS) for additional information and handling instructions.		*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Shee (SDS) for additional information and handling instructions.
	APPLICATIONS	
APPLICATIONS	Diseas Nata: Ontinnal a	tilutions should be determined by each laboratory for each application. Concrete Protocole are evolution in the Technical Information section on our website

Please Note: Optimal dilutions should be determined by each laboratory for each application. General Protocols are available in the Technical Information section on our website.			
Intracellular Staining by Flow Cytometry	Titration recommended for optimal concentration with starting range of 0.1-1 μg/1 million cells. Sample used for this experiment was HepG2 human hepatocellular carcinoma cells.		

PREPARATION AND STORAGE	
Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze.
	 12 months from date of receipt, 2 to 8 °C as supplied.

BACKGROUND

β-III Tubulin, also known as tubulin β-4, is regarded as a neuron-specific marker. The expression of β-III Tubulin has been suggested to be one of the earliest markers to signal commitment in primitive neuroepithelium.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Rev. 8/23/2023 Page 1 of 1



Global bio-techne.com info@bio-techne.com techsupport@bio-techne.com TEL +1 612 379 2956 USA TEL 800 343 7475 Canada TEL 855 668 8722 China TEL +86 (21) 52380373 Europe | Middle East | Africa TEL +44 (0)1235 529449