RD SYSTEMS a biotechne brand

Human IRAK4 Alexa Fluor® 405-conjugated Antibody

Monoclonal Rat IgG₁ Clone # 413613 Catalog Number: IC39191V 100 µg

| DESCRIPTION | | | |
|--------------------|--|--|--|
| Species Reactivity | Human | | |
| Specificity | Detects human IRAK4 in direct ELISAs. | | |
| Source | Monoclonal Rat IgG ₁ Clone # 413613 | | |
| Purification | Protein A or G purified from hybridoma culture supernatant | | |
| Immunogen | <i>E. coli-</i> derived recombinant human IRAK4 Met1-Ser460 Accession # Q9NWZ3 | | |
| Conjugate | Alexa Fluor 405 Excitation Wavelength: 405 nm Emission Wavelength: 421 nm | | |
| Formulation | Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. | | |

*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

| APPLICATIONS | | | | |
|---|---------------------------------|--|--|--|
| Please Note: Optimal dilutions should be determined by each laboratory for each application. General Protocols are available in the Technical Information section on our website. | | | | |
| | Recommended Concentration | Sample | | |
| Intracellular Staining by Flow Cytometry | 0.25-1 µg/10 ⁶ cells | K562 and Jurkat human cell line fixed with Flow Cytometry Fixation Buffer (Catalog # FC004) and permeabilized with Flow Cytometry Permeabilization/Wash Buffer I (Catalog # FC005) | | |

| PREPARATION AND STORAGE | | |
|-------------------------|---|--|
| Shipping | The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below. | |
| Stability & Storage | Protect from light. Do not freeze. | |
| | 12 months from date of receipt, 2 to 8 °C as supplied. | |

BACKGROUND

IL-1 receptor-associated kinases (IRAKs) are serine/threonine kinases that help mediate signaling from Toll-like receptor (TLR) and IL-1 receptor family members. Four human IRAKs have been identified: IRAK1, IRAK2, IRAK-M, and IRAK4. Upon TLR ligand challenge, IRAK4 knockout mice exhibit severely impaired NF- κ B activation. IRAK4 mutations have been described in patients with recurrent bacterial infections and poor inflammatory responses.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Rev. 1/15/2020 Page 1 of 1



Global bio-techne.com info@bio-techne.com techsupport@bio-techne.com TEL +1 612 379 2956 USA TEL 800 343 7475 **Canada** TEL 855 668 8722 **China** TEL +86 (21) 52380373 **Europe | Middle East | Africa** TEL +44 (0)1235 529449