

Catalog Number: LUH206

Pack Size: 100 Tests



#### **SPECIFICATIONS AND USE**


Recommended Sample Types Microparticle Region	<ul> <li>Cell culture supernates, serum, EDTA plasma, and heparin plasma.</li> <li>Region-32</li> </ul>
Components	• Microparticle Concentrate (Part 891062) is supplied as a 100X concentrated stock (0.075 mL) with preservatives.
	• Biotin-Antibody Concentrate (Part 892622) is supplied as a 100X concentrated stock solution (0.075 mL) with preservatives.
Other Supplies Required	• Luminex Performance Assay Human Base Kit A (Catalog Number LUH000) or Luminex Performance Assay Human Obesity Base Kit (Catalog Number LOB000).
Storage	<ul> <li>Store the unopened kit at 2-8 °C. Do not use past the expiration date on the label.</li> <li>Avoid freezing microparticles.</li> <li>Protect microparticles from light.</li> </ul>
	Defense the environments Dage (it is east for the Lympiney Device means Access are and

Instructions for Use

• Refer to the appropriate Base Kit insert for the Luminex Performance Assay procedure.

# **TYPICAL DATA**

This human IL-6 standard curve is provided only for demonstration. A standard curve must be generated each time an assay is run, utilizing values from the Standard Value Card included in the Base Kit.



Standard	pg/mL	MFI	Average	Corrected	
Blank	0	18 20	19		
1	2400	4875	5112	5094	
		5348		5071	
2	800	2910	3063	3045	
Z	800	3216	2002	3045	
3	267	1540	1570	1552	
C		1599	1370	1332	
4	89	575	606	588	
4		636	000	300	
5	30	256	266	248	
2	50	275	200	248	
6	10	104	107	89	
	10	109	107	09	
7	3	48	49	31	
/	2	50	49	51	

## **PERFORMANCE CHARACTERISTICS**

#### All data were collected with assays run as a multiplex.

**Sensitivity** - The Minimum Detectable Dose (MDD) was determined by adding two standard deviations to the mean MFI of twenty zero standard replicates and calculating the corresponding concentration. Forty-three assays were evaluated, and the MDD of human IL-6 ranged from 0.10-1.11 pg/mL. The mean MDD was 0.36 pg/mL. Intra-assay Precision (precision within an assay) - Three samples of known concentration were tested twenty times on one plate to assess precision within an assay.

Inter-assay Precision (precision between assays) - Three samples of known concentration were tested in twenty-five separate assays to assess precision between assays.

	Intra-assay Precision				Inter-assay Precision			
Sample	1	2	3		1	2	3	
n	20	20	20		25	25	25	
Mean (pg/mL)	30.2	180	903		34	189	900	
Standard Deviation	1.43	7.78	41.9		3.00	13.0	53.0	
% CV	4.7	4.3	4.6		8.8	6.9	5.9	

**Recovery and Linearity** - Samples containing and/or spiked with high concentrations of IL-6 were evaluated for recovery and were serially diluted to evaluate assay linearity.

Recovery					Linearity					
Sample Type	Average % Recovery	Range (%)				Cell culture supernates	Serum	EDTA Plasma	Heparin Plasma	
Cell culture 96 supernates	06	91-100		1:2	Average % of Expected	101	106	103	99	
	90			1:2	Range (%)	87-117	104-109	82-124	86-110	
Serum	108	95-118		1:4	Average % of Expected	103	112	122	101	
	108				Range (%)	86-127	103-117	116-129	85-114	
EDTA plasma	102	95-116	1:	1.0	Average % of Expected	102	115	126	100	
	102			1:8	Range (%)	82-130	97-126	121-131	87-110	
Heparin plasma	109	96-129								

**Specificity** - This assay recognizes natural and recombinant human IL-6. The assay was tested for cross-reactivity and interference with the following factors. Less than 0.5% cross-reactivity and interference was observed.

Recombinant human:			Recombinant mouse:	Recombinant rat:	Recombinant porcine:	Recombinant human multiplex partners	Recombinant human multiplex partners
6Ckine	IL-3 Ra	MCP-3	G-CSF	GM-CSF	GM-CSF	Panel A:	Obesity Panel:
CNTF	IL-4 R	MCP-4	GM-CSF	IFN-γ	IL-1α	ENA-78	Adiponectin
β-ECGF	IL-5 Ra	M-CSF	IFN-γ	IL-1α	IL-1β	FGF basic	CRP
FGF acidic	IL-6 R	TNF RI	IL-1α	IL-1β	IL-2	G-CSF	Factor D
FGF-4	IL-10 R	TNF-β	IL-1β	IL-2	IL-4	GM-CSF	IL-10
FGF-5	IL-3	VEGF <sub>121</sub>	IL-1ra	IL-4	IL-6	IFN-γ	Leptin
FGF-6	IL-7	VEGF165	IL-2	IL-6	IL-8	IL-1α	MCP-1
FGF-9	IL-9	VEGF-D	IL-4	IL-10	IL-10	IL-1β	Resistin
FGF-10	IL-11		IL-5	TNF-α		IL-1ra	Serpin E1
FGF-18	IL-12/IL-23 p40		IL-6			IL-2	TNF-a
GCP-2	IL-12 p70		IL-10			IL-4	
GROα	IL-13		IL-17			IL-5	
GR0β	IL-15		MIP-1a			IL-8	
GR0γ	IL-16		MIP-1β			IL-10	
I-309	IL-17		RANTES			IL-17	
IGF-I	IL-18		TNF-a			MCP-1	
IGF-II	LIF		Тро			MIP-1a	
IL-1 RI	LIF R					MIP-1β	
IL-1 RII	MIP-18					RANTES	
IL-2 Ra	MIP-3a					TNF-a	
IL-2 Rβ	ΜΙΡ-3β					Тро	
IL-2 Rγ	MCP-2					VEGF	

## **TECHNICAL HINTS**

- Protect the microparticles and streptavidin-PE from light at all times.
- Refer to the Base Kit Standard Value Card for reconstitution volume and values of the reconstituted standard.
- Diluted microparticles cannot be stored. Make a fresh dilution of microparticles each time the assay is run.
- The use of a vacuum manifold device made to accommodate a microplate is necessary for washing. Adjust the vacuum to between 15 and 40 cm Hg.
- Discrepancies may exist in values obtained for the same analyte utilizing different technologies.

Luminex Performance Assays afford the user the benefit of multi-analyte analysis of biomarkers in a complex sample. For each sample type, a single, multipurpose diluent is used to optimize recovery, linearity, and reproducibility. Such a multipurpose diluent may not optimize any single analyte to the same degree that a unique diluent selected for analysis of that analyte can optimize conditions. Therefore, some performance characteristics may be more variable than those for assays designed specifically for single analyte analysis.

All trademarks and registered trademarks are the property of their respective owners.