

Magnetic Luminex® Performance Assay Human IL-10 Kit

Catalog Number: LUHM217
Pack Size: 100 Tests

Recommended Sample Types

• Cell culture supernates, serum, EDTA plasma, and heparin plasma

Microparticle Region

• Region-29

Components

• Microparticle Concentrate (Part 894438) is supplied as a 100X concentrated stock (0.075 mL) with preservatives.

• Biotin-Antibody Concentrate (Part 892624) is supplied as a 100X concentrated stock

solution (0.075 mL) with preservatives.

Other Supplies Required

• Magnetic Luminex® Performance Assay Human Base Kit A (R&D Systems®,

Catalog # LUHM000).

or

Magnetic Luminex® Performance Assay Human Obesity Base Kit (R&D Systems®,

Catalog # LOBM000).

Storage

• Store the unopened kit at 2-8 °C. Do not use past the expiration date on the label.

· Avoid freezing microparticles.


· Protect microparticles from light.

Instructions for Use

• Refer to the Base Kit insert for the Magnetic Luminex® Performance Assay procedure.

TYPICAL DATA

This human IL-10 standard curve is provided only for demonstration. A standard curve must be generated each time an assay is run, utilizing values from the Standard Value Card included in the Base Kit.



Standard	(pg/mL)	MFI	Average	Corrected	
Blank	0	22	22	_	
		22			
1	3200	20,050	20,065	20,043	
		20,079			
2	1067	9300	9447	9425	
		9594			
3	356	3871	3889	3867	
		3907			
4	119	1671	1696	1674	
		1721			
5	40	734	736	714	
		738			
6	13	315	318	296	
		320			
7	4	157	159	137	
		161			
-				-	

PERFORMANCE CHARACTERISTICS

All data were collected with assays run as a multiplex.

Data obtained with polystyrene and magnetic beads were equivalent.

Sensitivity - The Minimum Detectable Dose (MDD) was determined by adding two standard deviations to the MFI of twenty zero standard replicates and calculating the corresponding concentration.

Forty-three assays were evaluated, and the MDD of human IL-10 ranged from 0.07-0.30 pg/mL. The mean MDD was 0.13 pg/mL.

bio-techne [®]	752780.1	4/17	FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.	
bio-techne.com info@bio-techne.com	North America	Europe • Middle East • Africa	China info.cn@bio-techne.com	Rest of World bio-techne.com/find-us/distributors
techsupport@bio-techne.com	TEL 800 343 7475	TEL +44 (0)1235 529449	TEL +86 (21) 52380373	TEL +1 612 379 2956

PRECISION

Intra-assay Precision (precision within an assay) - Three samples of known concentration were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (precision between assays) - Three samples of known concentration were tested in twenty-five separate assays to assess inter-assay precision.

	Intra-Assay Precision			Inter-Assay Precision		
Sample	1	2	3	1	2	3
n	20	20	20	25	25	25
Mean (pg/mL)	10.4	66	320	12	73	333
Standard deviation	0.54	3.4	20.7	1.0	5.0	24
CV (%)	5.2	5.2	6.4	10.1	7.4	7.3

RECOVERY

Samples containing and/or spiked with human IL-10 were evaluated for recovery.

Sample Type	Average % Recovery	Range	
Cell culture supernates	105	95-117%	
Serum	109	90-122%	
EDTA plasma	97	86-106%	
Heparin plasma	97	88-102%	

LINEARITY

Samples containing and/or spiked with human IL-10 were serially diluted to evaluate assay linearity.

		Cell culture supernates	Serum	EDTA plasma	Heparin plasma
1:2	Average % of Expected	101	99	100	96
	Range (%)	90-111	92-109	92-110	92-101
1:4	Average % of Expected	95	95	100	92
	Range (%)	83-119	87-102	90-111	88-98
1:8	Average % of Expected	96	91	96	89
	Range (%)	76-114	88-95	84-107	82-96

CORRELATION

This assay has been correlated to the Quantikine® ELISA Kit with a slope of 0.9-1.1 and an R² value greater than 0.9.

SPECIFICITY

Note: Refer to the base kit insert for a complete list of analytes tested for cross-reactivity and interference.

This assay recognizes natural and recombinant human IL-10.

TECHNICAL HINTS

- Protect the microparticles and streptavidin-PE from light at all times.
- Refer to the Base Kit Standard Value Card for reconstitution volume and values of the reconstituted standard.
- Diluted microparticles cannot be stored. Make a fresh dilution of microparticles each time the assay is run.
- The use of a magnetic device made to accommodate a microplate is necessary for washing.
- Discrepancies may exist in values obtained for the same analyte utilizing different technologies.

Magnetic Luminex® Performance Assays afford the user the benefit of multi-analyte analysis of biomarkers in a complex sample. For each sample type, a single, multipurpose diluent is used to optimize recovery, linearity, and reproducibility. Such a multipurpose diluent may not optimize any single analyte to the same degree that a unique diluent selected for analysis of that analyte can optimize conditions. Therefore, some performance characteristics may be more variable than those for assays designed specifically for single analyte analysis.

752780.1 LUHM210 PAGE 2 of 2