

Specifications:

Gene:	mGzmk
Accession:	NP_032222
Insert size:	805bp
Concentration:	10µg at 0.2µg/µL

mGranzyme K cDNA Plasmid

Gzmk granzyme K [*Mus musculus* (house mouse)]

Summary:


Granzyme K is a member of a group of related serine proteases from the cytoplasmic granules of cytotoxic lymphocytes. Cytolytic T lymphocytes (CTL) and natural killer (NK) cells share the remarkable ability to recognize, bind, and lyse specific target cells. They are thought to protect their host by lysing cells bearing on their surface 'nonself' antigens, usually peptides or proteins resulting from infection by intracellular pathogens. Granzyme K lacks consensus sequences for N-glycosylation present in other granzymes.

Description

This shuttle vector contains the complete ORF for the gene of interest, along with a Kozak consensus sequence for optimal translation initiation. It is inserted NotI to AscI. The gene insert is flanked with convenient multiple cloning sites which can be used to easily cut and transfer the gene cassette into your desired expression vector.

Preparation and Storage

Formulation	cDNA is provided in 10 mM Tris-Cl, pH 8.5
Shipping	Ships at ambient temperature
Stability	1 year from date of receipt when stored at -20°C to -80°C
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.



> RDC2029 Plasmid DNA Sequence

```

1  tcgcgcgctt  cggatgatgac  ggtgaaaaac  tctgacacat  gcagctcccc  gagacggtca  cagcttgctc  gtaagcggat  gccgggagca  gacaagcccc
101  tcaggggcgc  tcagcgggtg  ttggcgggtg  tcggggctgg  ctttaactatg  cggcatcaga  gcagattgta  ctgagagtgc  accatatgcg  gtgtgaaata
201  ccgcacacat  gcgtaaggag  aaaataccgc  atcaggcgcc  attcgccatt  caggctcgcc  aactgttggg  aagggcgatc  ggtcggggcc  tcttcgctat
301  tacgccagct  ggcgaaagg  ggatgtgctg  caaggcgatt  aagttgggta  acgccagggt  tttcccagtc  acgacgttgt  aaaacgacgg  ccagtgaatt
401  ggagacgtgt  taacaagctt  ggatccgata  tcgctagcgc  ggccgcacc  atgaggtttt  cttcatgggc  tctggtttcc  ctggtggctg  gcgtttatat
501  gtcttcagag  tgtttccata  ctgaaattat  tggaggggag  gaagtccagc  cgcatccag  gccatttatg  gogtccatcc  agtaccgag  caagcatatt
601  tgtggaggag  tctgtatcca  cccacagtgg  gtgctaacag  ccgcccactg  ctactcttgg  ttcccagag  gccactctcc  caccgtgggt  ttaggagcac
701  attctctttc  caagaatgag  cccatgaagc  agacatttga  aattaaaaag  ttcatcccat  tctcagcact  tcagtcgggt  tccgcatcgc  atgacatcat
801  gctgataaag  cttgcgactg  ctgcagaact  aaacaagaat  gtccaaactg  ttcaactggg  atccaaaaac  tatcttagag  atgggaccaa  atgccagggtg
901  actggctggg  gaaccaccaa  gccagatctg  ttaaccgcct  ctgataccct  gagagaagtc  actgttacca  tcataagtag  aaaacgctgt  aacagccaaa
1001  gctactacaa  ccacaaacct  gttataacca  aggacatgat  atgtgcagga  gatgccagag  gtcaaaagga  ttctgcaag  ggtgactotg  gtggcccttt
1101  gatctgcaaa  ggcattcttc  atgcccctag  ctctcagggc  tataaatgtg  gcctcgccaa  aaagcctgga  atctatacgc  tattaactaa  gaaataccag
1201  acctggatca  aaagcaagct  tgcgccatca  cgtgcacatt  aaaggcgccg  cagtatactc  tagagtcgac  acccggggaa  ttctcgagc  gctcgtctct
1301  agcttggcgt  aatcatggtc  atagctgttt  cctgtgtgaa  attgttatcc  gctcacaaat  ccacacaaca  tacgagccgg  aagcataaag  tgtaaagcct
1401  ggggtgccta  atgagtgagc  taactcacat  taattgcgtt  gcgctcactg  cccgctttcc  agtcgggaaa  cctgtcgtgc  cagctgcatt  aatgaatcgg
1501  ccaacgcgcg  gggagaggcg  gtttgcgtat  ttggcgctct  tccgcttctc  cgctcactga  ctgcgtgcgc  tcggctcgtc  ggtcgcggcg  agcggtatca
1601  gctcactcaa  aggcggtaat  acggttatcc  acagaaacag  gggataacgc  aggaaaagac  atgtgagcaa  aaggccagca  aaagccagg  aaccgtaaaa
1701  aggcgcgctt  gctggcgctt  ttccatagcg  tcgccccccc  tgacgagcat  cacaaaaatc  gacgctcaag  tcagaggtgg  cgaaaccoga  caggactata
1801  aagataccag  gcgtttcccc  ctggaagctc  cctcgtgcgc  tctcctgttc  cgacctgccc  gcttaccgga  tacctgtccg  cctttctccc  ttccggaagc
1901  gtggcgcttt  ctcaaatgct  acgtctagtg  cgtgtgaggt  cgttcgctcc  cgttcgctcc  aagctgggct  gtgtgcagca  acccccgttt  cagcccgacc
2001  gctgcgcctt  atccggtaac  tatcgtcttg  agtccaaccc  ggtaagacac  gacttatcgc  cactggcagc  agccactggt  aacaggatta  gcagagcgag
2101  gtatgtaggc  ggtgctacag  agttcttgaa  gtgtggcct  aactacggct  acactagaag  gacagtattt  ggtatctgcg  cctctgtgaa  gccagtacc
2201  ttccgaaaaa  gagttggtag  ctcttgatcc  ggcaaaaaa  ccacgcgtgg  tagcggtggt  tttttgttt  gcaagcagca  gattacgcgc  agaaaaaag
2301  gatctcaaga  agatcctttg  atcttttcta  cggggtctga  cgctcagtg  aacgaaaact  cacgttaagg  gattttggtc  atgagattat  caaaaaggat
2401  cttcacctag  atccttttaa  attaaaaatg  aagttttaaa  tcaatctaaa  gtatatatga  gtaaaacttg  tctgacagtt  accaatgctt  aatcagttag
2501  gcacctatct  cagcgatctg  tctatttctg  tcatccatag  ttgcctgact  ccccgctcgt  tagataacta  cgatacggga  gggcttacca  tctggcccca
2601  gtgctgcaat  gataccgcga  gacccacgct  caccggctcc  agatttatca  gcaataaacc  agccagcccg  aagggccgag  cgcagaagtg  gtcctgcaac
2701  tttatccgcc  tccatccagt  ctattaatg  ttgcgggaa  gctagagtaa  gtggttcgcc  agttaatagt  ttgcccacag  ttgttgccat  tgctacaggc
2801  atcgtggtgt  cagcgtcgtc  gtttggtag  gcttcattca  gctccggttc  ccaacgatca  aggcgagtta  catgatcccc  catgttgtgc  aaaaaagcgg
2901  ttagctcctt  cggtcctccg  atcgttgtca  gaagtaagtt  ggccgcagtg  ttatcactca  tgggttatgg  agcactgcat  aattctctta  ctgtcatgcc
3001  atccgtaaga  tgcctttctg  tgactggtga  gtactcaacc  aagtcattct  gagaatagtg  tatgcggcga  tatgcggcga  cagagttgct  cttgcccgcg
3101  gataataccg  cgccacatag  cagaacttta  aaagtgtcta  tcattggaaa  acgttctctg  gggcgaaaac  tctcaaggat  cttaccgctg  ttgagatcca
3201  gttcgatgta  acccactcgt  gcacccaact  gatcttcagc  atcttttact  ttcaccagcg  tttctgggtg  agcaaaaaa  ggaaggcaaa  atgccgcaaa
3301  aaagggaata  agggcgacac  ggaatgttg  aatactcata  ctcttccttt  ttcaatatta  ttgaagcatt  tatcagggtt  attgtctcat  gagcgatac
3401  atatttgaat  gtatttagaa  aaataaacia  atagggttcc  gcgcacatt  tccccgaaa  gtgccacctg  acgtctaaga  aaccattatt  atcatgacat
3501  taacctataa  aaataggcgt  atcacgagcg  cctttcgtc

```

> RDC2029 Translated Insert Sequence

```

1  mrfsswalvs  lvagvymsse  cfhteiiigr  evqphsrpfm  asiqrskhi  cggvlihpqw  vltaahcysw  fprghsptvv  lgahslskne  pmkqtfeikk
101  fipfslrlsq  sashdimlik  lrtaaelnkn  vqlhlgskn  ylrdgtkcv  tgwttkpdl  ltasdtlrev  tvtiisrkrc  nsqsyynhkp  vitkdmicag
201  dargqkdsck  gdsqplick  gifhalvsqg  ykcgiaakpg  iytlltkkyq  twiksklaps  rah

```