

Specifications:

Gene:	rVtcn1
Accession:	NP_001019415
Insert size:	862bp
Concentration:	10µg at 0.2µg/µL

Description

This shuttle vector contains the complete ORF for the gene of interest, along with a Kozak consensus sequence for optimal translation initiation. It is inserted NotI to AscI. The gene insert is flanked with convenient multiple cloning sites which can be used to easily cut and transfer the gene cassette into your desired expression vector.

Preparation and Storage

Formulation	cDNA is provided in 10 mM Tris-Cl, pH 8.5
Shipping	Ships at ambient temperature
Stability	1 year from date of receipt when stored at -20°C to -80°C
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.


rB7-H4 cDNA Plasmid

Vtcn1 V-set domain containing T cell activation inhibitor 1
[*Rattus norvegicus* (Norway rat)]

Also known as: RGD1311204

Summary:

B7-H4 belongs to the B7 costimulatory protein family. It negatively regulates T-cell-mediated immune response by inhibiting T-cell activation, proliferation, cytokine production and development of cytotoxicity. High levels of B7-H4 expression have been correlated with tumor progression. It may serve as a predictive marker for renal cell carcinoma.



> RDC2380 Plasmid DNA Sequence

```

1   tcgcgcgcttt  cggatgatgac  ggtgaaaaacc  tctgacacat  gcagctcccc  gagacggtca  cagcttgtct  gtaagcggat  gccgggagca  gacaagcccg
101  tcagggcgcg  tcagcgggtg  ttggcgggtg  tcggggctgg  ctttaactatg  cggcatcaga  gcagattgta  ctgagagtgc  accatatgcg  gtgtgaaata
201  ccgcacacgat  gcgtaaggag  aaaataccgc  atcaggcgcc  attcgccatt  caggctgcgc  aactgttggg  aagggcgatc  ggtgcgggcc  tcttcgctat
301  tacgccagct  ggcgaaagg  ggatgtgtg  caaggcgatt  aagttgggta  acgccagggt  tttcccagtc  acgacgttgt  aaaacgacgg  ccagtgaatt
401  ggagacgtgt  taacaagctt  ggatccgata  tcgctagcgc  ggccgcaccc  atggcttcc  tggggcagat  catcttttgg  agtattatta  acgtcatcat
501  catctctgct  ggagccatcg  tactcatcat  tggcttttgg  atttcaggca  agcaactcat  caccgtcaca  accttccact  cagctgaaa  tattggagag
601  gacgggaccc  tgagctgcac  gtttgaacct  gacatcaaac  tcaacggcat  cgtcatccag  tggctgaagg  aagcatcaa  aggtttggtc  catgagttca
701  aagagggcaa  agacgacctc  tcacagcagc  acgagatggt  tcgaggcgcg  acagcaggtg  ttgcagatca  agtgggtggt  ggcaatgctt  cctcgagact
801  gaaaaatgtg  cagctcactg  atgctggcac  ctacacatgt  tataccaca  cctcaaaagg  caaagggaat  gctaaccttg  agtataagac  tggagccttc
901  agtatgccag  agattaatgt  ggactataat  gccagttcgg  agagtctacg  ctgtgaggct  cctcgggtgt  tccccagcc  cacagtggcc  tgggcatctc
1001 aagttgacca  aggagccaac  ttctcagaag  tctccaacac  tagcttcgag  ttgaactctg  agaattgtac  catgaaggtc  gtatctgtgc  tctacaatgt
1101 cacaatcaac  aacacatact  cctgtatgat  cgaaaaagat  attgccaaag  ccacagggga  catcaaaagt  acagattcag  aggtcaaaag  gcgtagtcag
1201 ctggagctgc  tgaactctgg  gccttctcgg  tgtgtttctt  ctgtctctgc  cgtctctagc  gcaactctat  ctctctctg  ttgctctgat  ctaagataaa
1301 ggccgcgcag  tatactctag  agtcgacacc  cggggaattc  ctcgagcgct  cgtctctagc  ttggcgtaat  catggtcata  gctgtttcct  gtgtgaaatt
1401 gttatccgct  cacaattcca  cacaacatac  gagccggaag  cataaagtgt  aaagcctggg  gtgcctaatt  agtgagctaa  ctacatttaa  ttgcgttgcg
1501 ctcactgccc  gctttccagt  cgggaaacct  gtcgtccagc  ctgcattaat  gaatcgccca  acgcgcgggg  agaggcggtt  tgcgtattgg  gcgctcttcc
1601 gcttccctgc  tcaactgactc  gtcgctcgct  gtcgttcggc  tgcggtcagc  ggtatcagct  cactcaaaag  cggtaatcag  gttatccaca  gaatcgagg
1701 ataacgcagg  aaagaacatg  tgagcaaaag  gccagcaaaa  ggccagggaac  cgtaaaaagg  ccgctgttgt  ggcgttttcc  cataggctcc  gcccccctga
1801 cagagcatcac  aaaaactcgac  gctcaagtca  gaggtggcga  aacccgacag  gactataaaa  ataccaggcg  ataccaggcg  ttccccctg  gaagctccct  ctgcgctct
1901 ctgtttccga  cccctgcgct  tacggatata  ctgtccgctc  ttctcccttc  agtgaagcgt  gcgctttctc  aatgtctcac  gctgagtgat  ctgacttcgg
2001 tgtaggtcgt  tcgctccaag  ctgggtctgt  tgcacgaacc  ccccgttcag  cccgaccgct  gcgccttato  cggttaactat  cgtcttgagt  ccaaccoggt
2101 aagacacgac  ttatcgccac  tggcagcagc  cactggtaac  aggattagca  gagcgaggta  tgtaggcggt  gctacagagt  tcttgaagt  gtggcctaac
2201 tacggctaca  ctagaaggac  agtatttgg  tgcgaagcc  tgctgaagcc  agttaccttc  ggaaaaagag  ttggtagctc  ttgtccggc  aaacaaacca
2301 ccgctggtag  cgggtgtttt  tttgtttgca  agcagcagat  tacgcgcaga  aaaaaaggat  ctcaagaaga  tcctttgatc  ttttctacgg  ggtctgacgc
2401 tcagtggaa  gaaaactcac  gttacgggat  ttgtgctatg  agattatcaa  aaaggatctt  cactagatc  cttttaaatt  aaaaatgaag  ttttaataca
2501 atctaaagta  tatatgagta  aacttgggtc  gacagttaac  aatgcttaat  cagttaggca  cctatctcag  cgatctgtct  atttcgttca  tccatagttg
2601 cctgactccc  cgtcgtgtag  ataactacga  tacgggaggg  cttaccatct  ggccccagtg  ctgcaatgat  accgcgagac  ccacgtctac  cggctccaga
2701 tttatcagca  ataaaccagc  cagccggaag  ggccgagcgc  agaagtggtc  ctgcaacttt  atccgcctcc  atccagctca  ttaattgttg  ccgggaagct
2801 agagtaagta  gttcgccagt  taatagtttg  cgcaacggtg  ttgccattgc  tacaggcatc  gtggtgtcac  gctcgctggt  tggtagggt  tcattcagct
2901 ccggttccca  acgatcaagg  cgagttacat  gatccccat  gttgtgcaaa  aaagcggtta  gctccttcgg  tcctccgac  gttgtcagaa  gtaagtggc
3001 cgcagtgtaa  tcaactcatg  ttatggcagc  actgcataat  tctcttactg  tcatgcatc  cgtaagatgc  ttttctgtga  ctggtagta  ctcaaccaag
3101 tcattctgag  aatagtgtat  gcggcgaccg  agttgctctt  gcccgcgctc  aataccggat  aataccgcgc  cacatagcag  aactttaaaa  gtgctcatca
3201 ttgaaaaaac  ttcttcgggg  cgaaaactct  caaggatctt  accgctgttg  agatccagtt  cgatgtaacc  cactcgtgca  cccaactgat  ctccagcatc
3301 ttttactttc  accagcgttt  ctgggtgagc  aaaaacagga  aggcaaaatg  cgcacaaaaa  gggaataaag  gcgacacgga  aatgttgaat  actcatactc
3401 ttcctttttc  aatattattg  aagcatttat  cagggttatt  gtctcatgag  cggatacata  tttgaatgta  tttagaaaaa  taacaaaata  ggggttccgc
3501 gcacatttcc  ccgaaaagtg  ccacctgacg  tctaagaaac  cattattatc  atgacattaa  cctataaaaa  taggcgtatc  acgaggccct  ttctgc

```

> RDC2380 Translated Insert Sequence

```

1   maslgqiifw  siinviila  gaivliigfg  isgkhfitvt  tftsagnige  dgtlscftfp  diklngiviq  wlkegikglv  hefkegkddl  sqghemfrgr
101  tavfadqvvv  gnasrlrknv  qltdagtytc  yihtskgkgn  anleyktgaf  smpeinvdyn  asseslrcea  prwfpqptva  wasqvdqgan  fsevsntsfe
201  lnsenvtmkv  vsvlynvtn  ntyscmiend  iakatgdikv  tdsevkrrsq  lelinsgpp  cvssvsaagw  allslscclm  lr

```