

MATERIAL DATA SHEET**Di-Ubiquitin/Ub₂, 48-DCA-linked, Agarose, human recombinant**
Cat. # UCD-202

Dichloroacetone (DCA) forms of di-ubiquitin (Ub₂) are chemically linked via a DCA linker using cysteine introduced into the ubiquitin sequence at the appropriate residues to mimic a native isopeptide linkage. Reports have demonstrated that the DCA linkage is stable and cannot be processed by deubiquitinases. Additionally, the orientation of the two ubiquitins in DCA-linked Ub₂ is preserved in relation to di-ubiquitin containing a native isopeptide linkage. This makes DCA-linked Ub₂ useful for characterizing polyubiquitin-binding proteins in the absence of added deubiquitinase inhibitors.

Product Information

Quantity:	100 µl resin supplied as a 50% slurry in pH 5.5 20 mM sodium acetate, 0.01% sodium azide.
Stock:	Non-hydrolyzable Ub ₂ (K48) chains coupled to beads at X nmol/ml (X mg/ml resin).

Use and Storage

Use:	Equilibrate resin by washing with 10 volumes of desired start buffer. Binding and elution of material is dependent on individual experimental conditions and requirements.
Storage:	Polyubiquitin-agarose can be re-used if properly maintained. After use, clean resin with a wash cycle of 5 volumes 100 mM HEPES pH 8.0, 500 mM NaCl followed by 5 volumes 100 mM NaOAc pH 5.5, 500 mM NaCl. Repeat twice, then rinse resin with a low salt buffer. Store resin at 4°C in aqueous buffer (acetate, etc.) containing 0.01% sodium azide. DO NOT FREEZE.

Literature

References:	Carlile, C., <i>et al.</i> (2009) <i>J Biol Chem</i> , 284 : 29326-29334 Russell N. and Wilkinson K. (2004) <i>Biochemistry</i> 43 : 4844-4854 Yin L., <i>et al.</i> (2000) <i>Biochemistry</i> 39 : 10001-10010
--------------------	---

For Laboratory Research Use Only, Not For Use in Human

Rev: 08/02/2012

840 Memorial Drive, Cambridge, MA 02139 Phone: 617-241-7072 FAX: 617-492-3565
www.bostonbiochem.com

The contents of this datasheet (unless otherwise noted) are Copyright © 2012 Boston Biochem, Inc. All rights reserved. Duplication in whole or in part is strictly prohibited without the express written consent of Boston Biochem, Inc. "Boston Biochem" is a Trademark of Boston Biochem, Inc., 840 Memorial Drive, Cambridge, MA 02139.