

## DESCRIPTION

<b>Source</b>	Mouse myeloma cell line, NS0-derived mouse PDGF R beta protein		
	Mouse PDGF R $\beta$ (Leu32-Lys530) Accession # P05622	IEGRMD	Human IgG <sub>1</sub> (Pro100-Lys330)
	N-terminus		C-terminus
<b>N-terminal Sequence</b>	Leu32		
<b>Analysis</b>			
<b>Structure / Form</b>	Disulfide-linked homodimer		
<b>Predicted Molecular Mass</b>	82.6 kDa (monomer)		

## SPECIFICATIONS

<b>SDS-PAGE</b>	124 kDa, reducing conditions
<b>Activity</b>	Measured by its ability to inhibit the biological activity of PDGF-BB using NR6R-3T3 mouse fibroblast cells. Raines, E.W. <i>et al.</i> (1985) <i>Methods Enzymol.</i> <b>109</b> :749. The ED <sub>50</sub> for this effect is 6-48 ng/mL in the presence of 4 ng/mL recombinant human PDGF-BB.
<b>Endotoxin Level</b>	<0.10 EU per 1 $\mu$ g of the protein by the LAL method.
<b>Purity</b>	>90%, by SDS-PAGE under reducing conditions and visualized by silver stain.
<b>Formulation</b>	Lyophilized from a 0.2 $\mu$ m filtered solution in PBS. See Certificate of Analysis for details.

## PREPARATION AND STORAGE

<b>Reconstitution</b>	Reconstitute at 100 $\mu$ g/mL in sterile PBS.
<b>Shipping</b>	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.
<b>Stability &amp; Storage</b>	<p><b>Use a manual defrost freezer and avoid repeated freeze-thaw cycles.</b></p> <ul style="list-style-type: none"> <li>• 12 months from date of receipt, -20 to -70 °C as supplied.</li> <li>• 1 month, 2 to 8 °C under sterile conditions after reconstitution.</li> <li>• 3 months, -20 to -70 °C under sterile conditions after reconstitution.</li> </ul>

## BACKGROUND

The platelet-derived growth factor (PDGF) family consists of proteins derived from four genes (PDGF-A, -B, -C, and -D) that form disulfide-linked homodimers (PDGF-AA, -BB, -CC, and -DD) and a heterodimer (PDGF-AB) (1, 2). These proteins regulate diverse cellular functions by binding to and inducing the homo- or hetero-dimerization of two receptors (PDGF R $\alpha$  and R $\beta$ ). Whereas  $\alpha/\alpha$  homo-dimerization is induced by PDGF-AA, -BB, -CC, and -AB,  $\alpha/\beta$  hetero-dimerization is induced by PDGF-AB, -BB, -CC, and -DD, and  $\beta/\beta$  homo-dimerization is induced only by PDGF-BB, and -DD (1 - 4). Both PDGF R $\alpha$  and R $\beta$  are members of the class III subfamily of receptor tyrosine kinases (RTK) that also includes the receptors for M-CSF, SCF and Flt3-ligand. All class III RTKs are characterized by the presence of five immunoglobulin-like domains in their extracellular region and a split kinase domain in their intracellular region. Ligand-induced receptor dimerization results in autophosphorylation in trans resulting in the activation of several intracellular signaling pathways that can lead to cell proliferation, cell survival, cytoskeletal rearrangement, and cell migration. Many cell types, including fibroblasts and smooth muscle cells, express both the  $\alpha$  and  $\beta$  receptors. Others have only the  $\alpha$  receptors (oligodendrocyte progenitor cells, mesothelial cells, liver sinusoidal endothelial cells, astrocytes, platelets and megakaryocytes) or only the  $\beta$  receptors (myoblasts, capillary endothelial cells, pericytes, T cells, myeloid hematopoietic cells and macrophages). A soluble PDGF R $\alpha$  has been detected in normal human plasma and serum as well as in the conditioned medium of the human osteosarcoma cell line MG-63 (5). Both the recombinant mouse and human soluble PDGF R $\alpha$  bind PDGF with high affinity and are potent PDGF antagonists.

### References:

1. Betsholtz, C. *et al.* (2001) *BioEssays* **23**:494.
2. Ostman, A. and A.H. Heldin (2001) *Advances in Cancer Research* **80**:1.
3. Gilbertson, D. *et al.* (2001) *J. Biol. Chem.* **276**:27406.
4. LaRochells, W.J. *et al.* (2001) *Nature Cell Biol.* **3**:517.
5. Tiesman, J. and C.E. Hart (1993) *J. Biol. Chem.* **5**:9621.