

**DESCRIPTION**

<b>Source</b>	Human embryonic kidney cell, HEK293-derived human ERMAP protein		
	Human ERMAP (His30-Ala155) Accession # Q96PL5.1	IEGRMD	Human IgG <sub>1</sub> (Pro100-Lys330)
	N-terminus		C-terminus
<b>N-terminal Sequence Analysis</b>	His30 & Gly32		
<b>Structure / Form</b>	Disulfide-linked homodimer		
<b>Predicted Molecular Mass</b>	40 kDa		



**SPECIFICATIONS**

<b>SDS-PAGE</b>	40-55 kDa, under reducing conditions
<b>Activity</b>	Measured by its ability to inhibit anti-CD3 antibody induced IL-2 or IFN-gamma secretion by human T cells. ED <sub>50</sub> for this effect is 0.75-7.5 µg/mL.
<b>Endotoxin Level</b>	<0.10 EU per 1 µg of the protein by the LAL method.
<b>Purity</b>	>95%, by SDS-PAGE visualized with Silver Staining and quantitative densitometry by Coomassie® Blue Staining.
<b>Formulation</b>	Lyophilized from a 0.2 µm filtered solution in PBS. See Certificate of Analysis for details.

**PREPARATION AND STORAGE**

<b>Reconstitution</b>	Reconstitute at 200 µg/mL in PBS.
<b>Shipping</b>	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.
<b>Stability &amp; Storage</b>	<p><b>Use a manual defrost freezer and avoid repeated freeze-thaw cycles.</b></p> <ul style="list-style-type: none"> <li>• 12 months from date of receipt, -20 to -70 °C as supplied.</li> <li>• 1 month, 2 to 8 °C under sterile conditions after reconstitution.</li> <li>• 3 months, -20 to -70 °C under sterile conditions after reconstitution.</li> </ul>

**DATA**



**BACKGROUND**

Erythroid membrane-associated protein (ERMAP), also called Scianna blood group antigen, is a member of the butrophilin (BTN) and butrophilin-like (BTNL) family within the immunoglobulin superfamily. ERMAP consists of an extracellular domain (ECD) with one V-type Ig-like (IgV) fold, a transmembrane segment, and a cytoplasmic domain with a B30.2 motif (1). Within the ECD, human ERMAP shares 33% and 31% amino acid sequence identity with mouse and rat ERMAP, respectively. Interestingly, human ERMAP contains a single IgV-like domain in the ECD, while mouse and rat ERMAP contain both an IgV-like and IgC-like domain. However, the absence of the IgC-like domain in human ERMAP does not seem to affect its functions (2, 3). ERMAP is an adhesion receptor molecule highly expressed in erythroid tissues, on the cell surface of resting and activated antigen-presenting cells (APCs) and in some tumor tissues (1, 4, 5). Fusion proteins (ERMAP-Ig) of both human and mouse ERMAP inhibit T cell functions in vitro and administration of the fusion protein ameliorates autoimmune diseases, including experimental autoimmune encephalomyelitis and type 1 diabetes, in mice (5). ERMAP acts as a novel inhibitory molecule for T cells and macrophages (5).

**References:**

1. Xu, H. *et al.* (2001) *Genomics* **76**:2.
2. Compte, E.W. *et al.* (2004) *Eur. J. Immunol.* **34**:2089.
3. Su, Y.Y. *et al.* (2001) *Blood Cell Mol. Dis.* **27**:938
4. Abeler-Dorner, L. *et al.* (2012) *Trends Immunol.* **33**:34.
5. Su, M. *et al.* (2020) *Cell. Mol. Immunol.* doi: 10.1038/s41423-020-0494-8.