

DESCRIPTION

Source Mouse myeloma cell line, NS0-derived
His20-Val223, with a C-terminal 10-His tag
Accession # P02743

N-terminal Sequence Analysis His20

Structure / Form Noncovalently-linked homopentamer

Predicted Molecular Mass 25 kDa (monomer)

SPECIFICATIONS

SDS-PAGE 28-30 kDa, reducing conditions

Activity Measured by its binding ability in a functional ELISA.
When Recombinant Human Fcγ RIIA/CD32a (Catalog # 1330-CD) is immobilized at 2 µg/mL (100 µL/well), the concentration of Recombinant Human Pentraxin 2/SAP that produces 50% of the optimal binding response is 0.175-1 µg/mL.

Endotoxin Level <0.01 EU per 1 µg of the protein by the LAL method.

Purity >97%, by SDS-PAGE visualized with Silver Staining and quantitative densitometry by Coomassie® Blue Staining.

Formulation Lyophilized from a 0.2 µm filtered solution in PBS. See Certificate of Analysis for details.

PREPARATION AND STORAGE

Reconstitution Reconstitute at 100 µg/mL in sterile PBS.

Shipping The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Stability & Storage **Use a manual defrost freezer and avoid repeated freeze-thaw cycles.**

- 12 months from date of receipt, -20 to -70 °C as supplied.
- 1 month, 2 to 8 °C under sterile conditions after reconstitution.
- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

Pentraxin 2 (PTX2), also known as Serum Amyloid P Component (SAP), is a secreted serum glycoprotein that is a universal non-fibrillar component of amyloid deposits. These extracellular deposits of insoluble protein fibrils are the result of protein misfolding and can lead to tissue damage and disease (1, 2). PTX2 belongs to the pentaxin superfamily, whose members have the characteristic pentagonal discoid arrangement of five non-covalently bound subunits. Pentaxins bind to a variety of molecules in a calcium-dependent lectin-like manner through a pattern-recognition-binding site (1, 4, 5). Two subfamilies of pentaxins, the classical or short pentaxin subfamily that includes the serum C-reactive protein (CRP) and PTX2, and the fusion or long pentaxin subfamily whose members contain pentaxin-related carboxyl-terminal halves, are known (1).

PTX2 and CRP share approximately 50% amino acid sequence identity (2, 5). They are produced and secreted by liver hepatocytes and circulates in plasma. Mouse PTX2 is a major acute-phase protein whose plasma concentrations increase dramatically during an acute phase response (2). In human where CRP is the major acute-phase protein, the plasma concentration of human PTX2 remains relatively constant in response to tissue-damage (2, 5). The gene for PTX2 has been localized to human chromosome 1 where it is closely linked to the gene for CRP.

PTX2 associates ubiquitously with all amyloid deposits that are implicated in a diverse range of diseases including Alzheimer's and prion diseases, type 2 diabetes and various systemic amyloidoses (3, 6, 7). As a non-fibrillar component, PTX2 regulates the solubility of amyloid fibrils and protects them from degradation by proteolytic enzymes and phagocytic cells. In addition to its role in the pathogenesis of amyloidoses, PTX2 also has an important physiological function in innate immunity (8). It is an opsonin that interacts with all three types of human Fcγ receptors that mediate phagocytosis by human polymorphonuclear leukocytes. It has been proposed that PTX2 may function as an opsonin for a variety of ligands including autoantigens, apoptotic cells, chromatin, DNA, and micro-organisms.

References:

1. Goodman, A. *et al.* (1996) Cytokine Growth Factor Rev. **7**:191.
2. Steel, D. and A. Whitehead (1994) Immunol. Today **15**:81.
3. Hirschfield, G.M. and P.N. Hawkins (2003) Int. J. Biochem. Cell Biol. **35**:1608.
4. Emsley, J. *et al.* (1994) Nature **367**:338.
5. Mantzouranis, E. *et al.* (1985) J. Biol. Chem. **260**:7752.
6. Botto, M. *et al.* (1997) Nature Medicine **3**:855.
7. Pepys, M. *et al.* (2002) Nature **417**:254.
8. Bharadwaj, D. *et al.* (2001) J. Immunol. **166**:6735.