

**DESCRIPTION**

<b>Source</b>	<i>E. coli</i> -derived human IL-2 protein Ala21-Thr153, with an N-terminal Met Accession # P60568
<b>N-terminal Sequence Analysis</b>	Met
<b>Predicted Molecular Mass</b>	15 kDa




**SPECIFICATIONS**

<b>SDS-PAGE</b>	13 kDa, under reducing conditions.
<b>Activity</b>	Measured in a cell proliferation assay using CTLL-2 mouse cytotoxic T cells. Gearing, A.J.H. and C.B. Bird (1987) in <i>Lymphokines and Interferons, A Practical Approach</i> . Clemens, M.J. <i>et al.</i> (eds): IRL Press. 295. The ED <sub>50</sub> for this effect is 0.05-0.25 ng/mL.
<b>Endotoxin Level</b>	<0.10 EU per 1 µg of the protein by the LAL method.
<b>Purity</b>	>97%, by SDS-PAGE under reducing conditions and visualized by silver stain.
<b>Formulation</b>	Lyophilized from a 0.2 µm filtered solution in Acetonitrile and TFA. See Certificate of Analysis for details.

**PREPARATION AND STORAGE**

<b>Reconstitution</b>	Reconstitute at 100 µg/mL in sterile 100 mM Acetic Acid. Alternatively, reconstitute at 100 µg/mL in sterile deionized water and use within 24 hours, store at 2 to 8 °C.
<b>Shipping</b>	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
<b>Stability &amp; Storage</b>	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> <li>• 12 months from date of receipt, -20 to -70 °C as supplied.</li> <li>• 1 month, 2 to 8 °C under sterile conditions after reconstitution.</li> <li>• 3 months, -20 to -70 °C under sterile conditions after reconstitution.</li> </ul>

**DATA**



**BACKGROUND**

Recombinant Interleukin-2 (IL-2) is expressed in *E. coli* and has been engineered to contain the serine for cysteine substitution found in Proleukin® (aldesleukin). Recombinant IL-2 is widely used in cell culture for the expansion of T cells. IL-2 is expressed by CD4+ and CD8+ T cells,  $\gamma\delta$  T cells, B cells, dendritic cells, and eosinophils (1-3). Mature human IL-2 shares 56% and 66% amino acid (aa) sequence identity with mouse and rat IL-2, respectively. Human and mouse IL-2 exhibit cross-species activity (4). The receptor for IL-2 consists of three subunits that are present on the cell surface in varying preformed complexes (5-7). The 55 kDa IL-2 R $\alpha$  is specific for IL-2 and binds with low affinity. The 75 kDa IL-2 R $\beta$ , which is also a component of the IL-15 receptor, binds IL-2 with intermediate affinity. The 64 kDa common gamma chain  $\gamma$ c/IL-2 R $\gamma$ , which is shared with the receptors for IL-4, -7, -9, -15, and -21, does not independently interact with IL-2. Upon ligand binding, signal transduction is performed by both IL-2 R $\beta$  and  $\gamma$ c.

IL-2 is best known for its autocrine and paracrine activity on T cells. It drives resting T cells to proliferate and induces IL-2 and IL-2 R $\alpha$  synthesis (1, 2). It contributes to T cell homeostasis by promoting the Fas-induced death of naïve CD4+ T cells but not activated CD4+ memory lymphocytes (8). IL-2 plays a central role in the expansion and maintenance of regulatory T cells, although it inhibits the development of Th17 polarized cells (9-11). Thus, IL-2 may be a key cytokine in the natural suppression of autoimmunity (12, 13).

IL-2 expression and concentration can have either immunostimulatory effects at high doses or immunosuppressive effects at low doses due to its preferential binding to different receptor subunits expressed by various immune cell types. This has led to the generation of recombinant IL-2 variants aimed at modifying IL-2 receptor binding for increased antitumor efficacy (14, 15). These variants are typically used in combination with immune checkpoint inhibitors instead of as a monotherapy (14). IL-2 can be genetically engineered to express in NK cells for CAR T cell therapies, and in combination with other cytokines like IL-15, can increase cell viability and proliferation (16). In addition to adoptive cell transfer and checkpoint blockade inhibitors, cancer vaccines that boost immune responses have been combined with IL-2 treatment with promising results in recent studies (15).

In cell culture, IL-2 is a frequently used cytokine for the proliferation, differentiation, and increased antibody secretion of B cells as they transform into plasma cells *in vitro* (17). IL-2 is also a classically used cytokine for the expansion of NK cells, early differentiated T cells and effector memory Treg cells for adoptive cell transfer cancer immunotherapy (16, 18). GMP IL-2 is a commonly used supplement for the expansion of these cell types for cellular therapies.

**References:**

1. Ma, A. *et al.* (2006) *Annu. Rev. Immunol.* **24**:657.
2. Gaffen, S.L. and K.D. Liu (2004) *Cytokine* **28**:109.
3. Taniguchi, T. *et al.* (1983) *Nature* **302**:305.
4. Mosmann, T.R. *et al.* (1987) *J. Immunol.* **138**:1813.
5. Liparoto, S.F. *et al.* (2002) *Biochemistry* **41**:2543.
6. Wang, X. *et al.* (2005) *Science* **310**:1159.
7. Bodnar, A. *et al.* (2008) *Immunol. Lett.* **116**:117.
8. Jaleco, S. *et al.* (2003) *J. Immunol.* **171**:61.
9. Malek, T.R. (2003) *J. Leukoc. Biol.* **74**:961.
10. Laurence, A. *et al.* (2007) *Immunity* **26**:371.
11. Kryczek, I. *et al.* (2007) *J. Immunol.* **178**:6730.
12. Afzali, B. *et al.* (2007) *Clin. Exp. Immunol.* **148**:32.
13. Fehervari, Z. *et al.* (2006) *Trends Immunol.* **27**:109.
14. Xue, D. *et al.* (2021) *Antibody Therapeutics.* **4**(2): 123-133.
15. Wolfarth, A.A. *et al.* (2022) *Immune Netw.* **22**(1): e5.
16. Koehl, U. *et al.* (2015) *Oncoimmunology.* **5**(4).
17. Marsman, C. *et al.* (2022) *Front. In Immunol.* **13**(815449).
18. Chamucero-Millares, J.A. *et al.* (2021) *Cellular Immunology.* **360**(104257).