

DESCRIPTION

Source	<i>Spodoptera frugiperda</i> , Sf 21 (stably transfected)-derived Gln19-Ile144 Accession # P15248
N-terminal Sequence Analysis	No results obtained. Gln19 inferred from enzymatic pyroglutamate treatment revealing Gly20
Predicted Molecular Mass	14 kDa


SPECIFICATIONS

SDS-PAGE	12-24 kDa, reducing conditions
Activity	Measured in a cell proliferation assay using MO7e human megakaryocytic leukemic cells. Avanzi, G. <i>et al.</i> (1988) Br. J. Haematol. 69 :359. The ED ₅₀ for this effect is 0.03-0.3 ng/mL.
Endotoxin Level	<0.10 EU per 1 µg of the protein by the LAL method.
Purity	>95%, by SDS-PAGE visualized with Silver Staining and quantitative densitometry by Coomassie® Blue Staining.
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with BSA as a carrier protein. See Certificate of Analysis for details.

PREPARATION AND STORAGE

Reconstitution	Reconstitute at 100 µg/mL in PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -70 °C as supplied. ● 1 month, 2 to 8 °C under sterile conditions after reconstitution. ● 3 months, -20 to -70 °C under sterile conditions after reconstitution.

DATA



BACKGROUND

Interleukin-9 (IL-9), also known as P40 and MEA (mast cell growth-enhancing activity), is a 30-40 kDa glycosylated member of a cytokine family that includes Interleukins-2, -4, -7, -15, and -21. These proteins utilize heteromeric receptors containing the Common gamma chain (γ_c) in addition to ligand-specific subunits. IL-9 interacts selectively with IL-9 R which then associates with γ_c to form the functional receptor complex. IL-9 contributes to allergic inflammation, autoimmunity-induced inflammation, parasite clearance from the GI tract, and Treg-mediated immune suppression (1, 2). It enhances the expansion and recruitment of mast cells and eosinophils as well as the production of IgE and Th2 cytokines (3-6). It is required for anaphylactic responses to ingested allergens but not to systemic allergens (7). IL-9 plays multiple roles in the development and function of subsets within the CD4⁺ T cell lineage (8). It is expressed by activated Th9, Th17, Treg, and Th2 cells (3, 9-12). IL-9 acts as an autocrine growth and activation factor for Th17, Treg, and mast cells (3, 11, 13). It also can inhibit immune responses by enhancing the suppressive properties of Treg and by recruiting immune-suppressive mast cells to sites of inflammation (11, 12). Mature human IL-9 shares 57% amino acid sequence identity with mouse and rat IL-9 (14, 15).

References:

1. Noelle, R.J. and E.C. Nowak (2010) *Nat. Rev. Immunol.* **10**:683.
2. Goswami, R. and M.H. Kaplan (2011) *J. Immunol.* **186**:3283.
3. Nowak, E.C. *et al.* (2009) *J. Exp. Med.* **206**:1653.
4. Townsend, M.J. *et al.* (2000) *Immunity* **13**:573.
5. Leech, M.D. and R.K. Grencis (2006) *J. Immunol.* **176**:2505.
6. Fawaz, L.M. *et al.* (2007) *J. Allergy Clin. Immunol.* **120**:1208.
7. Osterfeld, H. *et al.* (2010) *J. Allergy Clin. Immunol.* **125**:469.
8. Jabeen, R. and M.H. Kaplan (2012) *Curr. Opin. Immunol.* **24**:303.
9. Tan, C. *et al.* (2010) *J. Immunol.* **185**:6795.
10. Jager, A. *et al.* (2009) *J. Immunol.* **183**:7169.
11. Elyaman, W. *et al.* (2009) *Proc. Natl. Acad. Sci. USA* **106**:12885.
12. Lu, L.-F. *et al.* (2006) *Nature* **442**:997.
13. Stassen, M. *et al.* (2000) *J. Immunol.* **164**:5549.
14. Yang, Y.-C. *et al.* (1989) *Blood* **74**:1880.
15. Renaud, J.-C. *et al.* (1990) *Cytokine* **2**:9.