

DESCRIPTION

Source	<i>Spodoptera frugiperda</i> , Sf 21 (baculovirus)-derived human HGF protein Gln32-Ser728 (Asp384Asn, Asp416Asn, Leu622Ser, His645Arg, Val678Ile) Accession # P14210.2
N-terminal Sequence Analysis	α chain: No results obtained; Gln32 predicted β chain: Val495
Structure / Form	Disulfide-linked heterodimer
Predicted Molecular Mass	53.7 kDa (α chain), 26 kDa (β chain)


SPECIFICATIONS

SDS-PAGE	60-70 kDa, non-reducing conditions 60 kDa, and 33-34 kDa, reducing conditions
Activity	Measured by its ability to induce IL-11 secretion by Saos-2 human osteosarcoma cells. Hjertner, O. <i>et al.</i> (1999) <i>Blood</i> 94 :3883. The ED ₅₀ for this effect is 0.1-0.5 ng/mL. The specific activity of Recombinant Human HGF is approximately 1.5 x 10 ³ IU/μg, which is calibrated against Recombinant Human HGF WHO International Standard (NIBSC code: 96/564).
Endotoxin Level	<1.0 EU per 1 μg of the protein by the LAL method.
Purity	>95%, by SDS-PAGE visualized with Silver Staining and quantitative densitometry by Coomassie® Blue Staining.
Formulation	Lyophilized from a 0.2 μm filtered solution in NaH ₂ PO ₄ and NaCl with BSA as a carrier protein. See Certificate of Analysis for details.

PREPARATION AND STORAGE

Reconstitution	Reconstitute at 100 μg/mL in sterile PBS containing at least 0.1% human or bovine serum albumin.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> • 12 months from date of receipt, -20 to -70 °C as supplied. • 1 month, 2 to 8 °C under sterile conditions after reconstitution. • 3 months, -20 to -70 °C under sterile conditions after reconstitution.

DATA



BACKGROUND

HGF, also known as scatter factor and hepatopoietin A, is a pleiotropic protein in the plasminogen subfamily of S1 peptidases. It is a multidomain molecule that includes an N-terminal PAN/APPLE-like domain, four Kringle domains, and a serine proteinase-like domain that has no detectable protease activity (1-5). Human HGF is secreted as an inactive 728 amino acid (aa) single chain propeptide. It is cleaved after the fourth Kringle domain by a serine protease to form bioactive disulfide-linked HGF with a 60 kDa α and 30 kDa β chain. Alternate splicing generates human HGF isoforms that lack the proteinase-like domain and different numbers of the Kringle domains. Human HGF shares 91%-94% aa sequence identity with bovine, canine, feline, mouse, and rat HGF. HGF binds heparan-sulfate proteoglycans and the widely expressed receptor tyrosine kinase, HGF R/c-MET (6, 7). HGF-dependent c-MET activation is implicated in the development of many human cancers (8). HGF regulates epithelial morphogenesis by inducing cell scattering and branching tubulogenesis (9, 10). HGF induces the up-regulation of integrin $\alpha 2 \beta 1$ in epithelial cells by a selective increase in $\alpha 2$ gene transcription (11). This integrin serves as a collagen I receptor, and its blockade disrupts epithelial cell branching tubulogenesis (11, 12). HGF can also alter epithelium morphology by the induction of nectin-1 α ectodomain shedding, an adhesion protein component of adherens junctions (13). In the thyroid, HGF induces the proliferation, motility, and loss of differentiation markers of thyrocytes and inhibits TSH-stimulated iodine uptake (14). HGF promotes the motility of cardiac stem cells in damaged myocardium (15).

References:

1. Karihaloo, A. *et al.* (2005) *Nephron Exp. Nephrol.* **100**:e40.
2. Hammond, D.E. *et al.* (2004) *Curr. Top. Microbiol. Immunol.* **286**:21.
3. Rosario, M. and W. Birchmeier (2004) *Dev. Cell* **7**:3.
4. Lesk, A.M. and W.D. Fordham (1996) *J. Mol. Biol.* **258**:501.
5. Nakamura, T. *et al.* (1989) *Nature* **342**:440.
6. Mizuno, K., *et al.* (1994) *J. Biol. Chem.* **269**:1131.
7. Gheradi, E. *et al.* (2003) *Proc. Natl. Acad. Sci.* **100**:12039.
8. Corso, S. *et al.* (2005) *Trends Mol. Med.* **11**:284.
9. Maeshima, A. *et al.* (2000) *Kid. Int.* **58**:1511.
10. Montesano, R. *et al.* (1991) *Cell* **67**:901.
11. Chiu, S-J. *et al.* (2002) *J. Biomed. Sci.* **9**:261.
12. Saelman, E.U.M. *et al.* (1995) *J. Cell Sci.* **108**:3531.
13. Tanaka, Y. *et al.* (2002) *Biochem. Biophys. Res. Commun.* **299**:472.
14. Mineo, R. *et al.* (1994) *Endocrinology* **145**:4355.
15. Urbanek, K. *et al.* (2005) *Circ. Res.* **97**:663.