

DESCRIPTION

Source Chinese Hamster Ovary cell line, CHO-derived
Ser19-Lys352
Accession # NP_149122

N-terminal Sequence Analysis Ser19

Predicted Molecular Mass 37.4 kDa

SPECIFICATIONS

SDS-PAGE 40 kDa, reducing conditions

Activity Measured by its ability to induce alkaline phosphatase production by MC3T3-E1 mouse preosteoblast cells. The ED₅₀ for this effect is typically 1.5-7.5 ng/mL.

Measured by its ability to induce Topflash reporter activity in HEK293T human embryonic kidney cells. Protein concentrations should be titrated based on cell type and if appropriate, passage number of the cell line. **Optimal concentrations should be determined by each laboratory for each application.**

Endotoxin Level <0.10 EU per 1 µg of the protein by the LAL method.

Purity >90%, by SDS-PAGE with silver staining.

Formulation Lyophilized from a 0.2 µm filtered solution in PBS, NaCl, EDTA and CHAPS. See Certificate of Analysis for details.

PREPARATION AND STORAGE

Reconstitution Reconstitute at 200 µg/mL in sterile PBS.


Shipping The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Stability & Storage Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

- 12 months from date of receipt, -20 to -70 °C as supplied.
- 1 month, 2 to 8 °C under sterile conditions after reconstitution.
- 3 months, -20 to -70 °C under sterile conditions after reconstitution.


DATA

SDS-PAGE



1 µg/lane of Recombinant Human Wnt-3a (High Purity) was resolved with SDS-PAGE under reducing (R) conditions and visualized by silver staining, showing a single band at 40 kDa.

Bioactivity



Recombinant Human Wnt-3a (High Purity, Catalog # 5036-WNP/CF) induces alkaline phosphatase production by the MC3T3-E1 mouse preosteoblast cell line. The ED₅₀ for this effect is typically 1.5-7.5 ng/mL.

BACKGROUND

Wnt-3a is one of 19 vertebrate members of the Wingless-type MMTV integration site (Wnt) family of highly conserved cysteine-rich secreted glycoproteins important for normal developmental processes (1). Wnts bind to the cell surface Frizzled family receptors in conjunction with low-density lipoprotein receptor-related protein family receptors (LRP5 or 6) resulting in the stabilization of intracellular β -catenin levels (2). As intracellular β -catenin levels rise, β -catenin binds to TCF/LEF transcription factors leading to expression of Wnt target genes (3). Endo-IWR 1 (Catalog # 3532, # PSM1324) is a cell-permeant small molecule inhibitor of Axin turnover that suppresses Wnt signal transduction by stabilizing the β -catenin destruction complex (4). Wnt-3a is a 44 kDa secreted hydrophobic glycoprotein containing a conserved pattern of 24 cysteine residues (5). Wnt-3a has two N-linked glycosylation sites (Asn 87, Asn 298), and Ser 209 is modified with palmitoleic acid (6). Glycosylation and acylation are essential for efficient Wnt secretion and biological activity, respectively (6, 7). Human Wnt-3a shares 96% amino acid (aa) identity with mouse mouse, bovine and canine Wnt-3a, and 89%, 86% and 84% aa identity with chicken, Xenopus and zebrafish Wnt-3a, respectively. It also shares 87% aa identity with Wnt3. During embryonic development, Wnt-3a is necessary for proper development of the hippocampus, anterior-posterior patterning, somite development, and tailbud formation (9-12). Wnt-3a also promotes self-renewal of hematopoietic stem cells, neural stem cells, and embryonic stem cells (13, 14).

References:

1. Willert, K. and Nusse, R. (2012) Cold Spring Harb. Perspect. Biol. **4**:a007864.
2. MacDonald, B.T. and X. He (2012) Cold Spring Harb. Perspect. Biol. **4**:a007880.
3. Korinek, V. *et al.* (1997) Science **275**:1784.
4. Chen, B. *et al.* (2009) Nat. Chem. Biol. **5**:100.
5. Smolich, B.D. *et al.* (1993) Mol. Biol. Cell **4**:1267.
6. Takada, R. *et al.* (2006) Dev. Cell **11**:791.
7. Komekado, H. (2007) Genes Cells **12**:521.
8. Dunty Jr. W. C. *et al.* (2008) Development **135**:85.
9. Ikeya, M. and S. Takada (2001) Mech. Dev. **103**:27.
10. Lee, S. M. *et al.* (2000) Development **127**:457.
11. Takada, S. *et al.* (1994) Genes Dev. **8**:174.
12. Willert, K. *et al.* (2003) Nature **423**:6938.
13. Kalani, M.Y. *et al.* (2008) Proc. Natl. Acad. Sci. USA **105**:16970.
14. Ten Berge, D. *et al.* (2011) Nat. Cell Biol. **13**:1070.