

Product Specifications

Amount	◆ 1 mg
Form	◆ Lyophilized
Purity	◆ > 98% by SDS-PAGE
Buffer	◆ 10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS) with 5% trehalose
Preservative	◆ 0.05% (w/v) sodium azide
Reconstitution	◆ Reconstitute with sterile PBS. If 0.5 mL of PBS is used, the antibody concentration will be 2 mg/mL.
Storage	◆ 2° - 8° C
Source	◆ Normal rat serum
Note	◆ Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts of water during disposal.