

Product Specifications

Amount	◆ 1 mg
Form	◆ Lyophilized
Purity	◆ > 90% by SDS-PAGE
Buffer	◆ PBS, pH 7.4 with 5% trehalose
Endotoxin	◆ < 0.10 EU per 1 µg of the antibody by the LAL method.
Preservative	◆ None
Reconstitution	◆ Reconstitute with sterile PBS. If 0.5 mL of PBS is used, the antibody concentration will be 2 mg/mL.
Storage	◆ 2° - 8° C
Source	◆ Normal rat serum